



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Manual Retail Prescription Form

MARCH 2021  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Azalea Self Propel Tall

### Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 24" Self Propulsion with User Manoeuvred Tyre Brakes and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm Width Adjustment)
- Seat Height 500mm
- Depth Adjustable Seat 480mm-550mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +25°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Ducktail Pad Armrests with Autolock
- Pearl Silver Frame Colour
- Chassis Extender with Solid Black Front Castors
- Quick Release Self Propel Wheels, Solid Tyres and Aluminium Handrims
- Pelvic Belt, Adjustable both sides
- Anti Tippers



**Standard Base price: £2,352**

Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
BNF9003	Rea Azalea Self Propel Tall	£2,352
BNF8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
<b>PROPULSION*</b>		
BNF0025	Self Propel 24" wheels with user brakes (NC 1076)	<input checked="" type="checkbox"/>
BNF0030	Self Propel 24" wheels with assistant hub brakes (NC 1074)	£129
<b>SEAT DIMENSIONS*</b>		
BNF0036	Tall (+ 50 mm In seat height, depth and chassis length)	<input checked="" type="checkbox"/>

SEAT WIDTH*			
BNF0040		390 mm/15" (Small chassis)	○
BNF0050		440 mm/17" (Medium chassis)	○
BNF0060		490 mm/19" (Large chassis)	○
SEAT WIDTH EXTENSION*			
BNF0065		+ 50 mm seat width pre-adjustment on armrests and legrests	○
BNF0085		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required	○
SEAT DEPTH*			
BNF0105		480-550 mm / 19" -22" (Azalea Tall - Incl Chassis Extension)	◆
SEAT HEIGHT*			
BNF0125		500 mm / 20" (Azalea Tall)	◆
SEAT ANGLE*			
BNF0000		Standard Seat Angle +1° to +25°	◆
BNF0127		Negative Seat Angle -1 To +23° (NC 0360)	○
SEAT TYPE*			
BNF0130		Depth Adjustable	◆
SEAT CUSHION*			
BNF0165		Cushion Floshape seat upholstery (NC 1111)	£52
BNF0170		Cushion Vicair Multifunctional (Black) (NC 0201, 0202)	£559
BNF0180		Without Seat Cushion, Seat Plate only (NC 0201, 0202)	○
SEAT UPHOLSTERY*			
BNF0201		Multi-stretch Black PU coated fabric	◆
BNF0202		Black Spacer TR36	○
BNF1111		No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion	○
BACKREST WIDTH EXTENSION*			
BNF0205		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required (NC 0221, 0223, 0225, 0228)	○
BNF0207		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required	○
BACKREST*			
BNF0210		Flex 3 Backrest, Fixed Plate (Width Adjustable) (NC 0233, 0234, 0240, 0343, 0771, 0905)	◆
BNF0220		Flex 3 Backrest Tension (Width Adjustable) (NC 0233, 0234, 0343, 0771, 0925)	£94
BNF0221		Matrx Elite (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£437
BNF0223		Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£485
BNF0225		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£437
BNF0228		Matrx Elite TR (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0820, 0880, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£361
BNF0230		Backrest Tubes Without Backrest (NC 0233, 0234, 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0345, 0740, 0771, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	○
MATRX BACKREST HEIGHT			
BNF0233		Matrx Backrest Height 410 mm / 16"	○
BNF0234		Matrx Backrest Height 510 mm / 20"	○

<b>BACKREST CUSHION*</b>				
BNF0240		Cover for tension adjustable upholstery	(NC 2222)	○
BNF0250		Cushion, Laguna (Lateral Support)	(NC 0905, 2222)	○
BNF0260		Cushion, Mistral 2 (Waist Support)	(NC 0322, 0905, 2222)	£133
BNF0270		Cushion, Passad 2 (Shoulder Support)	(NC 0905, 2222)	◆
BNF0280		Cushion, Shoulder High (Multi-stretch Black PU coated fabric)	(NC 0322, 0905, 2222)	£138
BNF0290		Cushion, Vicair Multifunctional (Black)	(NC 0321, 0322, 0905)	£537
BNF0300		Without Backrest Cover/ Cushion	(NC 0321, 0322, 0905)	○
<b>BACKREST UPHOLSTERY*</b>				
BNF0321		Multi-stretch Black PU coated fabric		◆
BNF0322		Black Spacer TR36		○
BNF2222		No upholstery, select with Matrx, Multifunctional and Without backrest options		○
<b>PUSH HANDLES*</b>				
BNF0335		Push Handles connected,no head/neckrest attachment		◆
BNF0343		Push Handles connected,with Rea head/neckrest attachment	(NC 0771)	£66
BNF0345		Push Bar, angle adjustable		£94
<b>SEAT TILT AND BACKREST REGULATION*</b>				
BNF0350		Assistant operated gas piston seat tilt and backrest recline	(NC 5555)	◆
BNF0360		Electrical operated seat tilt and backrest recline	(NC 1263)	£673
<b>TILT LOCK*</b>				
BNF1263		Tilt Lock (locking of pre-set angle)		◆
BNF5555		No Tilt Lock, Select with 0360 Electric operated seat and backrest		○
<b>CENTRAL LEGREST</b>				
BNF0370		Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444	(NC 0189, 0376, 0378, 0380, 0390, 0410, 0430, 0440, 0460, 0475, 0476, 0477, 0720, 0930, 0940, 0941, 0942, 1290)	£279
<b>LEGREST HANGER ATTACHMENT*</b>				
BNF0375		Standard position for legrest mounting		◆
SWO0189		Standard position with 0-30° Legrest Abduction function		£92
BNF0376		High position 50mm higher for leg length compensation		○
BNF0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest	(NC 0390, 0410, 0440, 0460)	○
<b>LEGREST HANGER LH*</b>				
BNF0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH	(NC 0475, 0480, 0490, 0942, 0945)	◆
BNF0390		Legrest Hanger fixed 80 deg, swing in/swing out LH		○
BNF0410		Legrest Amputee LH	(NC 0475, 0476, 0477, 0480, 0490, 0510, 0520, 0930, 0941, 0942, 0945)	£126
BNF0420		Without Legrest, LH	(NC 0475, 0476, 0477, 0480, 0490, 0510, 0520, 0930, 0941, 0942, 0945)	○
<b>LEGREST HANGER RH*</b>				
BNF0430		Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support RH	(NC 0475, 0540, 0550, 0942, 0945)	◆
BNF0440		Legrest Hanger fixed 80 deg, swing in/swing out RH		○
BNF0460		Legrest Amputee RH	(NC 0475, 0476, 0477, 0540, 0550, 0570, 0580, 0940, 0941, 0942, 0945)	£126
BNF0470		Without Legrest, RH	(NC 0475, 0476, 0477, 0540, 0550, 0570, 0580, 0940, 0941, 0942, 0945)	○

FOOTREST LH*			
BNF0475		Footrest one piece, depth and angle adjustable short tube	(NC 0477, 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940) £109
BNF0476		Footrest one piece, depth and angle adjustable long tube	(NC 0477, 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940) £109
BNF0477		Footbox, angle adjustable, padded with short tube	(NC 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940, 0941, 0942, 0945) £501
BNF0480		Footrest lower part, fixed with short tube and heel strap, LH	(NC 0930, 0941) ○
BNF0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH	○
BNF0510		Footrest lower part, fixed with long tube and heel strap, LH.	(NC 0930, 0941) ○
BNF0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH	○
BNF3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○
FOOTREST RH*			
BNF0540		Footrest lower part, fixed with short tube and heel strap, RH	(NC 0940, 0941) ○
BNF0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH	○
BNF0570		Footrest lower part, fixed with long tube and heel strap, RH	(NC 0940, 0941) ○
BNF0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH	○
BNF4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BNF0600		Armrest Ducktail Pad LH, height adjustable with autolock	(NC 1232) ◆
BNF0607		Armrest Reinforced Wide Pad, LH, height adjustable with autolock	(NC 1230, 1231, 1232, 1250, 1260) £82
BNF0610		Armrest Padded LH, height adjustable with autolock	(NC 1230, 1231, 1232, 1250, 1260) £59
BNF0615		Armrest Hemiplegic LH, height adjustable with autolock	(NC 0695, 1230, 1231, 1232, 1250, 1260) £165
BNF0620		Without Armrest LH Select when choosing BNA1232 Half Tray LH, Narrow Armrest is included	(NC 0695, 1230, 1231, 1250, 1260) ○
ARMREST RH*			
BNF0630		Armrest Ducktail Pad RH, height adjustable with autolock	(NC 1233) ◆
BNF0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock	(NC 1230, 1231, 1233, 1250, 1260) £82
BNF0640		Armrest Padded RH, height adjustable with autolock.	(NC 1230, 1231, 1233, 1250, 1260) £59
BNF0645		Armrest Hemiplegic, RH, height adjustable with autolock	(NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260) £165
BNF0650		Without Armrest RH Select when choosing BNF1233 Half Tray RH, Narrow Armrest is included	(NC 0695, 1230, 1231, 1240, 1245, 1250, 1260) ○
ARMREST OPTIONS. Side support cover 0680, must be selected if 0690 Foam for side support is required			
BNF0680		Side Support Cover, removable with tools (wipeable cover for the sideguard)	(NC 0695) £47
BNF0690		Foam for side support cover (Internal cushioned pad, which is only possible in combination with 0680)	(NC 0695) ○
BNF0695		Pad for sideguard, removable with hook and loop fastening	£47

SEAT OPTIONS			
BNF0710		Incontinence Cover	£37
BNF0720		Pommel Abduction Cushion (NC 0722)	£109
BNF0722		2 cm Increased seat height plate	£59
BNF0725		Tilt Scale for seat tilt	£41
BACKREST OPTIONS			
BNF0730		Headrest, Sideway Adjustable (NC 0750, 0770)	◆
BNF0750		Headrest with cheek support, sideway adjustable (NC 0740, 0770, 0780)	£99
BNF0740		Cover for headrest	£28
BNF0770		Neckrest (NC 0771)	○
BNF0771		Headrest Matrix, 250 mm Soft Pad (NC 0780, 0790)	£171
BNF0780		Cover For Neckrest	£20
BNF0782		Alternative headrest/neckrest mounting bracket	£21
BNF0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	○
BNF0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included	○
BNF0905		Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery	£56
BNF0910		Lumbar Support Cushion	£28
BNF0915		Tilt Scale for backrest recline	£41
BNF0920		Wedges (One Pair)	£9
BNF0925		Privacy Attachment (Material connecting backrest and seat)	£31
TRUNK SUPPORT LH			
BNF0800		Trunk Support Multi Adjustable with standard arm LH (NC 0900)	○
BNF0820		Trunk Support Multi Adjustable with extended height arm LH (NC 0900)	○
BNF0880		Trunk Support Swingaway LH (NC 0900)	£97
BNF0895		Trunk Support Multi Adjustable Swingaway LH (NC 0900)	£104
TRUNK SUPPORT RH			
BNF0810		Trunk Support Multi Adjustable with standard arm RH (NC 0900)	○
BNF0830		Trunk Support Multi Adjustable with extended height arm RH (NC 0900)	○
BNF0890		Trunk Support Swingaway RH (NC 0900)	£97
BNF0896		Trunk Support Multi Adjustable Swingaway RH (NC 0900)	£104
LEGREST OPTIONS			
BNF0930		Footplate Extender LH (NC 0941)	£31
BNF0940		Footplate Extender RH (NC 0941)	£31
BNF0941		Cover Footplate, Black Plush Fabric	£49
BNF0942		Calf Strap Padded	£67
BNF0945		Calf Strap (for One piece footrest)	○
CHASSIS COLOUR*			
BNF0944		Silver Pearl	◆
BNF0946		Azurite (Dark Blue)	○
BNF0947		Sand (Beige)	○
BNF0948		Electric Green	○
BNF0949		Happy Red	○

<b>CHASSIS SAFETY*</b>			
BNF0950	<input type="checkbox"/>	Chassis extender, castors extend 50mm forward for added stability	◆
<b>CASTOR ATTACHMENT*</b>			
BNF0960	<input type="checkbox"/>	Fixed castor attachment	◆
<b>FRONT CASTORS*</b>			
BNF1015	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	£19
BNF1025	<input type="checkbox"/>	200x45 mm (8" X 1 3/4") Black Solid tyre	◆
<b>REAR WHEEL SIZE*</b>			
BNF1062	<input type="checkbox"/>	24" / 600 mm Self Propel rear wheel with QR attachment	◆
<b>REAR WHEEL TYPE*</b>			
BNF1074	<input type="checkbox"/>	Standard self propel with colourless hub, clear anodized rim and Stainless steel spokes	◆
BNF1076	<input type="checkbox"/>	Hub Brake 22" or 24" Self Propel with colourless hub, clear anodized rim and Stainless steel spokes	○
<b>REAR WHEEL TYRE*</b>			
BNF1082	<input type="checkbox"/>	Pneumatic Profile Black 1 3/8" tyres	○
BNF1084	<input type="checkbox"/>	Pneumatic Puncture proof / flatfree 1 3/8" Black tyres	£57
BNF1092	<input type="checkbox"/>	Solid Black tyres	◆
<b>HANDRIM*</b>			
BNF1100	<input type="checkbox"/>	Aluminium handrim	◆
BNF1120	<input type="checkbox"/>	Plastic coated handrim	£45
<b>SPOKEGUARD</b>			
BNF1160	<input type="checkbox"/>	Frost spokeguard	£153
<b>BRAKE*</b>			
BNF1170	<input type="checkbox"/>	User manoeuvred tyre brakes	◆
BNF1180	<input type="checkbox"/>	One Arm Brake LH	£116
BNF1185	<input type="checkbox"/>	One Arm Brake RH	£116
<b>PELVIC BELT*</b>			
BNF1200	<input type="checkbox"/>	Pelvic belt	◆
<b>TRANSPORTATION BRACKETS*</b>			
BNF1225	<input type="checkbox"/>	Transportation Brackets, mounted to the front of the wheelchair	◆
<b>SAFETY OPTIONS</b>			
BNF1220	<input type="checkbox"/>	Reflector Kit (wheels and chassis)	£20

**OTHER OPTIONS**

BNF1230		Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	£127
BNF1231		Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	£127
BNF1232		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	£95
BNF1233		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	£95
BNF1240		Swing away table without elbow pads (NC 1235, 1250, 1260)	£140
BNF1245		Swing away table with elbow pads (NC 1235, 1260)	£140
BNF1235		Table tray padded cushion	£47
BNF1237		Half table tray padded cushion	£43
BNF1250		Locking device for table tray	£13
BNF1260		Mounting kit for table (if ordered separately)	£7
BNF1270		Pump	£20
BNF1280		Cane Holder	£46
BNF1290		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	£50
BNF1300		Harness attachment , must be selected in combination with Head/Neckrest	£79
BNF1310		Attachment for positioning belt	£29

**TOTAL PRICE**

**Non Conforming Options:**

BNF0201	Multi-stretch Black PU coated fabric
BNF0202	Black Spacer TR36
BNF1111	No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion
BNF0221	Matrx Elite (includes hardware mounted to frame, backrest supplied in box)
BNF0223	Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box)
BNF0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BNF0228	Matrx Elite TR (includes hardware mounted to frame, backrest supplied in box)
BNF0233	Matrx Backrest Height 410 mm / 16"
BNF0234	Matrx Backrest Height 510 mm / 20"
BNF0240	Cover for tension adjustable upholstery
BNF0250	Cushion, Laguna (Lateral Support)
BNF0260	Cushion, Mistral 2 (Waist Support)
BNF0270	Cushion, Passad 2 (Shoulder Support)
BNF0280	Cushion, Shoulder High (Multi-stretch Black PU coated fabric)
BNF0290	Cushion, Vicair Multifunctional (Black)
BNF0321	Multi-stretch Black PU coated fabric
BNF0322	Black Spacer TR36
BNF2222	No upholstery, select with Matrx, Multifunctional and Without backrest options
BNF0335	Push Handles connected,no head/neckrest attachment
BNF0343	Push Handles connected,with Rea head/neckrest attachment
BNF0345	Push Bar, angle adjustable
BNF0360	Electrical operated seat tilt and backrest recline
BNF1263	Tilt Lock (locking of pre-set angle)
BNF5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
SWO0189	Standard position with 0-30° Legrest Abduction function
BNF0376	High position 50mm higher for leg length compensation
BNF0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BNF0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH
BNF0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNF0410	Legrest Amputee LH
BNF0430	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support RH
BNF0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNF0460	Legrest Amputee RH
BNF0475	Footrest one piece, depth and angle adjustable short tube
BNF0476	Footrest one piece, depth and angle adjustable long tube
BNF0477	Footbox, angle adjustable, padded with short tube
BNF0480	Footrest lower part, fixed with short tube and heel strap, LH
BNF0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BNF0510	Footrest lower part, fixed with long tube and heel strap, LH.
BNF0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH
BNF0540	Footrest lower part, fixed with short tube and heel strap, RH
BNF0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BNF0570	Footrest lower part, fixed with long tube and heel strap, RH



BNF0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BNF0695	Pad for sideguard, removable with hook and loop fastening
BNF0720	Pommel Abduction Cushion
BNF0722	2 cm Increased seat height plate
BNF0750	Headrest with cheek support, sideway adjustable
BNF0740	Cover for headrest
BNF0770	Neckrest
BNF0771	Headrest Matrix, 250 mm Soft Pad
BNF0780	Cover For Neckrest
BNF0790	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BNF0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNF0905	Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery
BNF0910	Lumbar Support Cushion
BNF0920	Wedges (One Pair)
BNF0925	Privacy Attachment (Material connecting backrest and seat)
BNF0800	Trunk Support Multi Adjustable with standard arm LH
BNF0820	Trunk Support Multi Adjustable with extended height arm LH
BNF0880	Trunk Support Swingaway LH
BNF0895	Trunk Support Multi Adjustable Swingaway LH
BNF0810	Trunk Support Multi Adjustable with standard arm RH
BNF0830	Trunk Support Multi Adjustable with extended height arm RH
BNF0890	Trunk Support Swingaway RH
BNF0896	Trunk Support Multi Adjustable Swingaway RH
BNF0930	Footplate Extender LH
BNF0940	Footplate Extender RH
BNF0941	Cover Footplate, Black Plush Fabric
BNF0942	Calf Strap Padded
BNF0945	Calf Strap (for One piece footrest)
BNF1074	Standard self propel with colourless hub, clear anodized rim and Stainless steel spokes
BNF1076	Hub Brake 22" or 24" Self Propel with colourless hub, clear anodized rim and Stainless steel spokes
BNF1230	Fixed Table Tray without elbow pads includes mounting kit
BNF1231	Fixed Table Tray with elbow Pads includes mounting kit
BNF1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNF1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNF1240	Swing away table without elbow pads
BNF1245	Swing away table with elbow pads
BNF1235	Table tray padded cushion
BNF1237	Half table tray padded cushion
BNF1250	Locking device for table tray
BNF1260	Mounting kit for table (if ordered separately)
BNF1290	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):