



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

MARCH 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Azalea Transit

Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Choice of 2 Seat Heights 400mm or 450mm
- Depth Adjustable Seat 430mm-500mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +25°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Wide Pad Armrests with Autolock
- Pearl Silver Frame Colour, Solid Black Front Castors
- User Manoeuvred Tyre Brakes
- Pelvic Belt, Adjustable both sides
- Anti Tippers
- Front Transportation Brackets



Standard Base price: £2,190

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
BNE9002	Rea Azalea Transit	£2,190
BNE8000	SPECIAL REQUEST	POA
PROPULSION*		
BNE0020	Transit 16" wheels with hub brakes	<input checked="" type="checkbox"/>
SEAT WIDTH*		
BNE0040	390 mm/15" (Small chassis)	<input type="checkbox"/>
BNE0050	440 mm/17" (Medium chassis)	<input type="checkbox"/>
BNE0060	490 mm/19" (Large chassis) (NC 0206)	<input type="checkbox"/>

SEAT WIDTH EXTENSION*			
BNE0065		+ 50 mm pre-adjustment on armrests and legrests	○
BNE0075		+ 100 mm pre-adjustment on armrests and legrests (NC 0110)	○
BNE0085		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required (NC 0206)	○
SEAT DEPTH*			
BNE0100		430-500 mm / 17" - 20"	◆
SEAT HEIGHT*			
BNE0110		400 mm / 16" Seat Height (NC 0370, 1015, 1025)	○
BNE0120		450 mm / 18" Seat Height	○
SEAT ANGLE*			
BNE0127		Negative Seat Angle -1 To +23° (NC 0360)	○
BNE0000		Standard Seat Angle +1° to +25°	◆
SEAT TYPE*			
BNE0130		Depth Adjustable	◆
SEAT CUSHION*			
BNE0165		Cushion Floshape seat upholstery (NC 1111)	£52
SWO0190		Cushion Floshape seat upholstery Shorter depth 450mm deep. Only use with seat depth in shortest position (NC 0202, 1111)	£52
BNE0170		Cushion Vicair Multifunctional (Black) (NC 0201, 0202)	£559
BNE0180		Without Seat Cushion, Seat Plate only (NC 0201, 0202)	○
SEAT UPHOLSTERY*			
BNE0201		Multi-stretch Black PU coated fabric	◆
BNE0202		Black Spacer TR36	○
BNE1111		No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion	○
BACKREST WIDTH EXTENSION*			
BNE0205		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required (NC 0221, 0223, 0225, 0228)	○
BNE0206		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required (NC 0220, 0221, 0223, 0225, 0228)	○
BNE0207		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required	○
BACKREST*			
BNE0210		Flex 3 Backrest, Fixed Plate (Width Adjustable) (NC 0233, 0234, 0240, 0343, 0771, 0905)	◆
BNE0220		Flex 3 Backrest Tension (Width Adjustable) (NC 0233, 0234, 0343, 0771, 0925)	£94
BNE0221		Matrx Elite (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£437
BNE0223		Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£485
BNE0225		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£437
BNE0228		Matrx Elite TR (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0820, 0880, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£361
BNE0230		Backrest Tubes Without Backrest (NC 0233, 0234, 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0345, 0740, 0771, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	○
MATRX BACKREST HEIGHT			
BNE0233		Matrx Backrest Height 410 mm / 16"	○
BNE0234		Matrx Backrest Height 510 mm / 20"	○

BACKREST CUSHION*				
BNE0240		Cover for tension adjustable upholstery	(NC 2222)	○
BNE0250		Cushion, Laguna (Lateral Support)	(NC 0905, 2222)	○
BNE0260		Cushion, Mistral 2 (Waist Support)	(NC 0322, 0905, 2222)	£133
BNE0270		Cushion, Passad 2 (Shoulder Support)	(NC 0905, 2222)	◆
BNE0280		Cushion, Shoulder High. Higher Cushion to suit taller clients	(NC 0322, 0905, 2222)	£138
BNE0290		Cushion, Vicair Multifunctional (Black)	(NC 0321, 0322, 0905)	£537
BNE0300		Without Backrest Cover/ Cushion	(NC 0321, 0322, 0905)	○
BACKREST UPHOLSTERY*				
BNE0321		Multi-stretch Black PU coated fabric		◆
BNE0322		Black Spacer TR36		○
BNE2222		No upholstery, select with Matrix, Multifunctional and Without backrest options		○
PUSH HANDLES*				
BNE0335		Push Handles connected,no head/neckrest attachment		◆
BNE0343		Push Handles connected,with Rea head/neckrest attachment	(NC 0771)	£66
BNE0345		Push Bar, Angle Adjustable		£94
SEAT TILT AND BACKREST REGULATION*				
BNE0350		Assistant operated gas piston seat tilt and backrest recline	(NC 5555)	◆
BNE0360		Electrical operated seat tilt and backrest recline	(NC 1263)	£673
TILT LOCK*				
BNE1263		Tilt Lock (locking of pre-set angle)		◆
BNE5555		No Tilt Lock, Select with 0360 Electric operated seat and backrest		○
CENTRAL LEGREST				
BNE0370		Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444	(NC 0189, 0376, 0378, 0380, 0390, 0400, 0410, 0430, 0440, 0450, 0460, 0475, 0476, 0477, 0720, 0930, 0940, 0941, 0942, 1290)	£279
LEGREST HANGER ATTACHMENT*				
BNE0375		Standard position for legrest mounting		◆
SWO0189		Standard position with 0-30° Legrest Abduction function		£92
BNE0376		High position 50mm higher for leg length compensation		○
BNE0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest only	(NC 0390, 0400, 0410, 0440, 0450, 0460)	○
LEGREST HANGER LH*				
BNE0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH	(NC 0475, 0480, 0490, 0942, 0945)	◆
BNE0390		Legrest Hanger fixed 80 deg, swing in/swing out LH		○
BNE0400		Legrest Hanger fixed 90 deg, swing in/swing out LH	(NC 0480, 0520, 0930)	○
BNE0410		Legrest Amputee LH	(NC 0475, 0476, 0477, 0480, 0490, 0520, 0930, 0941, 0942, 0945)	£126
BNE0420		Without Legrest, LH	(NC 0475, 0476, 0477, 0480, 0490, 0520, 0930, 0941, 0942, 0945)	○
LEGREST HANGER RH*				
BNE0430		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support RH	(NC 0475, 0540, 0550, 0942, 0945)	◆
BNE0440		Legrest Hanger fixed 80 deg, swing in/swing out RH		○
BNE0450		Legrest Hanger fixed 90 deg, swing in/swing out RH	(NC 0540, 0580, 0940)	○
BNE0460		Legrest Amputee RH	(NC 0475, 0476, 0477, 0540, 0550, 0580, 0940, 0941, 0942, 0945)	£126
BNE0470		Without Legrest, RH	(NC 0475, 0476, 0477, 0540, 0550, 0580, 0940, 0941, 0942, 0945)	○

FOOTREST LH*			
BNE0475		Footrest one piece, depth and angle adjustable short tube (NC 0477, 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940)	£109
BNE0476		Footrest one piece, depth and angle adjustable long tube (NC 0477, 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940)	£109
BNE0477		Footbox, angle adjustable, padded with short tube (NC 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940, 0941, 0942, 0945)	£501
BNE0480		Footrest lower part, fixed with short tube and heel strap, LH (NC 0930, 0941)	○
BNE0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH	○
BNE0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH	○
BNE3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○
FOOTREST RH*			
BNE0540		Footrest lower part, fixed with short tube and heel strap, RH (NC 0940, 0941)	○
BNE0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH	○
BNE0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH	○
BNE4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BNE0600		Armrest Ducktail Pad LH, height adjustable with autolock (NC 1232)	○
BNE0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only (NC 1232)	◆
BNE0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	£82
BNE0610		Armrest Padded LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	£59
BNE0615		Armrest Hemiplegic LH, height adjustable with autolock (NC 0695, 1230, 1231, 1232, 1250, 1260)	£165
BNE0620		Without Armrest LH Select when choosing BNE1232 Half Tray LH, Narrow Armrest is included (NC 0695, 1230, 1231, 1250, 1260)	○
ARMREST RH*			
BNE0630		Armrest Ducktail Pad RH, height adjustable with autolock (NC 1233)	○
BNE0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only (NC 1233)	◆
BNE0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock (NC 1230, 1231, 1233, 1250, 1260)	£82
BNE0640		Armrest Padded RH, height adjustable with autolock. (NC 1230, 1231, 1233, 1250, 1260)	£59
BNE0645		Armrest Hemiplegic, RH, height adjustable with autolock (NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260)	£165
BNE0650		Without Armrest RH Select when choosing BNE1233 Half Tray RH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST OPTIONS. Side support cover 0680, must be selected if 0690 Foam for side support is required			
BNE0680		Side Support Cover, removable with tools (Wipeable cover for the sideguard) (NC 0695)	£47
BNE0690		Foam for side support cover (Internal cushioned pad, which is only possible in combination with 0680) (NC 0695)	○
BNE0695		Pad for sideguard, removable with hook and loop fastening	£47

SEAT OPTIONS			
BNE0710		Incontinence Cover	£37
BNE0720		Pommel Abduction Cushion (NC 0722)	£109
BNE0722		2 cm Increased seat height plate	£59
BNE0725		Tilt Scale for seat tilt	£41
BACKREST OPTIONS			
BNE0730		Headrest, Sideway Adjustable (NC 0750, 0770)	◆
BNE0750		Headrest with cheek support, sideway adjustable (NC 0740, 0770, 0780)	£99
BNE0740		Cover for headrest	£28
BNE0770		Neckrest (NC 0771)	○
BNE0771		Headrest Matrix, 250 mm Soft Pad (NC 0780, 0790)	£171
BNE0780		Cover For Neckrest	£20
BNE0782		Alternative headrest/neckrest mounting bracket	£21
BNE0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	○
BNE0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included	○
BNE0905		Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery	£56
BNE0910		Lumbar Support Cushion	£28
BNE0915		Tilt Scale for backrest recline	£41
BNE0920		Wedges (One Pair)	£9
BNE0925		Privacy Attachment (Material connecting backrest and seat)	£31
TRUNK SUPPORT LH			
BNE0800		Trunk Support Multi Adjustable with standard arm LH (NC 0900)	○
BNE0820		Trunk Support Multi Adjustable with extended height arm LH (NC 0900)	○
BNE0880		Trunk Support Swingaway LH (NC 0900)	£97
BNE0895		Trunk Support Multi Adjustable Swingaway LH (NC 0900)	£104
TRUNK SUPPORT RH			
BNE0810		Trunk Support Multi Adjustable with standard arm RH (NC 0900)	○
BNE0830		Trunk Support Multi Adjustable with extended height arm RH (NC 0900)	○
BNE0890		Trunk Support Swingaway RH (NC 0900)	£97
BNE0896		Trunk Support Multi Adjustable Swingaway RH (NC 0900)	£104
LEGREST OPTIONS			
BNE0930		Footplate Extender LH (NC 0941)	£31
BNE0940		Footplate Extender RH (NC 0941)	£31
BNE0941		Cover Footplate, Black Plush Fabric	£49
BNE0942		Calf Strap Padded	£67
BNE0945		Calf Strap (for One piece footrest)	○
CHASSIS COLOUR*			
BNE0944		Silver Pearl	◆
BNE0946		Azurite (Dark Blue)	○
BNE0947		Sand (Beige)	○
BNE0948		Electric Green	○
BNE0949		Happy Red	○

CHASSIS SAFETY			
BNE0950	<input type="checkbox"/>	Chassis extender, castors extend 50mm forward for added stability	○
CASTOR ATTACHMENT*			
BNE0960	<input type="checkbox"/>	Fixed castor attachment	◆
FRONT CASTORS*			
BNE0992	<input type="checkbox"/>	150x30 mm (6" x 1 1/4") Medium Hard Black tyre	○
BNE1015	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	○
BNE1025	<input type="checkbox"/>	200x45 mm (8" X 1 3/4") Black Solid tyre	◆
REAR WHEEL SIZE*			
BNE1047	<input type="checkbox"/>	16" / 405 mm Transit rear wheel with fixed attachment	◆
REAR WHEEL TYPE*			
BNE1076	<input type="checkbox"/>	Hub Brake 16" Transit with Black Plastic Hub and Rim	◆
REAR WHEEL TYRE*			
BNE1092	<input type="checkbox"/>	Solid Black tyres	◆
BRAKE*			
BNE1170	<input type="checkbox"/>	User manoeuvred tyre brakes	◆
SAFETY OPTIONS			
BNE1220	<input type="checkbox"/>	Reflector Kit (wheels and chassis)	£20
PELVIC BELT*			
BNE1200	<input type="checkbox"/>	Pelvic belt	◆
TRANSPORTATION BRACKETS*			
BNE1225	<input type="checkbox"/>	Transportation Brackets, mounted to the front of the wheelchair	◆
OTHER OPTIONS			
BNE1230	<input type="checkbox"/>	Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	£127
BNE1231	<input type="checkbox"/>	Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	£127
BNE1232	<input type="checkbox"/>	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	£95
BNE1233	<input type="checkbox"/>	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	£95
BNE1240	<input type="checkbox"/>	Swing away table without elbow pads (NC 1235, 1250, 1260)	£140
BNE1245	<input type="checkbox"/>	Swing away table with elbow pads (NC 1235, 1260)	£140
BNE1235	<input type="checkbox"/>	Table tray padded cushion	£47
BNE1237	<input type="checkbox"/>	Half table tray padded cushion	£43
BNE1250	<input type="checkbox"/>	Locking device for table tray	£13
BNE1260	<input type="checkbox"/>	Mounting kit for table (if ordered separately)	£7
BNE1270	<input type="checkbox"/>	Pump	£20
BNE1280	<input type="checkbox"/>	Cane Holder	£46
BNE1290	<input type="checkbox"/>	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	£50
BNE1300	<input type="checkbox"/>	Harness attachment , must be selected in combination with Head/Neckrest	£79
BNE1305	<input type="checkbox"/>	Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles	£171
BNE1310	<input type="checkbox"/>	Attachment for positioning belt	£29

TOTAL PRICE	
--------------------	--

Non Conforming Options:

BNE0110	400 mm / 16" Seat Height
BNE0201	Multi-stretch Black PU coated fabric
BNE0202	Black Spacer TR36
BNE1111	No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion
BNE0206	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BNE0220	Flex 3 Backrest Tension (Width Adjustable)
BNE0221	Matrx Elite (includes hardware mounted to frame, backrest supplied in box)
BNE0223	Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box)
BNE0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BNE0228	Matrx Elite TR (includes hardware mounted to frame, backrest supplied in box)
BNE0233	Matrx Backrest Height 410 mm / 16"
BNE0234	Matrx Backrest Height 510 mm / 20"
BNE0240	Cover for tension adjustable upholstery
BNE0250	Cushion, Laguna (Lateral Support)
BNE0260	Cushion, Mistral 2 (Waist Support)
BNE0270	Cushion, Passad 2 (Shoulder Support)
BNE0280	Cushion, Shoulder High. Higher Cushion to suit taller clients
BNE0290	Cushion, Vicair Multifunctional (Black)
BNE0321	Multi-stretch Black PU coated fabric
BNE0322	Black Spacer TR36
BNE2222	No upholstery, select with Matrx, Multifunctional and Without backrest options
BNE0335	Push Handles connected,no head/neckrest attachment
BNE0343	Push Handles connected,with Rea head/neckrest attachment
BNE0345	Push Bar, Angle Adjustable
BNE0360	Electrical operated seat tilt and backrest recline
BNE1263	Tilt Lock (locking of pre-set angle)
BNE5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
BNE0370	Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444
SWO0189	Standard position with 0-30° Legrest Abduction function
BNE0376	High position 50mm higher for leg length compensation
BNE0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest only
BNE0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH
BNE0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNE0400	Legrest Hanger fixed 90 deg, swing in/swing out LH
BNE0410	Legrest Amputee LH
BNE0430	Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support RH
BNE0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNE0450	Legrest Hanger fixed 90 deg, swing in/swing out RH
BNE0460	Legrest Amputee RH
BNE0475	Footrest one piece, depth and angle adjustable short tube
BNE0476	Footrest one piece, depth and angle adjustable long tube
BNE0477	Footbox, angle adjustable, padded with short tube
BNE0480	Footrest lower part, fixed with short tube and heel strap, LH

BNE0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BNE0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH
BNE0540	Footrest lower part, fixed with short tube and heel strap, RH
BNE0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BNE0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BNE0695	Pad for sideguard, removable with hook and loop fastening
BNE0720	Pommel Abduction Cushion
BNE0722	2 cm Increased seat height plate
BNE0750	Headrest with cheek support, sideways adjustable
BNE0740	Cover for headrest
BNE0770	Neckrest
BNE0771	Headrest Matrix, 250 mm Soft Pad
BNE0780	Cover For Neckrest
BNE0790	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BNE0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNE0905	Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery
BNE0910	Lumbar Support Cushion
BNE0920	Wedges (One Pair)
BNE0925	Privacy Attachment (Material connecting backrest and seat)
BNE0800	Trunk Support Multi Adjustable with standard arm LH
BNE0820	Trunk Support Multi Adjustable with extended height arm LH
BNE0880	Trunk Support Swingaway LH
BNE0895	Trunk Support Multi Adjustable Swingaway LH
BNE0810	Trunk Support Multi Adjustable with standard arm RH
BNE0830	Trunk Support Multi Adjustable with extended height arm RH
BNE0890	Trunk Support Swingaway RH
BNE0896	Trunk Support Multi Adjustable Swingaway RH
BNE0930	Footplate Extender LH
BNE0940	Footplate Extender RH
BNE0941	Cover Footplate, Black Plush Fabric
BNE0942	Calf Strap Padded
BNE0945	Calf Strap (for One piece footrest)
BNE1015	200x32 mm (8" x 1 1/4") Black Solid Tyre
BNE1025	200x45 mm (8" X 1 3/4") Black Solid tyre
BNE1230	Fixed Table Tray without elbow pads includes mounting kit
BNE1231	Fixed Table Tray with elbow Pads includes mounting kit
BNE1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNE1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNE1240	Swing away table without elbow pads
BNE1245	Swing away table with elbow pads
BNE1235	Table tray padded cushion
BNE1237	Half table tray padded cushion
BNE1250	Locking device for table tray
BNE1260	Mounting kit for table (if ordered separately)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):