



Manual Retail Prescription Form

Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MARCH 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Azalea Transit Tall

Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Seat Height 500mm
- Depth Adjustable Seat 430mm-500mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +25°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Wide Pad Armrests with Autolock
- Headrest, Sideway Adjustable, Pearl Silver Frame Colour
- Chassis Extender with Solid Black Front Castors
- User Manoeuvred Tyre Brakes, Pelvic Belt, Adjustable both sides
- Anti Tippers, Front Transportation Brackets



Standard Base price: £2,352

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
BNG9004	Rea Azalea Transit Tall	£2,352
BNG8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
PROPULSION*		
BNG0020	Transit 16" wheels with hub brakes	◆
SEAT DIMENSIONS*		
BNG0036	Tall	◆

SEAT WIDTH*			
BNG0040		390 mm/15" (Small chassis)	○
BNG0050		440 mm/17" (Medium chassis)	○
BNG0060		490 mm/19" (Large chassis) (NC 0206)	○
SEAT WIDTH EXTENSION*			
BNG0065		+ 50 mm pre-adjustment on armrests and legrests	○
BNG0075		+ 100 mm pre-adjustment on armrests and legrests	○
BNG0085		Without seat width extension (NC 0206)	○
SEAT DEPTH*			
BNG0105		480-550 mm / 19" -22" (Azalea Tall - Incl Chassis Extension)	◆
SEAT HEIGHT*			
BNG0125		500 mm / 20" (Azalea Tall)	◆
SEAT ANGLE*			
BNG0127		Negative Seat Angle -1 To +23° (NC 0360)	○
BNG0000		Standard Seat Angle +1° to +25°	◆
SEAT TYPE*			
BNG0130		Depth Adjustable	◆
SEAT CUSHION*			
BNG0165		Cushion Flosshape seat upholstery (NC 1111)	£52
BNG0170		Cushion Vicair Multifunctional (Black) (NC 0201, 0202)	£559
BNG0180		Without Seat Cushion, Seat Plate only (NC 0201, 0202)	○
SEAT UPHOLSTERY*			
BNG0201		Multi-stretch Black PU coated fabric	◆
BNG0202		Black Spacer TR36	○
BNG1111		No upholstery	○
BACKREST WIDTH EXTENSION*			
BNG0205		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required (NC 0221, 0223, 0225, 0228)	○
BNG0206		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required (NC 0220, 0221, 0223, 0225, 0228)	○
BNG0207		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required	○
BACKREST*			
BNG0210		Flex 3 Backrest, Fixed Plate (Width Adjustable) (NC 0233, 0234, 0240, 0343, 0771, 0905)	◆
BNG0220		Flex 3 Backrest Tension (Width Adjustable) (NC 0233, 0234, 0343, 0771, 0925)	£94
BNG0221		Matrx Elite (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£437
BNG0223		Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£485
BNG0225		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£437
BNG0228		Matrx Elite TR (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0800, 0820, 0880, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£361
BNG0230		Backrest Tubes Without Backrest (NC 0233, 0234, 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0345, 0740, 0771, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	○
MATRX BACKREST HEIGHT			
BNG0233		Matrx Backrest Height 410 mm / 16"	○
BNG0234		Matrx Backrest Height 510 mm / 20"	○

BACKREST CUSHION*				
BNG0240		Cover for tension adjustable upholstery	(NC 2222)	○
BNG0250		Cushion, Laguna (Lateral Support)	(NC 0905, 2222)	○
BNG0260		Cushion, Mistral 2 (Waist Support)	(NC 0322, 0905, 2222)	£133
BNG0270		Cushion, Passad 2 (Shoulder Support)	(NC 0905, 2222)	◆
BNG0280		Cushion, Shoulder High (Multi-stretch Black PU coated fabric)	(NC 0322, 0905, 2222)	£138
BNG0290		Cushion, Vicair Multifunctional (Black)	(NC 0321, 0322, 0905)	£537
BNG0300		Without Backrest Cover/ Cushion	(NC 0321, 0322, 0905)	○
BACKREST UPHOLSTERY*				
BNG0321		Multi-stretch PU coated fabric		◆
BNG0322		Black Spacer TR36		○
BNG2222		No upholstery, select with Matrx, Multifunctional and Without backrest options		○
PUSH HANDLES*				
BNG0335		Push Handles connected,no head/neckrest attachment		◆
BNG0343		Push Handles connected,with Rea head/neckrest attachment	(NC 0771)	£66
BNG0345		Push Bar, angle adjustable		£94
SEAT TILT AND BACKREST REGULATION*				
BNG0350		Assistant operated gas piston seat tilt and backrest recline	(NC 5555)	◆
BNG0360		Electrical operated seat tilt and backrest recline	(NC 1263)	£673
TILT LOCK*				
BNG1263		Tilt Lock (locking of pre-set angle)		◆
BNE5555		No Tilt Lock, Select with 0360 Electric operated seat and backrest		○
CENTRAL LEGREST				
BNG0370		Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444	(NC 0189, 0376, 0378, 0380, 0390, 0410, 0430, 0440, 0460, 0475, 0476, 0477, 0720, 0930, 0940, 0941, 0942, 1290)	£279
LEGREST HANGER ATTACHMENT*				
BNG0375		Standard position for legrest mounting		◆
SWO0189		Standard position with 0-30° Legrest Abduction function		£92
BNG0376		High position 50mm higher for leg length compensation		○
BNG0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest	(NC 0390, 0410, 0440, 0460)	○
LEGREST HANGER LH*				
BNG0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH	(NC 0475, 0480, 0490, 0942, 0945)	○
BNG0390		Legrest Hanger fixed 80 deg, swing in/swing out LH		◆
BNG0410		Legrest Amputee LH	(NC 0475, 0476, 0477, 0480, 0490, 0510, 0520, 0930, 0941, 0942, 0945)	£126
BNG0420		Without Legrest, LH	(NC 0475, 0476, 0477, 0480, 0490, 0510, 0520, 0930, 0941, 0942, 0945)	○
LEGREST HANGER RH*				
BNG0430		Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support RH	(NC 0475, 0540, 0550, 0942, 0945)	◆
BNG0440		Legrest Hanger fixed 80 deg, swing in/swing out RH		○
BNG0460		Legrest Amputee RH	(NC 0475, 0476, 0477, 0540, 0550, 0570, 0580, 0940, 0941, 0942, 0945)	£126
BNG0470		Without Legrest, RH	(NC 0475, 0476, 0477, 0540, 0550, 0570, 0580, 0940, 0941, 0942, 0945)	○

FOOTREST LH*			
BNG0475		Footrest one piece, depth and angle adjustable short tube	(NC 0477, 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940) £109
BNG0476		Footrest one piece, depth and angle adjustable long tube	(NC 0477, 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940) £109
BNG0477		Footbox, angle adjustable, padded with short tube	(NC 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940, 0941, 0942, 0945) £501
BNG0480		Footrest lower part, fixed with short tube and heel strap, LH	(NC 0930, 0941) ○
BNG0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH	○
BNG0510		Footrest lower part, fixed with long tube and heel strap, LH.	(NC 0930, 0941) ○
BNG0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH	○
BNG3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○
FOOTREST RH*			
BNG0540		Footrest lower part, fixed with short tube and heel strap, RH	(NC 0940, 0941) ○
BNG0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH	○
BNG0570		Footrest lower part, fixed with long tube and heel strap, RH	(NC 0940, 0941) ○
BNG0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH	○
BNG4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BNG0600		Armrest Ducktail Pad LH, height adjustable with autolock	(NC 1232) ○
BNG0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only	(NC 1232) ◆
BNG0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock	(NC 1230, 1231, 1232, 1250, 1260) £82
BNG0610		Armrest Padded LH, height adjustable with autolock	(NC 1230, 1231, 1232, 1250, 1260) £59
BNG0615		Armrest Hemiplegic LH, height adjustable with autolock	(NC 0695, 1230, 1231, 1232, 1240, 1245, 1250, 1260) £165
BNG0620		Without Armrest LH Select when choosing BNG1232 Half Tray LH, Narrow Armrest is included	(NC 0695, 1230, 1231, 1240, 1245, 1250, 1260) ○
ARMREST RH*			
BNG0630		Armrest Ducktail Pad RH, height adjustable with autolock	(NC 1233) ○
BNG0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only	(NC 1233) ◆
BNG0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock	(NC 1230, 1231, 1233, 1250, 1260) £82
BNG0640		Armrest Padded RH, height adjustable with autolock.	(NC 1230, 1231, 1233, 1250, 1260) £59
BNG0645		Armrest Hemiplegic, RH, height adjustable with autolock	(NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260) £165
BNG0650		Without Armrest RH Select when choosing BNG1233 Half Tray RH, Narrow Armrest is included	(NC 0695, 1230, 1231, 1240, 1245, 1250, 1260) ○
ARMREST OPTIONS. Side support cover 0680, must be selected if 0690 Foam for side support is required			
BNG0680		Side Support Cover, removable with tools (wipeable cover for the sideguard)	(NC 0695) £47
BNG0690		Foam for side support cover (Internal cushioned pad, which is only possible in combination with 0680)	(NC 0695) ○
BNG0695		Pad for sideguard, removable with hook and loop fastening	£47

SEAT OPTIONS			
BNG0710		Incontinence Cover	£37
BNG0720		Pommel Abduction Cushion	£109
BNG0725		Tilt Scale for seat tilt	£41
BACKREST OPTIONS			
BNG0730		Headrest, Sideway Adjustable (NC 0750, 0770)	◆
BNG0750		Headrest with cheek support, sideway adjustable (NC 0740, 0770, 0780)	£99
BNG0740		Cover for headrest	£28
BNG0770		Neckrest (NC 0771)	○
BNG0771		Headrest Matrix, 250 mm Soft Pad (NC 0780, 0790)	£171
BNG0780		Cover For Neckrest	£20
BNG0782		Alternative headrest/neckrest mounting bracket	£21
BNG0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	○
BNG0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included	◆
BNG0905		Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery	£56
BNG0910		Lumbar Support Cushion	£28
BNG0915		Tilt Scale for backrest recline	£41
BNG0920		Wedges (One Pair)	£9
BNG0925		Privacy Attachment (Material connecting backrest and seat)	£31
TRUNK SUPPORT LH			
BNG0800		Trunk Support Multi Adjustable with standard arm LH (NC 0900)	○
BNG0820		Trunk Support Multi Adjustable with extended height arm LH (NC 0900)	○
BNG0880		Trunk Support Swingaway LH (NC 0900)	£97
BNG0895		Trunk Support Multi Adjustable Swingaway LH (NC 0900)	£97
TRUNK SUPPORT RH			
BNG0810		Trunk Support Multi Adjustable with standard arm RH (NC 0900)	○
BNG0830		Trunk Support Multi Adjustable with extended height arm RH (NC 0900)	○
BNG0890		Trunk Support Swingaway RH (NC 0900)	£97
BNG0896		Trunk Support Multi Adjustable Swingaway RH (NC 0900)	£97
LEGREST OPTIONS			
BNG0930		Footplate Extender LH (NC 0941)	£31
BNG0940		Footplate Extender RH (NC 0941)	£31
BNG0941		Cover Footplate, Black Plush Fabric	£49
BNG0942		Calf Strap Padded	£67
BNG0945		Calf Strap (for One piece footrest)	○
CHASSIS COLOUR*			
BNG0944		Silver Pearl	◆
BNG0946		Azurite (Dark Blue)	○
BNG0947		Sand (Beige)	○
BNG0948		Electric Green	○
BNG0949		Happy Red	○
CHASSIS SAFETY*			
BNG0950		Chassis extender, castors extend 50mm forward for added stability	◆

CASTOR ATTACHMENT*			
BNG0960	<input type="checkbox"/>	Fixed castor attachment	◆
FRONT CASTORS*			
BNG1015	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	£19
BNG1025	<input type="checkbox"/>	200x45 mm (8" X 1 3/4") Black Solid tyre	◆
REAR WHEEL SIZE*			
BNG1047	<input type="checkbox"/>	16" / 405 mm Transit rear wheel with fixed attachment	◆
REAR WHEEL TYPE*			
BNG1076	<input type="checkbox"/>	Hub Brake 16" Transit with Black Plastic Hub and Rim	○
REAR WHEEL TYRE*			
BNG1092	<input type="checkbox"/>	Solid Black tyres	◆
BRAKE			
BNG1170	<input type="checkbox"/>	User manoeuvred tyre brakes	◆
SAFETY OPTIONS			
BNG1220	<input type="checkbox"/>	Reflector Kit (wheels and chassis)	£20
PELVIC BELT*			
BNG1200	<input type="checkbox"/>	Pelvic belt	◆
TRANSPORTATION BRACKETS*			
BNG1225	<input type="checkbox"/>	Transportation Brackets, mounted to the front of the wheelchair	◆
OTHER OPTIONS			
BNG1230	<input type="checkbox"/>	Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	£127
BNG1231	<input type="checkbox"/>	Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	£127
BNG1232	<input type="checkbox"/>	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	£95
BNG1233	<input type="checkbox"/>	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	£95
BNG1240	<input type="checkbox"/>	Swing away table without elbow pads (NC 1235, 1250, 1260)	£140
BNG1245	<input type="checkbox"/>	Swing away table with elbow pads (NC 1235, 1260)	£140
BNG1235	<input type="checkbox"/>	Table tray padded cushion	£47
BNG1237	<input type="checkbox"/>	Half table tray padded cushion	£43
BNG1250	<input type="checkbox"/>	Locking device for table tray	£13
BNG1260	<input type="checkbox"/>	Mounting kit for table (if ordered separately)	£7
BNG1270	<input type="checkbox"/>	Pump	£20
BNG1280	<input type="checkbox"/>	Cane Holder	£46
BNG1290	<input type="checkbox"/>	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	£50
BNG1300	<input type="checkbox"/>	Harness attachment , must be selected in combination with Head/Neckrest	£79
BNG1310	<input type="checkbox"/>	Attachment for positioning belt	£29

TOTAL PRICE	
--------------------	--

Non Conforming Options:

BNG0201	Multi-stretch Black PU coated fabric
BNG0202	Black Spacer TR36
BNG1111	No upholstery
BNG0206	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BNG0220	Flex 3 Backrest Tension (Width Adjustable)
BNG0221	Matrx Elite (includes hardware mounted to frame, backrest supplied in box)
BNG0223	Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box)
BNG0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BNG0228	Matrx Elite TR (includes hardware mounted to frame, backrest supplied in box)
BNG0233	Matrx Backrest Height 410 mm / 16"
BNG0234	Matrx Backrest Height 510 mm / 20"
BNG0240	Cover for tension adjustable upholstery
BNG0250	Cushion, Laguna (Lateral Support)
BNG0260	Cushion, Mistral 2 (Waist Support)
BNG0270	Cushion, Passad 2 (Shoulder Support)
BNG0280	Cushion, Shoulder High (Multi-stretch Black PU coated fabric)
BNG0290	Cushion, Vicair Multifunctional (Black)
BNG0321	Multi-stretch PU coated fabric
BNG0322	Black Spacer TR36
BNG2222	No upholstery, select with Matrx, Multifunctional and Without backrest options
BNG0335	Push Handles connected,no head/neckrest attachment
BNG0343	Push Handles connected,with Rea head/neckrest attachment
BNG0345	Push Bar, angle adjustable
BNG0360	Electrical operated seat tilt and backrest recline
BNG1263	Tilt Lock (locking of pre-set angle)
BNE5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
SWO0189	Standard position with 0-30° Legrest Abduction function
BNG0376	High position 50mm higher for leg length compensation
BNG0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BNG0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH
BNG0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNG0410	Legrest Amputee LH
BNG0430	Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support RH
BNG0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNG0460	Legrest Amputee RH
BNG0475	Footrest one piece, depth and angle adjustable short tube
BNG0476	Footrest one piece, depth and angle adjustable long tube
BNG0477	Footbox, angle adjustable, padded with short tube
BNG0480	Footrest lower part, fixed with short tube and heel strap, LH
BNG0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BNG0510	Footrest lower part, fixed with long tube and heel strap, LH.
BNG0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH
BNG0540	Footrest lower part, fixed with short tube and heel strap, RH

BNG0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BNG0570	Footrest lower part, fixed with long tube and heel strap, RH
BNG0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BNG0695	Pad for sideguard, removable with hook and loop fastening
BNG0720	Pommel Abduction Cushion
BNG0750	Headrest with cheek support, sideway adjustable
BNG0740	Cover for headrest
BNG0770	Neckrest
BNG0771	Headrest Matrix, 250 mm Soft Pad
BNG0780	Cover For Neckrest
BNG0790	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BNG0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNG0905	Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery
BNG0910	Lumbar Support Cushion
BNG0920	Wedges (One Pair)
BNG0925	Privacy Attachment (Material connecting backrest and seat)
BNG0800	Trunk Support Multi Adjustable with standard arm LH
BNG0820	Trunk Support Multi Adjustable with extended height arm LH
BNG0880	Trunk Support Swingaway LH
BNG0895	Trunk Support Multi Adjustable Swingaway LH
BNG0810	Trunk Support Multi Adjustable with standard arm RH
BNG0830	Trunk Support Multi Adjustable with extended height arm RH
BNG0890	Trunk Support Swingaway RH
BNG0896	Trunk Support Multi Adjustable Swingaway RH
BNG0930	Footplate Extender LH
BNG0940	Footplate Extender RH
BNG0941	Cover Footplate, Black Plush Fabric
BNG0942	Calf Strap Padded
BNG0945	Calf Strap (for One piece footrest)
BNG1230	Fixed Table Tray without elbow pads includes mounting kit
BNG1231	Fixed Table Tray with elbow Pads includes mounting kit
BNG1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNG1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNG1240	Swing away table without elbow pads
BNG1245	Swing away table with elbow pads
BNG1235	Table tray padded cushion
BNG1237	Half table tray padded cushion
BNG1250	Locking device for table tray
BNG1260	Mounting kit for table (if ordered separately)
BNG1290	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):