



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

MARCH 2021  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Dahlia 45

### Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Choice of Seat Heights 380mm - 500mm (15" - 20")
- Depth Adjustable Seat 430mm - 500mm (17" - 20")
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Black Multi-stretch PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +45° and Backrest Recline 30° with Tilt Lock
- Fixed 80° Legrests Swing In Swing Out. Height, Depth and Angle Adjustable Footplates
- Height Adjustable Ducktail Pad Armrests with Autolock
- Black Frame Colour, Solid Black Front Castors
- User Manoeuvred Tyre Brakes, Pelvic Belt, Adj both sides
- Anti Tippers, Front Transportation Brackets



**Standard Base price: £2,390**

Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
BXA0015	INVACARE DAHLIA 45 °	£2,390
BXA8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
<b>PROPULSION*</b>		
BXA0020	Transit 16" wheels with hub brakes	<input checked="" type="checkbox"/>

SEAT WIDTH*			
BXA0050		390 mm/15" (Small chassis)	○
BXA0060		440 mm/17" (Medium chassis)	○
BXA0070		490 mm/19" (Large chassis)	○
SEAT WIDTH EXTENSION*			
BXA0080		+ 50 mm seat width pre-adjustment on armrests and legrests	○
BXA0090		+ 100 mm seat width pre-adjustment on armrests and legrests	○
BXA0100		Without seat width extension, set at 0mm adjustment (NC 0420, 0430) but can be adjusted wider if required	○
SEAT DEPTH*			
BXA0110		430-500 mm / 17" - 20"	◆
SEAT HEIGHT*			
BXA0130		380 mm / 15" (NC 0630, 1415, 1435, 1455, 1475)	○
BXA0170		480 mm / 19" (NC 1400)	○
BXA0180		500 mm / 20" (NC 0680, 0730, 1400, 1415, 1435)	○
SEAT TILT ANGLE*			
BXA0200		Full Tilt: Positive Seat Angle 45°	◆
BACKREST RECLINE ANGLE*			
BXA0207		Full Backrest Recline: Angle 30°	◆
SEAT CUSHION*			
BXA0230		Cushion Flosshape (NC 1111)	£48
SWO0190		Cushion Flosshape seat upholstery Shorter depth 450mm deep. Only use with seat depth in shortest position (NC 0305, 1111)	£48
BXA0270		Cushion Vicair Multifunctional (Black) (NC 0300, 0305)	£559
BXA0280		Without Seat Cushion, Seat Plate only (NC 0300, 0305)	◆
SEAT UPHOLSTERY*			
BXA0300		Black Multi-stretch PU coated fabric	◆
BXA0305		Black Spacer TR36	○
BXA1111		No upholstery, select with Matrix, Multifunctional and Without backrest options	○
BACKREST TUBES*			
BXA0310		Bent Back Tubes (20° from 380 mm above the seatplate - Matrix only) (NC 0390, 0400)	○
BXA0320		Straight Back Tubes (Flex 3 and Matrix)	◆
BACKREST RECLINE PIVOT POINT*			
BXA0330		Plus 50mm width extension, outrigs backrest pivot points by 5cm overall to allow a wider backrest to be fitted to a smaller chassis (NC 0430, 0595)	£78
BXA0340		Without Backrest Pivot Point Width Extension	◆
BACKREST*			
BXA0350		Matrix Elite (includes hardware mounted to frame, backrest supplied in box) (NC 0420, 0430, 0450, 0460, 0470, 0480, 0490, 0500, 0530, 0535, 1100, 1110, 1120, 1130, 1200, 1210, 1220, 1225, 1230, 1240, 1250)	£437
BXA0360		Matrix PB Elite Deep (includes hardware mounted to frame, backrest supplied in box) (NC 0420, 0430, 0450, 0460, 0470, 0480, 0490, 0500, 0530, 0535, 1100, 1110, 1120, 1130, 1200, 1210, 1220, 1225, 1230, 1240, 1250)	£485
BXA0365		Matrix Elite TR (includes hardware mounted to frame, backrest supplied in box) (NC 0420, 0430, 0450, 0460, 0470, 0480, 0490, 0500, 0530, 0535, 1100, 1110, 1120, 1130, 1200, 1210, 1220, 1225, 1230, 1240, 1250)	£361
BXA0390		Flex 3 Backrest, Fixed Plate (Width Adjustable) (NC 0413, 0417, 0450, 0570, 1072, 1225)	◆
BXA0400		Flex 3 Backrest Tension (Width Adjustable) (NC 0413, 0417, 0430, 0570, 1072, 1250)	£94
BXA0410		Backrest Tubes Without Backrest (NC 0413, 0417, 0420, 0430, 0450, 0460, 0470, 0480, 0490, 0500, 0530, 0535, 1072, 1100, 1110, 1120, 1130, 1200, 1210, 1220, 1225, 1230, 1240, 1250)	○

<b>MATRIX BACKREST HEIGHT</b>			
BXA0413		Matrix Backrest Height 410 mm / 16"	○
BXA0417		Matrix Backrest Height 510 mm / 20"	○
<b>BACKREST WIDTH EXTENSION*</b>			
BXA0420		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required	○
BXA0430		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required	○
BXA0440		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required	○
<b>BACKREST CUSHION*</b>			
BXA0450		Cover for tension adjustable upholstery (NC 2222)	○
BXA0460		Cushion, Laguna (Lateral Support) (NC 1225, 2222)	○
BXA0470		Cushion, Mistral 2 (Waist Support) (NC 0535, 1225, 2222)	£133
BXA0480		Cushion, Passad 2 (Shoulder Support) (NC 1225, 2222)	◆
BXA0490		Cushion, Shoulder High. Higher Cushion to suit taller clients (NC 0535, 1225, 2222)	£138
BXA0500		Cushion, Vicair Multifunctional (Black) (NC 0530, 0535, 1225)	£537
BXA0510		Without Backrest Cover/ Cushion (NC 0530, 0535, 1225)	○
<b>BACKREST UPHOLSTERY*</b>			
BXA0530		Black Multi-stretch PU coated fabric	◆
BXA0535		Black Spacer TR36	○
BXA2222		No upholstery, select with Matrix, Multifunctional and Without backrest options	○
<b>PUSH HANDLES*</b>			
BXA0540		Push Bar, Angle Adjustable	£94
BXA0560		Push Handles connected, no head/neckrest attachment	◆
BXA0570		Push Handles connected, with head/neckrest attachment (NC 1072)	£66
<b>SEAT TILT*</b>			
BXA0580		Assistant Operated gas piston seat tilt (NC 5555)	◆
BXA0595		Electrical Operated seat and backrest tilt (NC 0600, 0630)	£673
<b>BACKREST REGULATION*</b>			
BXA0600		Assistant Operated gas piston for backrest recline	◆
BXA5555		Electric backrest regulation, select with option 0595	○
<b>CENTRAL LEGREST</b>			
BXA0630		Central Legrest, includes calfpad and footrest. Also select options 0700, 0750, 3333, 4444 (NC 0641, 0650, 0655, 0660, 0670, 0680, 0690, 0710, 0720, 0730, 0740, 0760, 0770, 1000, 1270, 1280, 1300, 1310, 1800)	£279
<b>LEGREST HANGER ATTACHMENT*</b>			
BXA0640		Standard position for legrest mounting	◆
SWO0641		Standard position with 0-30° Legrest Abduction function	£92
BXA0650		High position 50mm higher for leg length compensation	○
BXA0655		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest (NC 0670, 0680, 0690, 0720, 0730, 0740)	○
<b>LEGREST HANGER LH*</b>			
BXA0660		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH (NC 0760, 0790, 1290)	£92
BXA0670		Legrest Hanger fixed 80 deg, swing in/swing out LH	◆
BXA0680		Legrest Hanger fixed 90 deg, swing in/swing out LH (NC 0810, 1270, 1455, 1475)	○
BXA0690		Legrest Amputee LH (NC 0760, 0770, 0790, 0810, 1270, 1300, 1310)	£126
BXA0700		Without Legrest, LH (NC 0760, 0770, 0790, 0810, 1270, 1300, 1310)	○

LEGREST HANGER RH*				
BXA0710		Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support RH	(NC 0760, 0830, 1290)	£92
BXA0720		Legrest Hanger fixed 80 deg, swing in/swing out RH		◆
BXA0730		Legrest Hanger fixed 90 deg, swing in/swing out RH	(NC 0850, 1280, 1455, 1475)	○
BXA0740		Legrest Amputee RH	(NC 0760, 0770, 0830, 0850, 1280, 1300, 1310)	£126
BXA0750		Without Legrest, RH	(NC 0760, 0770, 0830, 0850, 1280, 1300, 1310)	○
FOOTREST LH*				
BXA0760		Footrest one piece, depth and angle adjustable short tube	(NC 0790, 0810, 0830, 0850, 1270, 1280)	£109
BXA0770		Footrest one piece, depth and angle adjustable long tube	(NC 0790, 0810, 0830, 0850, 1270, 1280)	£109
BXA0790		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH		○
BXA0810		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH		○
BXA3333		Without footrest, select with options Central legrest, Amputee legrest or Without legrest		○
FOOTREST RH*				
BXA0830		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH		○
BXA0850		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH		○
BXA4444		Without footrest, select with options Central legrest, Amputee legrest, One Piece or Without legrest		○
ARMREST LH*				
BXA0870		Armrest Ducktail Pad LH, Height Adjustable from 200 mm - 320mm with Autolock	(NC 1716)	◆
BXA0875		Armrest Extended Height, Ducktail pad LH, Adjusts from 255 mm - 355 mm with Autolock	(NC 1716)	○
BXA0880		Armrest Wide Pad LH, Height Adjustable from 200 mm - 320mm with autolock.	(NC 1716)	○
BXA0890		Armrest Hemiplegic LH, Height Adjustable from 200 mm - 320 mm with autolock	(NC 0980, 1710, 1715, 1716, 1720, 1730, 1735, 1740, 1750)	£165
BXA0910		Without Armrest LH Select when choosing BXA1716 Half Tray LH, Narrow Armrest is included	(NC 0980, 1710, 1715, 1720, 1730, 1735, 1740, 1750)	○
ARMREST RH*				
BXA0930		Armrest Ducktail Pad RH, Height Adjustable from 200 mm - 320mm with Autolock	(NC 1717)	◆
BXA0935		Armrest Extended Height, Ducktail pad RH, Adjusts from 255 mm - 355 mm with Autolock	(NC 1717)	○
BXA0940		Armrest Wide Pad RH, Height Adjustable from 200 mm - 320mm with autolock.	(NC 1717)	○
BXA0950		Armrest Hemiplegic RH, Height Adjustable from 200 mm - 320 mm with autolock	(NC 0980, 1710, 1715, 1717, 1720, 1730, 1735, 1740, 1750)	£165
BXA0970		Without Armrest RH Select when choosing BXA1717 Half Tray RH, Narrow Armrest is included	(NC 0980, 1710, 1715, 1720, 1730, 1735, 1740, 1750)	○
ARMREST OPTIONS				
BXA0980		Pad for sideguard, removable with hook and loop fastening		£47

SEAT OPTIONS			
BXA0990		Incontinence Cover	£37
BXA1000		Pommel Abduction Cushion (NC 1020)	£109
BXA1010		Seat Plate Extension Plates. Width, Forward, Rear or Assymetric (Supplied as a kit of 2 metal plates unmounted)	£39
BXA1020		2 cm Increased seat height plate. Allows an increased 20mm of seat height by adding a foam plate to the seat board with hook and loop	£59
BXA1030		Tilt Scale for seat tilt	£41
BACKREST OPTIONS			
BXA1040		Headrest, Sideway Adjustable (NC 1070, 1072, 1080)	○
BXA1060		Headrest with cheek support, sideway adjustable (NC 1050, 1070, 1080)	£99
BXA1050		Cover for headrest (NC 1070, 1072, 1080)	£28
BXA1070		Neckrest	◆
BXA1072		Headrest Matrix, 250 mm Soft Pad (NC 1080, 1090)	£171
BXA1080		Cover For Neckrest	£20
BXA1085		Alternative headrest/neckrest mounting bracket	£21
BXA1090		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	£45
BXA1220		Trunk Support Attachment. If trunk supports are chosen, the attachment bracket is already included	○
BXA1225		Lateral Positioning Pads. Attaches with hook and loop between the straps and the cover of the tension adjustable upholstery	£56
BXA1230		Lumbar Support Cushion	£28
BXA1240		Wedges (One Pair)	£9
BXA1250		Privacy Attachment (Material connected between backrest and seat)	£31
BXA1260		Tilt Scale for backrest recline	£41
TRUNK SUPPORTS LH			
BXA1100		Trunk Support Multi Adjustable with standard arm LH (NC 1220)	○
BXA1120		Trunk Support Multi Adjustable with extended height arm LH (NC 1220)	○
BXA1200		Trunk Support Multi Adjustable Swingaway LH (NC 1220)	£104
TRUNK SUPPORTS RH			
BXA1110		Trunk Support Multi Adjustable with standard arm RH (NC 1220)	○
BXA1130		Trunk Support Multi Adjustable with extended height arm RH (NC 1220)	○
BXA1210		Trunk Support Multi Adjustable Swingaway RH (NC 1220)	£104
LEGREST OPTIONS			
BXA1270		Footplate Extender LH (NC 1300)	£31
BXA1280		Footplate Extender RH (NC 1300)	£31
BXA1290		Calf Strap (for One piece footrest)	○
BXA1300		Cover Footplate, Black Plush Fabric	£49
BXA1310		Calf Strap Padded	£67
CHASSIS COLOUR*			
BXA1365		Black Frame	◆
CASTOR ATTACHMENT*			
BXA1370		Fixed	◆

<b>FRONT CASTORS*</b>			
BXA1400	<input type="checkbox"/>	125x27 mm (5" x 1") Medium Hard Black Tyre	○
BXA1415	<input type="checkbox"/>	140x45 mm (5 1/2" x 1 3/4") Soft Black Tyre	£19
BXA1435	<input type="checkbox"/>	150X30 mm (6" x 1 1/4") Medium Hard Black Tyre	○
BXA1455	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	○
BXA1475	<input type="checkbox"/>	200x45 mm (8" x 1 3/4") Black Solid Tyre	◆
<b>REAR WHEEL SIZE*</b>			
BXA1505	<input type="checkbox"/>	16" / 405 mm Transit rear wheel with fixed attachment	◆
<b>REAR WHEEL TYPE*</b>			
BXA1546	<input type="checkbox"/>	Hub Brake 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes	◆
<b>REAR WHEEL TYRE*</b>			
BXA1585	<input type="checkbox"/>	Solid Black Tyres	◆
<b>BRAKE*</b>			
BXA1680	<input type="checkbox"/>	User manoeuvred tyre brakes	◆
<b>PELVIC BELT*</b>			
BXA1690	<input type="checkbox"/>	Pelvic belt	◆
<b>SAFETY OPTIONS</b>			
BXA1700	<input type="checkbox"/>	Reflector Kit (wheels and chassis)	£20
<b>OTHER OPTIONS</b>			
BXA1710	<input type="checkbox"/>	Table Tray without elbow pads includes mounting kit (NC 1725, 1750)	£127
BXA1715	<input type="checkbox"/>	Table Tray with elbow Pads includes mounting kit (NC 1725, 1750)	£127
BXA1716	<input type="checkbox"/>	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1720, 1750)	£95
BXA1717	<input type="checkbox"/>	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1720, 1750)	£95
BXA1730	<input type="checkbox"/>	Swing away table without elbow pads (NC 1720, 1725, 1740, 1750)	£140
BXA1735	<input type="checkbox"/>	Swing away table with elbow pads (NC 1720, 1725, 1740, 1750)	£140
BXA1720	<input type="checkbox"/>	Table tray padded cushion (NC 1750)	£47
BXA1725	<input type="checkbox"/>	Half table tray padded cushion (NC 1750)	£43
BXA1740	<input type="checkbox"/>	Locking device for table tray (NC 1750)	£13
BXA1750	<input type="checkbox"/>	Mounting kit for table (if ordered separately)	£7
BXA1780	<input type="checkbox"/>	Pump	£20
BXA1790	<input type="checkbox"/>	Cane Holder	£46
BXA1800	<input type="checkbox"/>	Clinic Kit	£50
BXA1810	<input type="checkbox"/>	Harness Attachment Bracket Seat (Pair) on seat	£21
BXA1820	<input type="checkbox"/>	Harness Attachment Bracket Upper	£21
BXA1830	<input type="checkbox"/>	Kerb Lever (NC 1860)	£39
BXA1860	<input type="checkbox"/>	Accessories Tray to suit Ventilators etc. Positioned at the rear of the frame between the rear wheels	£248
BXA1890	<input type="checkbox"/>	Restrictor Clip kit, Plastic clips which attach to the gas struts to limit the seat tilt or backrest recline	£6
<b>TRANSPORTATION BRACKETS*</b>			
BXA1850	<input type="checkbox"/>	Transportation Brackets	◆

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

BXA0300	Black Multi-stretch PU coated fabric
BXA0305	Black Spacer TR36
BXA1111	No upholstery, select with Matrix, Multifunctional and Without backrest options
BXA0390	Flex 3 Backrest, Fixed Plate (Width Adjustable)
BXA0400	Flex 3 Backrest Tension (Width Adjustable)
BXA0413	Matrix Backrest Height 410 mm / 16"
BXA0417	Matrix Backrest Height 510 mm / 20"
BXA0420	+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required
BXA0430	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BXA0450	Cover for tension adjustable upholstery
BXA0460	Cushion, Laguna (Lateral Support)
BXA0470	Cushion, Mistral 2 (Waist Support)
BXA0480	Cushion, Passad 2 (Shoulder Support)
BXA0490	Cushion, Shoulder High. Higher Cushion to suit taller clients
BXA0500	Cushion, Vicair Multifunctional (Black)
BXA0530	Black Multi-stretch PU coated fabric
BXA0535	Black Spacer TR36
BXA2222	No upholstery, select with Matrix, Multifunctional and Without backrest options
BXA0570	Push Handles connected, with head/neckrest attachment
BXA0595	Electrical Operated seat and backrest tilt
BXA0600	Assistant Operated gas piston for backrest recline
BXA5555	Electric backrest regulation, select with option 0595
BXA0630	Central Legrest, includes calfpad and footrest. Also select options 0700, 0750, 3333, 4444
SWO0641	Standard position with 0-30° Legrest Abduction function
BXA0650	High position 50mm higher for leg length compensation
BXA0655	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BXA0660	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH
BXA0670	Legrest Hanger fixed 80 deg, swing in/swing out LH
BXA0680	Legrest Hanger fixed 90 deg, swing in/swing out LH
BXA0690	Legrest Amputee LH
BXA0710	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support RH
BXA0720	Legrest Hanger fixed 80 deg, swing in/swing out RH
BXA0730	Legrest Hanger fixed 90 deg, swing in/swing out RH
BXA0740	Legrest Amputee RH
BXA0760	Footrest one piece, depth and angle adjustable short tube
BXA0770	Footrest one piece, depth and angle adjustable long tube
BXA0790	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BXA0810	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH
BXA0830	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BXA0850	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BXA0980	Pad for sideguard, removable with hook and loop fastening
BXA1000	Pommel Abduction Cushion
BXA1020	2 cm Increased seat height plate. Allows an increased 20mm of seat height by adding a foam plate to the seat board with hook and loop

BXA1050	Cover for headrest
BXA1070	Neckrest
BXA1072	Headrest Matrix, 250 mm Soft Pad
BXA1080	Cover For Neckrest
BXA1090	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BXA1220	Trunk Support Attachment. If trunk supports are chosen, the attachment bracket is already included
BXA1225	Lateral Positioning Pads. Attaches with hook and loop between the straps and the cover of the tension adjustable upholstery
BXA1230	Lumbar Support Cushion
BXA1240	Wedges (One Pair)
BXA1250	Privacy Attachment (Material connected between backrest and seat)
BXA1100	Trunk Support Multi Adjustable with standard arm LH
BXA1120	Trunk Support Multi Adjustable with extended height arm LH
BXA1200	Trunk Support Multi Adjustable Swingaway LH
BXA1110	Trunk Support Multi Adjustable with standard arm RH
BXA1130	Trunk Support Multi Adjustable with extended height arm RH
BXA1210	Trunk Support Multi Adjustable Swingaway RH
BXA1270	Footplate Extender LH
BXA1280	Footplate Extender RH
BXA1290	Calf Strap (for One piece footrest)
BXA1300	Cover Footplate, Black Plush Fabric
BXA1310	Calf Strap Padded
BXA1400	125x27 mm (5" x 1") Medium Hard Black Tyre
BXA1415	140x45 mm (5 1/2" x 1 3/4") Soft Black Tyre
BXA1435	150X30 mm (6" x 1 1/4") Medium Hard Black Tyre
BXA1455	200x32 mm (8" x 1 1/4") Black Solid Tyre
BXA1475	200x45 mm (8" x 1 3/4") Black Solid Tyre
BXA1710	Table Tray without elbow pads includes mounting kit
BXA1715	Table Tray with elbow Pads includes mounting kit
BXA1716	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BXA1717	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BXA1730	Swing away table without elbow pads
BXA1735	Swing away table with elbow pads
BXA1720	Table tray padded cushion
BXA1725	Half table tray padded cushion
BXA1740	Locking device for table tray
BXA1750	Mounting kit for table (if ordered separately)
BXA1800	Clinic Kit
BXA1860	Accessories Tray to suit Ventilators etc. Positioned at the rear of the frame between the rear wheels





For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):