



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

MARCH 2021  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Dahlia 45 Base

### Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Choice of Seat Heights 380mm - 500mm (15" - 20")
- Depth Adjustable Seat 430mm - 500mm (17" - 20")
- Wooden Seat Board
- Backrest Tubes without Backrest Plate
- Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +45° and Backrest Recline 30° with Tilt Lock
- Black Frame Colour, Solid Black Front Castors
- User Manoeuvred Tyre Brakes, Pelvic Belt, Adj both sides
- Anti Tippers, Front Transportation Brackets



**Standard Base price: £1,421**

Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
BXC0019	Dahlia 45 ° Without seating	£1,421
BXC8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
<b>PROPULSION*</b>		
BXC0020	Transit 16" wheels with hub brakes	◆

<b>SEAT WIDTH*</b>			
BXC0050	<input type="checkbox"/>	390 mm/15" (Small chassis)	<input type="radio"/>
BXC0060	<input type="checkbox"/>	440 mm/17" (Medium chassis)	<input type="radio"/>
BXC0070	<input type="checkbox"/>	490 mm/19" (Large chassis)	<input type="radio"/>
<b>SEAT WIDTH EXTENSION*</b>			
BXC0080	<input type="checkbox"/>	+ 50 mm seat width pre-adjustment on armrests and legrests	<input type="radio"/>
BXC0090	<input type="checkbox"/>	+ 100 mm seat width pre-adjustment on armrests and legrests	<input type="radio"/>
BXC0100	<input type="checkbox"/>	Without seat width extension, set at 0mm adjustment but can be adjusted wider if required	<input type="radio"/>
<b>SEAT DEPTH*</b>			
BXC0110	<input type="checkbox"/>	430-500 mm / 17" - 20"	<input checked="" type="radio"/>
<b>SEAT HEIGHT*</b>			
BXC0130	<input type="checkbox"/>	380 mm / 15"	(NC 0630, 1415, 1435, 1455, 1475) <input type="radio"/>
BXC0170	<input type="checkbox"/>	480 mm / 19"	(NC 1400) <input type="radio"/>
BXC0180	<input type="checkbox"/>	500 mm / 20"	(NC 0680, 0730, 1400, 1415, 1435) <input type="radio"/>
<b>SEAT TILT ANGLE*</b>			
BXC0200	<input type="checkbox"/>	Full Tilt: Positive Seat Angle 45°	<input checked="" type="radio"/>
<b>BACKREST RECLINE ANGLE*</b>			
BXC0207	<input type="checkbox"/>	Full Backrest Recline: Angle 30°	<input checked="" type="radio"/>
<b>SEAT DEPTH RANGE*</b>			
BXC0210	<input type="checkbox"/>	Depth Adjustable 430-500 mm / 17" - 20"	<input checked="" type="radio"/>
<b>SEAT CUSHION*</b>			
BXC0280	<input type="checkbox"/>	Without Seat Cushion, Seat Plate only	<input checked="" type="radio"/>
<b>BACKREST TUBES*</b>			
BXC0310	<input type="checkbox"/>	Bent Back Tubes (20° from 380 mm above the seatplate - Matrx only)	<input type="radio"/>
BXC0320	<input type="checkbox"/>	Straight Back Tubes (Flex 3 and Matrx)	<input checked="" type="radio"/>
<b>BACKREST RECLINE PIVOT POINT*</b>			
BXC0330	<input type="checkbox"/>	Plus 50mm width extension, outrigs backrest pivot points by 5cm overall to allow a wider backrest to be fitted to a smaller chassis	(NC 0595) <input type="radio"/>
BXC0340	<input type="checkbox"/>	Without Backrest Pivot Point Width Extension	<input type="radio"/>
<b>BACKREST*</b>			
BXC0410	<input type="checkbox"/>	Backrest Tubes Without Backrest	<input checked="" type="radio"/>
<b>PUSH HANDLES*</b>			
BXC0540	<input type="checkbox"/>	Push Bar, Angle Adjustable	£94
BXC0560	<input type="checkbox"/>	Push Handles connected, no head/neckrest attachment	<input checked="" type="radio"/>
<b>SEAT TILT*</b>			
BXC0580	<input type="checkbox"/>	Assistant Operated gas piston seat tilt	<input checked="" type="radio"/>
BXC0595	<input type="checkbox"/>	Electrical Operated seat and backrest tilt	(NC 0600, 0620, 0630) £673
<b>BACKREST REGULATION*</b>			
BXC0600	<input type="checkbox"/>	Assistant Operated gas piston for backrest recline	<input checked="" type="radio"/>
BXA5555	<input type="checkbox"/>	Electric backrest regulation, select with option 0595	<input type="radio"/>
BXC0620	<input type="checkbox"/>	Fixed Backrest Angle (Barrel Screw Adjustment)	<input type="radio"/>
<b>CENTRAL LEGREST</b>			
BXC0630	<input type="checkbox"/>	Central Legrest, includes calfpad and footrest. Also select options 0700, 0750, 3333, 4444	(NC 0641, 0650, 0655, 0660, 0670, 0680, 0690, 0710, 0720, 0730, 0740, 0760, 0770, 1000, 1270, 1280, 1300, 1310, 1800) £279

LEGREST HANGER ATTACHMENT*			
BXC0640		Standard position for legrest mounting	◆
SWO0641		Standard position with 0-30° Legrest Abduction function	£92
BXC0650		High position 50mm higher for leg length compensation	○
BXC0655		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest (NC 0670, 0680, 0690, 0720, 0730, 0740)	○
LEGREST HANGER LH*			
BXC0660		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH (NC 0760, 0790, 1290)	£92
BXC0670		Legrest Hanger fixed 80 deg, swing in/swing out LH	£52
BXC0680		Legrest Hanger fixed 90 deg, swing in/swing out LH (NC 0810, 1270, 1455, 1475)	£52
BXC0690		Legrest Amputee LH (NC 0760, 0770, 0790, 0810, 1270, 1300, 1310)	£126
BXC0700		Without Legrest, LH (NC 0760, 0770, 0790, 0810, 1270, 1300, 1310)	○
LEGREST HANGER RH*			
BXC0710		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support RH (NC 0760, 0830, 1290)	£92
BXC0720		Legrest Hanger fixed 80 deg, swing in/swing out RH	£52
BXC0730		Legrest Hanger fixed 90 deg, swing in/swing out RH (NC 0850, 1280, 1455, 1475)	£52
BXC0740		Legrest Amputee RH (NC 0760, 0770, 0830, 0850, 1280, 1300, 1310)	£126
BXC0750		Without Legrest, RH (NC 0760, 0770, 0830, 0850, 1280, 1300, 1310)	○
FOOTREST LH*			
BXC0760		Footrest one piece, depth and angle adjustable short tube (NC 0790, 0810, 0830, 0850, 1270, 1280)	£184
BXC0770		Footrest one piece, depth and angle adjustable long tube (NC 0790, 0810, 0830, 0850, 1270, 1280)	£184
BXC0790		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH	○
BXC0810		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH	○
BXC3333		Without footrest, select with options Central legrest, Amputee legrest or Without legrest	○
FOOTREST RH*			
BXC0830		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH	○
BXC0850		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH	○
BXC4444		Without footrest, select with options Central legrest, Amputee legrest, One Piece or Without legrest	○
ARMREST LH*			
BXC0870		Armrest Ducktail Pad LH, Height Adjustable from 200 mm - 320mm with Autolock	£77
BXC0875		Armrest Extended Height, Ducktail pad LH, Adjusts from 255 mm - 355 mm with Autolock	£77
BXC0880		Armrest Wide Pad LH, Height Adjustable from 200 mm - 320mm with autolock.	£77
BXC0890		Armrest Hemiplegic LH, Height Adjustable from 200 mm - 320 mm with autolock (NC 0980, 1710, 1715, 1720, 1730, 1735, 1740, 1750)	£184
BXC0910		Without Armrest LH (NC 0980, 1710, 1715, 1720, 1730, 1735, 1740, 1750)	○
ARMREST RH*			
BXC0930		Armrest Ducktail Pad RH, Height Adjustable from 200 mm - 320mm with Autolock	£77
BXC0935		Armrest Extended Height, Ducktail pad RH, Adjusts from 255 mm - 355 mm with Autolock	£77
BXC0940		Armrest Wide Pad RH, Height Adjustable from 200 mm - 320mm with autolock.	£77
BXC0950		Armrest Hemiplegic RH, Height Adjustable from 200 mm - 320 mm with autolock (NC 0980, 1710, 1715, 1720, 1730, 1735, 1740, 1750)	£184
BXC0970		Without Armrest RH (NC 0980, 1710, 1715, 1720, 1730, 1735, 1740, 1750)	○

<b>ARMREST OPTIONS</b>			
BXC0980	<input type="checkbox"/>	Pad for sideguard, removable with hook and loop fastening	£47
<b>SEAT OPTIONS</b>			
BXC0990	<input type="checkbox"/>	Incontinence Cover	£37
BXC1000	<input type="checkbox"/>	Pommel Abduction Cushion (NC 1020)	£109
BXC1010	<input type="checkbox"/>	Seat Plate Extension Plates. Width, Forward, Rear or Assymmetric (Supplied as a kit of 2 metal plates unmounted)	£39
BXC1020	<input type="checkbox"/>	2 cm Increased seat height plate. Allows an increased 20mm of seat height by adding a foam plate to the seat board with hook and loop	£59
BXC1030	<input type="checkbox"/>	Tilt Scale for seat tilt	£41
<b>BACKREST OPTIONS</b>			
BXC1260	<input type="checkbox"/>	Tilt Scale for backrest recline	£41
<b>LEGREST OPTIONS</b>			
BXC1270	<input type="checkbox"/>	Footplate Extender LH (NC 1300)	£31
BXC1280	<input type="checkbox"/>	Footplate Extender RH (NC 1300)	£31
BXC1290	<input type="checkbox"/>	Calf Strap (for One piece footrest)	○
BXC1300	<input type="checkbox"/>	Cover Footplate, Black Plush Fabric	£49
BXC1310	<input type="checkbox"/>	Calf Strap Padded	£67
<b>CHASSIS COLOUR*</b>			
BXC1365	<input type="checkbox"/>	Black Frame	◆
<b>CASTOR ATTACHMENT*</b>			
BXC1370	<input type="checkbox"/>	Fixed	○
<b>FRONT CASTORS*</b>			
BXC1400	<input type="checkbox"/>	125x27 mm (5" x 1") Medium Hard Black Tyre	○
BXC1415	<input type="checkbox"/>	140x45 mm (5 1/2" x 1 3/4") Soft Black Tyre	○
BXC1435	<input type="checkbox"/>	150X30 mm (6" x 1 1/4") Medium Hard Black Tyre	○
BXC1455	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	○
BXC1475	<input type="checkbox"/>	200x45 mm (8" x 1 3/4") Black Solid Tyre	◆
<b>REAR WHEEL SIZE*</b>			
BXC1505	<input type="checkbox"/>	16" / 405 mm Transit rear wheel with fixed attachment	◆
<b>REAR WHEEL TYPE*</b>			
BXC1546	<input type="checkbox"/>	Hub Brake 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes	◆
<b>REAR WHEEL TYRE*</b>			
BXC1585	<input type="checkbox"/>	Solid Black Tyres	◆
<b>BRAKE*</b>			
BXC1680	<input type="checkbox"/>	User manoeuvred tyre brakes	◆
<b>PELVIC BELT*</b>			
BXC1690	<input type="checkbox"/>	Pelvic belt	◆
<b>SAFETY OPTIONS</b>			
BXC1700	<input type="checkbox"/>	Reflector Kit (wheels and chassis)	£19

**OTHER OPTIONS**

BXC1710		Table Tray without elbow pads includes mounting kit	(NC 1750)	£127
BXC1715		Table Tray with elbow Pads includes mounting kit	(NC 1750)	£127
BXC1730		Swing away table without elbow pads	(NC 1720, 1740, 1750)	£140
BXC1735		Swing away table with elbow pads	(NC 1720, 1740, 1750)	£140
BXC1720		Table tray padded cushion	(NC 1750)	£47
BXC1740		Locking device for table tray	(NC 1750)	£13
BXC1750		Mounting kit for table (if ordered separately)		£7
BXC1780		Pump		£20
BXC1790		Cane Holder		£46
BXC1800		Clinic Kit		£50
BXC1810		Harness Attachment Bracket Seat (Pair) on seat		£21
BXC1820		Harness Attachment Bracket Upper		£21
BXC1830		Kerb Lever	(NC 1860)	£39
BXC1860		Accessories Tray to suit Ventilators etc. Positioned at the rear of the frame between the rear wheels		£248
BXC1890		Restrictor Clip kit, Plastic clips which attach to the gas struts to limit the seat tilt or backrest recline		£6

**TRANSPORTATION BRACKETS\***

BXC1850		Transportation Brackets		◆
---------	--	-------------------------	--	---

**TOTAL PRICE**

**Non Conforming Options:**

BXC0595	Electrical Operated seat and backrest tilt
BXC0600	Assistant Operated gas piston for backrest recline
BXC0620	Fixed Backrest Angle (Barrel Screw Adjustment)
BXC0630	Central Legrest, includes calfpad and footrest. Also select options 0700, 0750, 3333, 4444
SWO0641	Standard position with 0-30° Legrest Abduction function
BXC0650	High position 50mm higher for leg length compensation
BXC0655	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BXC0660	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH
BXC0670	Legrest Hanger fixed 80 deg, swing in/swing out LH
BXC0680	Legrest Hanger fixed 90 deg, swing in/swing out LH
BXC0690	Legrest Amputee LH
BXC0710	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support RH
BXC0720	Legrest Hanger fixed 80 deg, swing in/swing out RH
BXC0730	Legrest Hanger fixed 90 deg, swing in/swing out RH
BXC0740	Legrest Amputee RH
BXC0760	Footrest one piece, depth and angle adjustable short tube
BXC0770	Footrest one piece, depth and angle adjustable long tube
BXC0790	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BXC0810	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH
BXC0830	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BXC0850	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BXC0980	Pad for sideguard, removable with hook and loop fastening
BXC1000	Pommel Abduction Cushion
BXC1020	2 cm Increased seat height plate. Allows an increased 20mm of seat height by adding a foam plate to the seat board with hook and loop
BXC1270	Footplate Extender LH
BXC1280	Footplate Extender RH
BXC1290	Calf Strap (for One piece footrest)
BXC1300	Cover Footplate, Black Plush Fabric
BXC1310	Calf Strap Padded
BXC1400	125x27 mm (5" x 1") Medium Hard Black Tyre
BXC1415	140x45 mm (5 1/2" x 1 3/4") Soft Black Tyre
BXC1435	150X30 mm (6" x 1 1/4") Medium Hard Black Tyre
BXC1455	200x32 mm (8" x 1 1/4") Black Solid Tyre
BXC1475	200x45 mm (8" x 1 3/4") Black Solid Tyre
BXC1710	Table Tray without elbow pads includes mounting kit
BXC1715	Table Tray with elbow Pads includes mounting kit
BXC1730	Swing away table without elbow pads
BXC1735	Swing away table with elbow pads
BXC1720	Table tray padded cushion
BXC1740	Locking device for table tray
BXC1750	Mounting kit for table (if ordered separately)
BXC1800	Clinic Kit
BXC1860	Accessories Tray to suit Ventilators etc. Positioned at the rear of the frame between the rear wheels





For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):