



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

MARCH 2021  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Storm 4 Modulite

### Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO 7176 -Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable Seat Widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat Depths 410mm - 510mm Stepless Adjustment (16" - 20")
- Tension Adjustable Backrest
- Standard T-Armrests, detachable with Short armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F with Fixed Footplate
- LinX Control REM110 with standard bracket
- Choice of 6 Colours
- Puncture Proof Castors
- Puncture Proof Drive Wheels
- 73.5 Ah Batteries
- 6 km/h (4mph) Motors
- 8A (200 - 250V) Charger
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers
- Rear wheel suspension

**Standard Base price: £6,302**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ATB0000	Storm 4 Modulite	£6,302
ATZ0000	Special Storm 4 Modulite (Please provide details in special section)	POA
<b>SEATING OPTIONS*</b>		
ATB2100	Modulite Seating Unit (Flex3 seating unit)	◆
ATB2111	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆

SEAT WIDTH*			
ATBXXXX		Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	○
ATB2121		Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2701, 2703, 2705)	○
ATB2122		Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest) (NC 2701, 2703, 2705)	○
ATB2123		Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2702, 2704, 2706)	○
ATB2124		Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest) (NC 2702, 2704, 2706)	○
SEAT DEPTH*			
ATB2132		410mm (16") (adjustable 410mm - 510mm stepless adjustment) (NC 2196, 2701, 2702, 2703, 2704, 2705, 2706)	○
ATB2134		460mm (18") (adjustable 410mm - 510mm stepless adjustment)	○
ATB2136		510mm (20") (adjustable 410mm - 510mm stepless adjustment) (NC 2143)	○
SEAT OPTION			
ATBSWO		"H" Frame Interface Kit	£426
SEAT CUSHION*			
ATB2142		Matrx Flo-tech Lite	£107
ATB2143		Matrx Flo-tech Contour Visco (NC 2136, 2155)	£107
ATB2147		Matrx Libra 90mm thick (NC 2155)	£372
ASB2140		Without seat cushion	○
BACKREST SUPPORT*			
ATB2162		Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	◆
ATB2163		Backposts (Short) (NC 2172, 2173)	○
PUSH HANDLES FOR BACKPOSTS*			
ATB2107		Backposts WITH push handles	◆
ATB2108		Backposts WITHOUT push handles, Push handles will be cut down on the vertical section. (NC 2220, 2219, 2227, 2228, 2245)	£34
BACKREST UPHOLSTERY / CUSHION*			
ATB2171		Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only (NC 2152, 2155, 1496, 2229, 2230, 2247, 2322, 2323, 2324, 2325, 2326, 2327)	◆
ATB2179		Flex3 Tension Adjustable (540mm or 480mm) Black multi stretch PU coated fabric only (NC 2151, 2155, 2229, 2230, 2247, 2322, 2323, 2324, 2325, 2326, 2327)	£131
ATB2172		Flex3 adjustable Backrest Plates (without Backrest Cushion) (NC 2163, 2128)	£213
ATB2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion (NC 2163, 2128)	£346
ATB2174		Matrx Posture Back, height 400mm (Mesh fabric) (NC 1496)	£472
ATB2175		Matrx Posture Back, height 500mm (Mesh fabric) (NC 1496)	£472
ATB2176		Matrx Elite Full Back, height 500mm (Mesh fabric) (NC 1496, 2121)	£472
ATBS177		Matrx Elite TR Back, height 500mm (Mesh fabric) (NC 1496, 2121)	£472
ATB2178		Matrx Elite Deep Back, height 500mm (Mesh fabric) (NC 1496)	£524
ATB2170		Without Backrest Upholstery (NC 1496, 2151, 2152, 2155, 2230, 2247)	○
UPHOLSTREY TYPE*			
ATB2151		Black fabric (NC 2179, 2170)	◆
ATB2152		Black multi stretch PU coated fabric (NC 2171, 2170)	○
ATB2155		Black Spacer fabric (NC 2171, 2179, 2143, 2147)	○

<b>ARMREST (PAIR)*</b>			
ATB2191		T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	◆
ATB2181		T-Armrest Extended height, Range of adjustment from seat plate 295mm - 365mm.	£34
ATB2183		In-rigged Armrest (reduces seat width by 25mm each side)	○
ATB2193		Flip-up Armrest (height and depth adjustable) (height adjustable from 230mm - 360mm)	£55
ATB2192		Following Armrest (height adjustable from 230mm - 300mm)	£157
ATB2194		Following Armrest (height adjustable form 300mm - 360mm)	£157
<b>ARMPAD (PAIR)*</b>			
ATB2195		Short Armpad (325mm)	◆
ATB2196		Long Armpad (375mm) (NC 2132)	○
ATB2197		Long Ergonomic Armpad (397mm) (NC 2132, 2230, 2247)	○
ATB2199		Soft Armpad (REA Azalea Wide Pad) (NC 2132, 2179)	£82
<b>HIP SUPPORTS LH (Not available with 410mm seat depth aNd MaNual CeNtre legrests)</b>			
ATB2184		SMALL 140mm x 100mm, for flip up and following armrests only (NC 2701, 2702, 2703, 2704, 2705, 2706)	£215
ATB2186		MEDIUM 230mm x 100mm, for flip up and following armrests only (NC 2701, 2702, 2703, 2704, 2705, 2706)	£215
ATB2188		LARGE 330mm x 100mm, for flip up and following armrests only (NC 2701, 2702, 2703, 2704, 2705, 2706)	£215
<b>HIP SUPPORTS RH (Not available with 410mm seat depth aNd MaNual CeNtre legrests)</b>			
ATB2185		SMALL 140mm x 100mm, for flip up and following armrests only (NC 2701, 2702, 2703, 2704, 2705, 2706)	£215
ATB2187		MEDIUM 230mm x 100mm, for flip up and following armrests only (NC 2701, 2702, 2703, 2704, 2705, 2706)	£215
ATB2189		LARGE 330mm x 100mm, for flip up and following armrests only (NC 2701, 2702, 2703, 2704, 2705, 2706)	£215
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>			
ATB2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
<b>HEADREST / NECKREST - Adequate headrest optioN must be selected*</b>			
ATB1350		neckrest: height and depth adjustable (NC 2230, 2247)	£69
ATB1432		Headrest, height, depth and sideways adjustable (NC 2230, 2247)	£193
ATB1427		Headrest with cheek support: height and depth adjustable (NC 2230, 2247)	£193
ATB2323		Elan Standard Pad 150mm (6") wide with standard cover complete with Elan Headrest Hardware (NC 2170, 2230, 2247)	£280
ATB2324		Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware (NC 2170, 2230, 2247)	£280
ATB2322		Elan Standard Pad 355mm (14") wide with standard cover complete with Elan Headrest Hardware (NC 2170, 2230, 2247)	£280
ATB2325		Elan Occipital Pad 230mm x 180mm (9" x 7") with standard cover complete with Elan Headrest Hardware (NC 2170, 2230, 2247)	£280
ATB2326		Elan Occipital Pad 305mm x 200mm (12" x 8") with standard cover complete with Elan Headrest Hardware (NC 2170, 2230, 2247)	£280
ATB2327		Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") with standard cover complete with Elan Headrest Hardware (NC 2170, 2230, 2247)	£377
ATB1496		non Invacare Headrest but profile mount fitting Flex3 hard back (NC 2170, 2171)	○
ATB1495		Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, no Headrest - no mounting hardware. (NC 2108)	○

<b>LATERAL SUPPORTS</b>			
ATB1446		Rigid pad (Azalea type, multi stretch PU coated fabric) (NC 2170, 2171)	£143
ATB1445		Multi-adjustable pad (Azalea type, multi stretch PU coated fabric) (NC 2170, 2171)	£143
<b>SEAT HEIGHT AND SEAT RISER*</b>			
ATB1023		Fixed Seat Height 450mm (17.75")	○
ATB1025		Fixed Seat Height 480mm (19")	○
ATB1030		Electrical Seat Riser (400mm - 650mm) (16" - 25.5")	£1,036
ATB1040		Electrical Seat Riser (440mm - 690mm) incl. adapter (17.25" - 27")	£1,101
<b>SEAT TILT*</b>			
ATB1002		Manual Seat Tilt by spindle	◆
ATB1014		Electrical Tilt 0° to 25°	£514
<b>BACKREST RECLINER*</b>			
ATB0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
ATB0977		Powered Recliner (0° - 30°) Foldable	£431
<b>LEGREST FIXATION POSITION*</b>			
ATB0590		Position for seat cushion up to 70mm thick	◆
<b>FOOTREST HANGER TYPE (PAIR)*</b>			
ATB0503		Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	○
ATB0533		Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm (NC 0512)	£199
ATB0532		ADM Manual : Angle adjustable with quick release lever 0° - 80°, length 290mm - 460mm (NC 0512, 0522)	£195
ATB0531		ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm (NC 0512, 0522)	£571
ATB0501		Without Hanger (NC 0512, 0504, 0513, 0522)	○
<b>FOOTPLATE TYPE*</b>			
ATB0512		Fixed Footplate (NC 0501, 0533, 0532, 0531)	○
ATB0504		Angle Adjustable Footplate (NC 0501)	£62
ATB0513		Angle and depth adjustable Footplate (Adapt) (NC 0501)	£101
ATB0522		Angle Adjustable Footboard (NC 1244, 0501, 0532, 0531)	£130
ATBXXXX		Without footplates for option 0501 without hanger (NC 0503, 0533, 0532, 0531)	○
<b>CONTROLS*</b>			
ATB2201		REM110 Drive only (NC 0977, 0531, 1437, 2209, 2210, 2230, 2247, 2257, 2233, 2240, 2270, 2219, 2245)	◆
ATB2203		REM211 Drive & Actuators (NC 1437, 2209, 2210, 2230, 2247, 2257, 2233, 2240, 2270)	○
ATB2204		REM216 Drive, Actuators & Lights (NC 2209, 2210, 2230, 2247, 2257, 2233, 2240, 2270)	£129
ATB2205		REM400, touch-screen Interface, built in Bluetooth HID (NC 2209, 2210)	£291
ATB2206		REM500, without joystick for alternative access, built in Bluetooth HID (NC 2227, 2228, 2270)	£663

CONTROL BRACKETRY*				
ATB1207		Standard Bracket RH	(NC 1404, 2206)	○
ATB1208		Standard Bracket LH	(NC 1405, 2206)	○
ATB1707		Swing Away RH Bracket	(NC 1404, 2206)	£113
ATB1708		Swing Away LH Bracket	(NC 1405, 2206)	£113
ATB2213		Centre mount swing away to RH nucleus Bracket	(NC 1404, 1405)	£426
ATB2214		Centre mount swing away to LH nucleus Bracket	(NC 1404, 1405)	£426
ATB2227		Bracket RH mounted for Attendant	(NC 2206, 1404, 1405, 1287, 2128)	£54
ATB2228		Bracket LH mounted for Attendant	(NC 2206, 1404, 1405, 1286, 2128)	£54
ATZ2239		Individual mounting option per (special) - Requires additional information / specification		POA
ATB2209		REM500 Swing away RH multi-adjustable Bracket	(NC 2201, 2203, 2204, 2205, 1404)	○
ATB2210		REM500 Swing away LH multi-adjustable Bracket	(NC 2201, 2203, 2204, 2205)	○
REM400 TOGGLE OPTION				
ATB2270		Assignable Toggle Switch Module for REM400, (as default set to left On / OFF, right speed dial)	(NC 2201, 2203, 2204)	£269
CHIN CONTROL				
ATB2229		Chin Control (Requires Primary Remote REM400 or REM500) (Only in combination with Headrest 2125) (Please select chin control bracket below)	(NC 2201, 2203, 2204, 2171, 2179, 2128)	£1,702
CHIN CONTROL BRACKETRY				
ATB1289		Chin Control Swing Away Bracket RH Manual	(NC 2201, 2203, 2204, 2171, 2179, 2322)	○
ATB1288		Chin Control Swing Away Bracket LH Manual	(NC 2201, 2203, 2204, 2171, 2179, 2322)	○
ATB1287		Chin Control Swing Away Bracket RH Powered	(NC 2201, 2203, 2204, 2171, 2179, 1404, 2227, 1415, 1312)	£2,144
ATB1286		Chin Control Swing Away Bracket LH Powered	(NC 2201, 2203, 2204, 2171, 2179, 2228, 1415, 1339, 1363)	£2,144
DUAL CONTROL				
ATB2220		REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)	(NC 1221, 1222, 2227, 2228, 2245, 2128, 2230, 2247, 2257, 2233, 2240)	£337
ATB2219		Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)	(NC 1221, 1222, 2227, 2228, 2245, 2128, 2201)	£337
DUAL CONTROL BRACKETRY				
ATB1207		Dual Control Standard Bracket RH	(NC 1221, 1222, 2227, 2228, 2245)	○
ATB1208		Dual Control Standard Bracket LH	(NC 1221, 1222, 2227, 2228, 2245)	○
INVACARE LINX IDC DUAL CONTROL				
ATB2245		IDC Dual control (Please select IDC dual control bracket below)	(NC 2220, 2219, 2227, 2228, 2201)	£641
INVACARE LINX IDC DUAL CONTROL BRACKETRY				
ATB1221		IDC Dual control Detachable Bracket RH	(NC 2220, 2219, 2227, 2228)	○
ATB1222		IDC Dual control Detachable Bracket LH	(NC 2220, 2219, 2227, 2228)	○
LINX COMPACT REMOTE				
ATB2231		Proportional Compact Remote DLX-CR400 (Please also select Proportional Bracket)	(NC 2201, 2203, 2204)	£706
ATB2232		Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry)	(NC 2201, 2203, 2204)	£943
LINX COMPACT REMOTE BRACKET				
ATB2207		Proportional Compact RH Mount Bracket	(NC 2201, 2203, 2204)	○
ATB2208		Proportional Compact LH Mount Bracket	(NC 2201, 2203, 2204)	○
ATB2213		Proportional Compact Centre Mount swing away to RH (nucleus)	(NC 2201, 2203, 2204)	£426
ATB2214		Proportional Compact Centre Mount swing away to LH (nucleus)	(NC 2201, 2203, 2204, 1404, 1405)	£426
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2201, 2203, 2204)	POA

LINX HEAD CONTROL			
ATB2230		Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M4	(NC 2171, 2179, 2197, 2201, 2203, 2204, 2220, 2229, 2233, 2323, 2324, 2322, 2325, 2326, 2327, 1350, 1432, 1427, 2257, 2240) £3,652
ATB2247		Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings	(NC 2171, 2179, 2197, 2201, 2203, 2204, 2220, 2229, 2233, 2323, 2324, 2322, 2325, 2326, 2327, 1350, 1432, 1427, 2257, 2240) £4,277
FOUR-SWITCH PROXIMITY ARRAY REMOTE			
ATZ2257		Four-Switch Proximity Array (ASL106)	(NC 2201, 2203, 2204, 2220, 2229, 2230, 2247, 2233, 2240, 2128) £5,656
FOUR-SWITCH PROXIMITY ARRAY BRACKETRY			
ATZ2275		Centre Mount, swing-away to RH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)	(NC 1404, 1405, 1415) £426
ATZ2276		Centre Mount, swing-away to LH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)	(NC 1404, 1405, 1339) £426
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	POA
WIRELESS ACCESSORIES FOR ATOM			
ATB2249		Wireless Triple Switch Receiver (ASL557-3) (Sub D plug)	£291
ATB2250		Wireless Mouse Emulator (ASL558)	£698
LINX MICRO EXTREMITY CONTROL			
ATB2233		MEC Proportional Joystick (Micro Extremity Control) (Please select LinX MEC Bracketry)	(NC 2201, 2203, 2204, 2220, 2229, 2230, 2247, 2257, 2240, 2128) £3,546
LINX MICRO EXTREMITY CONTROL BRACKET			
ATB2235		MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please also select MEC Bracketry)	(NC 1404, 1415) ○
ATB2236		MEC Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1405, 1339) ○
ATB2213		Centre Mount, swing away to RH, nucleus	(NC 1404, 1405, 1415) £426
ATB2214		Centre Mount, swing away to LH, nucleus	(NC 1404, 1405, 1339) £426
ATB2241		Centre Mount, swing away to RH, nucleus (Inc Tray ASL603, 100m x 150mm)	(NC 1404, 1405, 1415) £426
ATB2242		Centre Mount, swing away to LH, nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405, 1339) £426
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	POA
LINX HEAVY DUTY REMOTE			
ATZ2240		Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LinX Heavy Duty Bracketry)	(NC 2201, 2203, 2204, 2220, 2229, 2230, 2247, 2257, 2233, 2128) £1,977
LINX HEAVY DUTY REMOTE BRACKET			
ATZ2235		Tray Mount RH, (ASL603, 100mm x 150mm)	(NC 1404, 1415) ○
ATZ2236		Tray Mount LH, (ASL603, 100mm x 150mm)	(NC 1339) ○
ATZ2213		Centre Mount, swing away to RH, nucleus	(NC 1404, 1405, 1415) £426
ATZ2214		Centre Mount, swing away to LH, nucleus	(NC 1404, 1405, 1339) £426
ATZ2241		Centre Mount, swing away to RH, nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405, 1415) £426
ATZ2242		Centre Mount, swing away to LH, nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405, 1339) £426
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	POA
LINX 10-WAY SWITCH			
ATB2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch)	£431
LINX 10-WAY SWITCH BRACKET			
ATB2217		10-Way Switch Bracket RH, Mounted for User	(NC 2209, 2213, 2275, 2235, 2241) ○
ATB2218		10-Way Switch Bracket LH, Mounted for User	(NC 2210, 2214, 2276, 2236, 2242) ○
ATB2227		10-Way Switch Bracket RH, Mounted for Attendant	(NC 1287) £54
ATB2228		10-Way Switch Bracket LH, Mounted for Attendant	(NC 1286) £54

<b>LINX ACCESSORY</b>			
ATB2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	£86
<b>LINX USB CHARGER PORT</b>			
ATB2253		USB Charger Port DLX-USB02-A (NC 2701, 2702, 2703, 2704, 2705, 2706)	£136
<b>LINX USB CHARGER PORT BRACKET</b>			
ATB2207		USB Charger Port RH Mount	○
ATB2208		USB Charger Port LH Mount	○
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
ATB2300		To be enclosed only	○
ATB2301		To be connected and programmed (special)	POA
<b>BUDDY BUTTONS</b>			
ATB1685		Red Buddy Button	£75
ATB1686		Yellow Buddy Button	£75
ATB1687		Green Buddy Button	£75
ATB1690		Blue Buddy Button	£75
<b>EGG SWITCH</b>			
ATB2254		Egg Switch, black	£75
<b>PIKO BUTTONS</b>			
ATB1711		Red Piko Button	£135
ATB1712		Yellow Piko Button	£135
ATB1713		Green Piko Button	£135
ATB1714		Orange Piko Button	£135
ATB1715		Blue Piko Button	£135
ATB1716		Black Piko Button	£135
<b>TOGGLES / BUTTONS</b>			
ATZ2261		Stereo Toggle - direct control for Tilt only	£229
ATZ2262		Stereo Button - direct control for Tilt only	£220
ATZ2263		Four-Way Toggle - direct control for Tilt & Recline	£287
ATZ2264		Four-Way Button - direct control for Tilt & Recline	£306
ATZ2207		Toggle / Button RH Mount	○
ATZ2208		Toggle / Button LH Mount	○
ATZ2225		Without Mounting	○
<b>BODYPOINT JOYSTICK KNOBS</b>			
ATB2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£61
ATB2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£61
ATB2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£61
ATB2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£61
<b>REMOTE STOP SWITCH</b>			
ATZ2258		Remote Stop Switch (ASL504)	£896

SHROUD COLOUR*				
ATB0202		Blue		○
ATB0228		Black		○
ATB0220		Champagne		○
ATB0204		Transparent Cherry Red		○
ATB0222		White		○
ATB0216		Olive green		○
ATZ0299		Special shroud colour (RAL code required)		POA
FRONT WHEEL*				
ATB0851		3.00 - 6 Grey Pneumatic		○
ATB0853		3.00 - 6 Grey Puncture Proof (sealed)		○
ATB0855		3.00 - 6 Grey Puncture Proof (solid inlay)		£92
REAR WHEEL*				
ATB0771		3.00 - 8 Grey Pneumatic		○
ATB0773		3.00 - 8 Grey Puncture Proof (sealed)		○
ATB0750		3.00 - 8 Grey Puncture Proof (solid inlay)		£92
MOTOR TYPE*				
ATB1277		6 km/h (4mph)		◆
ATB1908		10 km/h (6mph) (includes LED lights system with indicators)	(NC 1437)	£442
ATB1909		12 km/h (8mph) (includes LED light system with indicators)	(NC 1437)	£442
BATTERY*				
ATB1364		70 Ah Gel		◆
SUSPENSION*				
ATB1240		Suspension low (up to 80kg user weight)		○
ATB1241		Suspension middle (from 70kg up to 100kg user weight)		○
ATB1242		Suspension hard (from 90kg up to 150kg user weight)		○
G TRAC				
ATB1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)		£703
CHARGER*				
ATB1343		OptiCharge Charger 8A (200 - 250V)		◆
ATB1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster		£129
LIGHTS*				
ATB1437		LED Light system with indicators	(NC 1908, 1909)	£409
ATB1254		Without Light system		◆
ANTI TIPPER*				
ATB1205		Anti-tipper (standard version)		○
KERB RIDER				
ATB1244		Central Kerbrider (with parking position)	(NC 0522)	£181
TRANSPORTATION KIT*				
ATB1450		Preparation for DAHL-Docking station, (Adequate headrest is compulsory, see options below)	(NC 1244)	£349
ATB1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)		◆



FIXATION BELT*				
ATB1329		Lap Belt with Auto Lock	(NC 2701, 2702, 2703, 2704, 2705, 2706)	○
ATB2113		Lap belt with auto lock (flexible mount in side rails)	(NC 2701, 2702, 2703, 2704, 2705, 2706)	◆
ATB1332		H-Belt		£65
ATB2701		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£81
ATB2702		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£81
ATB2703		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£74
ATB2704		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£74
ATB2705		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£74
ATB2706		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£74
REAR VIEW MIRROR				
ATB1415		Rear view Mirror RH - Must be selected with 10 km/h Motor	(NC 2209, 2213, 1287, 2275, 2235, 2241)	£36
ATB1339		Rear view Mirror LH - Must be selected with 10 km/h Motor	(NC 2210, 2214, 1286, 2276, 2236, 2242)	£36
BAG HOOKS				
ATB2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	(NC 2172, 2173, 2227, 2228, 2220, 2219, 2229, 2245, 2257, 2233, 2240)	£66
LUGGAGE RACK				
ATB1394		Luggage Rack		£281
CANE HOLDER				
ATB1312		Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	(NC 1287)	£36
ATB1363		Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	(NC 1286)	£36
TRAY				
ATB1404		Tray Swing Away RH	(NC 1207, 1707, 2213, 2214, 2209, 1287, 2235, 2236, 2241, 2242, 2207, 2275, 2276)	£299
ATB1405		Tray Swing Away LH	(NC 1208, 1708, 2213, 2214, 2210, 1286, 2235, 2236, 2241, 2242, 2208, 2275, 2276)	£299
FENDER				
ATB1416		Fender for front Wheels		£103
ATB1417		Fender for rear Wheels		£103

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ATB2121	Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest)
ATB2122	Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)
ATB2132	410mm (16") (adjustable 410mm - 510mm stepless adjustment)
ATB2136	510mm (20") (adjustable 410mm - 510mm stepless adjustment)
ATB2143	Matrx Flo-tech Contour Visco
ATB2147	Matrx Libra 90mm thick
ATB2163	Backposts (Short)
ATB2108	Backposts WITHOUT push handles, Push handles will be cut down on the vertical section.
ATB2171	Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only
ATB2179	Flex3 Tension Adjustable (540mm or 480mm) Black multi stretch PU coated fabric only
ATB2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ATB2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
ATB2170	Without Backrest Upholstery
ATB2151	Black fabric
ATB2152	Black multi stretch PU coated fabric
ATB2155	Black Spacer fabric
ATB2196	Long Armpad (375mm)
ATB2197	Long Ergonomic Armpad (397mm)
ATB2184	SMALL 140mm x 100mm, for flip up and following armrests only
ATB2186	MEDIUM 230mm x 100mm, for flip up and following armrests only
ATB2188	LARGE 330mm x 100mm, for flip up and following armrests only
ATB2185	SMALL 140mm x 100mm, for flip up and following armrests only
ATB2187	MEDIUM 230mm x 100mm, for flip up and following armrests only
ATB2189	LARGE 330mm x 100mm, for flip up and following armrests only
ATB1350	neckrest: height and depth adjustable
ATB1432	Headrest, height, depth and sideways adjustable
ATB1427	Headrest with cheek support: height and depth adjustable
ATB2323	Elan Standard Pad 150mm (6") wide with standard cover complete with Elan Headrest Hardware
ATB2324	Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware
ATB2322	Elan Standard Pad 355mm (14") wide with standard cover complete with Elan Headrest Hardware
ATB2325	Elan Occipital Pad 230mm x 180mm (9" x 7") with standard cover complete with Elan Headrest Hardware
ATB2326	Elan Occipital Pad 305mm x 200mm (12" x 8") with standard cover complete with Elan Headrest Hardware
ATB2327	Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") with standard cover complete with Elan Headrest Hardware
ATB1496	non Invacare Headrest but profile mount fitting Flex3 hard back
ATB0977	Powered Recliner (0° - 30°) Foldable
ATB0503	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm
ATB0533	Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ATB0532	ADM Manual : Angle adjustable with quick release lever 0° - 80°, length 290mm - 460mm
ATB0531	ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm
ATB0501	Without Hanger
ATB0512	Fixed Footplate
ATB0504	Angle Adjustable Footplate
ATB0513	Angle and depth adjustable Footplate (Adapt)

ATB0522	Angle Adjustable Footboard
ATB2201	REM110 Drive only
ATB2203	REM211 Drive & Actuators
ATB2204	REM216 Drive, Actuators & Lights
ATB2205	REM400, touch-screen Interface, built in Bluetooth HID
ATB2206	REM500, without joystick for alternative access, built in Bluetooth HID
ATB1207	Standard Bracket RH
ATB1208	Standard Bracket LH
ATB1707	Swing Away RH Bracket
ATB1708	Swing Away LH Bracket
ATB2213	Centre mount swing away to RH nucleus Bracket
ATB2214	Centre mount swing away to LH nucleus Bracket
ATB2227	Bracket RH mounted for Attendant
ATB2228	Bracket LH mounted for Attendant
ATB2209	REM500 Swing away RH multi-adjustable Bracket
ATB2210	REM500 Swing away LH multi-adjustable Bracket
ATB2270	Assignable Toggle Switch Module for REM400, (as default set to left On / OFF, right speed dial)
ATB2229	Chin Control (Requires Primary Remote REM400 or REM500) (Only in combination with Headrest 2125) (Please select chin control bracket below)
ATB1287	Chin Control Swing Away Bracket RH Powered
ATB1286	Chin Control Swing Away Bracket LH Powered
ATB2220	REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)
ATB2219	Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)
ATB1207	Dual Control Standard Bracket RH
ATB1208	Dual Control Standard Bracket LH
ATB2245	IDC Dual control (Please select IDC dual control bracket below)
ATB1221	IDC Dual control Detachable Bracket RH
ATB1222	IDC Dual control Detachable Bracket LH
ATB2207	Proportional Compact RH Mount Bracket
ATB2208	Proportional Compact LH Mount Bracket
ATB2213	Proportional Compact Centre Mount swing away to RH (nucleus)
ATB2214	Proportional Compact Centre Mount swing away to LH (nucleus)
ATB2230	Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M4
ATB2247	Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings
ATZ2257	Four-Switch Proximity Array (ASL106)
ATZ2275	Centre Mount, swing-away to RH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)
ATZ2276	Centre Mount, swing-away to LH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)
ATB2233	MEC Proportional Joystick (Micro Extremity Control) (Please select LinX MEC Bracketry)
ATB2235	MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please also select MEC Bracketry)
ATB2236	MEC Tray Mount LH (ASL603, 100mm x 150mm)
ATB2213	Centre Mount, swing away to RH, nucleus
ATB2214	Centre Mount, swing away to LH, nucleus
ATB2241	Centre Mount, swing away to RH, nucleus (Inc Tray ASL603, 100m x 150mm)
ATB2242	Centre Mount, swing away to LH, nucleus (Inc Tray ASL603, 100mm x 150mm)
ATZ2240	Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LinX Heavy Duty Bracketry)

ATZ2235	Tray Mount RH, (ASL603, 100mm x 150mm)
ATZ2236	Tray Mount LH, (ASL603, 100mm x 150mm)
ATZ2213	Centre Mount, swing away to RH, nucleus
ATZ2214	Centre Mount, swing away to LH, nucleus
ATZ2241	Centre Mount, swing away to RH, nucleus (Inc Tray ASL603, 100mm x 150mm)
ATZ2242	Centre Mount, swing away to LH, nucleus (Inc Tray ASL603, 100mm x 150mm)
ATB2227	10-Way Switch Bracket RH, Mounted for Attendant
ATB2228	10-Way Switch Bracket LH, Mounted for Attendant
ATB2253	USB Charger Port DLX-USB02-A
ATB2207	USB Charger Port RH Mount
ATB2208	USB Charger Port LH Mount
ATZ2207	Toggle / Button RH Mount
ATZ2208	Toggle / Button LH Mount
ATB1908	10 km/h (6mph) (includes LED lights system with indicators)
ATB1909	12 km/h (8mph) (includes LED light system with indicators)
ATB1437	LED Light system with indicators
ATB1244	Central Kerbrider (with parking position)
ATB1329	Lap Belt with Auto Lock
ATB2113	Lap belt with auto lock (flexible mount in side rails)
ATB2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ATB2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ATB2703	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ATB2704	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ATB2705	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ATB2706	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ATB1415	Rear view Mirror RH - Must be selected with 10 km/h Motor
ATB1339	Rear view Mirror LH - Must be selected with 10 km/h Motor
ATB2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
ATB1312	Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)
ATB1363	Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)
ATB1404	Tray Swing Away RH
ATB1405	Tray Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):