



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Powerchairs Retail Prescription Form

MARCH 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm4 X-plore Recaro

Standard Features

- Max user weight 150kg (23.5 stone)
- 6 km/h, 10 km/h or 12 km/h
- Recaro seating N Joy Mechanical, Including Headrest
- Seat Width 490mm - 530mm (360mm seat)
- Seat Depth 460mm - 510mm
- Fixed armrests adjustable in height 290mm - 380mm
- Manually settable seat Tilt
- Angle adjustable legrest hangers (VARI F) with Fixed footplates
- LiNX REM110 Control with Standard Bracket
- Choice of 6 colours
- Puncture proof (sealed) drive wheels and front castors
- 73.5 Ah GEL batteries
- 8A (200 - 250V) Charger
- Lapbelt (adjustable both sides)

Standard Base price: £7,664



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ATE0017	Storm 4 X-plore-Recaro	£7,664
ATZ0017	Storm 4 X-plore-RecaroSpecial Basic Model (Please provide details in special section)	POA
SEATING TYPE*		
ATE1152	Recaro Seat Unit	◆
SEAT WIDTH*		
AFT0370	490mm - 530mm	◆
SEAT DEPTH*		
ATE0329	460mm - 510mm	◆

Customer Ref:

Order Number:

SEAT TYPE*			
ATE0906		Seat N-Joy (Mechanical) upholstery fabric (seat depth 130mm front / 70mm back) Always without skirt guard	◆
ATE0927		Seat Mobility (Powered) in combination with True Track Plus: Upholstery fabric (seat depth 130mm front / 70mm back)	£478
ARMRESTS*			
ATE0636		Fixed Armrest 290mm - 380mm	◆
ATE0690		Without Armrests	○
ARMPADS*			
ATE0647		Short Armpad (325mm)	◆
ATE0649		Long Armpad (375mm)	○
HEADREST*			
ATE1355		Recaro Headrest	◆
SEAT HEIGHT AND SEAT RISER*			
ATE1023		Fixed Seat Height 450mm	◆
ATE1025		Fixed Seat Height 480mm	○
ATE1040		Electrical Seat Riser (440mm - 680mm) including adapter (mandatory for X-plore)	£1,101
SEAT TILT*			
ATE1002		Manual Seat Tilt by spindle	◆
ATE1014		Electrical Tilt from 0° to 25°	£514
PREPARED FOR ELECTRICAL ADJUSTMENTS			
ATE1436		Prepared for electrical Legrest & electrical adjustments	○
LEGREST FIXATION POSITION*			
ATE0591		Position for seat cushion over 70mm thick	◆
FOOTREST HANGER TYPE (PAIR)*			
ATE0503		Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	◆
ATE0533		Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	£200
ATE0532		ADM Manual : Angle adjustable with quick release lever 0° - 80°, length 290mm - 460mm	£181
ATE0531		ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm	£571
FOOTPLATE TYPE*			
ATE0512		Fixed Footplate	◆
ATE0504		Angle Adjustable Footplate	£63
ATE0513		Angle and depth adjustable Footplate (Adapt)	£101
CONTROLS*			
ATE2201		REM110 Drive only	◆
ATE2203		REM211 Drive & Actuators	○
ATE2204		REM216 Drive, Actuators & Lights	£130
ATE2205		REM400, touch-screen Interface, built in Bluetooth HID	£290
ATE2206		REM500, without joystick for alternative access, built in Bluetooth HID	£662

CONTROL BRACKETRY*			
AFT1207		Standard Bracket RH	○
ATE1208		Standard Bracket LH	○
ATE1707		Swing Away RH	£113
ATE1708		Swing Away LH	£113
ATE2213		Centre mount swing away to RH Nucleus Bracket	£425
ATE2214		Centre mount swing away to LH Nucleus Bracket	£425
ATE2227		Bracket RH mounted for Attendant	○
ATE2228		Bracket LH mounted for Attendant	○
ATZ2239		Individual mounting option per (special) - Requires additional information / specification	POA
ATB2209		REM500 Swing away RH multi-adjustable Bracket	○
ATB2210		REM500 Swing away LH multi-adjustable Bracket	○
REM400 TOGGLE OPTION			
ATB2270		Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)	£270
DUAL CONTROL			
ATE2220		REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)	£338
ATE2219		Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)	£338
DUAL CONTROL BRACKET			
ATE1207		Dual Control Standard Bracket RH	○
ATE1208		Dual Control Standard Bracket LH	○
INVACARE IDC DUAL CONTROL			
ATE1645		IDC Dual control (Please select IDC dual control bracket below)	£641
INVACARE IDC DUAL CONTROL BRACKET			
ATE1221		IDC Dual control Detachable Bracket RH	○
ATE1222		IDC Dual control Detachable Bracket LH	○
LINX COMPACT REMOTE			
ATE2231		Proportional Compact Remote DLX-CR400 (Please also select Proportional Bracket)	£705
ATE2232		Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry)	£943
LINX COMPACT REMOTE BRACKET			
ATE2207		Proportional Compact RH Mount Bracket	○
ATE2208		Proportional Compact LH Mount Bracket	○
ATE2213		Proportional Compact Centre Mount swing away to RH (Nucleus)	£425
ATE2214		Proportional Compact Centre Mount swing away to LH (Nucleus)	£425
ATE2215		Individual mounting option per request (Special) Requires additional information / specification	POA
CHIN CONTROL			
ATE2229		Chin Control (Requires Primary Remote REM400 or REM500) (Only in combination with Headrest 2125) (Please select chin control bracket below)	£1,702

CHIN CONTROL BRACKET			
ATE1289		Chin Control Swing Away Bracket RH Manual	○
ATE1288		Chin Control Swing Away Bracket LH Manual	○
ATE1287		Chin Control Swing Away Bracket RH Powered	£2,144
ATE1286		Chin Control Swing Away Bracket LH Powered	£2,144
FOUR-SWITCH PROXIMITY ARRAY REMOTE			
ATZ2257		Four-Switch Proximity Array (ASL106)	£5,655
FOUR-SWITCH PROXIMITY ARRAY BRACKETRY			
ATZ2275		Centre Mount, swing-away to RH, Nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)	£425
ATZ2276		Centre Mount, swing-away to LH, Nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)	£425
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	POA
WIRELESS ACCESSORIES FOR ATOM			
ZTE2249		Wireless Triple Switch Receiver (ASL557-3) (Sub D plug)	£290
ATE2250		Wireless Mouse Emulator (ASL558)	£698
LINX HEAVY DUTY REMOTE			
ATZ2240		Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LiNX Heavy Duty Bracketry)	£1,977
LINX HEAVY DUTY REMOTE BRACKET			
ATZ2235		Tray Mount RH, (ASL603, 100mm x 150mm)	○
ATZ2236		Tray Mount LH, (ASL603, 100mm x 150mm)	○
ATZ2213		Centre Mount, swing away to RH, Nucleus	£425
ATZ2214		Centre Mount, swing away to LH, Nucleus	£425
ATZ2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	£425
ATZ2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	£425
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	POA
LINX 10-WAY SWITCH			
ATE2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch)	£431
LINX 10-WAY SWITCH BRACKET			
ATE2217		10-Way Switch Bracket RH, Mounted for User	○
ATZ2218		10-Way Switch Bracket LH, Mounted for User	○
ATZ2227		10-Way Switch Bracket RH, Mounted for Attendant	£54
ATZ2228		10-Way Switch Bracket LH, Mounted for Attendant	£54
LINX ACCESSORY			
ATE2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	£86
FURTHER LINX COMPONENTS / ACCESSORIES			
ATE2300		To be enclosed only	○
ATZ2301		To be connected and programmed (special)	POA

BUDDY BUTTONS		
ATE1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
ATE1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
ATE1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
ATE1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
EGG SWITCH		
ATE2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
PIKO BUTTONS		
ATE1711		Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ATE1712		Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ATE1713		Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ATE1714		Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ATE1715		Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ATE1716		Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
TOGGLES / BUTTONS		
ATZ2261		Stereo Toggle - direct control for Tilt only £230
ATZ2262		Stereo Button - direct control for Tilt only £219
ATZ2263		Four-Way Toggle - direct control for Tilt & Recline £286
ATZ2264		Four-Way Button - direct control for Tilt & Recline £306
ATZ2207		Toggle / Button RH Mount ○
ATZ2208		Toggle / Button LH Mount ○
ATZ2225		Without Mounting ○
BODYPOINT JOYSTICK KNOBS		
ATE2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A £62
ATE2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A £62
ATE2723		Domed Rubber, Joystick Handle, Bodypoint PC107A £62
ATE2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A £62
REMOTE STOP SWITCH		
ATZ2258		Remote Stop Switch (ASL504) £897
SHROUD COLOUR*		
ATE0202		Blue ○
ATE0228		Black ○
ATE0220		Champagne ○
ATE0204		Transparent Cherry Red ○
ATE0222		White ○
ATE0216		Olive green ○
ATE0299		Special Shroud Colour (RAL code required) POA

CASTORS*			
ATE0851	<input type="checkbox"/>	3.00 - 6 Grey Pneumatic	○
ATE0853	<input type="checkbox"/>	3.00 - 6 Grey Puncture Proof sealed	○
ATE0855	<input type="checkbox"/>	3.00 - 6 Grey Puncture Proof (solid inlay)	£92
DRIVE WHEEL*			
ATE0771	<input type="checkbox"/>	3.00 - 8 Grey Pneumatic	○
ATE0773	<input type="checkbox"/>	3.00 - 8 Grey Puncture Proof(sealed)	○
ATE0750	<input type="checkbox"/>	3.00 - 8 Grey Puncture Proof (solid inlay)	£92
MOTOR TYPE*			
ATE1277	<input type="checkbox"/>	6 km/h (4 mph)	○
ATE1908	<input type="checkbox"/>	10 km/h (6mph) (includes LED lights system with indicators)	£442
ATE1909	<input type="checkbox"/>	12 km/h (8mph) (includes LED light system with indicators) (G-Trac recommended)	£442
BATTERY*			
ATE1364	<input type="checkbox"/>	73.5 Ah (GEL)	◆
REAR SUSPENSION*			
ATE1233	<input type="checkbox"/>	Hitec spring X-plore	◆
FRONT SUSPENSION*			
ATE1238	<input type="checkbox"/>	Front suspension (suggested for up to 100kg)	○
ATE1239	<input type="checkbox"/>	Front suspension (suggested for above 100kg)	○
CHARGER*			
ATE1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
AFT1313	<input type="checkbox"/>	OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£130
LIGHTS*			
ATE1437	<input type="checkbox"/>	LED Light system with indicators (included with 10 km/h & 12 km/h motor)	○
ATE1254	<input type="checkbox"/>	Without Light system	○
ANTI-TIPPERS*			
ATE1205	<input type="checkbox"/>	Anti-tipper (standard version)	◆
KERB-RIDER			
ATE1244	<input type="checkbox"/>	Kerb climber assembly	£181
TRANSPORTATION KIT*			
ATE1450	<input type="checkbox"/>	Preparation for DAHL-Docking station, (Adequate headrest is compulsory, see options below)	£349
ATE1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory).	◆
LAP BELT*			
ATE1329	<input type="checkbox"/>	Lapbelt (adjustable both sides) fitted to base of backrest	○
ATE1332	<input type="checkbox"/>	H-Belt	£65
REAR VIEW MIRROR			
ATE1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h & 12 km/h Motor	£36
ATE1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h & 12 km/h Motor	£36
LUGGAGE RACK			
ATE1394	<input type="checkbox"/>	Luggage Rack	£281

CANE HOLDER		
ATE1312		Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote) £36
ATE1363		Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote) £36
FENDER		
ATE1416		Fender for front castors £103
ATE1417		Fender for rear / drive wheels £103
TRAY		
ATE1404		Swing away RH, Select opposite side to user remote £300
ATE1405		Swing away LH ,Select opposite side to user remote £300

TOTAL PRICE	
--------------------	--



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions
Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):