



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

MARCH 2021  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Storm 4 X-plore

### Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO7176-Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable seat Widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat Depth 410mm - 510mm Stepless adjustment (16" - 20")
- Tension Adjustable Backrest
- Standard T-Armrests, detachable with Short Armpads
- Manual settable seat Tilt 0° - 20°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F hangers with Fixed Footplate
- LiNX control REM110 with standard bracket
- Choice of 6 colours
- Puncture proof castors
- Puncture proof Drive Wheels
- 73.5 Ah GEL batteries
- 6 km/h (4mph) Motors
- 8A (200 - 250V) Charger
- 4 Point-tie down system
- Lapbelt with autolock (adjustable both sides)
- Wheel lock / parking brake
- Anti-Tippers
- X-plore suspension 50mm all round

**Standard Base price: £6,458**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ATE0017	Storm 4 X-plore	£6,458
ATZ0017	Special Storm 4 X-plore	POA
<b>SEATING TYPE*</b>		
ATE2100	Invacare Modulite 'Flex3'	◆
<b>SEAT FRAME*</b>		
ATE2111	Telescopic Seat Frame with metal seat plate	◆

<b>SEAT WIDTH*</b>			
ATEXXXX	<input type="checkbox"/>	Seat width 330mm (13")	○
ATE2121	<input type="checkbox"/>	Set at 380mm	○
ATE2122	<input type="checkbox"/>	Set at 430mm	○
ATE2123	<input type="checkbox"/>	Set at 480mm	○
ATE2124	<input type="checkbox"/>	Set at 530mm	○
<b>SEAT DEPTH*</b>			
ATE2132	<input type="checkbox"/>	410mm (16") (NC 2196, 2197)	○
ATE2134	<input type="checkbox"/>	460mm (18")	○
ATE2136	<input type="checkbox"/>	510mm (20")	○
<b>SEAT OPTION</b>			
ATESWO	<input type="checkbox"/>	"H" Frame Interface Kit	£425
<b>SEAT CUSHION*</b>			
ATE2142	<input type="checkbox"/>	Matrx Flo-tech Lite	£107
ATE2143	<input type="checkbox"/>	Matrx Flo-tech Contour Visco	£107
ATE2147	<input type="checkbox"/>	Matrx Libra 90mm thick	£371
ATE2140	<input type="checkbox"/>	Without seat cushion	◆
<b>BACKREST SUPPORT*</b>			
ATE2162	<input type="checkbox"/>	Backposts 540mm (Tall)	○
ATE2163	<input type="checkbox"/>	Backposts 480mm (Short) (NC 2173)	○
<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ATE2107	<input type="checkbox"/>	Backposts WITH push handles	◆
ATE2108	<input type="checkbox"/>	Backposts WITHOUT push handles (NC 1332, 2227, 2228, 2220, 2219, 2245)	£34
<b>BACKREST UPHOLSTREY / CUSHION*</b>			
ATE2171	<input type="checkbox"/>	Basic Tension Adjustable Backrest (NC 1445, 1446, 1496, 2155, 2230, 2247)	◆
ATE2179	<input type="checkbox"/>	Flex3 Tension Adjustable (540mm or 480mm) (NC 2151, 2192, 2193, 2194, 2199, 1446, 2155, 2229, 2230, 2247, 2128)	£132
ATE2172	<input type="checkbox"/>	Flex3 adjustable Backrest Plates (without Backrest Cushion) (NC 2163, 2128)	£213
ATE2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion (NC 2163, 2128)	£346
ATE2174	<input type="checkbox"/>	Matrx Posture Back, height 400mm (Mesh fabric) (NC 1446, 1496, 2124)	£472
ATE2175	<input type="checkbox"/>	Matrx Posture Back, height 500mm (Mesh fabric) (NC 1446, 1496)	£472
ATE2176	<input type="checkbox"/>	Matrx Elite Back, height 500mm (Mesh fabric) (NC 1446, 1496)	£472
ATE2177	<input type="checkbox"/>	Matrx Elite TR Back, height 500mm (Mesh fabric) (NC 1446, 1496)	£472
ATE2178	<input type="checkbox"/>	Matrx Elite Deep Back, height 500mm (Mesh fabric) (NC 1446, 1496, 2128)	£524
ATE2170	<input type="checkbox"/>	Without Backrest Upholstery (NC 1446, 1496, 2230, 2247)	○
<b>UPHOLSTERY COLOUR*</b>			
ATE2151	<input type="checkbox"/>	Black fabric (NC 2170)	○
ATE2152	<input type="checkbox"/>	Black multi stretch PU coated fabric (NC 2170)	○
ATE2155	<input type="checkbox"/>	Black Spacer fabric (NC 2171, 2179, 2143, 2147)	○
ATZ1002	<input type="checkbox"/>	Without backrest upholstery (for option 2170) (NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178)	○
<b>ARMREST*</b>			
ATE2191	<input type="checkbox"/>	T-Armrest Standard	◆
AFT2181	<input type="checkbox"/>	T-Armrest Extended height	£34
ATE2183	<input type="checkbox"/>	In-rigged Armrests	○
ATE2193	<input type="checkbox"/>	Flip-up Armrest	£55
ATE2192	<input type="checkbox"/>	Following Armrest	£158
ATE2194	<input type="checkbox"/>	Following Armrest	£158

<b>ARMPAD*</b>			
ATE2195		Short Armpad (325mm)	◆
ATE2196		Long Armpad (375mm)	○
ATE2197		Long Ergonomic Armpad (397mm)	(NC 2193) ○
ATE2199		Soft Armpad	£82
<b>HIP SUPPORTS LH (Not available with 410mm seat depth and Manual Centre legrests)</b>			
ATE2184		SMALL 140mm x 100mm	£215
ATE2186		MEDIUM 230mm x 100mm	£215
ATE2188		LARGE 330mm x 100mm	£215
<b>HIP SUPPORTS RH (Not available with 410mm seat depth and Manual Centre legrests)</b>			
ATE2185		SMALL 140mm x 100mm	£215
ATE2187		MEDIUM 230mm x 100mm	£215
ATE2189		LARGE 330mm x 100mm	£215
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>			
ATE2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
ATE1350		Neckrest: height and depth Adjustable	(NC 2247) £69
ATE1432		Headrest: height, depth and sideways Adjustable	(NC 2247) £193
ATE1427		Headrest with cheek support: height and depth Adjustable	(NC 2247) £193
ATE2323		Elan Standard Pad 150mm (6")	(NC 2247) £280
ATE2324		Elan Standard Pad 355mm (10")	(NC 2247) £280
ATE2322		Elan Standard Pad 355mm (14")	(NC 2247) £280
ATE2325		Elan Occipital Pad 230mm x 180mm (9" x 7")	(NC 2247) £280
ATE2326		Elan Occipital Pad 305mm x 200mm (12" x 8")	(NC 2247) £280
ATE2327		Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10")	(NC 2247) £377
ATE1496		All Bracketry as option 1432 without Headrest cushion	(NC 2171, 2174, 2175, 2176, 2177, 2178) ○
ATE1495		Push bar	○
<b>LATERALS</b>			
ATE1446		Rigid pad (Azalea type, multi stretch PU coated fabric)	(NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194, 2179) £143
ATE1445		Multi-adjustable pad (Azalea type, multi stretch PU coated fabric)	(NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194) £143
<b>SEAT HEIGHT SETTING &amp; SEAT RISER*</b>			
ATE1023		Fixed Seat Height 450mm	◆
ATE1025		Fixed Seat Height 480mm	○
ATE1040		Electrical Seat Riser (440mm - 690mm)	(NC 2201) £1,101
<b>SEAT TILT*</b>			
ATE1002		Manual Seat Tilt by spindle	○
ATE1014		Electrical Tilt from 0° to 25°	£514
<b>BACKREST RECLINER*</b>			
ATE0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
ATE0977		Powered Recliner (0° - 30°) Foldable	(NC 2201) £431
<b>PREPARED FOR ELECTRIC ADJUSTMENT</b>			
ATE1436		Prepared for electrical Legrests & electrical adjustments	(NC 0531, 2201) ○

FOOTREST HANGER TYPE (PAIR)*			
ATE0503		Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	◆
ATE0533		Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	£200
ATE0532		ADM Manual : Angle adjustable with quick release lever, 0° - 80°, length 290mm - 460mm	£195
ATE0531		ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm	(NC 0522, 1436, 2201) £571
ATE0501		Without Hanger	○
FOOTPLATE TYPE*			
ATE0512		Fixed Footplate	◆
ATE0504		Angle Adjustable Footplate	£63
ATE0513		Angle and depth adjustable Footplate (Adapt)	£101
ATE0522		Angle Adjustable Footboard	(NC 1244) £131
ATZXXXX		Without footplates for option 0501 without hanger	(NC 0503, 0533, 0532, 0531) ○
CONTROLS*			
ATE2201		REM110 Drive only	(NC 1437, 1436, 1040, 1014, 0977, 0531, 2231, 2232, 2229, 2230, 2247, 2233, 2270, 2219, 2245, 1908, 1909) ◆
ATE2203		REM211 Drive & Actuators	(NC 1437, 2231, 2232, 2229, 2230, 2247, 2233, 2270, 1908, 1909) ○
ATE2204		REM216 Drive, Actuators & Lights	(NC 2232, 2229, 2230, 2247, 2233, 2270) £130
ATE2205		REM400, touch-screen Interface, built in Bluetooth HID	£290
ATE2206		REM500, without joystick for alternative access, built in Bluetooth HID	(NC 2270) £662
CONTROL BRACKETRY*			
ATE1207		Standard Bracket RH	(NC 2206) ○
ATE1208		Standard Bracket LH	(NC 2206) ○
ATE1707		Swing Away RH Bracket	(NC 2206) £113
ATE1708		Swing Away LH Bracket	(NC 2206) £113
ATE2213		Centre mount swing away to RH Nucleus Bracket	(NC 1405) £425
ATE2214		Centre mount swing away to LH Nucleus Bracket	(NC 1405) £425
ATE2227		Bracket RH mounted for Attendant	(NC 1404, 1405, 1287, 2128) £54
ATE2228		Bracket LH mounted for Attendant	(NC 1404, 1405, 1286, 2128) £54
ATZ2239		Individual mounting option per (special) - Requires additional information / specification	POA
ATE2209		REM500 Swing away RH multi-adjustable Bracket	(NC 2203, 2204, 2205, 1404) ○
ATE2210		REM500 Swing away LH multi-adjustable Bracket	(NC 2203, 2204, 2205) ○
REM400 TOGGLE OPTION			
ATE2270		Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)	(NC 2203, 2204) £270
DUAL CONTROL			
ATE2220		REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)	(NC 1222, 2227, 2228, 2245, 2128, 2230, 2247, 2257, 2233, 2240) £338
ATE2219		Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)	(NC 1222, 2227, 2228, 2245, 2128, 2201) £338
DUAL CONTROL BRACKET			
ATE1207		Dual Control Standard Bracket RH	○
ATE1208		Dual Control Standard Bracket LH	○
INVACARE IDC DUAL CONTROL			
ATE2245		Invacare IDC Dual control (Please select AQE1221 OR AQE1222 for bracket option below)	£641
INVACARE IDC DUAL CONTROL BRACKET			
ATE1221		IDC Dual Detachable Bracket RH	○
ATE1222		IDC Dual Detachable Bracket LH	○

LINX COMPACT REMOTE				
ATE2231		Proportional Compact Remote DLX-CR400 (Please also select Proportional Bracket)	(NC 2203, 2204)	£705
ATE2232		Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry)	(NC 2203, 2204)	£943
LINX COMPACT REMOTE BRACKET				
ATE2207		Proportional Compact RH Mount Bracket	(NC 2203, 2204)	○
ATE2208		Proportional Compact LH Mount Bracket	(NC 2203, 2204)	○
ATE2213		Proportional Compact Centre Mount swing away to RH (Nucleus)	(NC 2203, 2204)	£425
ATE2214		Proportional Compact Centre Mount swing away to LH (Nucleus)	(NC 2203, 2204, 1404, 1405)	£425
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2203, 2204)	POA
CHIN CONTROL				
ATE2229		Chin Control (Requires Primary Remote REM400 or REM500) (Only in combination with Headrest 2125) (Please select chin control bracket below)	(NC 2203, 2204, 2171, 2179, 2128)	£1,702
CHIN CONTROL BRACKET				
ATE1289		Chin Control Swing Away Bracket RH Manual		○
ATE1288		Chin Control Swing Away Bracket LH Manual		○
ATE1287		Chin Control Swing Away Bracket RH Powered		£2,144
ATE1286		Chin Control Swing Away Bracket LH Powered		£2,144
LINX HEAD CONTROL				
ATE2230		Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M4	(NC 2179, 2197, 2201, 2203, 2204, 2220, 2229, 2233, 2323, 2324, 2322, 2325, 2326, 2327, 1350, 1432, 1427, 2257, 2240)	£3,652
ATE2247		Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings	(NC 2179, 2197, 2201, 2203, 2204, 2220, 2229, 2233, 2323, 2324, 2322, 2325, 2326, 2327, 1350, 1432, 1427, 2257, 2240)	£4,277
FOUR-SWITCH PROXIMITY ARRAY REMOTE				
ATZ2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2220, 2229, 2230, 2247, 2233, 2240, 2128)	£5,655
FOUR-SWITCH PROXIMITY ARRAY BRACKETRY				
ATZ2275		Centre Mount, swing-away to RH, Nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)	(NC 1405, 1415)	£425
ATZ2276		Centre Mount, swing-away to LH, Nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)	(NC 1405, 1339)	£425
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification		POA
WIRELESS ACCESSORIES FOR ATOM				
ATE2249		Wireless Triple Switch Receiver (ASL557-3) (Sub D plug)		£290
ATE2250		Wireless Mouse Emulator (ASL558)		£698
LINX MICRO EXTREMITY CONTROL				
ATE2233		MEC Proportional Joystick (Micro Extremity Control) (Please select LiNX MEC Bracketry)	(NC 2203, 2204, 2220, 2229, 2230, 2247, 2257, 2240, 2128)	£3,546
LINX MICRO EXTREMITY CONTROL BRACKET				
ATE2235		MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please also select MEC Bracketry)	(NC 1415)	○
ATE2236		MEC Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1339)	○
ATE2213		Centre Mount, swing away to RH, Nucleus	(NC 1405, 1415)	£425
ATE2214		Centre Mount, swing away to LH, Nucleus	(NC 1405, 1339)	£425
ATE2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100m x 150mm)	(NC 1405, 1415)	£425
ATE2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1405, 1339)	£425
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification		POA

<b>LINX HEAVY DUTY REMOTE</b>			
ATZ2240		Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LiNX Heavy Duty Bracketry)	(NC 2203, 2204, 2220, 2229, 2230, 2247, 2257, 2233, 2128) £1,977
<b>LINX HEAVY DUTY REMOTE BRACKET</b>			
ATZ2235		Tray Mount RH, (ASL603, 100mm x 150mm)	(NC 1415) ○
ATZ2236		Tray Mount LH, (ASL603, 100mm x 150mm)	(NC 1339) ○
ATZ2213		Centre Mount, swing away to RH, Nucleus	(NC 1405, 1415) £425
ATZ2214		Centre Mount, swing away to LH, Nucleus	(NC 1405, 1339) £425
ATZ2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1405, 1415) £425
ATZ2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1405, 1339) £425
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	POA
<b>LINX 10-WAY SWITCH</b>			
ATE2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch)	£431
<b>LINX 10-WAY SWITCH BRACKET</b>			
ATE2217		10-Way Switch Bracket RH, Mounted for User	(NC 2213, 2275, 2235, 2241) ○
ATE2218		10-Way Switch Bracket LH, Mounted for User	(NC 2214, 2276, 2236, 2242) ○
ATE2227		10-Way Switch Bracket RH, Mounted for Attendant	£54
ATE2228		10-Way Switch Bracket LH, Mounted for Attendant	£54
<b>LINX ACCESSORY</b>			
ATE2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	£86
<b>LINX USB CHARGER PORT</b>			
ATE2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) £136
<b>LINX USB CHARGER PORT BRACKET</b>			
ATE2207		USB Charger Port RH Mount	○
ATE2208		USB Charger Port LH Mount	○
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
ATE2300		To be enclosed only	○
ATE2301		To be connected and programmed (special)	POA
<b>BUDDY BUTTONS</b>			
ATE1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£75
ATE1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£75
ATE1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£75
ATE1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£75
<b>EGG SWITCH</b>			
ATE2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	£75



PIKO BUTTONS		
ATE1711		Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ATE1712		Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ATE1713		Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ATE1714		Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ATE1715		Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ATE1716		Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
TOGGLES / BUTTONS		
ATZ2261		Stereo Toggle - direct control for Tilt only £230
ATZ2262		Stereo Button - direct control for Tilt only £219
ATZ2263		Four-Way Toggle - direct control for Tilt & Recline £286
ATZ2264		Four-Way Button - direct control for Tilt & Recline £306
ATZ2207		Toggle / Button RH Mount ○
ATZ2208		Toggle / Button LH Mount ○
ATZ2225		Without Mounting ○
BODYPOINT JOYSTICK KNOBS		
ATE2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A £62
ATE2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A £62
ATE2723		Domed Rubber, Joystick Handle, Bodypoint PC107A £62
ATE2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A £62
REMOTE STOP SWITCH		
ATZ2258		Remote Stop Switch (ASL504) £897
SHROUD COLOUR*		
ATE0202		Blue ○
ATE0228		Black ○
ATE0220		Champagne ○
ATE0204		Transparent Cherry Red ○
ATE0222		White ○
ATE0216		Olive green ○
ATE0299		Special Shroud Colour (RAL code required) POA
CASTOR WHEEL*		
ATE0851		3.00 - 6 Grey Pneumatic ○
ATE0853		3.00 - 6 Grey Puncture Proof (sealed) ○
ATE0855		3.00 - 6 Grey Puncture Proof (solid inlay) £92
DRIVE WHEEL*		
ATE0771		3.00 - 8 Grey Pneumatic ○
ATE0773		3.00 - 8 Grey Puncture Proof (sealed) ○
ATE0750		3.00 - 8 Grey Puncture Proof (solid inlay) £92
MOTOR TYPE*		
ATE1277		6 km/h (4mph) ○
ATE1908		10 km/h (6mph) includes LED light system with indicators £442
ATE1909		12 km/h (8mph) includes LED light system with indicators £442

<b>BATTERY*</b>			
ATE1364	<input type="checkbox"/>	73.5 Ah GEL	◆
<b>REAR SUSPENSION*</b>			
ATE1233	<input type="checkbox"/>	Hitec spring X-plore	○
<b>FRONT SUSPENSION*</b>			
ATE1238	<input type="checkbox"/>	Front suspension (suggested for up to 100kg)	○
ATE1239	<input type="checkbox"/>	Front suspension (suggested for above 100kg)	○
<b>G TRAC</b>			
ATE1215	<input type="checkbox"/>	G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£703
<b>CHARGER*</b>			
ATE1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
ATE1313	<input type="checkbox"/>	OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£130
<b>LIGHTS*</b>			
ATE1437	<input type="checkbox"/>	LED Light system with indicators	£410
ATZ1006	<input type="checkbox"/>	Lights Included with Motor Package 10 & 12kmh	○
ATE1254	<input type="checkbox"/>	Without Light system (NC 1909)	○
<b>ANTI-TIPPERS*</b>			
ATE1205	<input type="checkbox"/>	Anti-tipper (standard version)	◆
<b>KERB CLIMBER</b>			
ATE1244	<input type="checkbox"/>	Kerb climber assembly	£181
<b>TRANSPORTATION*</b>			
ATE1450	<input type="checkbox"/>	Preparation for DAHL-Docking station, (Adequate headrest is compulsory, see options below)	£349
ATE1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory).	◆
<b>FIXATION BELT*</b>			
ATE1329	<input type="checkbox"/>	Lapbelt (adjustable both sides) fitted to base of backrest	○
ATE2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	◆
ATE1332	<input type="checkbox"/>	H-Belt	£65
ATE2701	<input type="checkbox"/>	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£80
ATE2702	<input type="checkbox"/>	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£80
ATE2703	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£74
ATE2704	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£74
ATE2705	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£74
ATE2706	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£74
<b>REAR VIEW MIRROR</b>			
ATE1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	£36
ATE1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	£36
<b>BAG HOOKS</b>			
ATE2128	<input type="checkbox"/>	Bag Hooks (Max load up to 6kg = impact on max user weight) (NC 2172, 2173, 2227, 2228, 2220, 2219, 2229, 2245, 2257, 2233, 2240)	£66



<b>LUGGAGE RACK</b>		
ATE1394		Luggage Rack £281
<b>CANE HOLDER</b>		
ATE1312		Cane holder RH (In combination with Attendant bracket or Dual control the Can Holder has to be mounted on the opposite side of the remote) £36
ATE1363		Cane holder LH (In combination with Attendant bracket or Dual control the Can Holder has to be mounted on the opposite side of the remote) £36
<b>TRAY</b>		
ATE1404		Swing away RH, Select opposite side to user remote £300
ATE1405		Swing away LH ,Select opposite side to user remote £300
<b>FENDER</b>		
ATE1416		Fender for front Wheels £103
ATE1417		Fender for rear Wheels £103

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ATE2121	Set at 380mm
ATE2122	Set at 430mm
ATE2124	Set at 530mm
ATE2132	410mm (16")
ATE2143	Matrx Flo-tech Contour Visco
ATE2147	Matrx Libra 90mm thick
ATE2163	Backposts 480mm (Short)
ATE2171	Basic Tension Adjustable Backrest
ATE2179	Flex3 Tension Adjustable (540mm or 480mm)
ATE2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ATE2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
ATE2174	Matrx Posture Back, height 400mm (Mesh fabric)
ATE2175	Matrx Posture Back, height 500mm (Mesh fabric)
ATE2176	Matrx Elite Back, height 500mm (Mesh fabric)
ATE2177	Matrx Elite TR Back, height 500mm (Mesh fabric)
ATE2178	Matrx Elite Deep Back, height 500mm (Mesh fabric)
ATE2170	Without Backrest Upholstery
ATE2151	Black fabric
ATE2155	Black Spacer fabric
ATE2193	Flip-up Armrest
ATE2192	Following Armrest
ATE2194	Following Armrest
ATE2196	Long Armpad (375mm)
ATE2197	Long Ergonomic Armpad (397mm)
ATE2199	Soft Armpad
ATE2184	SMALL 140mm x 100mm
ATE2186	MEDIUM 230mm x 100mm
ATE2188	LARGE 330mm x 100mm
ATE2185	SMALL 140mm x 100mm
ATE2187	MEDIUM 230mm x 100mm
ATE2189	LARGE 330mm x 100mm
ATE1350	Neckrest: height and depth Adjustable
ATE1432	Headrest: height, depth and sideways Adjustable
ATE1427	Headrest with cheek support: height and depth Adjustable
ATE2323	Elan Standard Pad 150mm (6")
ATE2324	Elan Standard Pad 355mm (10")
ATE2322	Elan Standard Pad 355mm (14")
ATE2325	Elan Occipital Pad 230mm x 180mm (9" x 7")
ATE2326	Elan Occipital Pad 305mm x 200mm (12" x 8")
ATE2327	Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10")
ATE1496	All Bracketry as option 1432 without Headrest cushion
ATE1446	Rigid pad (Azalea type, multi stretch PU coated fabric)
ATE1445	Multi-adjustable pad (Azalea type, multi stretch PU coated fabric)

ATE1040	Electrical Seat Riser (440mm - 690mm)
ATE1014	Electrical Tilt from 0° to 25°
ATE0977	Powered Recliner (0° - 30°) Foldable
ATE1436	Prepared for electrical Legrests & electrical adjustments
ATE0503	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm
ATE0533	Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ATE0532	ADM Manual : Angle adjustable with quick release lever, 0° - 80°, length 290mm - 460mm
ATE0531	ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm
ATE0522	Angle Adjustable Footboard
ATE2201	REM110 Drive only
ATE2203	REM211 Drive & Actuators
ATE2204	REM216 Drive, Actuators & Lights
ATE2205	REM400, touch-screen Interface, built in Bluetooth HID
ATE2206	REM500, without joystick for alternative access, built in Bluetooth HID
ATE2213	Centre mount swing away to RH Nucleus Bracket
ATE2214	Centre mount swing away to LH Nucleus Bracket
ATE2227	Bracket RH mounted for Attendant
ATE2228	Bracket LH mounted for Attendant
ATE2270	Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)
ATE2220	REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)
ATE2219	Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)
ATE2245	Invacare IDC Dual control (Please select AQE1221 OR AQE1222 for bracket option below)
ATE1222	IDC Dual Detachable Bracket LH
ATE2231	Proportional Compact Remote DLX-CR400 (Please also select Proportional Bracket)
ATE2232	Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry)
ATE2213	Proportional Compact Centre Mount swing away to RH (Nucleus)
ATE2214	Proportional Compact Centre Mount swing away to LH (Nucleus)
ATE2229	Chin Control (Requires Primary Remote REM400 or REM500) (Only in combination with Headrest 2125) (Please select chin control bracket below)
ATE1287	Chin Control Swing Away Bracket RH Powered
ATE1286	Chin Control Swing Away Bracket LH Powered
ATE2230	Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M4
ATE2247	Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings
ATZ2257	Four-Switch Proximity Array (ASL106)
ATZ2275	Centre Mount, swing-away to RH, Nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)
ATZ2276	Centre Mount, swing-away to LH, Nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)
ATE2233	MEC Proportional Joystick (Micro Extremity Control) (Please select LiNX MEC Bracketry)
ATE2235	MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please also select MEC Bracketry)
ATE2236	MEC Tray Mount LH (ASL603, 100mm x 150mm)
ATE2213	Centre Mount, swing away to RH, Nucleus
ATE2214	Centre Mount, swing away to LH, Nucleus
ATE2241	Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100m x 150mm)
ATE2242	Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ATZ2240	Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LiNX Heavy Duty Bracketry)
ATZ2235	Tray Mount RH, (ASL603, 100mm x 150mm)

ATZ2236	Tray Mount LH, (ASL603, 100mm x 150mm)
ATZ2213	Centre Mount, swing away to RH, Nucleus
ATZ2214	Centre Mount, swing away to LH, Nucleus
ATZ2241	Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ATZ2242	Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ATE2227	10-Way Switch Bracket RH, Mounted for Attendant
ATE2228	10-Way Switch Bracket LH, Mounted for Attendant
ATE2253	USB Charger Port DLX-USB02-A
ATE1908	10 km/h (6mph) includes LED light system with indicators
ATE1909	12 km/h (8mph) includes LED light system with indicators
ATE1437	LED Light system with indicators
ATE1244	Kerb climber assembly
ATE1329	Lapbelt (adjustable both sides) fitted to base of backrest
ATE2113	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion
ATE1332	H-Belt
ATE2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ATE2703	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ATE2704	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ATE2705	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ATE2706	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ATE1415	Rear view Mirror RH - Must be selected with 10 km/h Motor
ATE1339	Rear view Mirror LH - Must be selected with 10 km/h Motor
ATE2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
ATE1404	Swing away RH, Select opposite side to user remote
ATE1405	Swing away LH ,Select opposite side to user remote



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):