



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Powerchairs Retail Prescription Form

MARCH 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 Low Rider

Standard Features

- Max user weight 160kg (25.2 Stone)
- Crash tested ISO7176-Part 19
- Compact, manoeuvrable, CWD TDX SP2
- Telescopic adjustable Seat Widths from 330mm - 530mm (13" - 21") with an additional 50mm (2") by arms
- Seat Depths 410mm - 560mm Stepless Adjustment (16" - 22")
- Tension Adjustable backrest
- Standard T-Armrests, detachable with Short Armpad
- Seat lifter with powered Tilt (300mm, 30°) CoG Shift
- Choice of seat height either 420mm or 403mm
- Manual settable Recliner 0° - 30° , in 3.75° steps
- Standard 80° Hanger with fixed footplate
- LiNX Control REM211 with standard bracket
- Selection of 10 Stunning Shroud and Rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Technology
- 6 km/h 4 Pole Motors
- 8A (200 - 25V) Charge
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £7,115



Q'STRAIT
 QLK-150
 Compatible

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
ARQ0071	TDX SP2 Low Rider 73ah Batteries 420mm seat to ground height	○
ARQ0072	TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height	○
ARZ0071	Special TDX SP2 Low Rider 73ah Batteries 420mm seat to ground height	POA
ARZ0072	Special TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height	POA

User Weight*			
ARQ1918	<input type="checkbox"/>	User weight less 100 kg	<input type="radio"/>
ARQ1919	<input type="checkbox"/>	User weight 101 to 125 kg	<input type="radio"/>
ARQ1921	<input type="checkbox"/>	User weight 126 to 136 kg	<input type="radio"/>
ARQ1922	<input type="checkbox"/>	User weight 137 to 150 kg	<input type="radio"/>
ARQ1923	<input type="checkbox"/>	User weight 151 to 160 kg (NC 2192, 2193, 2194, 1908)	<input type="radio"/>
SEATING TYPE*			
ARQ2100	<input type="checkbox"/>	Modulite Seating Unit	£7,115
SEAT FRAME*			
ARQ2111	<input type="checkbox"/>	Telescopic Seat Frame with metal seat plate	◆
SEAT WIDTH*			
ARQ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19)) (NC 2703, 2705)	<input type="radio"/>
ARQ2123	<input type="checkbox"/>	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21")) (NC 2704, 2706)	<input type="radio"/>
ARQ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23")) (NC 2704, 2706)	<input type="radio"/>
SEAT DEPTH RAIL LENGTH*			
ARQ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm)	<input type="radio"/>
ARQ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm) (NC 2134)	<input type="radio"/>
SEAT DEPTH*			
ARQ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm) (NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706)	<input type="radio"/>
ARQ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm)	<input type="radio"/>
ARQ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm)	<input type="radio"/>
ARQ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm)	<input type="radio"/>
SEAT OPTION			
ARQSWO	<input type="checkbox"/>	"H" Frame Interface Kit	£425
SEAT CUSHION*			
ARQ2142	<input type="checkbox"/>	Matrx Flo-tech Lite	£107
ARQ2143	<input type="checkbox"/>	Matrx Flo-tech Contour Visco (NC 2124, 2138, 2155)	£107
ARQ2147	<input type="checkbox"/>	Matrx Libra 90mm thick	£372
ARQ2140	<input type="checkbox"/>	Without seat cushion	<input type="radio"/>
BACKREST SUPPORT*			
ARQ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	<input type="radio"/>
ARQ2163	<input type="checkbox"/>	Backposts 480mm (Short) (NC 2173)	<input type="radio"/>
PUSH HANDLES FOR BACKPOSTS*			
ARQ2107	<input type="checkbox"/>	Backposts WITH push handles	<input type="radio"/>
ARQ2108	<input type="checkbox"/>	Backposts WITHOUT push handles (NC 2220, 2219, 2227, 2228, 2245)	£34

BACKREST UPHOLSTERY / CUSHION*			
ARQ2171		Basic Tension Adjustable Backrest	(NC 1445, 1446, 1496, 2229, 2155, 2230, 2247, 2322, 2323, 2324, 2325, 2326, 2327) ○
ARQ2179		Flex3 Tension Adjustable (540mm or 480mm) Black multi stretch PU coated fabric only	(NC 2155, 2229, 2230, 2247, 2322, 2324, 2325, 2326, 2327) £132
ARQ2172		Flex3 Adjustable Backrest Plates (without Backrest Cushion)	(NC 2128) £213
ARQ2173		Flex3 Adjustable Backrest Plates with Laguna Backrest cushion	(NC 2128) £346
ARQ2174		Matrx Posture Back, height 400mm (Mesh fabric)	(NC 1445, 1496, 2124) £472
ARQ2175		Matrx Posture Back, height 500mm (Mesh fabric)	(NC 1445, 1496, 2124) £472
ARQ2176		Matrx Elite Back, height 500mm (Mesh fabric)	(NC 1445, 1496, 2124) £472
ARQ2177		Matrx Elite TR Back, height 500mm (Mesh fabric)	(NC 1445, 1496) £472
ARQ2178		Matrx Elite Deep Back, height 500mm (Mesh fabric)	(NC 1446, 1496) £524
ARQ2170		Without Backrest Upholstery	(NC 1446, 1496, 2151, 2152) ○
UPHOLSTERY TYPE*			
ARQ2151		Black fabric	(NC 2170) ○
ARQ2152		Black multi stretch PU coated fabric	(NC 2170) ○
ARQ2155		Black Spacer fabric	(NC 2147, 2171, 2179) ○
ARMREST*			
ARQ2191		T-Armrest Standard, Range of adjustment	○
ARQ2181		T-Armrest Extended height	£34
ARQ2183		In-rigged Armrests	£24
ARQ2193		Flip-up Armrest (from 230mm - 360mm)	(NC 1446) £55
ARQ2192		Following Armrest (from 230mm - 300mm)	(NC 1446) £158
ARQ2194		Following Armrest (from 300mm - 360mm)	(NC 1446) £158
ARMPAD*			
ARQ2195		Short Armpad (325mm)	○
ARQ2196		Long Armpad (375mm)	○
ARQ2197		Long Ergonomic Armpad (397mm)	(NC 2230, 2247) ○
ARQ2199		Soft Arpad (REA Azalea Wide Pad)	£82
HIP SUPPORTS LH			
ARQ2184		SMALL 140mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £215
ARQ2186		MEDIUM 230mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £215
ARQ2188		LARGE 330mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £215
HIP SUPPORTS RH			
ARQ2185		SMALL 140mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £215
ARQ2187		MEDIUM 230mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £215
ARQ2189		LARGE 330mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £215
MOUNTING HARDWARE FOR HIP SUPPORTS			
ARQ2165		Mounting hardware for hip supports (8.1" / 206mm high)	○

HEADREST / NECKREST*			
ARQ1350		Neckrest height and depth adjustable, Max User weight 136kg (NC 2247)	£69
ARQ1432		Headrest, height, depth and sideways adjustable, Max User weight 136kg (NC 2247)	○
ARQ1427		Headrest with cheek support, height and depth adjustable, Max User weight 136kg (NC 2247)	£193
ARQ2323		Elan Standard Pad 150mm (6") wide (NC 2247)	£280
ARQ2324		Elan Standard Pad 355mm (10") wide (NC 2247)	£280
ARQ2322		Elan Standard Pad 355mm (14") wide (NC 2247)	£280
ARQ2325		Elan Occipital Pad 230mm x 180mm (9" x 7") (NC 2247)	£280
ARQ2326		Elan Occipital Pad 305mm x 200mm (12" x 8") (NC 2247)	£280
ARQ2327		Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") (NC 2247)	£377
ARQ1496		All Bracketry as above option without Headrest cushion (NC 2171)	○
ARQ1495		Push bar dia. 22mm	○
LATERAL SUPPORT			
ARQ1446		Rigid pad (Azalea type, multi stretch PU coated fabric) (NC 2171)	£143
ARQ1445		Multi-adjustable pad (Azalea type, multi stretch PU coated fabric) (NC 2171)	£143
SEAT LIFTER / SEAT TILT*			
ARQ1094		Seat lifter with Electrical Tilt	○
SEAT TILT OPTION			
ARQ2054		5° Fixed Anterior Tilt Assist mount	£344
BACKREST RECLINER*			
ARQ0975		Manual Settable Recliner	○
ARQ0977		Powered Recliner (0° - 30°) Foldable	£431
FOOTREST HANGER TYPE*			
ARQ0556		Centre Mounted: Manual : Hanger, centre mounted (NC 0504, 0513)	£318
ARQ0589		Centre Mounted : LNX Powered: Centre mounted legrest (NC 0504, 0513)	£1,400
ARQ0526		Side Mounted : Standard 80° hangers	○
ARQ0503		Side Mounted : Vari F : Angle adjustable	£102
ARQ0533		Side Mounted : Vari A : Angle adjustable (NC 0559)	£200
ARQ0532		Side Mounted : ADM Manual (NC 0559)	£189
ARQ0531		Side Mounted : ADE Powered (NC 0559)	£571
ARQ0501		Without Hanger (NC 0504, 0513, 0559)	○
FOOTPLATE TYPE*			
ARQ0512		Fixed Footplate (NC 0532, 0531, 0556, 0589, 0501)	○
ARQ0504		Angle Adjustable Footplate (NC 0589, 0501)	£61
ARQ0513		Angle and depth adjustable Footplate (Adapt) (NC 0589, 0501)	£101
ARQ0559		Angle Adjustable Footplate for Centre mounted legrest (NC 0503, 0533, 0532, 0531, 0501)	○
ARQXXXX		Without footplates f (NC 0589, 0526, 0503, 0533, 0532, 0531)	POA
CONTROLS*			
ARQ2203		REM211 Drive & Actuators (NC 2209, 2210, 2229, 2270)	○
ARQ2204		REM216 Drive, Actuators & Lights (NC 2210, 2229, 2270)	£130
ARQ2205		REM400, touch-screen Interface, built in Bluetooth HID (NC 2210)	£290
ARQ2206		REM500 (NC 2228, 2270)	£662

CONTROL BRACKETRY*				
ARQ1207		Standard Bracket RH	(NC 2206)	○
ARQ1208		Standard Bracket LH	(NC 2206)	○
ARQ1707		Swing Away RH Bracket		£113
ARQ1708		Swing Away LH Bracket	(NC 2206)	£113
ARQ2213		Centre mount swing away to RH Nucleus		£425
ARQ2214		Centre mount swing away to LH Nucleus	(NC 1405)	£425
ARQ2227		Bracket RH mounted for Attendant	(NC 1287)	£54
ARQ2228		Bracket LH mounted for Attendant	(NC 1286)	£54
ARZ2239		Individual mounting option per (special)		POA
ARQ2209		REM500 Swing away RH multi-adjustable	(NC 2203, 2204, 2205)	○
ARQ2210		REM500 Swing away LH multi-adjustable	(NC 2203, 2204, 2205)	○
REM400 TOGGLE OPTION				
ARQ2270		Assignable Toggle Switch Module for REM400,	(NC 2203, 2204)	£270
CHIN CONTROL				
ARQ2229		Proportional Chin Joystick	(NC 2203, 2204, 2171, 2179, 2128)	£1,702
CHIN CONTROL BRACKETRY				
ARQ1289		Chin Control Swing Away Bracket RH Manual	(NC 2203, 2204, 2171, 2179, 2322)	○
ARQ1288		Chin Control Swing Away Bracket LH Manual	(NC 2203, 2204, 2171, 2179, 2322)	○
ARQ1287		Chin Control Swing Away Bracket RH Powered	(NC 2203, 2204, 2171, 2179, 2227, 1404, 1415, 1312)	£2,144
ARQ1286		Chin Control Swing Away Bracket LH Powered	(NC 2203, 2204, 2171, 2179, 2228, 1405, 1339, 1363)	£2,144
DUAL CONTROL				
ARQ2220		REM050 Dual Control	(NC 1222, 2227, 2228, 2245, 2128)	£242
ARQ2219		Compact Dual Control (DLX-ACU200)	(NC 1222, 2227, 2228, 2245, 2128)	£338
DUAL CONTROL BRACKETRY				
ARQ1207		Dual Control Standard Bracket RH		○
ARQ1208		Dual Control Standard Bracket LH		○
INVACARE LiNX IDC DUAL CONTROL				
ARQ2245		IDC Dual control	(NC 2220, 2227, 2228)	£641
INVACARE LiNX IDC DUAL CONTROL BRACKETRY				
ARQ1221		IDC Dual Control Detachable Bracket RH	(NC 2220, 2227, 2228, 2128)	○
ARQ1222		IDC Dual Control Detachable Bracket LH	(NC 2220, 2227, 2228)	○
LINX COMPACT REMOTE				
ARQ2231		Proportional Compact Remote DLX-CR400	(NC 2203, 2204)	£702
ARQ2232		Proportional Compact Remote Low Force DLX-CR400LF	(NC 2203, 2204)	£943
LINX COMPACT REMOTE BRACKET				
ARQ2207		Proportional Compact RH Mount Bracket	(NC 2203, 2204)	○
ARQ2208		Proportional Compact LH Mount Bracket	(NC 2203, 2204)	○
ARQ2213		Proportional Compact Centre Mount swing away to RH (Nucleus)	(NC 2203, 2204)	£425
ARQ2214		Proportional Compact Centre Mount swing away to LH (Nucleus)	(NC 2203, 2204)	£425
ARZ2215		Individual mounting option per request (Special)	(NC 2203, 2204)	POA
LINX HEAD CONTROL				
ARQ2230		Head Control with ultrasonic proximity sensors	(NC 2203, 2204, 2323, 2324, 2322, 2325, 2326, 2327)	£3,652
ARQ2247		Head Control with ultrasonic proximity sensors	(NC 2203, 2204, 2323, 2324, 2322, 2325, 2326, 2327)	£4,277
LINX MICRO EXTREMITY CONTROL				
ARZ2233		MEC Proportional Joystick	(NC 2203, 2204)	£3,546

LINX MICRO EXTREMITY CONTROL BRACKET				
ARZ2235		MEC Tray Mount RH	(NC 2203, 2204)	○
ARZ2236		MEC Tray Mount LH	(NC 2203, 2204)	○
ARZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	£425
ARZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	£425
ARZ2241		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	£425
ARZ2242		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	£425
ARZ2215		Individual mounting option per request (Special)	(NC 2203, 2204)	POA
LINX HEAVY DUTY REMOTE				
ARZ2240		Proportional Heavy Duty Remote Joystick with mode switch	(NC 2203, 2204)	£1,977
LINX HEAVY DUTY REMOTE BRACKET				
ARZ2235		Tray Mount RH, (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2236		Tray Mount LH, (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	£425
ARZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	£425
ARZ2241		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	£425
ARZ2242		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	£425
ARZ2215		Individual mounting option per request (Special)	(NC 2203, 2204)	POA
LINX 10-WAY SWITCH				
ARQ2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		£431
LINX 10-WAY SWITCH BRACKET				
ARQ2217		10-Way Switch Bracket RH, Mounted for User		○
ARQ2218		10-Way Switch Bracket LH, Mounted for User		○
ARQ2227		10-Way Switch Bracket RH, Mounted for Attendant		£54
ARQ2228		10-Way Switch Bracket LH, Mounted for Attendant		£54
LINX ACCESSORY				
ARQ2252		Loom providing jack sockets for Drive and Mode		£86
LINX USB CHARGER PORT				
ARQ2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706)	£136
LINX USB CHARGER PORT BRACKET				
ARQ2207		USB Charger Port RH Mount		○
ARQ2208		USB Charger Port LH Mount		○
FURTHER LINX COMPONENTS / ACCESSORIES				
ARQ2300		To be enclosed only		○
ARZ2301		To be connected and programmed (special)		POA
BUDDY BUTTONS				
ARQ1685		Red Buddy Button		£75
ARQ1686		Yellow Buddy Button		£75
ARQ1687		Green Buddy Button		£75
ARQ1690		Blue Buddy Button		£75
EGG SWITCH				
ARQ2254		Egg Switch, black		£75

PIKO BUTTONS			
ARQ1711		Red Piko Button	£135
ARQ1712		Yellow Piko Button	£135
ARQ1713		Green Piko Button	£135
ARQ1714		Orange Piko Button	£135
ARQ1715		Blue Piko Button	£135
ARQ1716		Black Piko Button	£135
TOGGLES / BUTTONS			
ARZ2261		Stereo Toggle - direct control for Tilt only	£230
ARZ2262		Stereo Button - direct control for Tilt only	£219
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	£286
ARZ2264		Four-Way Button - direct control for Tilt & Recline	£306
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
ARQ2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£62
ARQ2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£62
ARQ2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£62
ARQ2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£62
REMOTE STOP SWITCH			
ARZ2258		Remote Stop Switch (ASL504)	£897
SHROUD COLOUR*			
ARQ0228		Wet Black	○
ARQ2601		SpaceStation Silver	○
ARQ2602		Island Blue	○
ARQ2603		Grape JellyBean (Purple)	○
ARQ2604		Monster Green	○
ARQ2605		Tangy Orange	○
ARQ2606		WhiteOut	○
ARQ2607		Rockstar Red	○
ARQ2608		Lights Out Black - Matt	○
ARQ2609		Invacare Blue	○
RIM COLOUR*			
ARQ2611		SpaceStation Silver	○
ARQ2612		Island Blue	○
ARQ2613		Grape JellyBean (Purple)	○
ARQ2614		Monster Green	○
ARQ2615		Tangy Orange	○
ARQ2616		WhiteOut	○
ARQ2617		Rockstar Red	○
ARQ2618		Lights Out Black - Matt	○
ARQ2619		Invacare Blue	○
ARQ2620		Wet Black	○
ARQ2621		Without Rim insert	○

FORKS*			
ARQ2631		Single sided fork (New) includes Black castor	£301
ARQ2632		Double sided fork	○
CENTRE WHEELS (DRIVING WHEELS)*			
ARQ0704		355mm x 76mm (14" x 3") Black Tyres, Pneumatic	○
ARQ0733		355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	○
MOTOR TYPE*			
ARQ1294		3 km/h (2mph)	○
ARQ1282		6 km/h (4mph) (4 pole)	○
ARQ1284		8 km/h (5mph) (4 pole)	£442
ARQ1908		10 km/h (6mph) (includes light system and indicators) G-Trac recommended	(NC 2203, 1437, 1254) £442
ARQ1909		12 km/h (8mph) (includes light system and indicators) G-Trac recommended	(NC 2203, 1437, 1254) £442
G TRACK			
ARQ1215		G-Trac	£703
CHARGER*			
ARQ1343		OptiCharge Charger 8A (200 - 250V)	○
ARQ1313		OptiCharge Charger 10A (200 - 250V)	£130
LIGHTS*			
ARQ1437		Light system with Indicators & LED Headlights	£410
ARQ1254		Without Light system	○
ARQ1006		Lights Included with Motor Package 10 & 12kmh	(NC 1282) ○
TRANSPORTATION KIT			
ARQ1450		Preparation for DAHL Docking station:	£499
ARQ1494		4-Point tie-down	○
FIXATION BELT*			
ARQ1329		Lapbelt (adjustable both sides)	(NC 2702, 2703, 2704, 2705, 2706) ○
ARQ2113		Lapbelt (adjustable both sides)	(NC 2702, 2703, 2704, 2705, 2706) ○
ARQ1332		H-Belt	£63
ARQ2701		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £80
ARQ2702		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £80
ARQ2703		Bodypoint 2-Point Hip Belt, Centre pull - push button	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £74
ARQ2704		Bodypoint 2-Point Hip Belt, Centre pull - push button	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £74
ARQ2705		Bodypoint 2-Point Hip Belt, Dual pull Side release	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £74
ARQ2706		Bodypoint 2-Point Hip Belt, Dual pull Side release	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £74
REAR-VIEW MIRROR			
ARQ1415		Rear View Mirror RH	(NC 2213, 1287) £36
ARQ1339		Rear View Mirror LH	(NC 2214, 1286) £36
BAG HOOKS			
ARQ2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	(NC 2173, 2227, 2228, 2220, 2219, 2229, 2245, 2233, 2240) £66
CANE HOLDER			
ARQ1312		Cane holder RH	£36
ARQ1363		Cane holder LH	£36

TRAY				
------	--	--	--	--

ARQ1404		Tray Swing Away RH	(NC 1707, 2213, 2214, 2209, 1287)	£289
ARQ1405		Tray Swing Away LH	(NC 1708, 2213, 2214, 2210, 1286)	£289

TOTAL PRICE				
--------------------	--	--	--	--

Non Conforming Options:

ARQ2122	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19))
ARQ2123	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21"))
ARQ2124	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23"))
ARQ2132	410mm (16") (410mm - 510mm)
ARQ2134	460mm (18") (410mm - 510mm)
ARQ2138	560mm (22") (510mm - 560mm)
ARQ2147	Matrx Libra 90mm thick
ARQ2171	Basic Tension Adjustable Backrest
ARQ2179	Flex3 Tension Adjustable (540mm or 480mm) Black multi stretch PU coated fabric only
ARQ2173	Flex3 Adjustable Backrest Plates with Laguna Backrest cushion
ARQ2170	Without Backrest Upholstery
ARQ2151	Black fabric
ARQ2152	Black multi stretch PU coated fabric
ARQ2155	Black Spacer fabric
ARQ2193	Flip-up Armrest (from 230mm - 360mm)
ARQ2192	Following Armrest (from 230mm - 300mm)
ARQ2194	Following Armrest (from 300mm - 360mm)
ARQ2196	Long Armpad (375mm)
ARQ2197	Long Ergonomic Armpad (397mm)
ARQ2184	SMALL 140mm x 100mm
ARQ2186	MEDIUM 230mm x 100mm
ARQ2188	LARGE 330mm x 100mm
ARQ2185	SMALL 140mm x 100mm
ARQ2187	MEDIUM 230mm x 100mm
ARQ2189	LARGE 330mm x 100mm
ARQ2323	Elan Standard Pad 150mm (6") wide
ARQ2324	Elan Standard Pad 355mm (10") wide
ARQ2322	Elan Standard Pad 355mm (14") wide
ARQ2325	Elan Occipital Pad 230mm x 180mm (9" x 7")
ARQ2326	Elan Occipital Pad 305mm x 200mm (12" x 8")
ARQ2327	Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10")
ARQ1496	All Bracketry as above option without Headrest cushion
ARQ1446	Rigid pad (Azalea type, multi stretch PU coated fabric)
ARQ1445	Multi-adjustable pad (Azalea type, multi stretch PU coated fabric)
ARQ0556	Centre Mounted: Manual : Hanger, centre mounted
ARQ0589	Centre Mounted : LNX Powered: Centre mounted legrest
ARQ0526	Side Mounted : Standard 80° hangers
ARQ0503	Side Mounted : Vari F : Angle adjustable
ARQ0533	Side Mounted : Vari A : Angle adjustable
ARQ0532	Side Mounted : ADM Manual
ARQ0531	Side Mounted : ADE Powered
ARQ0501	Without Hanger
ARQ0504	Angle Adjustable Footplate

ARQ0513	Angle and depth adjustable Footplate (Adapt)
ARQ0559	Angle Adjustable Footplate for Centre mounted legrest
ARQ2203	REM211 Drive & Actuators
ARQ2204	REM216 Drive, Actuators & Lights
ARQ2205	REM400, touch-screen Interface, built in Bluetooth HID
ARQ2206	REM500
ARQ1707	Swing Away RH Bracket
ARQ1708	Swing Away LH Bracket
ARQ2213	Centre mount swing away to RH Nucleus
ARQ2214	Centre mount swing away to LH Nucleus
ARQ2227	Bracket RH mounted for Attendant
ARQ2228	Bracket LH mounted for Attendant
ARQ2209	REM500 Swing away RH multi-adjustable
ARQ2210	REM500 Swing away LH multi-adjustable
ARQ2270	Assignable Toggle Switch Module for REM400,
ARQ2229	Proportional Chin Joystick
ARQ1287	Chin Control Swing Away Bracket RH Powered
ARQ1286	Chin Control Swing Away Bracket LH Powered
ARQ2220	REM050 Dual Control
ARQ2219	Compact Dual Control (DLX-ACU200)
ARQ2245	IDC Dual control
ARQ1222	IDC Dual Control Detachable Bracket LH
ARQ2213	Proportional Compact Centre Mount swing away to RH (Nucleus)
ARQ2214	Proportional Compact Centre Mount swing away to LH (Nucleus)
ARQ2230	Head Control with ultrasonic proximity sensors
ARQ2247	Head Control with ultrasonic proximity sensors
ARZ2233	MEC Proportional Joystick
ARZ2213	Centre Mount, swing away to RH, Nucleus
ARZ2214	Centre Mount, swing away to LH, Nucleus
ARZ2240	Proportional Heavy Duty Remote Joystick with mode switch
ARZ2213	Centre Mount, swing away to RH, Nucleus
ARZ2214	Centre Mount, swing away to LH, Nucleus
ARQ2227	10-Way Switch Bracket RH, Mounted for Attendant
ARQ2228	10-Way Switch Bracket LH, Mounted for Attendant
ARQ2253	USB Charger Port DLX-USB02-A
ARQ1282	6 km/h (4mph) (4 pole)
ARQ1908	10 km/h (6mph) (includes light system and indicators) G-Trac recommended
ARQ1437	Light system with Indicators & LED Headlights
ARQ1254	Without Light system
ARQ1329	Lapbelt (adjustable both sides)
ARQ2113	Lapbelt (adjustable both sides)
ARQ2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARQ2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARQ2703	Bodypoint 2-Point Hip Belt, Centre pull - push button

ARQ2704	Bodypoint 2-Point Hip Belt, Centre pull - push button
ARQ2705	Bodypoint 2-Point Hip Belt, Dual pull Side release
ARQ2706	Bodypoint 2-Point Hip Belt, Dual pull Side release
ARQ1415	Rear View Mirror RH
ARQ1339	Rear View Mirror LH
ARQ2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
ARQ1312	Cane holder RH
ARQ1363	Cane holder LH
ARQ1404	Tray Swing Away RH
ARQ1405	Tray Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions
Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):