



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MARCH 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 Max

Standard Features



- Max user weight 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Compact, manoeuvrable, CWD TDX SP2
- Seat Widths from 490mm - 610mm
- Seat Depths 560mm - 630mm
- Tension Adjustable backrest
- Armrest with reinforced depth adjustable bracket with Short armpad
- Manual settable seat Tilt
- Manual settable backrest Recline
- VARI F hangers with fixed footplate
- LiNX Control REM110 with standard bracket
- Selection of 10 Stunning Shroud and Rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3")
Black Tyres
- Suspension and Sure Step Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 73.5 Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £5,225



Q'STRAIT
QLK-150
Compatible

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ARH0000	TDX SP2 Max	£5,225
ARZ0000	Special TDX SP2 Max (Please provide details in special section)	POA
SEATING TYPE*		
ARH1145	MAX Seating system	◆

SEAT WIDTH SETTING*			
ARH0348	<input type="checkbox"/>	Seat Width of Seat Frame 490mm (+100mm adjustment in arms)	<input type="radio"/>
ARH0378	<input type="checkbox"/>	Seat Width of Seat Frame 550mm (+100mm adjustment in arms)	<input type="radio"/>
ARH0380	<input type="checkbox"/>	Seat Width of Seat Frame 610mm (+100mm adjustment in arms) (NC 0533)	<input type="radio"/>
SEAT DEPTH SETTING SETTABLE IN 3 STEPS*			
ARH0396	<input type="checkbox"/>	Seat Depth 560mm	<input type="radio"/>
ARH0397	<input type="checkbox"/>	Seat Depth 595mm	<input type="radio"/>
ARH0398	<input type="checkbox"/>	Seat Depth 630mm	<input type="radio"/>
SEAT CUSHION*			
ARH1484	<input type="checkbox"/>	MAX Seat Cushion	<input type="radio"/>
ARH1114	<input type="checkbox"/>	Without Seat Cushion (10mm seat board only)	<input type="radio"/>
ADAPTOR TO REDUCE THE SEAT DEPTH BY 130MM			
ARH0399	<input type="checkbox"/>	Adaptor to shorten the seat depth by 130mm (Required if seat width 490mm is selected in combination with seat depth 595mm or 630mm)	<input type="radio"/>
BACKREST UPHOLSTERY*			
ARH1106	<input type="checkbox"/>	Tension Adjustable	<input checked="" type="checkbox"/>
UPHOLSTERY COLOUR*			
ARH0278	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input checked="" type="checkbox"/>
ARMREST INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*			
ARH0635	<input type="checkbox"/>	Fixed Armrest, 250mm - 340mm tool adjustable in height, detachable	<input checked="" type="checkbox"/>
ARH0636	<input type="checkbox"/>	Fixed Armrest, 290mm - 380mm tool adjustable in height, detachable	<input type="radio"/>
ARMPADS*			
ARH0647	<input type="checkbox"/>	Short armpad (325mm)	<input checked="" type="checkbox"/>
ARH0649	<input type="checkbox"/>	Long armpad (375mm)	<input type="radio"/>
PUSH BAR*			
ARH1307	<input type="checkbox"/>	Push bar	<input checked="" type="checkbox"/>
HEADREST - Adequate headrest option must be selected*			
ARH1432	<input type="checkbox"/>	Headrest: height, depth and sideways adjustable (Azalea type)	<input checked="" type="checkbox"/>
ARH1427	<input type="checkbox"/>	Headrest with cheek support: height and depth adjustable (Azalea type)	£193
ARH1495	<input type="checkbox"/>	Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware.	<input type="radio"/>
SEAT TILT*			
ARH1091	<input type="checkbox"/>	Manually settable Seat Tilt (Max user weight 180kg)	<input checked="" type="checkbox"/>
ARH1095	<input type="checkbox"/>	Electric heavy duty basic Tilt (0° - 20°) select either 470mm or 495mm seat to ground height. (Max user weight 180kg)	£373
SEAT HEIGHT SETTING*			
ARH1024	<input type="checkbox"/>	Fixed seat height 470mm (18.5")	<input type="radio"/>
ARH1046	<input type="checkbox"/>	Fixed seat height 495mm (19.5")	<input type="radio"/>
BACKREST RECLINER*			
ARH0958	<input type="checkbox"/>	Manually settable recliner by spindle (0° - 25°)	<input checked="" type="checkbox"/>
ARH0947	<input type="checkbox"/>	Electrical recliner (0° - 25°)	£431

FOOTREST HANGER TYPE*			
ARH0556		Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate)	(NC 0504, 0513, 0599, 0522) £318
ARH0589		Centre Mounted : LNX Powered: Centre mounted legrest, length compensating (including calf pad and angle adjustable footplate)	(NC 0504, 0513, 0599, 0522) £1,400
ARH0564		Reinforced hangers, fixed angle (not detachable without tools) incl Calf pads	○
ARH0503		Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	◆
ARH0533		Side Mounted : Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	(NC 0559) £200
ARH0532		Side Mounted : ADM Manual : Angle adjustable with quick release lever, 0° - 80°, length 290mm - 460mm	(NC 0559) £189
ARH0531		Side Mounted : ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm	(NC 0513, 0599, 0559) £571
FOOTPLATE TYPE*			
ARH0512		Fixed Footplate	(NC 0533, 0532, 0556, 0589) ○
ARH0504		Angle Adjustable Footplate	(NC 0589) £63
ARH0513		Angle and depth adjustable Footplate (Adapt)	(NC 0589) £101
ARH0599		Angle adjustable flip-back, heavy-duty footboard (Recommended for SW 550mm)	(NC 0531, 0532, 0533, 0556, 0589) £285
ARH0522		Angle Adjustable Footboard (for max user weight up to 180 kg) (Recommended for SW 550 mm)	(NC 0531, 0532, 0556, 0589) £192
ARH0559		Angle Adjustable Footplate for Centre mounted legrest	(NC 0503, 0533, 0532, 0531) ○
CONTROLS*			
ARH2201		REM110 Drive only	(NC 1437, 1095, 0947, 0589, 0531, 2219, 2245) ◆
ARH2203		REM211 Drive & Actuators	(NC 1437) ○
ARH2204		REM216 Drive, Actuators & Lights	£130
CONTROL BRACKETRY*			
ARH1207		Standard Bracket RH	○
ARH1208		Standard Bracket LH	○
ARH1707		Swing Away RH Bracket	£113
ARH1708		Swing Away LH Bracket	£113
ARH2227		Bracket RH mounted for Attendant	£54
ARH2228		Bracket LH mounted for Attendant	£54
ARH1289		Tray / Centre control Swing Away Bracket RH	(NC 2203, 2204) ○
ARH1288		Tray / Centre control Swing Away Bracket LH	(NC 2203, 2204) ○
DUAL CONTROL			
ARH2220		REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)	(NC 1221, 1222) £338
ARH2219		Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)	(NC 1222, 2201) £338
DUAL CONTROL BRACKETRY			
ARH1207		Dual Control Standard Bracket RH	○
ARH1208		Dual Control Standard Bracket LH	○
INVACARE LINX IDC DUAL CONTROL			
ARH2245		Invacare IDC Dual control (Please select IDC Dual control bracket below)	(NC 2219, 2227, 2228) £641
INVACARE LINX IDC DUAL CONTROL BRACKETRY			
ARH1221		IDC Dual control Detachable Bracket RH	(NC 2219, 2227, 2228) ○
ARH1222		IDC Dual control Detachable Bracket LH	(NC 2219, 2227, 2228) ○

LINX 10-WAY SWITCH		
ARH2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch) (Please select 10 way switch bracket below) £431
LINX 10-WAY SWITCH BRACKET		
ARH2217		10-Way Switch Bracket RH, Mounted for User ○
ARH2218		10-Way Switch Bracket LH, Mounted for User ○
ARH2227		10-Way Switch Bracket RH, Mounted for Attendant £54
ARH2228		10-Way Switch Bracket LH, Mounted for Attendant £54
LINX ACCESSORY		
ARH2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch) £86
FURTHER LINX COMPONENTS / ACCESSORIES		
ARH2300		To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected) ○
ARZ2301		To be connected and programmed (special) POA
BUDDY BUTTONS		
ARH1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
ARH1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
ARH1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
ARH1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
EGG SWITCH		
ARH2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
PIKO BUTTONS		
ARH1711		Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ARH1712		Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ARH1713		Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ARH1714		Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ARH1715		Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
ARH1716		Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
TOGGLES / BUTTONS		
ARZ2261		Stereo Toggle - direct control for Tilt only £230
ARZ2262		Stereo Button - direct control for Tilt only £219
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline £286
ARZ2264		Four-Way Button - direct control for Tilt & Recline £306
ARZ2207		Toggle / Button RH Mount ○
ARZ2208		Toggle / Button LH Mount ○
ARZ2225		Without Mounting ○

BODYPOINT JOYSTICK KNOBS			
ARH2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£62
ARH2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£62
ARH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£62
ARH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£62
REMOTE STOP SWITCH			
ARZ2258		Remote Stop Switch (ASL504)	£897
SHROUD COLOUR*			
ARH0228		Wet Black	○
ARH2601		SpaceStation Silver	○
ARH2602		Island Blue	○
ARH2603		Grape JellyBean (Purple)	○
ARH2604		Monster Green	○
ARH2605		Tangy Orange	○
ARH2606		WhiteOut	○
ARH2607		Rockstar Red	○
ARH2608		Lights Out Black - Matt	○
ARH2609		Invacare Blue	○
RIM COLOUR*			
ARH2611		SpaceStation Silver	○
ARH2612		Island Blue	○
ARH2613		Grape JellyBean (Purple)	○
ARH2614		Monster Green	○
ARH2615		Tangy Orange	○
ARH2616		WhiteOut	○
ARH2617		Rockstar Red	○
ARH2618		Lights Out Black - Matt	○
ARH2619		Invacare Blue	○
ARH2620		Wet Black	○
ARH2621		Without Rim insert	○
FORKS*			
ARH2631		Single sided fork (New) includes Black castor	£301
ARH2632		Double sided fork	○
CENTRE WHEELS (DRIVING WHEELS)*			
ARH0704		355mm x 76mm (14" x 3") Black Tyres, Pneumatic	○
ARH0733		355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	○
MOTOR TYPE*			
ARH1294		3 km/h (2mph)	○
ARH1282		6 km/h (4mph) (4 pole)	◆
ARH1284		8 km/h (5mph) (4 pole)	£442
BATTERY*			
ARH1260		73.5 Ah	◆

G TRAC		
ARH1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain) £703
CHARGER*		
ARH1343		OptiCharge Charger 8A (200 - 250V) ◆
ARH1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster £130
LIGHTS*		
ARH1437		Light system with Indicators & LED Headlights £410
ARH1254		Without light system ◆
TRANSPORTATION KIT		
ARH1450		Preparation for DAHL Docking station: Docking station for car implementation, please order docking station separately (Adequate headrest is compulsory, see options below) £499
ARH1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) ◆
FIXATION BELT*		
ARH1329		MAX Lapbelt (adjustable both sides) fitted to base of backrest ◆
REAR VIEW MIRROR		
ARH1415		Rear-view Mirror RH £36
ARH1339		Rear-view Mirror LH £36
CANE HOLDER		
ARH1312		Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote) £36
ARH1363		Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote) £36

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARH1095	Electric heavy duty basic Tilt (0° - 20°) select either 470mm or 495mm seat to ground height. (Max user weight 180kg)
ARH0947	Electrical recliner (0° - 25°)
ARH0556	Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate)
ARH0589	Centre Mounted : LNX Powered: Centre mounted legrest, length compensating (including calf pad and angle adjustable footplate)
ARH0503	Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm
ARH0533	Side Mounted : Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARH0532	Side Mounted : ADM Manual : Angle adjustable with quick release lever, 0° - 80°, length 290mm - 460mm
ARH0531	Side Mounted : ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm
ARH0504	Angle Adjustable Footplate
ARH0513	Angle and depth adjustable Footplate (Adapt)
ARH0599	Angle adjustable flip-back, heavy-duty footboard (Recommended for SW 550mm)
ARH0522	Angle Adjustable Footboard (for max user weight up to 180 kg) (Recommended for SW 550 mm)
ARH0559	Angle Adjustable Footplate for Centre mounted legrest
ARH2201	REM110 Drive only
ARH2203	REM211 Drive & Actuators
ARH2204	REM216 Drive, Actuators & Lights
ARH2227	Bracket RH mounted for Attendant
ARH2228	Bracket LH mounted for Attendant
ARH2219	Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)
ARH2245	Invacare IDC Dual control (Please select IDC Dual control bracket below)
ARH1221	IDC Dual control Detachable Bracket RH
ARH1222	IDC Dual control Detachable Bracket LH
ARH2227	10-Way Switch Bracket RH, Mounted for Attendant
ARH2228	10-Way Switch Bracket LH, Mounted for Attendant
ARH1437	Light system with Indicators & LED Headlights



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):