



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MARCH 2021
Customer Ref:

Customer Ref*, Contact No*, Delivery Address*, Account No*, Order No*, Order Qty*, Order date*

TDX SP2 NB Express Spec 3-4

Standard Features

- Max user weight 150kg (23.5 Stone)
Crash tested ISO 7176 -Part 19
Telescopic adjustable seat width 430mm (17") (adjustable)
440mm - 480mm)
Seat depth set at 460mm (adjustable 410mm - 510mm)
Flex3 adjustable backrest, Matrix Flo-tech Lite cushion
Standard T-armrests, Detachable with long armpad
Powered Tilt (0°-20°)
Manual settable backrest Recline (0°- 30° in 3.75° steps)
Standard 80° hangers with fixed Footplates
LiNX REM110 control unit, Standard RH Bracket (can be swapped to LH)
6" Black castors, 14" solid Black drive wheels
Sure Step and Stability Lock
Matt Black body shrouds, Black rims
6 km/h (4mph) motors and 50Ah Batteries, 8Ah charger
4 Point tie down system, lapbelt
Express shipment! - orders will be shipped same day if ordered before 2pm, next day if ordered after 2pm!

Standard Base price: £6,616



Q'STRAIT QLK-150 Compatible

Please make your selection with a cross in the box [] Free of charge [O] standard feature [◆] POA special request

Customer weight = KG

Table with 2 columns: Code, Description, Price. Includes rows for TDX SP2 NB EXPRESS STOCK with codes 1650310 and 1650312.

TOTAL PRICE