



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

JUNE 2023  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Action 2NG Stock Spec 1

### Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Lap belt (adjustable both sides)
- Anti-tippers (24" wheels only)
- Anthracite Black
- Choice of 5 seat widths 15" - 19"
- Depth adjustable 16" to 17"
- Black padded nylon seat upholstery
- Half folding 7° backrest 17" (430mm)
- Swing in/out 80° detachable hanger with adjustable heel strap
- Removable, lightweight and height adjustable full arm (200mm > 260mm)
- Black 24" x 1.3/8" quick release rear wheel solid tyre or
- 12" X 2" fixed axle, black solid tyre wheel forward position
- Black 8" x 1.1/4" (200mm x 25mm) solid front castor aluminium fork
- Push to lock brakes
- NET PRICE ONLY



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>A2NG SELF PROPEL</b>		
BQ1529465-AF	A2NG 15" W X 16" to 17" D	£195
BQ1529466-AF	A2NG 16" W X 16" to 17" D	£195
BQ1529469-AF	A2NG 17" W X 16" to 17" D	£195
BQ1529470-AF	A2NG 18" W X 16" to 17" D	£195
BQ1529485-AF	A2NG 19" W X 16" to 17" D	£195

**A2NG TRANSIT**

BQ1529486- AF		A2NG 15" W X 16" to 17" D	£195
BQ1529487- AF		A2NG 16" W X 16" to 17" D	£195
BQ1529490- AF		A2NG 17" W X 16" to 17" D	£195
BQ1529491- AF		A2NG 18" W X 16" to 17" D	£195
BQ1529492- AF		A2NG 19" W X 16" to 17" D	£195

**TOTAL PRICE**