



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JUNE 2023  
 Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## Action 3NG DHP

### Standard Features

- Max user weight 125kg
- Crash tested ISO 7176-Part 19
- Dual Hand Rim Self Propulsion (DHP)
- Choice of 6 seat widths (15" - 20") 380mm to 505mm
- Choice of 2 seat depths (16" or 18") 400mm or 450mm
- Choice of 3 seat heights (17" - 19") 435mm to 485mm
- Black Nylon Upholstery
- Fixed 7° (17") 430mm backrest
- Lumbar support padded backrest with shaped cover
- 80° Swing in/Swing out legrest hangers
- Flip up footpates
- Swing away and detachable desk armrests
- Quick release axle
- Choice of 8 standard colours (with 7 special colour options)
- Push or Pull to lock brakes (configuration dependant)
- Transportation kit and lap belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

**Standard Base price: £985**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

| Code                 | Description                |           |                       |
|----------------------|----------------------------|-----------|-----------------------|
| <b>BASE PRODUCT*</b> |                            |           |                       |
| CEH0006              | ACTION 3NG DHR             | (NC 1084) | <input type="radio"/> |
| CEHSPEC              | SPECIAL                    |           | POA                   |
| <b>PROPULSION*</b>   |                            |           |                       |
| CEH0140              | DUAL HANDRIM RH PROPULSION | (NC 1330) | £985                  |
| CEH0150              | DUAL HANDRIM LH PROPULSION | (NC 1320) | £985                  |

| SEAT WIDTH*          |  |   |      |
|----------------------|--|---|------|
| CEH0180              |  | 15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm) (NC 0852)   | ○    |
| CEH0190              |  | 16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)   | ○    |
| CEH0200              |  | 17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)   | ◆    |
| CEH0210              |  | 18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)   | ○    |
| CEH0220              |  | 19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)   | ○    |
| CEH0230              |  | 20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)   | ○    |
| SEAT DEPTH*          |  |   |      |
| CEH0240              |  | 16" / 400 MM SEAT DEPTH (NC 0879)   | ○    |
| CEH0250              |  | 18" / 450 MM SEAT DEPTH (NC 0878)   | ◆    |
| CEH0252              |  | 16" TO 18" / 400 TO 450 MM ADJUSTABLE SEAT DEPTH (NC 0330, 0350)  | £21  |
| SEAT HEIGHT*         |  |   |      |
| CEH0270              |  | 19" / 485 MM SEAT HEIGHT  | ○    |
| CEH0280              |  | 18" / 460 MM SEAT HEIGHT  | ◆    |
| CEH0290              |  | 17" / 435 MM SEAT HEIGHT (NC 1010, 1030)  | ○    |
| SEAT UPHOLSTERY*     |  |   |      |
| CEH0330              |  | NYLON W/O HOOK AND LOOP STRIPS  | ◆    |
| CEH0340              |  | NYLON WITH HOOK AND LOOP STRIPS   | ○    |
| CEH0345              |  | TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM   | ○    |
| CEH0350              |  | PADDED NYLON  | £30  |
| BACKREST*            |  |   |      |
| CEH0380              |  | FIXED 0° BACKREST 2 HEIGHTS 16"/ 400 MM or 17" / 430 MM (NC 0471, 0477, 0945)   | ○    |
| CEH0390              |  | FIXED 0° BACKREST 20" / 510 MM HEIGHT (NC 0471, 0477, 0945)   | ○    |
| CEH0395              |  | FIXED 7° BACKREST HT ADJ 15" / 380 MM TO 18" / 460 MM FIX PUSH HANDLES (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)        | £67  |
| CEH0396              |  | FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 TO 18" / 460 MM HT ADJ PUSH HANDLES (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940) | £80  |
| CEH0410              |  | FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES (NC 0471, 0477, 0945)                                      | ◆    |
| CEH0420              |  | FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES (NC 0471, 0477, 0945)  | ○    |
| CEH0425              |  | FIXED 7° BACKREST 17" / 430 MM HEIGHT HT ADJ PUSHHANDLE (NC 0471, 0477, 0945)   | £43  |
| CEH0432              |  | ANGLE ADJUSTABLE AND FOLD DOWN (-15° to + 15°) 17" (NC 0465, 0471, 0474, 0478, 0480, 0481, 0882)                                  | £86  |
| CEH0455              |  | FOLDING 7° BACKREST 17" (NC 0471, 0477, 0945)   | £62  |
| CEH0460              |  | FOLDING 7° BACKREST 20" (NC 0471, 0477, 0945)   | £62  |
| BACKREST UPHOLSTERY* |  |   |      |
| CEH0465              |  | WITHOUT BACK UPHOLSTERY (NC 0882)   | ○    |
| CEH0471              |  | STANDARD NYLON SLING (JEMIMA) (NC 0882, 0905, 0906, 0907)   | ○    |
| CEH0474              |  | LUMBAR SUPPORT BI MATERIAL PADDED COVER   | ○    |
| CEH0477              |  | CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER  | £55  |
| CEH0478              |  | CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE (NC 0940, 0945)   | £62  |
| CEH0480              |  | CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE (NC 0940, 0945)  | £120 |
| CEH0481              |  | CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE (NC 0940, 0945)   | £120 |

| LEGREST HANGER RH* |  |   |   |      |
|--------------------|--|---|---|------|
| CEH0510            |  | 80° RH SWING IN/SWING OUT   | (NC 1111, 66666)                                    | ◆    |
| CEH0520            |  | ELEVATING LEGREST RH SWING IN/SWING OUT   | (NC 0611, 0930, 0935, 0936, 66666)                  | £116 |
| CEH0545            |  | AMPUTEE LEGREST RH  | (NC 0600, 0610, 0611, 0930, 0936, 1111)             | £133 |
| LEGREST HANGER LH* |  |   |   |      |
| CEH0560            |  | 80° LH SWING IN/SWING OUT   | (NC 2222)   | ◆    |
| CEH0570            |  | ELEVATING LEGREST LH SWING IN/SWING OUT   | (NC 0611, 0930, 0935, 0936, 55555)                  | £116 |
| CEH0595            |  | AMPUTEE LEGREST LH  | (NC 0611, 0620, 0630, 0930, 0936, 2222)             | £133 |
| FOOTREST RH*       |  |   |   |      |
| CEH0600            |  | FLIP UP RH FOOTPLATE  | (NC 0937, 0938)                                     | ◆    |
| CEH0610            |  | FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS. |   | £43  |
| SWO1111            |  | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)   |   | £43  |
| CEH0611            |  | FLIP UP ONE PIECE   | (NC 0620, 0630, 0935, 0936, 0937, 0938)             | £78  |
| CEE66666           |  | NO FOOTREST   |   | ○    |
| FOOTREST LH*       |  |   |   |      |
| CEF0620            |  | FLIP UP LH FOOTPLATE  | (NC 0937, 0938)                                     | ◆    |
| CEH0630            |  | FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS. |   | £43  |
| SWO2222            |  | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)   |   | £43  |
| CEE55555           |  | NO FOOTREST   |   | ○    |
| ARMREST RH*        |  |   |   |      |
| CEH0640            |  | REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED   | (NC 0760, 0770, 0780, 0790, 0795, 0796, 1480, 1483) | ○    |
| CEH0650            |  | REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED  | (NC 0760, 0770, 0780, 0790, 0795, 0796)             | ○    |
| CEH0660            |  | SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT   | (NC 0740, 0750, 0780, 0790, 0795, 0796, 1480, 1483) | ◆    |
| CEH0670            |  | SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT  | (NC 0740, 0750, 0780, 0790, 0795, 0796)             | ○    |
| CEH0680            |  | SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)   | (NC 0740, 0750, 0760, 0770, 0795, 0796, 1480, 1483) | £56  |
| CEH0690            |  | SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH  | (NC 0740, 0750, 0760, 0770, 0795, 0796, 1483)       | £56  |
| CEH0695            |  | REMOVABLE RH (N°4) A DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. UTOLOCK SAFETY DEVICE INCLUDED   | (NC 0740, 0750, 0760, 0770, 0780, 0790, 1480, 1483) | £60  |
| CEH0696            |  | REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED   | (NC 0740, 0750, 0760, 0770, 0780, 0790, 1483)       | £60  |
| CEH0732            |  | REMOVABLE RH (N°4) ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED  | (NC 1480, 1483)                                     | £59  |

| ARMREST LH* |  |   |     |
|-------------|--|---|-----|
| CEH0740     |  | REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1484)                                       | ○   |
| CEH0750     |  | REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED  | ○   |
| CEH0760     |  | SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT (NC 1480, 1484)   | ◆   |
| CEH0770     |  | SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT  | ○   |
| CEH0780     |  | SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) (NC 1480, 1484)   | £56 |
| CEH0790     |  | SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH (NC 1484)  | £56 |
| CEH0795     |  | REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. ) AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1484)          | £60 |
| CEH0796     |  | REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 1484) | £60 |
| CEH0832     |  | REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1484)        | £59 |

#### SEAT CUSHION

|         |  |  |      |
|---------|--|--|------|
| CEH0849 |  | BASIC 50 MM FOAM CUSHION - NYLON AND ANTI-SLIP COVER | £24  |
| CEH0869 |  | MATRX FLO-TECH LITE CUSHION - Low Risk               | £121 |
| CEH0871 |  | MATRX FLO-TECH CONTOUR VISCO CUSHION - High Risk     | £169 |
| CEH0873 |  | MATRX FLO-TECH IMAGE CUSHION - High Risk             | £288 |
| CEH0874 |  | MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk     | £398 |
| CEH0876 |  | MATRX LIBRA CUSHION - Very High Risk                 | £410 |

#### SEAT CUSHION DEPTH

|         |  |                       |   |
|---------|--|-----------------------|---|
| CEH0878 |  | FOR 400 MM SEAT DEPTH | ○ |
| CEH0879 |  | FOR 450 MM SEAT DEPTH | ○ |

#### OTHER SEATING OPTIONS

|         |  |  |      |
|---------|--|--|------|
| CEH0852 |  | SEAT RIGIDIZER   | £39  |
| CEH0882 |  | LATERAL POSITIONING PADS (PAIR)  | £61  |
| CEH0905 |  | NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0940, 0945)        | £131 |
| CEH0906 |  | HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0940, 0945)              | £159 |
| CEH0907 |  | HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0940, 0945) | £183 |
| CEH0930 |  | CALF STRAP   | ◆    |
| CEH0935 |  | LONGER FOOTREST TUBE (NC 0936)   | £30  |
| CEH0936 |  | SHORTER FOOTREST SUPPORT   | £42  |
| CEH0937 |  | BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED                                | £114 |
| CEH0938 |  | BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED                                | £114 |
| CEH0940 |  | BACKBRACE (NC 0945)  | £34  |
| CEH0945 |  | FOLDABLE BACKBRACE   | £38  |

| <b>COLOUR*</b>              |                          |  |                       |
|-----------------------------|--------------------------|--|-----------------------|
| CEH0961                     | <input type="checkbox"/> | SLATE GREY (RAL7024 FINE TEXTURE)        | <input type="radio"/> |
| CEH0964                     | <input type="checkbox"/> | PEARL WHITE                              | <input type="radio"/> |
| CEH0965                     | <input type="checkbox"/> | ELECTRIC GREEN                           | <input type="radio"/> |
| CEH0966                     | <input type="checkbox"/> | MATT BLACK                               | <input type="radio"/> |
| CEH0967                     | <input type="checkbox"/> | BLUE POP                                 | <input type="radio"/> |
| CEH0968                     | <input type="checkbox"/> | HYPNOTIC PURPLE                          | <input type="radio"/> |
| CEH0969                     | <input type="checkbox"/> | RED MATT                                 | <input type="radio"/> |
| CEH0972                     | <input type="checkbox"/> | WATER BLUE                               | <input type="radio"/> |
| CEH0980                     | <input type="checkbox"/> | ANTHRACITE BLACK                         | <input type="radio"/> |
| CEH0991                     | <input type="checkbox"/> | SPECIAL COLOUR ON DEMAND ONLY            | £169                  |
| SWO3333                     | <input type="checkbox"/> | DARK RED - SPECIAL                       | £169                  |
| SWO4444                     | <input type="checkbox"/> | PEARL GREY - SPECIAL                     | £169                  |
| SWO5555                     | <input type="checkbox"/> | INTENSE YELLOW                           | £169                  |
| SWO6666                     | <input type="checkbox"/> | VIBRANT ORANGE                           | £169                  |
| SWO7777                     | <input type="checkbox"/> | POLAR BLUE                               | £169                  |
| SWO8888                     | <input type="checkbox"/> | CANDY PINK                               | £169                  |
| <b>CASTOR HOUSING*</b>      |                          |  |                       |
| CEH0995                     | <input type="checkbox"/> | STANDARD                                 | ◆                     |
| CEH0996                     | <input type="checkbox"/> | LIGHT WEIGHT                             | £28                   |
| <b>FRONT CASTORS*</b>       |                          |  |                       |
| CEH1010                     | <input type="checkbox"/> | 8" X 1.1/4" BLACK SOLID TYRE             | ◆                     |
| CEH1030                     | <input type="checkbox"/> | 8" X 2" BLACK SOLID TYRE                 | £28                   |
| CEH1035                     | <input type="checkbox"/> | 6" X 1.1/4" BLACK SOLID TYRE MEDIUM HARD | £28                   |
| CEH1050                     | <input type="checkbox"/> | 6" X 1.1/4" BLACK SOLID TYRE             | £28                   |
| <b>REAR WHEEL SIZE*</b>     |                          |  |                       |
| CEH1068                     | <input type="checkbox"/> | 24"                                      | ◆                     |
| <b>REAR WHEELS*</b>         |                          |  |                       |
| CEH1075                     | <input type="checkbox"/> | ONE ARM DRIVE (NC 1084)                  | ◆                     |
| <b>TYRES*</b>               |                          |  |                       |
| CEH1082                     | <input type="checkbox"/> | SOLID BLACK 1.3/8" / 35MM                | £21                   |
| CEH1084                     | <input type="checkbox"/> | PNEUMATIC PROFILE BLACK 1.3/8" / 35MM    | ◆                     |
| <b>REAR WHEEL AXLE*</b>     |                          |  |                       |
| CEH1180                     | <input type="checkbox"/> | Q/R AXLE                                 | ◆                     |
| <b>REAR WHEEL POSITION*</b> |                          |  |                       |
| CEH1190                     | <input type="checkbox"/> | STANDARD REAR WHEELS POSITION            | <input type="radio"/> |
| CEH1200                     | <input type="checkbox"/> | PASSIVE REAR WHEELS POSITION             | ◆                     |
| <b>PUMP</b>                 |                          |  |                       |
| CEH1230                     | <input type="checkbox"/> | PUMP AND SUPPORT                         | £14                   |
| <b>CHASSIS OPTION</b>       |                          |  |                       |
| CEH1241                     | <input type="checkbox"/> | FRONT WHEELS EXTENSION                   | £22                   |
| <b>BRAKES*</b>              |                          |  |                       |
| CEH1270                     | <input type="checkbox"/> | PUSH TO LOCK BRAKES                      | ◆                     |
| CEH1275                     | <input type="checkbox"/> | PULL TO LOCK BRAKES                      | <input type="radio"/> |
| CEH1320                     | <input type="checkbox"/> | OAD WIRE CONTROLLED BRAKE RH (NC 1390)   | £45                   |
| CEH1330                     | <input type="checkbox"/> | OAD WIRE CONTROLLED BRAKE LH (NC 1390)   | £45                   |

| TRANSPORTATION KIT* |  |  |                               |
|---------------------|--|--|-------------------------------|
| CEH1345             |  | TRANSPORTATION BRACKET   | ◆                             |
| SAFETY OPTIONS*     |  |  |                               |
| CEH1360             |  | LAP BELT   | (NC 1361, 1362, 1363, 1364) ◆ |
| CEH1361             |  | BODYPOINT 2-PT POSTURE BELT M                                  | (NC 1362, 1363, 1364) £83     |
| CEH1362             |  | BODYPOINT 2-PT POSTURE BELT L                                  | (NC 1363, 1364) £83           |
| CEH1363             |  | BODYPOINT 4-PT POSTURE BELT M                                  | (NC 1364) £98                 |
| CEH1364             |  | BODYPOINT 4-PT POSTURE BELT L                                  | £98                           |
| CEH1370             |  | REAR ANTI-TIPPERS PAIR   | (NC 1375, 1440, 1441) £36     |
| CEH1375             |  | REAR ANTI-TIPPERS PAIR WITH TIPPING AID                        | (NC 1440, 1441) £72           |
| CEH1380             |  | REFLECTORS KIT   | £8                            |
| CEH1390             |  | BRAKE EXTENSION PAIR   | £14                           |
| OTHER OPTIONS       |  |  |                               |
| CEH1430             |  | CANE HOLDER  | £23                           |
| CEH1440             |  | TIPPING AID RH   | £14                           |
| CEH1441             |  | TIPPING AID LH   | £14                           |
| CEH1480             |  | ACRYLIC TRAY   | (NC 1483, 1484, 1485) £116    |
| CEH1483             |  | HALF TABLE TRAY RH   | (NC 1484) £89                 |
| CEH1484             |  | HALF TABLE TRAY LH   | £89                           |
| CEH1485             |  | COMFORT PADDING FOR HALF TABLE TRAY                            | £44                           |
| SWO10001            |  | ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618 | £107                          |

|                    |  |
|--------------------|--|
| <b>TOTAL PRICE</b> |  |
|--------------------|--|

**Non Conforming Options:**

|          |   |
|----------|---|
| CEH0330  | NYLON W/O HOOK AND LOOP STRIPS  |
| CEH0350  | PADDED NYLON  |
| CEH0465  | WITHOUT BACK UPHOLSTERY   |
| CEH0471  | STANDARD NYLON SLING (JEMIMA)   |
| CEH0474  | LUMBAR SUPPORT BI MATERIAL PADDED COVER   |
| CEH0477  | CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER  |
| CEH0478  | CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE   |
| CEH0480  | CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE  |
| CEH0481  | CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE   |
| CEH0600  | FLIP UP RH FOOTPLATE  |
| CEH0610  | FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS. |
| SWO1111  | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)   |
| CEH0611  | FLIP UP ONE PIECE   |
| CEE66666 | NO FOOTREST   |
| CEF0620  | FLIP UP LH FOOTPLATE  |
| CEH0630  | FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS. |
| SWO2222  | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)   |
| CEE55555 | NO FOOTREST   |
| CEH0740  | REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED   |
| CEH0750  | REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED  |
| CEH0760  | SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT   |
| CEH0770  | SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT  |
| CEH0780  | SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)   |
| CEH0790  | SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH  |
| CEH0795  | REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. ) AUTOLOCK SAFETY DEVICE INCLUDED  |
| CEH0796  | REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED   |
| CEH0878  | FOR 400 MM SEAT DEPTH   |
| CEH0879  | FOR 450 MM SEAT DEPTH   |
| CEH0852  | SEAT RIGIDIZER  |
| CEH0882  | LATERAL POSITIONING PADS (PAIR)   |
| CEH0905  | NECKREST BLACK MULTI STRETCH PU COATED FABRIC   |
| CEH0906  | HEADREST BLACK MULTI STRETCH PU COATED FABRIC   |
| CEH0907  | HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC  |
| CEH0930  | CALF STRAP  |
| CEH0935  | LONGER FOOTREST TUBE  |
| CEH0936  | SHORTER FOOTREST SUPPORT  |
| CEH0937  | BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED   |
| CEH0938  | BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED   |
| CEH0940  | BACKBRACE   |
| CEH0945  | FOLDABLE BACKBRACE  |

|         |   |
|---------|---|
| CEH1010 | 8" X 1.1/4" BLACK SOLID TYRE            |
| CEH1030 | 8" X 2" BLACK SOLID TYRE                |
| CEH1084 | PNEUMATIC PROFILE BLACK 1.3/8" / 35MM   |
| CEH1320 | OAD WIRE CONTROLLED BRAKE RH            |
| CEH1330 | OAD WIRE CONTROLLED BRAKE LH            |
| CEH1361 | BODYPOINT 2-PT POSTURE BELT M           |
| CEH1362 | BODYPOINT 2-PT POSTURE BELT L           |
| CEH1363 | BODYPOINT 4-PT POSTURE BELT M           |
| CEH1364 | BODYPOINT 4-PT POSTURE BELT L           |
| CEH1375 | REAR ANTI-TIPPERS PAIR WITH TIPPING AID |
| CEH1390 | BRAKE EXTENSION PAIR                    |
| CEH1440 | TIPPING AID RH                          |
| CEH1441 | TIPPING AID LH                          |
| CEH1480 | ACRYLIC TRAY                            |
| CEH1483 | HALF TABLE TRAY RH                      |
| CEH1484 | HALF TABLE TRAY LH                      |
| CEH1485 | COMFORT PADDING FOR HALF TABLE TRAY     |





For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):