



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JUNE 2023
Customer Ref:

Customer Ref\*, Contact No\*, Delivery Address\*, Account No\*, Order No\*, Order Qty\*, Order date\*

Action 3NG Self-Propel

Standard Features

- Max user weight 125kg (19.7 stone)
Crash tested ISO- 7176-Part 19
Self Propulsion with 22" or 24" Quick release rear wheels
Choice of 6 seat widths 380mm to 505mm (15" to 20")
Choice of 2 seat depths 400mm or 450mm (16" or 18")
Choice of 6 seat heights 385mm to 510mm (15" to 20")
Black nylon upholstery
Fixed 7° 430mm backrest
Lumbar support padded backrest with shaped cover
80° Swing in / Swing out legrest hangers
Flip up footplates
Swing away and detachable desk armrests
Choice of 8 colours (with 7 special colours)
Push or pull to lock brakes (configuration dependant)
Transportation kit and lap belt (adjustable both sides)
Adjustable heel straps
User tool kit

Standard Base price: £641



Please make your selection with a cross in the box [ ] Free of charge [ ] standard feature [ ] POA special request [ ]

Customer weight = KG

Table with 2 columns: Code, Description. Row: CEE0100 [ ] SELF PROPULSION £641

| <b>SEAT WIDTH*</b>      |  |   |     |
|-------------------------|--|---|-----|
| CEE44444                |  | 13" / 330 MM SEAT WIDTH - SELECT SEAT WIDTH<br>OPTION CEE0180 AND IN-RIGGED ARMREST<br>OPTION CEE0952 (lapbelt adjustment 420mm -<br>1520mm)                    | £41 |
| CEE0180                 |  | 15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm<br>- 1520mm) (NC 0852)  | ○   |
| CEE0190                 |  | 16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm<br>- 1520mm)  | ○   |
| CEE0200                 |  | 17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm<br>- 1520mm)  | ◆   |
| CEE0210                 |  | 18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm<br>- 1880mm)  | ○   |
| CEE0220                 |  | 19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm<br>- 1880mm)  | ○   |
| CEE0230                 |  | 20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm<br>- 1880mm)  | ○   |
| <b>SEAT DEPTH*</b>      |  |   |     |
| CEE0240                 |  | 16" / 400 MM SEAT DEPTH (NC 0879)   | ○   |
| CEE0250                 |  | 18" / 450 MM SEAT DEPTH (NC 0878)   | ◆   |
| CEE0252                 |  | 16" TO 18" / 400 TO 450 MM ADJUSTABLE SEAT<br>DEPTH (NC 0330, 0350)   | £21 |
| <b>SEAT HEIGHT*</b>     |  |   |     |
| CEE0260                 |  | 20" / 510 MM SEAT HEIGHT (NC 0511, 0561, 0738, 1035, 1050, 1054, 1056, 1066)  | ○   |
| CEE0270                 |  | 19" / 485 MM SEAT HEIGHT (NC 0511, 0561, 1054)  | ○   |
| CEE0280                 |  | 18" / 460 MM SEAT HEIGHT  | ◆   |
| CEE0290                 |  | 17" / 435 MM SEAT HEIGHT (NC 1010, 1030)  | ○   |
| CEE0300                 |  | 16" / 410 MM SEAT HEIGHT (NC 0935, 1010, 1030, 1056, 1370, 1440, 1441)  | ○   |
| CEE0310                 |  | 15" / 385 MM SEAT HEIGHT , (NC 0511, 0561, 0660, 0670, 0738, 0760, 0770, 0838, 0839, 0935,<br>1010, 1030, 1035, 1050, 1056, 1068, 1073, 1280, 1370, 1440, 1441) | ○   |
| <b>SEAT UPHOLSTERY*</b> |  |   |     |
| CEE0330                 |  | NYLON W/O HOOK AND LOOP STRIPS  | ◆   |
| CEE0340                 |  | NYLON WITH HOOK AND LOOP STRIPS   | ○   |
| CEE0345                 |  | TENSION ADJUSTABLE SEAT UPHOLSTERY &<br>DEPTH ADJUSTABLE -50MM  | ○   |
| CEE0350                 |  | PADDED NYLON  | £30 |

| <b>BACKREST - Anti-tippers should be selected with reclining backrests 0435, 0435, 0440, 0445 and 0450*</b> |  |   |  |      |
|---|--|---|--|------|
| CEE0380   |  | FIXED 0° BACKREST 2 HEIGHTS 16"/ 400 MM or 17" / 430 MM                       | (NC 0471, 0477, 0945)                                      | ○    |
| CEE0390   |  | FIXED 0° BACKREST 20" / 510 MM HEIGHT   | (NC 0471, 0477, 0945)                                      | ○    |
| CEE0395   |  | FIXED 7° BACKREST HT ADJ 15" / 380 MM TO 18" / 460 MM FIX PUSH HANDLES        | (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)        | £67  |
| CEE0396   |  | FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 TO 18" / 460 MM HT ADJ PUSH HANDLES | (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)        | £80  |
| CEE0410   |  | FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES        | (NC 0471, 0477, 0945)                                      | ◆    |
| CEE0420   |  | FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES                        | (NC 0471, 0477, 0945)                                      | ○    |
| CEE0425   |  | FIXED 7° BACKREST 17" / 430 MM HEIGHT HT ADJ PUSHHANDLE                       | (NC 0471, 0477, 0945)                                      | £43  |
| CEE0428   |  | FIXED 7° BACKREST ACTIVE 14" / 350 MM HT ADJ TO 16" / 400 MM FIX PUSH HANDLES | (NC 0471, 0477, 0480, 0481, 0905, 0906, 0907, 0940, 1280)  | ○    |
| CEE0432   |  | ANGLE ADJUSTABLE AND FOLD DOWN (15° to + 15°) 17"                             | (NC 0465, 0471, 0474, 0478, 0480, 0481, 0882, 0952)        | £86  |
| CEE0435   |  | RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 15" 16" 17"                      | (NC 0471, 0477, 0738, 0838, 0940, 0945, 1190)              | £175 |
| CEE0445   |  | RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 15" 16" 17"                   | (NC 0471, 0477, 0738, 0838, 0940, 0945, 1190)              | £233 |
| CEE0450   |  | RECLINING (GAS SPRING) 0/30° 3 HEIGHTS 18" 19" 20"                            | (NC 0471, 0477, 0738, 0838, 0940, 0945, 1190)              | £233 |
| CEE0455   |  | FOLDING 7° BACKREST 17"   | (NC 0471, 0477, 0945)                                      | £62  |
| CEE0460   |  | FOLDING 7° BACKREST 20"   | (NC 0471, 0477, 0945)                                      | £62  |
| <b>BACKREST UPHOLSTERY*</b>   |  |   |  |      |
| CEE0465   |  | WITHOUT BACK UPHOLSTERY   | (NC 0882)  | ○    |
| CEE0471   |  | STANDARD NYLON SLING (JEMIMA)   | (NC 0882, 0905, 0906, 0907)                                | ○    |
| CEE0474   |  | LUMBAR SUPPORT BI MATERIAL PADDED COVER                                       |  | ○    |
| CEE0477   |  | CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER                      |  | £55  |
| CEE0478   |  | CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE       | (NC 0940, 0945)  | £62  |
| CEE0480   |  | CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE            | (NC 0940, 0945)  | £120 |
| CEE0481   |  | CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE       | (NC 0940, 0945)  | £120 |
| <b>LEGREST HANGER RH*</b>   |  |   |  |      |
| CEE0510   |  | 80° RH SWING IN/SWING OUT   | (NC 1111, 66666)   | ◆    |
| CEE0511   |  | 90° RH SWING IN/SWING OUT   | (NC 0600, 0936, 1010, 1030, 1035, 1050, 1056, 1111, 66666) | ○    |
| CEE0520   |  | ELEVATING LEGREST RH SWING IN/SWING OUT                                       | (NC 0611, 0930, 0935, 0936, 66666)                         | £116 |
| CEE0545   |  | AMPUTEE LEGREST RH  | (NC 0600, 0610, 0611, 0930, 0936, 1111)                    | £133 |
| <b>LEGREST HANGER LH*</b>   |  |   |  |      |
| CEE0560   |  | 80° LH SWING IN/SWING OUT   | (NC 2222)  | ◆    |
| CEE0561   |  | 90° LH SWING IN/SWING OUT   | (NC 0620, 0936, 1010, 1030, 1035, 1050, 1056, 2222)        | ○    |
| CEE0570   |  | ELEVATING LEGREST LH SWING IN/SWING OUT                                       | (NC 0611, 0930, 0935, 0936, 55555)                         | £116 |
| CEE0595   |  | AMPUTEE LEGREST LH  | (NC 0611, 0620, 0630, 0930, 0936, 2222)                    | £133 |

| FOOTREST RH* |  |   |     |
|--------------|--|---|-----|
| CEE0600      |  | FLIP UP RH FOOTPLATE (NC 0937, 0938)  | ◆   |
| CEE0610      |  | FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS. | £43 |
| SWO1111      |  | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)   | £43 |
| CEE0611      |  | FLIP UP ONE PIECE (NC 0620, 0630, 0935, 0936, 0937, 0938)   | £78 |
| CEE66666     |  | NO FOOTREST   | ○   |
| FOOTREST LH* |  |   |     |
| CEE0620      |  | FLIP UP LH FOOTPLATE (NC 0937, 0938)  | ◆   |
| CEE0630      |  | FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS. | £43 |
| SWO2222      |  | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)   | £43 |
| CEE55555     |  | NO FOOTREST   | ○   |
| ARMREST RH*  |  |   |     |
| CEE0640      |  | REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 0760, 0770, 0780, 0790, 0795, 0796, 0838, 0839, 1480, 1483)   | ○   |
| CEE0650      |  | REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 0760, 0770, 0780, 0790, 0795, 0796, 0838, 0839)  | ○   |
| CEE0660      |  | SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT (NC 0740, 0750, 0780, 0790, 0795, 0796, 0838, 0839, 0952, 0955, 1480, 1483)   | ◆   |
| CEE0670      |  | SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT (NC 0740, 0750, 0780, 0790, 0795, 0796, 0838, 0839, 0952, 0955)  | ○   |
| CEE0680      |  | SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170 to 295 MM) (NC 0740, 0750, 0760, 0770, 0795, 0796, 0838, 0839, 0952, 0955, 1480, 1483)  | £56 |
| CEE0690      |  | SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170 to 295 MM) AND DEPTH (NC 0740, 0750, 0760, 0770, 0795, 0796, 0838, 0839, 0952, 0955, 1483)   | £56 |
| CEE0695      |  | REMOVABLE RH (N°4) DESK ADJ. IN HEIGHT (210345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0740, 0750, 0760, 0770, 0780, 0790, 0838, 0839, 0952, 0955, 1480, 1483)   | £60 |
| CEE0696      |  | REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0740, 0750, 0760, 0770, 0780, 0790, 0838, 0839, 0952, 0955, 1483)  | £60 |
| CEE0732      |  | REMOVABLE RH (N°4) ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED (NC 0838, 0839, 0952, 1480, 1483)  | £59 |
| CEE0738      |  | REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM) (NC 0740, 0750, 0760, 0770, 0780, 0790, 0795, 0796, 0832, 0839, 0952, 0955, 1066, 1200, 1201, 1480, 1483)   | £35 |

| ARMREST LH*        |  |   |      |
|--------------------|--|---|------|
| CEE0740            |  | REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1484)   | ○    |
| CEE0750            |  | REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED  | ○    |
| CEE0760            |  | SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT (NC 0952, 0955, 1480, 1484)   | ◆    |
| CEE0770            |  | SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT (NC 0952, 0955)  | ○    |
| CEE0780            |  | SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170 to 295 MM) (NC 0952, 0955, 1480, 1484)  | £56  |
| CEE0790            |  | SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170 to 295 MM) AND DEPTH (NC 0952, 0955, 1484)   | £56  |
| CEE0795            |  | REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0952, 0955, 1480, 1484)   | £60  |
| CEE0796            |  | REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0952, 0955, 1484)                                      | £60  |
| CEE0832            |  | REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION (NC 0952, 1480, 1484)   | £59  |
| CEE0838            |  | REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM) (NC 0952, 0955, 1066, 1200, 1201, 1484)                           | £35  |
| CEE0839            |  | TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) (NC 0952, 0955, 1200, 1201, 1345, 1430, 1480, 1484) | ○    |
| SEAT CUSHION       |  |   |      |
| CEE0849            |  | BASIC 50 MM FOAM CUSHION NYLON AND ANTISLIP COVER (NC 0952)   | £24  |
| CEE0869            |  | MATRX FLO-TECH LITE CUSHION - Low Risk (NC 0952)  | £121 |
| CEE0871            |  | MATRX FLO-TECH CONTOUR VISCO CUSHION - High Risk (NC 0952)  | £169 |
| CEE0873            |  | MATRX FLO-TECH IMAGE CUSHION - High Risk (NC 0952)  | £288 |
| CEE0874            |  | MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk (NC 0952)  | £398 |
| CEE0876            |  | MATRX LIBRA CUSHION - Very High Risk (NC 0952)  | £410 |
| SEAT CUSHION DEPTH |  |   |      |
| CEE0878            |  | FOR 400 MM SEAT DEPTH   | ○    |
| CEE0879            |  | FOR 450 MM SEAT DEPTH   | ○    |

| OTHER SEATING OPTIONS |  |  |      |
|-----------------------|--|--|------|
| CEE0852               |  | SEAT RIGIDIZER (NC 0952)   | £39  |
| CEE0882               |  | LATERAL POSITIONING PADS (PAIR)  | £61  |
| CEE0905               |  | NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0940, 0945)        | £131 |
| CEE0906               |  | HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0940, 0945)              | £159 |
| CEE0907               |  | HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0940, 0945) | £183 |
| SWO3333               |  | MULT ADJUSTABLE LATERALS   | £284 |
| CEE0930               |  | CALF STRAP   | ○    |
| CEE0935               |  | LONGER FOOTREST TUBE   | £30  |
| CEE0936               |  | SHORTER FOOTREST SUPPORT   | £42  |
| CEE0937               |  | BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED                                | £114 |
| CEE0938               |  | BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED                                | £114 |
| CEE0940               |  | BACKBRACE (NC 0945)  | £34  |
| CEE0945               |  | FOLDABLE BACKBRACE   | £38  |
| CEE0952               |  | INRIGGED ARMREST (NC 1480, 1483, 1484)   | £41  |
| CEE0955               |  | W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY                                  | ○    |

| COLOUR* |  |                                   |      |
|---------|--|-----------------------------------|------|
| CEE0961 |  | SLATE GREY (RAL7024 FINE TEXTURE) | ○    |
| CEE0964 |  | PEARL WHITE                       | ○    |
| CEE0965 |  | ELECTRIC GREEN                    | ○    |
| CEE0966 |  | MATT BLACK                        | ○    |
| CEE0967 |  | BLUE POP                          | ○    |
| CEE0968 |  | HYPNOTIC PURPLE                   | ○    |
| CEE0969 |  | RED MATT                          | ○    |
| CEE0972 |  | WATER BLUE                        | ○    |
| CEE0980 |  | ANTHRACITE BLACK                  | ○    |
| CEE0991 |  | SPECIAL COLOUR ON DEMAND ONLY     | £169 |
| SWO4444 |  | DARK RED SPECIAL                  | £169 |
| SWO5555 |  | PEARL GREY SPECIAL                | £169 |
| SWO6666 |  | INTENSE YELLOW                    | £169 |
| SWO7777 |  | VIBRANT ORANGE                    | £169 |
| SWO8888 |  | POLAR BLUE                        | £169 |
| SWO9999 |  | CANDY PINK                        | £169 |

| CASTOR HOUSING* |  |              |     |
|-----------------|--|--------------|-----|
| CEE0995         |  | STANDARD     | ◆   |
| CEE0996         |  | LIGHT WEIGHT | £28 |

| FRONT CASTORS* |  |  |     |
|----------------|--|--|-----|
| CEE1010        |  | 8" X 1.1/4" BLACK SOLID TYRE                 | ◆   |
| CEE1030        |  | 8" X 2" BLACK SOLID TYRE                     | £28 |
| CEE1035        |  | 6" X 1.1/4" BLACK SOLID TYRE MEDIUM HARD     | £28 |
| CEE1050        |  | 6" X 1.1/4" BLACK SOLID TYRE                 | £28 |
| CEE1054        |  | 5" X 1" BLACK SOLID TYRE                     | £28 |
| CEE1056        |  | 5.1/2" X 1.3/4" BLACK SOLID TYRE MEDIUM HARD | £28 |

| REAR WHEEL SIZE*     |  |   |   |      |
|----------------------|--|---|---|------|
| CEE1066              |  | 22"   | (NC 1092, 1260)                         | ○    |
| CEE1068              |  | 24"   | (NC 1261)                               | ◆    |
| REAR WHEELS*         |  |   |   |      |
| CEE1071              |  | STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES )  | (NC 1092, 1170, 1280)                   | ◆    |
| CEE1072              |  | SELF PROPEL REINFORCED REAR WHEEL                                     | (NC 1092, 1280)                         | £23  |
| CEE1073              |  | DRUM BRAKE SELFPROPEL 22 AND 24" or TRANSIT 12"                       | (NC 1092, 1201, 1270, 1275)             | ○    |
| CEE1076              |  | LIGHTWEIGHT (BLACK HUB AND RIM, 24 BLACK RADIAL SPOKES)               | (NC 1170, 1201, 1230, 1260, 1261, 1280) | £59  |
| TYRES*               |  |   |   |      |
| CEE1082              |  | SOLID BLACK 1.3/8" / 35MM   |   | ◆    |
| CEE1084              |  | PNEUMATIC PROFILE BLACK 1.3/8" / 35MM                                 |   | £21  |
| CEE1088              |  | PNEUMATIC LOW PROFILE BLACK 1" / 25MM                                 |   | ○    |
| CEE1092              |  | PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM |   | £39  |
| HANDRIM*             |  |   |   |      |
| CEE1140              |  | ALUMINIUM PREANODIZED HANDRIM   |   | ◆    |
| CEE1150              |  | HANDRIM COATED IN BLACK ANTISLIP PLASTIC                              |   | £34  |
| REAR WHEEL AXLE*     |  |   |   |      |
| CEE1170              |  | FIXED AXLE  |   | ○    |
| CEE1180              |  | Q/R AXLE  |   | ◆    |
| REAR WHEEL POSITION* |  |   |   |      |
| CEE1190              |  | STANDARD REAR WHEELS POSITION   |   | ◆    |
| CEE1200              |  | PASSIVE REAR WHEELS POSITION  |   | ○    |
| CEE1201              |  | DUAL AMPUTEE REAR WHEEL POSITION                                      | (NC 1275, 1280)                         | £72  |
| PUMP                 |  |   |   |      |
| CEE1230              |  | PUMP AND SUPPORT  |   | £14  |
| CHASSIS OPTION       |  |   |   |      |
| CEE1241              |  | FRONT WHEELS EXTENSION  |   | £22  |
| CEE1260              |  | SPOKE GUARDS PAIR 24" BLK   | (NC 1261)                               | £78  |
| CEE1261              |  | SPOKE GUARDS PAIR 560 AND 515 MM (22" AND 20") CLR                    |   | £98  |
| BRAKES*              |  |   |   |      |
| CEE1270              |  | PUSH TO LOCK BRAKES   |   | ◆    |
| CEE1275              |  | PULL TO LOCK BRAKES   |   | ○    |
| CEE1280              |  | ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)                  |   | £231 |
| TRANSPORTATION KIT*  |  |   |   |      |
| CEE1345              |  | TRANSPORTATION BRACKETS   |   | ◆    |
| SAFETY OPTIONS*      |  |   |   |      |
| CEE1360              |  | LAP BELT  | (NC 1361, 1362, 1363, 1364)             | ◆    |
| CEE1361              |  | BODYPOINT 2-PT POSTURE BELT M   | (NC 1362, 1363, 1364)                   | £83  |
| CEE1362              |  | BODYPOINT 2-PT POSTURE BELT L   | (NC 1363, 1364)                         | £83  |
| CEE1363              |  | BODYPOINT 4-PT POSTURE BELT M   | (NC 1364)                               | £98  |
| CEE1364              |  | BODYPOINT 4-PT POSTURE BELT L   |   | £98  |
| CEE1370              |  | REAR ANTITIPPERS PAIR   | (NC 1375, 1440, 1441)                   | £36  |
| CEE1375              |  | REAR ANTITIPPERS PAIR WITH TIPPING AID                                | (NC 1440, 1441)                         | £72  |
| CEE1380              |  | REFLECTORS KIT  |   | £8   |
| CEE1390              |  | BRAKE EXTENSION PAIR  |   | £14  |

| OTHER OPTIONS |  |  |      |
|---------------|--|--|------|
| CEE1430       |  | CANE HOLDER  | £23  |
| CEE1440       |  | TIPPING AID RH   | £14  |
| CEE1441       |  | TIPPING AID LH   | £14  |
| CEE1480       |  | ACRYLIC TRAY (NC 1483, 1484, 1485)                             | £116 |
| CEE1483       |  | HALF TABLE TRAY RH (NC 1484)                                   | £89  |
| CEE1484       |  | HALF TABLE TRAY LH   | £89  |
| CEE1485       |  | COMFORT PADDING FOR HALF TABLE TRAY                            | £44  |
| SWO10002      |  | ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618 | £107 |

|                    |  |
|--------------------|--|
| <b>TOTAL PRICE</b> |  |
|--------------------|--|



**Non Conforming Options:**

|          |   |
|----------|---|
| CEE0330  | NYLON W/O HOOK AND LOOP STRIPS  |
| CEE0350  | PADDED NYLON  |
| CEE0465  | WITHOUT BACK UPHOLSTERY   |
| CEE0471  | STANDARD NYLON SLING (JEMIMA)   |
| CEE0474  | LUMBAR SUPPORT BI MATERIAL PADDED COVER   |
| CEE0477  | CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER  |
| CEE0478  | CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE   |
| CEE0480  | CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE  |
| CEE0481  | CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE   |
| CEE0511  | 90° RH SWING IN/SWING OUT   |
| CEE0561  | 90° LH SWING IN/SWING OUT   |
| CEE0600  | FLIP UP RH FOOTPLATE  |
| CEE0610  | FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS. |
| SWO1111  | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)   |
| CEE0611  | FLIP UP ONE PIECE   |
| CEE66666 | NO FOOTREST   |
| CEE0620  | FLIP UP LH FOOTPLATE  |
| CEE0630  | FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS. |
| SWO2222  | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)   |
| CEE55555 | NO FOOTREST   |
| CEE0660  | SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT   |
| CEE0670  | SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT  |
| CEE0738  | REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)   |
| CEE0740  | REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED   |
| CEE0750  | REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED  |
| CEE0760  | SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT   |
| CEE0770  | SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT  |
| CEE0780  | SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170 to 295 MM)  |
| CEE0790  | SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170 to 295 MM) AND DEPTH   |
| CEE0795  | REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED   |
| CEE0796  | REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED  |
| CEE0832  | REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION   |
| CEE0838  | REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)   |
| CEE0839  | TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM)   |
| CEE0878  | FOR 400 MM SEAT DEPTH   |
| CEE0879  | FOR 450 MM SEAT DEPTH   |
| CEE0852  | SEAT RIGIDIZER  |
| CEE0882  | LATERAL POSITIONING PADS (PAIR)   |

|         |   |
|---------|---|
| CEE0905 | NECKREST BLACK MULTI STRETCH PU COATED FABRIC                         |
| CEE0906 | HEADREST BLACK MULTI STRETCH PU COATED FABRIC                         |
| CEE0907 | HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC      |
| CEE0930 | CALF STRAP  |
| CEE0935 | LONGER FOOTREST TUBE  |
| CEE0936 | SHORTER FOOTREST SUPPORT  |
| CEE0937 | BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED                     |
| CEE0938 | BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED                     |
| CEE0940 | BACKBRACE   |
| CEE0945 | FOLDABLE BACKBRACE  |
| CEE0952 | INRIGGED ARMREST  |
| CEE0955 | W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY                       |
| CEE1010 | 8" X 1.1/4" BLACK SOLID TYRE  |
| CEE1030 | 8" X 2" BLACK SOLID TYRE  |
| CEE1035 | 6" X 1.1/4" BLACK SOLID TYRE MEDIUM HARD                              |
| CEE1050 | 6" X 1.1/4" BLACK SOLID TYRE  |
| CEE1054 | 5" X 1" BLACK SOLID TYRE  |
| CEE1056 | 5.1/2" X 1.3/4" BLACK SOLID TYRE MEDIUM HARD                          |
| CEE1066 | 22"   |
| CEE1068 | 24"   |
| CEE1073 | DRUM BRAKE SELFPROPEL 22 AND 24" or TRANSIT 12"                       |
| CEE1092 | PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM |
| CEE1170 | FIXED AXLE  |
| CEE1190 | STANDARD REAR WHEELS POSITION   |
| CEE1200 | PASSIVE REAR WHEELS POSITION  |
| CEE1201 | DUAL AMPUTEE REAR WHEEL POSITION                                      |
| CEE1230 | PUMP AND SUPPORT  |
| CEE1260 | SPOKE GUARDS PAIR 24" BLK   |
| CEE1261 | SPOKE GUARDS PAIR 560 AND 515 MM (22" AND 20") CLR                    |
| CEE1270 | PUSH TO LOCK BRAKES   |
| CEE1275 | PULL TO LOCK BRAKES   |
| CEE1280 | ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)                  |
| CEE1345 | TRANSPORTATION BRACKETS   |
| CEE1361 | BODYPOINT 2-PT POSTURE BELT M   |
| CEE1362 | BODYPOINT 2-PT POSTURE BELT L   |
| CEE1363 | BODYPOINT 4-PT POSTURE BELT M   |
| CEE1364 | BODYPOINT 4-PT POSTURE BELT L   |
| CEE1370 | REAR ANTITIPPERS PAIR   |
| CEE1375 | REAR ANTITIPPERS PAIR WITH TIPPING AID                                |
| CEE1430 | CANE HOLDER   |
| CEE1440 | TIPPING AID RH  |
| CEE1441 | TIPPING AID LH  |
| CEE1480 | ACRYLIC TRAY  |
| CEE1483 | HALF TABLE TRAY RH  |

---

|         |                                     |
|---------|-------------------------------------|
| CEE1484 | HALF TABLE TRAY LH                  |
| CEE1485 | COMFORT PADDING FOR HALF TABLE TRAY |

---



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):