



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JUNE 2023
Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Action 5 Rigid

Standard Features

- Max user weight 130kg
- Crash tested ISO 7176 - Part 9
- Self Propulsion with 22" or 24" quick release rear wheels
- Choice of 7 seat widths 15" - 21" (380mm to 530mm)
- Choice of 7 seat depths 14" - 20" (350mm to 500mm)
- Choice of 8 front seat to floor heights 14" - 21" (360mm to 535mm)
- Choice of 7 rear seat to floor heights 14" - 20" (360mm to 510mm)
- Tension adjustable nylon seat upholstery
- Fixed backrest (90° to the seat rail)
- Fixed backrest, nylon and mesh padded cover
- Fixed front end or 80° swing in/swing out legrest hangers
- Flip up footplates
- sideguards / skirtguards
- Choice of 10 exciting colours and 2 component accents
- 2° Camber
- Push to lock brakes
- Transportation kit with lap belt
- Adjustable heel straps
- User tool kit



Standard Base price: £1,288

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

| Code | Description | | |
|----------------------|---------------------------------------|---|--------|
| BASE PRODUCT* | | | |
| SJE0016 | ACTION 5 RIGID | (NC 1510) | ○ |
| SJESPEC | SPECIAL | | POA |
| FRAME* | | | |
| SJE0200 | SWING AWAY FRAME | (NC 3230, 3440, 3450) | £1,288 |
| SJE0210 | FIXED FRONT END 80 DEGREE ANGLE FRAME | (NC 1111, 3010, 3020, 3030, 3040, 3050, 3110, 3120, 3130, 3140, 3150, 3333, 3400, 3410, 3420, 3430) | £1,288 |
| PROPULSION* | | | |
| SJE0300 | SELF PROPULSION | | ◆ |

| SEAT WIDTH* | | | |
|----------------------|--|---|---|
| SJE1030 | | 15" SW | ○ |
| SJE1040 | | 16" SW | ○ |
| SJE1050 | | 17" SW | ○ |
| SJE1060 | | 18" SW | ○ |
| SJE1070 | | 19" SW | ○ |
| SJE1075 | | 20" SW | (NC 1100, 1110, 1170) ○ |
| SJE1078 | | 21" SW | (NC 1100, 1110, 1170) ○ |
| SEAT DEPTH* | | | |
| SJE1100 | | 14" SD | (NC 1270, 3040, 3140, 7080, 7100, 7110, 7600, 7610) ○ |
| SJE1110 | | 15" SD | (NC 1270, 7600) ○ |
| SJE1120 | | 16" SD | (NC 1270) ○ |
| SJE1130 | | 17" SD | ○ |
| SJE1140 | | 18" SD | ○ |
| SJE1150 | | 19" SD | ○ |
| SJE1160 | | 20" SD | ○ |
| SJE1170 | | 14" - 16" ADJ SD | (NC 1510, 7600) ○ |
| SJE1180 | | 16" - 18" ADJ SD | (NC 1510) ○ |
| SJE1190 | | 18" - 20" ADJ SD | (NC 1510) ○ |
| FRONT SEAT HEIGHT* | | | |
| SJE1200 | | 14" FSH | (NC 1310, 1320, 1330, 1340, 1350, 1360, 3420, 3430, 3450, 7040, 7055, 7070, 7080, 7100, 7110) ○ |
| SJE1210 | | 15" FSH | (NC 1320, 1330, 1340, 1350, 1360, 3430, 7000, 7040, 7070, 7080, 7100, 7110) ○ |
| SJE1220 | | 16" FSH | (NC 1330, 1340, 1350, 1360, 7080, 7100, 7110) ○ |
| SJE1230 | | 17" FSH | (NC 1300, 1340, 1350, 1360, 7110) ○ |
| SJE1240 | | 18" FSH | (NC 1300, 1310, 1350, 1360) ○ |
| SJE1250 | | 19" FSH | (NC 1300, 1310, 1320, 1360) ○ |
| SJE1260 | | 20" FSH | (NC 1300, 1310, 1320, 1330, 7000, 7040) ○ |
| SJE1270 | | 21" FSH | (NC 1300, 1310, 1320, 1330, 1340, 3030, 3130, 7000, 7040, 7055, 7070, 7080, 7100) ○ |
| REAR SEAT HEIGHT* | | | |
| SJE1300 | | 14" RSH | (NC 7250, 7260) ○ |
| SJE1310 | | 15" RSH | ○ |
| SJE1320 | | 16" RSH | ○ |
| SJE1330 | | 17" RSH | ○ |
| SJE1340 | | 18" RSH | ○ |
| SJE1350 | | 19" RSH | (NC 7200) ○ |
| SJE1360 | | 20" RSH | (NC 7200, 7230, 7240) ○ |
| SEAT UPHOLSTERY* | | | |
| SJE1510 | | NYLON FIXED | ○ |
| SJE1520 | | TENSION AND DEPTH ADJUSTABLE | ○ |
| BACKREST ATTACHMENT* | | | |
| SJE2000 | | FIXED 90° | ○ |
| SJE2020 | | ANGLE ADJUSTABLE (-15° TO +15°) AND FOLDING | £100 |

| BACKREST FIXED HEIGHT - INCLUDING FIXED PUSH HANDLES | | | | |
|--|--|--|--|------|
| SJE2130 | | 12" BACK HT | (NC 2300, 2310, 2400, 2410, 2420, 2505, 8160, 8200) | ○ |
| SJE2140 | | 13" BACK HT | (NC 2300, 2310, 2400, 2410, 2420, 2505, 8160, 8200) | ○ |
| SJE2143 | | 14" BACK HT | (NC 2300, 2310, 2400, 2410, 2420, 2505, 8160, 8200) | ○ |
| SJE2150 | | 15" BACK HT | (NC 2300, 2310, 2400, 2410, 2420, 8200) | ○ |
| SJE2153 | | 16" BACK HT | (NC 2300, 2310, 2400, 2410, 2420) | ○ |
| SJE2160 | | 17" BACK HT | (NC 2300, 2310, 2400, 2410, 2420) | ○ |
| SJE2163 | | 18" BACK HT | (NC 2300, 2310, 2400, 2410, 2420) | ○ |
| SJE2170 | | 19" BACK HT | (NC 2300, 2310, 2400, 2410, 2420) | ○ |
| SJE2173 | | 20" BACK HT | (NC 2300, 2310, 2400, 2410, 2420) | ○ |
| BACKREST HEIGHT ADJUSTABLE - PUSH HANDLES MUST BE SELECTED IN THE SECTION BELOW | | | | |
| SJE2300 | | HT ADJ 14.5" - 16" | (NC 2505) | £25 |
| SJE2310 | | HT ADJ 16" - 19.5" | | £25 |
| PUSH HANDLES - MUST BE SELECTED IN COMBINATION WITH HEIGHT ADJUSTABLE BACKRESTS | | | | |
| SJE2400 | | W/O PUSH HANDLES | (NC 8160, 9100, 9200) | ○ |
| SJE2410 | | REAR SET | (NC 2505) | £49 |
| SJE2420 | | INTEGRATED | | £81 |
| BACKREST UPHOLSTERY* | | | | |
| SJE2500 | | W/O BACK UPH - DELIVERED WITH SIMPLE SLEEVE | | ○ |
| SJE2505 | | FIXED NYLON AND MESH | (NC 5600) | ◆ |
| SJE2510 | | CONTOUR ADJUSTABLE | | £74 |
| LEGREST HANGER RH* | | | | |
| SJE3000 | | W/O HANGER RH | (NC 5400) | ○ |
| SJE3010 | | 70° RH | (NC 1111, 3440, 3450, 22222) | ○ |
| SJE3020 | | 80° RH | (NC 1111, 3440, 3450, 22222) | ○ |
| SJE3030 | | 90° RH | (NC 1111, 3200, 3440, 3450, 7080, 7100, 7110, 22222) | ○ |
| SJE3040 | | ELEVATING LEGREST RH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GROUND IF UL IS MORE THAN 400MM | (NC 3230, 3400, 3440, 3450, 5400, 22222, 8120) | £116 |
| SJE3050 | | AMPUTEE LEGREST RH | (NC 1111, 3200, 3210, 3230, 5400, 8100, 8120) | £133 |
| LEGREST HANGER LH* | | | | |
| SJE3100 | | W/O HANGER LH | (NC 5400) | ○ |
| SJE3110 | | 70° LH | (NC 3333, 3440, 3450, 44444) | ○ |
| SJE3120 | | 80° LH | (NC 3333, 3440, 3450, 44444) | ○ |
| SJE3130 | | 90° LH | (NC 3300, 3333, 3440, 3450, 44444, 7080, 7100, 7110) | ○ |
| SJE3140 | | ELEVATING LEGREST LH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GROUND IF UL IS MORE THAN 400MM | (NC 3400, 3440, 3450, 44444, 5400, 8120) | £116 |
| SJE3150 | | AMPUTEE LEGREST LH | (NC 3300, 3310, 3333, 5400, 8100, 8120) | £133 |
| FOOTREST RH* | | | | |
| SJE3200 | | FLIP UP RH | (NC 3410) | ◆ |
| SJE3210 | | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS. | | £43 |
| SWO1111 | | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) | | £43 |
| SJE3230 | | ONE PIECE | (NC 3400, 3410, 3420, 3430) | £186 |
| SJE22222 | | NO FOOTPLATE REQUIRED | | ○ |

| FOOTREST LH* | | | |
|-----------------------|--|--|-----|
| SJE3300 | | FLIP UP LH (NC 3410) | ◆ |
| SJE3310 | | FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS. | £43 |
| SWO3333 | | FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) | £43 |
| SJE44444 | | NO FOOTPLATE REQUIRED | ○ |
| FOOTREST TUBE* | | | |
| SJE3400 | | SWING AWAY FRAME HIGH MOUNTED TUBE | ○ |
| SJE3410 | | SWING AWAY FRAME SHORT TUBE | ○ |
| SJE3420 | | SWING AWAY FRAME STANDARD TUBE | ○ |
| SJE3430 | | SWING AWAY FRAME LONG TUBE | ○ |
| SJE3440 | | FIXED FRAME STANDARD TUBE | ○ |
| SJE3450 | | FIXED FRAME LONG TUBE | ○ |
| MUDGUARDS/SIDEGUARDS* | | | |
| SJE3700 | | W/O SIDEGUARD/MUDGUARD | ○ |
| SJE3750 | | SIDEGUARDS (NC 4010, 4020, 4030, 4040, 4110, 4120, 4130, 4140) | ◆ |
| SJE3800 | | MUDGUARDS (NC 4010, 4020, 4030, 4040, 4110, 4120, 4130, 4140, 8160) | £74 |
| ARMREST RH* | | | |
| SJE4000 | | W/O ARMREST | ○ |
| SJE4010 | | T-ARM DESK ARMPAD RH (NC 7600, 7610) | £74 |
| SJE4020 | | T-ARM FULL ARMPAD RH | £74 |
| SJE4030 | | SWING AWAY DESK ARMPAD RH | £99 |
| SJE4040 | | SWING AWAY FULL ARMPAD RH | £99 |
| ARMREST LH* | | | |
| SJE4100 | | W/O ARMREST | ○ |
| SJE4110 | | T-ARM DESK ARMPAD LH (NC 7600, 7610) | £74 |
| SJE4120 | | T-ARM FULL ARMPAD LH | £74 |
| SJE4130 | | SWING AWAY DESK ARMPAD LH | £99 |
| SJE4140 | | SWING AWAY FULL ARMPAD LH | £99 |
| SEAT CUSHION* | | | |
| SJE55555 | | WITHOUT SEAT CUSHION | ○ |
| SJE5000 | | BASIC FOAM CUSHION | £21 |
| SEATING OPTIONS | | | |
| SJE5400 | | CALF STRAP | ○ |
| SJE5600 | | FIXED BACKBRACE | £21 |

| FRAME COLOUR* | | | |
|------------------|--|---|------|
| SJE6000 | | PEARL WHITE | ○ |
| SJE6010 | | ANTHRACITE BLACK | ○ |
| SJE6020 | | POP BLUE | ○ |
| SJE6030 | | DARK RED | ○ |
| SJE6040 | | ELECTRIC GREEN | ○ |
| SJE6050 | | HYPNOTIC PURPLE | ○ |
| SJE6060 | | PEARL GREY | ○ |
| SJE6070 | | MATT BLACK | ○ |
| SJE6080 | | INTENSE YELLOW | ○ |
| SJE6090 | | VIBRATING ORANGE | ○ |
| SJE6099 | | SPECIAL COLOUR | £169 |
| SWO10002 | | SLATE GREY - SPECIAL COLOUR 3 WEEK LEAD TIME | £169 |
| SWO10003 | | RED MATT - SPECIAL COLOUR 3 WEEK LEAD TIME | £169 |
| SWO10004 | | POLAR BLUE - SPECIAL COLOUR 3 WEEK LEAD TIME | £169 |
| SWO10005 | | CANDY PINK - SPECIAL COLOUR 3 WEEK LEAD TIME | £169 |
| ACCENT COLOUR* | | | |
| SJE6100 | | WITHOUT ACCENT COLOUR SILVER FORK AND REARWHEEL SUPPORT | ○ |
| SJE6110 | | BLUE ACCENT | ○ |
| SJE6200 | | PURPLE ACCENT | ○ |
| FRONT CASTORS* | | | |
| SJE7000 | | 4" X 1.1/4 SOFTROLL | £37 |
| SJE7040 | | 4" X 1.3/8" ALU SOFTROLL | £37 |
| SJE7055 | | 5" X 1" SOLID | ◆ |
| SJE7070 | | 5" X 1.3/8" ALU SOFTROLL | £37 |
| SJE7080 | | 6" X 1" SOLID | ○ |
| SJE7100 | | 6" X 1.3/8" ALU SOFTROLL | £37 |
| SJE7110 | | 7" X 1.1/4" SOLID | ○ |
| REAR WHEELS* | | | |
| SJE7200 | | 22" STANDARD, CROSS SPOKED WHEEL (NC 7320, 7410) | ○ |
| SJE7230 | | 24" STANDARD, CROSS SPOKED WHEEL (NC 7320) | ◆ |
| SJE7240 | | 24" LIGHTWEIGHT WHEEL (NC 7410, 8000, 8160) | £59 |
| SJE7250 | | 25" LIGHTWEIGHT WHEEL (NC 7302, 7304, 7315, 7410, 8000, 8160) | £59 |
| SJE7260 | | 25" SPINERGY WHEELS (NC 7302, 7304, 7315, 7410, 7420, 8000, 8160) | £561 |
| REAR WHEEL TYRE* | | | |
| SJE7302 | | SOLID BLACK 1.3/8" | ◆ |
| SJE7304 | | PNEUMATIC TYRE BLACK 1.3/8" | ○ |
| SJE7315 | | PNEUMATIC LOW PROFILE BLACK 1" | ○ |
| SJE7320 | | PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM | £70 |
| HANDRIM* | | | |
| SJE7400 | | ALUMINIUM HANDRIM | ◆ |
| SJE7410 | | STAINLESS STEEL | ○ |
| SJE7420 | | HANDRIM COATED IN BLACK ANTI-SLIP | £49 |

| REAR WHEEL AXLE* | | | |
|-----------------------------|--------------------------|---|------|
| SJE7500 | <input type="checkbox"/> | Q/R AXLE | ◆ |
| REAR WHEEL POSITION* | | | |
| SJE7600 | <input type="checkbox"/> | VERY ACTIVE | ○ |
| SJE7610 | <input type="checkbox"/> | ACTIVE | ○ |
| SJE7620 | <input type="checkbox"/> | MEDIUM ACTIVE | ○ |
| SJE7630 | <input type="checkbox"/> | STANDARD | ◆ |
| CAMBER ANGLE* | | | |
| SJE7700 | <input type="checkbox"/> | 0° CAMBER | ○ |
| SJE7710 | <input type="checkbox"/> | 1° CAMBER | ○ |
| SJE7720 | <input type="checkbox"/> | 2° CAMBER | ◆ |
| SJE7730 | <input type="checkbox"/> | 3° CAMBER | ○ |
| SJE7740 | <input type="checkbox"/> | 4° CAMBER | ○ |
| CHASSIS OPTION | | | |
| SJE8000 | <input type="checkbox"/> | SPOKE PROTECTORS | £81 |
| PUMP | | | |
| SJE8030 | <input type="checkbox"/> | PUMP - HIGH PRESSURE | £41 |
| BRAKES* | | | |
| SJE8100 | <input type="checkbox"/> | PUSH TO LOCK BRAKES | ◆ |
| SJE8110 | <input type="checkbox"/> | PULL TO LOCK BRAKES | ○ |
| SJE8120 | <input type="checkbox"/> | PUSH TO LOCK BRAKES WITH EXTENTION | £17 |
| SJE8130 | <input type="checkbox"/> | PULL TO LOCK BRAKES WITH EXTENTION | £17 |
| SJE8160 | <input type="checkbox"/> | ATTENDANT HUB BRAKES INCLUDES STANDARD USER PUSH TO LOCK BRAKE (OR PULL VERSION IF IT APPEARS AS THE ONLY CHOICE ABOVE) | £231 |
| TRANSPORTATION KIT | | | |
| SJE8200 | <input type="checkbox"/> | TRANSPORTATION BRACKETS | ◆ |
| SAFETY OPTIONS* | | | |
| SJE8300 | <input type="checkbox"/> | LAP BELT | ◆ |
| SJE8500 | <input type="checkbox"/> | REAR ANTI-TIPPER ACTIVE RH (NC 9200) | £64 |
| SJE8550 | <input type="checkbox"/> | REAR ANTI-TIPPER ACTIVE LH (NC 9100) | £64 |
| SJE8600 | <input type="checkbox"/> | REFLECTORS KIT (WHEELS) | £17 |
| ATTENDANT OPTION | | | |
| SJE9000 | <input type="checkbox"/> | CANE HOLDER | £23 |
| SJE9100 | <input type="checkbox"/> | TIPPING AID LH | £49 |
| SJE9200 | <input type="checkbox"/> | TIPPING AID RH | £49 |
| OTHER OPTIONS | | | |
| SWO10001 | <input type="checkbox"/> | ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618 | £107 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|----------|--|
| SJE1100 | 14" SD |
| SJE1110 | 15" SD |
| SJE1170 | 14" - 16" ADJ SD |
| SJE1270 | 21" FSH |
| SJE1300 | 14" RSH |
| SJE1310 | 15" RSH |
| SJE1320 | 16" RSH |
| SJE1330 | 17" RSH |
| SJE1340 | 18" RSH |
| SJE1350 | 19" RSH |
| SJE1360 | 20" RSH |
| SJE1510 | NYLON FIXED |
| SJE2300 | HT ADJ 14.5" - 16" |
| SJE2310 | HT ADJ 16" - 19.5" |
| SJE2400 | W/O PUSH HANDLES |
| SJE2410 | REAR SET |
| SJE2420 | INTEGRATED |
| SJE2505 | FIXED NYLON AND MESH |
| SJE3010 | 70° RH |
| SJE3020 | 80° RH |
| SJE3030 | 90° RH |
| SJE3040 | ELEVATING LEGREST RH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GROUND IF UL IS MORE THAN 400MM |
| SJE3050 | AMPUTEE LEGREST RH |
| SJE3110 | 70° LH |
| SJE3120 | 80° LH |
| SJE3130 | 90° LH |
| SJE3140 | ELEVATING LEGREST LH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GROUND IF UL IS MORE THAN 400MM |
| SJE3150 | AMPUTEE LEGREST LH |
| SJE3200 | FLIP UP RH |
| SJE3210 | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS. |
| SWO1111 | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) |
| SJE3230 | ONE PIECE |
| SJE22222 | NO FOOTPLATE REQUIRED |
| SJE3300 | FLIP UP LH |
| SJE3310 | FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS. |
| SWO3333 | FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) |
| SJE44444 | NO FOOTPLATE REQUIRED |
| SJE3400 | SWING AWAY FRAME HIGH MOUNTED TUBE |

| | |
|---------|---|
| SJE3410 | SWING AWAY FRAME SHORT TUBE |
| SJE3420 | SWING AWAY FRAME STANDARD TUBE |
| SJE3430 | SWING AWAY FRAME LONG TUBE |
| SJE3440 | FIXED FRAME STANDARD TUBE |
| SJE3450 | FIXED FRAME LONG TUBE |
| SJE4010 | T-ARM DESK ARMPAD RH |
| SJE4020 | T-ARM FULL ARMPAD RH |
| SJE4030 | SWING AWAY DESK ARMPAD RH |
| SJE4040 | SWING AWAY FULL ARMPAD RH |
| SJE4110 | T-ARM DESK ARMPAD LH |
| SJE4120 | T-ARM FULL ARMPAD LH |
| SJE4130 | SWING AWAY DESK ARMPAD LH |
| SJE4140 | SWING AWAY FULL ARMPAD LH |
| SJE5400 | CALF STRAP |
| SJE5600 | FIXED BACKBRACE |
| SJE7000 | 4" X 1.1/4 SOFTROLL |
| SJE7040 | 4" X 1.3/8" ALU SOFTROLL |
| SJE7055 | 5" X 1" SOLID |
| SJE7070 | 5" X 1.3/8" ALU SOFTROLL |
| SJE7080 | 6" X 1" SOLID |
| SJE7100 | 6" X 1.3/8" ALU SOFTROLL |
| SJE7110 | 7" X 1.1/4" SOLID |
| SJE7200 | 22" STANDARD, CROSS SPOKED WHEEL |
| SJE7230 | 24" STANDARD, CROSS SPOKED WHEEL |
| SJE7240 | 24" LIGHTWEIGHT WHEEL |
| SJE7250 | 25" LIGHTWEIGHT WHEEL |
| SJE7260 | 25" SPINERGY WHEELS |
| SJE7302 | SOLID BLACK 1.3/8" |
| SJE7304 | PNEUMATIC TYRE BLACK 1.3/8" |
| SJE7315 | PNEUMATIC LOW PROFILE BLACK 1" |
| SJE7320 | PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM |
| SJE7410 | STAINLESS STEEL |
| SJE7420 | HANDRIM COATED IN BLACK ANTI-SLIP |
| SJE7600 | VERY ACTIVE |
| SJE7610 | ACTIVE |
| SJE8000 | SPOKE PROTECTORS |
| SJE8100 | PUSH TO LOCK BRAKES |
| SJE8120 | PUSH TO LOCK BRAKES WITH EXTENTION |
| SJE8160 | ATTENDANT HUB BRAKES INCLUDES STANDARD USER PUSH TO LOCK BRAKE (OR PULL VERSION IF IT APPEARS AS THE ONLY CHOICE ABOVE) |
| SJE8200 | TRANSPORTATION BRACKETS |
| SJE9100 | TIPPING AID LH |
| SJE9200 | TIPPING AID RH |



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):