



Yes, you can:

Manual Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JUNE 2023

Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Clematis Pro

Standard Features

- Max user weight 135kg/21 Stone
- Crash tested ISO 7176 Part 19
- Self propulsion, 24" rear wheels
- Choice of 3 seat widths 390mm, 440mm and 490mm (+25mm width adj)
- Seat height 450mm
- Depth adjustable seat 420mm - 480mm
- Backrest plate, Laguna 2 cushion with height adjustable push handles
- Flo-shape seat cushion, Multi stretch PU coated upholstery TR26
- Assistant operated gas piston 25 Deg seat tilt and 30 Deg backrest recline
- Angle adjustable legrests with calf support
- Height adjustable standard pad armrests (including autolock)
- Pearl beige frame
- Quick release self propel wheels, solid tyres and aluminium handrims
- Anti tip devices
- Headrest, Pelvic belt, User manoeuvred tyre brakes



Standard Base price: £1,813

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
SNA0001	CLEMATIS PRO	£1,813
SNASPEC	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	○
PROPULSION*		
SNA0100	SELF PROPULSION (24" REAR WHEELS)	◆
FRAME*		
SNA0200	FIXED FRAME WITH 25MM WIDTH ADJUSTMENT EACH SIDE	◆
SEAT WIDTH*		
SNA0300	390 MM	○
SNA0310	440 MM	○
SNA0320	490 MM	○
SEAT DEPTH*		
SNA0360	420 - 480 MM	◆
SEAT HEIGHT*		
SNA0390	450 MM	◆
SEAT CUSHION*		
SNA0430	CUSHION FLO-SHAPE BLACK MULTI STRETCH PU COATED TR26	◆
BACKREST*		
SNA0500	RIGID BACKREST PLATE	◆
BACKREST CUSHION*		
SNA0530	CUSHION LAGUNA 2 BLACK MULTI STRETCH PU COATED FABRIC	◆
PUSH HANDLES*		
SNA0560	HEIGHT ADJUSTABLE PUSH HANDLES	◆
SNA0570	PUSH BAR ANGLE ADJUSTABLE	£102
SEAT TILT AND BACKREST RECLINE*		
SNA0600	ASSISTANT OPERATED GAS PISTON 30 DEG BACKREST RECLINE AND 25 DEG SEAT TILT	◆
TILT LOCK*		
SNA0630	TILT LOCK (LOCKING OF PRE-SET ANGLE)	◆
LEGREST HANGER ATTACHMENT*		
SNA0700	STANDARD LEGREST ATTACHMENT	◆
LEGREST HANGER LH*		
SNA0720	LEGREST ANGLE ADJUSTABLE, INCL CALF SUPPORT, LH	◆
LEGREST HANGER RH*		
SNA0740	LEGREST ANGLE ADJUSTABLE, INCL CALF SUPPORT, RH	◆
FOOTREST LH*		
SNA0770	FOOTREST LOWER PART, DEPTH, WIDTH AND ANGLE ADJUSTABLE WITH HEEL STRAP, LH	◆
FOOTREST RH*		
SNA0790	FOOTREST LOWER PART, DEPTH, WIDTH AND ANGLE ADJUSTABLE WITH HEEL STRAP, RH	◆

ARMREST LH*			
SNA0800	<input type="checkbox"/>	DUCKTAIL PAD, LH DEPTH ADJUSTABLE	○
SNA0810	<input type="checkbox"/>	CONCAVE ARMREST PAD, LH DEPTH ADJUSTABLE	◆
ARMREST RH*			
SNA0850	<input type="checkbox"/>	DUCKTAIL PAD, RH DEPTH ADJUSTABLE	○
SNA0860	<input type="checkbox"/>	CONCAVE ARMREST PAD, RH DEPTH ADJUSTABLE	◆
SEAT OPTION			
SNA0900	<input type="checkbox"/>	SIDE SUPPORT COVER	£51
SNA0910	<input type="checkbox"/>	POMMEL ABDUCTION CUSHION	£118
ARMREST SPACER*			
SNA0920	<input type="checkbox"/>	ARMREST SPACER, INCREASES WIDTH TO A TOTAL +25 MM PER SIDE	◆
BACKREST OPTIONS			
SNA1100	<input type="checkbox"/>	HEADREST ADJUSTABLE WITH HANDWHEEL	◆
SNA1110	<input type="checkbox"/>	HEADREST WITH CHEEK SUPPORT, SIDEWISE ADJUSTABLE (NC 1130)	£107
SNA1120	<input type="checkbox"/>	NECKREST FIXED (NC 1130)	£31
SNA1130	<input type="checkbox"/>	COVER FOR HEADREST	£44
SNA1140	<input type="checkbox"/>	COVER FOR NECKREST	£22
TRUNK SUPPORT LH			
SNA1200	<input type="checkbox"/>	TRUNK SUPPORT SMALL, LH	£44
SNA1210	<input type="checkbox"/>	TRUNK SUPPORT SWINGAWAY (LH)	£48
TRUNK SUPPORT RH			
SNA1220	<input type="checkbox"/>	TRUNK SUPPORT SMALL, RH	£44
SNA1230	<input type="checkbox"/>	TRUNK SUPPORT SWINGAWAY (RH)	£48
FOOTPLATE PAD			
SNA1300	<input type="checkbox"/>	FOOTPLATE PADDING	£54
CHASSIS COLOUR*			
SNA1400	<input type="checkbox"/>	PEARL BEIGE	◆
CASTOR ATTACHMENT*			
SNA1410	<input type="checkbox"/>	CASTOR ATTACHMENT QUICK RELEASE	◆
FRONT CASTORS*			
SNA1440	<input type="checkbox"/>	200X45 MM BLACK SOLID TYRE	◆
REAR WHEEL TYPE*			
SNA1510	<input type="checkbox"/>	24" (610 MM) REAR WHEELS QR ATTACHMENT	◆
SNA1525	<input type="checkbox"/>	24" (610 MM) DRUM BRAKE WHEELS QR ATTACHMENT	£133
REAR WHEEL TYRE*			
SNA1540	<input type="checkbox"/>	BLACK TYRES SOLID	◆
HANDRIM*			
SNA1550	<input type="checkbox"/>	HANDRIM ALUMINIUM	◆
SPOKEGUARD			
SNA1560	<input type="checkbox"/>	FROST	£166
BRAKE*			
SNA1600	<input type="checkbox"/>	OCCUPANT OPERATED BRAKES	◆
SAFETY OPTIONS*			
SNA1650	<input type="checkbox"/>	POSTURE BELT	◆
SNA1660	<input type="checkbox"/>	ANTI-TIP DEVICE INCLUDED IN THE CHASSIS	◆

TRANSPORTATION KIT*

SNA1680		TRANSPORTATION KIT (WITHOUT POSTURE BELT)	◆
---------	--	---	---

OTHER OPTIONS

SNA1700		TABLE TRAY INCL MOUNTING KIT, WITH ELBOW PADS (NC 1800)	£138
SNA1710		TABLE TRAY INCL MOUNTING KIT (NC 1800)	£138
SNA1740		TABLE TRAY CUSHION (TO BE ORDERED IN COMBINATION WITH THE FIXED & SWINGAWAY TABLE TRAY)	£51
SNA1790		LOCKING DEVICE FOR TABLE TRAY	£15
SNA1800		MOUNTING KIT TABLE (IF TABLE ORDERED SEPARATELY)	£9

TOTAL PRICE			
--------------------	--	--	--

Non Conforming Options:

SNA1130 COVER FOR HEADREST

SNA1800 MOUNTING KIT TABLE (IF TABLE ORDERED SEPARATELY)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):