



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JUNE 2023

Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Clematis Pro Stock

Standard Features

- Max user weight: 135kg (21.5 stone)
- Crash tested ISO 7176-Part 19
- 12" Transit wheels black solid tyres
- Choice of 3 widths 390mm, 440mm and 490mm with +25mm spacer on armrests
- Seat depth, stepless adjustable 420mm-480mm
- Seat height 450mm
- Floshape seat cushion with multi stretch PU coated fabric TR26
- Plastic rigid backrest plate with Swingaway trunk supports
- Laguna back cushion with multi stretch PU coated fabric TR26
- Height adjustable push handles
- Assistant operated gas strut 30 Deg backrest recline and 25 Deg seat tilt with tilt lock
- Angle adjustable legrests with depth and angle adjustable footplates and heelstraps
- Concave armrests with armrest spacer adds additional 25mm of seat space
- Headrest, adjustable with handwheel
- Pearl Beige frame colour
- Quick release castor attachment with 200mm x 45mm medium hard black castors
- Foot operated brakes, Pelvic belt, Flip up Anti tip devices
- Transportation kit



Standard Base price: £1,650

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
Rea Clematis Pro*		
DV1662394	W 390mm (Small chassis) x D 430-500mm (Depth adjustable)	£1,650
DV1662395	W 440mm (Medium chassis) x D 430-500mm (Depth adjustable)	£1,650
DV1662439	W 490mm (Large chassis) x D 430-500mm (Depth adjustable)	£1,650

TOTAL PRICE	
--------------------	--