



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JUNE 2023
Customer Ref:

Customer Ref*, Contact No*, Delivery Address*, Account No*, Order No*, Order Qty*, Order date*

Storm 4 Max
Standard Features

- Max User Weight 220KG (34.5 Stone)
Max Seating
Seat Width from 490mm - 610mm
Seat Depths 560mm - 630mm
multi stretch PU coated fabric Upholstery
Seat Height 480mm - 510mm
Manual Settable Tilt with Manual Settable Recline
Reinforced Armrests
Reinforced Legrests with Angle adjustable flip-back footboard
Headrest, sideways and depth adjustable
LiNX Control REM110 with standard bracket
Choice of 6 Colours
Puncture Proof Castors
Puncture Proof Drive Wheels
6 km/h (4mph) Motors
8A (200 - 250V) Charger
4 Point tie down system
73.5 Ah (Gel) batteries
Bolted on anti tippers

Standard Base price: £8,023



Please make your selection with a cross in the box [] Free of charge [O] standard feature [◆] POA special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include CHAIR TYPE* (AWA0000, AWZ0000), SEATING TYPE* (AWA1145), and BACK CUSHION* (AWA1106).

SEAT WIDTH*			
AWA0378	<input type="checkbox"/>	Seat Width 550mm	○
AWA0380	<input type="checkbox"/>	Seat Width 610mm (NC 0503)	○
SEAT CUSHION*			
AWA1459	<input type="checkbox"/>	Seat cushion Max (550mm x 635mm)	○
AWA1460	<input type="checkbox"/>	Seat cushion Max (610mm x 635mm)	○
SEAT DEPTH SETTABLE IN 3 STEPS*			
AWA0396	<input type="checkbox"/>	Seat Depth 560mm	○
AWA0397	<input type="checkbox"/>	Seat Depth 595mm	○
ADAPTER TO REDUCE SEAT DEPTH			
AWA0399	<input type="checkbox"/>	Adapter to shorten the seat depth by 130mm	£51
UPHOLSTERY TYPE*			
AWA0278	<input type="checkbox"/>	multi stretch PU coated fabric	◆
ARMRESTS INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*			
AWA0635	<input type="checkbox"/>	Fixed armrest 250mm - 340mm height adjustable removable	○
ARMPAD*			
AWA0647	<input type="checkbox"/>	Standard Armpad 325mm	○
AWA0649	<input type="checkbox"/>	Armpad extended plus 50mm	○
PUSH BAR*			
AWA1307	<input type="checkbox"/>	Push bar - height adjustable	○
HEADREST / NECKREST - Adequate headrest option must be selected*			
AWA1432	<input type="checkbox"/>	Headrest sideways and depth adjustable	◆
AWA1427	<input type="checkbox"/>	Headrest with cheek support sideways adjustable	£205
AWA1495	<input type="checkbox"/>	Equivalent non Invacare Headrest supporting the requirements of ISO7176-19:2001 Provided independently from this order	○
SEAT HEIGHT*			
AWA1025	<input type="checkbox"/>	Fixed seat height 480mm	○
AWA1006	<input type="checkbox"/>	Fixed seat height 510mm	○
SEAT TILT*			
AWA1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	○
AWA1014	<input type="checkbox"/>	Electrical Tilt 0° to 15°	£540
BACKREST RECLINER*			
AWA0958	<input type="checkbox"/>	Settable recliner by spindle (0° - 25°)	○
AWA0947	<input type="checkbox"/>	Electrical recliner (with high pivot, shearing compensation)	£477
FOOTREST HANGER TYPE (PAIR)*			
AWA0564	<input type="checkbox"/>	Reinforced hanger fixed angle (not detachable without tools) (NC 0522)	○
AWA0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm (NC 0599)	£107
CALF PADS			
AWA0583	<input type="checkbox"/>	Calf Pads adjustable	○
FOOTBOARD*			
AWA0599	<input type="checkbox"/>	Angle adjustable flip back heavy duty footboard (NC 0503)	◆
AWA0522	<input type="checkbox"/>	Angle adjustable footboard (for max user weight 200kg) (NC 0564)	○
CONTROLS*			
AWA2201	<input type="checkbox"/>	REM110 Drive only (NC 1437)	◆
AWA2203	<input type="checkbox"/>	REM211 Drive & Actuators (NC 1437)	○
AWA2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	£137

CONTROL BRACKETRY*			
AWA1207		Standard Bracket RH	○
AWA1208		Standard Bracket LH	○
AWA1707		Swing Away RH Bracket	£119
AWA1708		Swing Away LH Bracket	£119
AWA2227		Bracket RH mounted for Attendant (NC 2219, 2220)	£57
AWA2228		Bracket LH mounted for Attendant (NC 2219, 2220)	£57
DUAL CONTROL			
AWA2220		REM050 Dual Control (Please also select user control and bracket) (Please select dual control bracket) (NC 2227, 2228, 2245)	£355
AWA2219		Compact Dual Control (ACU200) (Please select dual control bracket) (NC 2227, 2228, 2245, 2201)	£355
DUAL CONTROL BRACKETEY			
AWA1207		Dual Control Standard Bracket RH (NC 2227, 2228, 2245)	○
AWA1208		Dual Control Standard Bracket LH	○
INVACARE LiNX IDC DUAL CONTROL			
AWA2245		Invacare LiNX IDC Dual control (Please select IDC dual control bracket below) (NC 2227, 2228, 2220, 2219, 2201)	£673
INVACARE LiNX IDC DUAL CONTROL BRACKETRY			
AWA1221		IDC Dual control Detachable Bracket RH (NC 2227, 2228, 2220, 2219)	○
AWA1222		LiNX IDC Dual control Detachable Bracket LH (NC 2227, 2228, 2220, 2219)	○
LINX 10-WAY SWITCH			
AWA2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch)	£453
LINX 10-WAY SWITCH BRACKET			
AWA2217		10-Way Switch Bracket RH, Mounted for User	○
AWA2218		10-Way Switch Bracket LH, Mounted for User	○
AWA2227		10-Way Switch Bracket RH, Mounted for Attendant	£57
AWA2228		10-Way Switch Bracket LH, Mounted for Attendant	£57
LINX ACCESSORY			
AWA2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	£90
FURTHER LiNX COMPONENTS / ACCESSORIES			
AWA2300		To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	○
BUDDY BUTTONS			
AWA1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
AWA1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
AWA1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
AWA1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
EGG SWITCH			
AWA2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79

PIKO BUTTONS			
AWA1711	<input type="checkbox"/>	Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
AWA1712	<input type="checkbox"/>	Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
AWA1713	<input type="checkbox"/>	Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
AWA1714	<input type="checkbox"/>	Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
AWA1715	<input type="checkbox"/>	Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
AWA1716	<input type="checkbox"/>	Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
TOGGLES / BUTTONS			
AWZ2261	<input type="checkbox"/>	Stereo Toggle - direct control for Tilt only	£242
AWZ2262	<input type="checkbox"/>	Stereo Button - direct control for Tilt only	£230
AWZ2263	<input type="checkbox"/>	Four-Way Toggle - direct control for Tilt & Recline	£300
AWZ2264	<input type="checkbox"/>	Four-Way Button - direct control for Tilt & Recline	£321
AWZ2207	<input type="checkbox"/>	Toggle / Button RH Mount	○
AWZ2208	<input type="checkbox"/>	Toggle / Button LH Mount	○
AWZ2225	<input type="checkbox"/>	Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
AWA2721	<input type="checkbox"/>	3" U-Shaped, Joystick Handle, Bodypoint PC101A	£65
AWA2722	<input type="checkbox"/>	4" U-Shaped, Joystick Handle, Bodypoint PC102A	£65
AWA2723	<input type="checkbox"/>	Domed Rubber, Joystick Handle, Bodypoint PC107A	£65
AWA2724	<input type="checkbox"/>	Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£65
REMOTE STOP SWITCH			
AWZ2258	<input type="checkbox"/>	Remote Stop Switch (ASL504)	£942
SHROUD COLOUR*			
AWA0202	<input type="checkbox"/>	Blue	○
AWA0228	<input type="checkbox"/>	Black	○
AWA0220	<input type="checkbox"/>	Champagne	○
AWA0204	<input type="checkbox"/>	Transparent Cherry Red	○
AWZ0299	<input type="checkbox"/>	Special Shroud Colour (RAL Code must be supplied)	POA
FRONT WHEEL*			
AWA0840	<input type="checkbox"/>	10" Grey Pneumatic	○
AWA0858	<input type="checkbox"/>	10" Grey Puncture Proof Sealed	○
AWA0842	<input type="checkbox"/>	10" Grey Puncture Proof Solid Inlay	◆
REAR WHEEL*			
AWA0771	<input type="checkbox"/>	14" Grey Pneumatic	○
AWA0773	<input type="checkbox"/>	14" Grey Puncture Proof Sealed	○
AWA0750	<input type="checkbox"/>	14" Grey Puncture Proof Solid Inlay	◆
MOTOR TYPE*			
AWA1277	<input type="checkbox"/>	6 km/h	◆
BATTERY*			
AWA1364	<input type="checkbox"/>	73.5 AH Gel	◆
SUSPENSION*			
AWA1203	<input type="checkbox"/>	Heavy Duty Suspension	◆

PERFORMANCE ACCESSORIES			
AWA1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£738
CHARGER*			
AWA1343		OptiCharge Charger 8A (200 - 250V)	◆
AWA1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£137
LIGHTS*			
AWA1254		No Lights	◆
AWA1437		LED Light system with indicators	£431
TRANSPORTATION KIT*			
AWA1494		4 Point Tie Down (Tested version for using wheelchair as a vehicle seat according to ISO7176-19:2001) Adequate headrest is compulsory	◆
FIXATION BELT*			
AWA1455		Fixation Belt Extra long with Auto Lock	◆
REAR VIEW MIRROR			
AWA1415		Rear view Mirror RH - Must be selected with 10 km/h Motor (NC 2217)	£38
AWA1339		Rear view Mirror LH - Must be selected with 10 km/h Motor (NC 2218)	£38
ACCESSORIES			
AWA1312		Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	£38
AWA1363		Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	£38
AWA1417		Fenders for rear wheels	£109

TOTAL PRICE	
--------------------	--

Non Conforming Options:

AWA0564	Reinforced hanger fixed angle (not detachable without tools)
AWA0503	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm
AWA0599	Angle adjustable flip back heavy duty footboard
AWA0522	Angle adjustable footboard (for max user weight 200kg)
AWA2201	REM110 Drive only
AWA2227	Bracket RH mounted for Attendant
AWA2228	Bracket LH mounted for Attendant
AWA2220	REM050 Dual Control (Please also select user control and bracket) (Please select dual control bracket)
AWA2219	Compact Dual Control (ACU200) (Please select dual control bracket)
AWA2245	Invacare LiNX IDC Dual control (Please select IDC dual control bracket below)
AWA2217	10-Way Switch Bracket RH, Mounted for User
AWA2218	10-Way Switch Bracket LH, Mounted for User
AWA2227	10-Way Switch Bracket RH, Mounted for Attendant
AWA2228	10-Way Switch Bracket LH, Mounted for Attendant
AWA1437	LED Light system with indicators



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):