



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JAN-MAR 2024
Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

AVIVA RX40 Modulite M27 707015 - Weekly Price £59.18

Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flo-tech lite cushion
- Flex3 Laguna backrest
- Standard T-Armrests, detachable with Long armpad
- Neckrest depth adjustable
- Seat lifter with powered tilt (300 mm, 30°) with CoG shift
- Manual settable backrest Recline
- Standard hanger 80° pair with Fixed Footplates
- LiNX control REM216 with swing away bracket
- 9" Castor, puncture proof solid tyre
- Solid inlay Dual density foam insert
- 12 km/h Ultra High Speed Package with LED Light system with Indicators, 30 Lux headlights & Mirror
- 73Ah batteries
- Anti tippers
- 6-Point tie-down system
- Lapbelt with autolock (adjustable both sides)



Standard Base price: £9,804

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

| Code | Description | |
|-----------------------|---|--------|
| BASE MODEL* | | |
| JBG001 | AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries | £9,804 |
| CHASSIS WIDTH* | | |
| JBL2810 | Chassis Narrow (560 mm/22") | ◆ |
| SEATING TYPE* | | |
| JBG2100 | Invacare® Modulite Seat Unit | ◆ |

| | | | |
|---|--------------------------|---|---|
| SEAT FRAME* | | | |
| JBG2111 | <input type="checkbox"/> | Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms | ◆ |
| SEAT WIDTH* | | | |
| JBG2121 | <input type="checkbox"/> | 380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests | ○ |
| JBG2122 | <input type="checkbox"/> | 430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests | ○ |
| JBG2123 | <input type="checkbox"/> | 480 mm Seat Plate | ○ |
| JBG2124 | <input type="checkbox"/> | 530 mm Seat Plate | ○ |
| SEAT DEPTH* | | | |
| JBG2132 | <input type="checkbox"/> | Seat Depth 410 mm | ○ |
| JBG2134 | <input type="checkbox"/> | Seat Depth 460 mm | ○ |
| JBG2136 | <input type="checkbox"/> | Seat Depth 510 mm | ○ |
| SEAT CUSHION* | | | |
| JBG2142 | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite | ◆ |
| BACKREST SUPPORT* | | | |
| JBG2162 | <input type="checkbox"/> | Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket | ◆ |
| PUSH HANDLES FOR BACKPOSTS | | | |
| JBG2108 | <input type="checkbox"/> | Backposts without push handles | ◆ |
| BACKREST UPHOLSTERY / CUSHION* | | | |
| JBG2173 | <input type="checkbox"/> | Flex3 adjustable Backrest Plates with Laguna Backrest cushion | ◆ |
| UPHOLSTERY TYPE* | | | |
| JBG2152 | <input type="checkbox"/> | Black multi stretch PU coated fabric Upholstery | ○ |
| JBG2155 | <input type="checkbox"/> | Black Spacer Fabric Upholstery | ○ |
| ARMREST (Pair)* | | | |
| JBG2191 | <input type="checkbox"/> | T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm | ◆ |
| ARMPAD (Pair)* | | | |
| JBG2196 | <input type="checkbox"/> | Long Armpad (375 mm) | ◆ |
| HEADREST / NECKREST - Adequate headrest option must be selected* | | | |
| JBG1350 | <input type="checkbox"/> | Neckrest depth adjustable | ◆ |
| SEAT LIFTER / SEAT TILT* | | | |
| JBG1094 | <input type="checkbox"/> | Seat lifter with powered tilt (300 mm, 30°) with CoG shift | ◆ |
| SEAT TO GROUND HEIGHT* | | | |
| JBG2857 | <input type="checkbox"/> | Seat-to-Floor Height 435 mm | ○ |
| JBG2858 | <input type="checkbox"/> | Seat-to-Floor Height 460 mm | ○ |
| JBG2859 | <input type="checkbox"/> | Seat-to-Floor Height 485 mm | ○ |
| BACKREST RECLINER* | | | |
| JBG0975 | <input type="checkbox"/> | Manual Settable Recliner (in 3.75° steps) 90° - 120° | ◆ |
| FOOTREST HANGER TYPE* | | | |
| JBG0526 | <input type="checkbox"/> | Standard hanger 80° pair Lower leg length: 350 - 450 mm | ◆ |
| FOOTPLATE TYPE* | | | |
| JBG0512 | <input type="checkbox"/> | Fixed Footplate pair | ◆ |
| CONTROLS* | | | |
| JBG2204 | <input type="checkbox"/> | DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights | ◆ |

| CONTROL BRACKETRY* | | | |
|----------------------------|--------------------------|---|----------------------------------|
| JBG1707 | <input type="checkbox"/> | Swing Away Bracket RH | <input type="radio"/> |
| JBG1708 | <input type="checkbox"/> | Swing Away Bracket LH | <input type="radio"/> |
| SHROUD COLOUR* | | | |
| JBG2871 | <input type="checkbox"/> | Mid range - Black/Anthracite Shroud | <input type="radio"/> |
| JBG2872 | <input type="checkbox"/> | Mid range - Silver Shroud | <input type="radio"/> |
| JBG2873 | <input type="checkbox"/> | Mid range - Bengal Red Shroud | <input type="radio"/> |
| JBG2874 | <input type="checkbox"/> | Mid range - Mercury Blue Shroud | <input type="radio"/> |
| JBG2875 | <input type="checkbox"/> | Mid range - Midnight Purple Shroud | <input type="radio"/> |
| RIM COLOUR* | | | |
| JBG2882 | <input type="checkbox"/> | Mid range - Black Rim Inserts | <input type="radio"/> |
| JBG2883 | <input type="checkbox"/> | Mid range - Silver Rim Inserts | <input type="radio"/> |
| JBG2884 | <input type="checkbox"/> | Mid range - Bengal Red Rim Inserts | <input type="radio"/> |
| JBG2885 | <input type="checkbox"/> | Mid range - Mercury Blue Rim Inserts | <input type="radio"/> |
| JBG2886 | <input type="checkbox"/> | Mid range - Midnight Purple Rim Inserts | <input type="radio"/> |
| FORKS* | | | |
| JBG2632 | <input type="checkbox"/> | Double sided fork Without castor rim insert (enlarges wide base to 656 mm) | <input checked="" type="radio"/> |
| FRONT WHEEL* | | | |
| JBG0836 | <input type="checkbox"/> | 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre | <input checked="" type="radio"/> |
| DRIVE WHEEL* | | | |
| JBG0705 | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert | <input checked="" type="radio"/> |
| C.T.C SUSPENSION* | | | |
| JBG2870 | <input type="checkbox"/> | CTC Suspension Control - Traction - Comfort | <input checked="" type="radio"/> |
| SUSPENSION SETTING* | | | |
| JBG1910 | <input type="checkbox"/> | Soft (suggested for up to 70 kg User Weight) (NC 1094) | <input type="radio"/> |
| JBG1912 | <input type="checkbox"/> | Medium (suggested for up to 90 kg User Weight) (NC 1094) | <input type="radio"/> |
| JBG1911 | <input type="checkbox"/> | Firm (suggested for up to 110 kg User Weight) | <input type="radio"/> |
| JBG1915 | <input type="checkbox"/> | Very firm (suggested for up to 136 kg User Weight) | <input type="radio"/> |
| MOTOR TYPE* | | | |
| JBG1909 | <input type="checkbox"/> | 12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights | <input checked="" type="radio"/> |
| BATTERY TYPE* | | | |
| JBG1260 | <input type="checkbox"/> | 73,5 Ah Gel | <input checked="" type="radio"/> |
| CHARGER* | | | |
| JBG1343 | <input type="checkbox"/> | Opti Charge Charger 8A (IP44) | <input checked="" type="radio"/> |
| ANTI-TIPPERS* | | | |
| JBG1205 | <input type="checkbox"/> | Anti-tipper | <input checked="" type="radio"/> |
| TRANSPORTATION KIT* | | | |
| JBG1401 | <input type="checkbox"/> | 6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory | <input checked="" type="radio"/> |
| LAP BELT* | | | |
| JBG2113 | <input type="checkbox"/> | Lap belt with Auto Lock Flexible mount in side rails | <input checked="" type="radio"/> |
| REAR VIEW MIRROR | | | |
| JBG1415 | <input type="checkbox"/> | Rear view mirror RH | <input type="radio"/> |
| JBG1339 | <input type="checkbox"/> | Rear view mirror LH | <input type="radio"/> |
| TOTAL PRICE | | | <input type="text"/> |

Non Conforming Options:

JBG1094 Seat lifter with powered tilt (300 mm, 30°) with CoG shift