



Yes, you can:

# Motability Retail Prescription Form

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

JAN-MAR2024  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## E-Fix & ACTION 3NG M82 707362 - Weekly Price £46.32

### Standard Features

- Action 3NG + E-Fix
- Max user weight 125kg (19.7 stone)
- Crash tested ISO- 7176-Part 19
- Seat widths 380mm to 505mm (15" -20")
- Seat depths 400mm or 450mm (16" or 18")
- Seat heights 460mm to 510mm (18" -20")
- Nylon with hook and loop strips seat upholstery
- Fixed 7°Backrest 2 heights 400mm/430mm
- Lumbar support padded backrest with shaped cover
- 80°swing in/swing out hangers
- Flip up footplates
- Removable full arms
- Basic 50mm foam cushion
- Adjustable calf strap
- Choice of 8 standard colours
- 8" Black solid front wheels, with 24" Black solid tyre spoked rim
- Aluminium Pre-anodized handrim
- Quick release axle
- Push to lock brakes
- Transportation kit (including buckle belt)
- **Alber Efix E35 1566293 (24" wheels)**
- Brackets
- Alber anti-tippers
- Adapter Kit
- Mobility Plus package
- Fitting at Invacare



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
CEE0003	ACTION 3NG SELF PROPEL	(NC 1195) ◆
<b>PROPULSION*</b>		
CEE0100	SELF PROPULSION	◆
CEE9259	E-Fix E36	£6,086

<b>SEAT WIDTH*</b>			
CEE0180	<input type="checkbox"/>	15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0190	<input type="checkbox"/>	16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0200	<input type="checkbox"/>	17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	◆
CEE0210	<input type="checkbox"/>	18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0220	<input type="checkbox"/>	19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0230	<input type="checkbox"/>	20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
<b>SEAT DEPTH*</b>			
CEE0240	<input type="checkbox"/>	16" / 400 MM SEAT DEPTH (NC 0879)	○
CEE0250	<input type="checkbox"/>	18" / 450 MM SEAT DEPTH (NC 0878)	◆
<b>SEAT HEIGHT*</b>			
CEE0260	<input type="checkbox"/>	20" / 510 MM SEAT HEIGHT	○
CEE0270	<input type="checkbox"/>	19" / 485 MM SEAT HEIGHT	○
CEE0280	<input type="checkbox"/>	18" / 460 MM SEAT HEIGHT	◆
<b>SEAT UPHOLSTERY*</b>			
CEE0330	<input type="checkbox"/>	NYLON W/O HOOK AND LOOP STRIPS	◆
<b>BACKREST*</b>			
CEE0410	<input type="checkbox"/>	FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES	◆
<b>BACKREST UPHOLSTERY*</b>			
CEE0474	<input type="checkbox"/>	LUMBAR SUPPORT BI MATERIAL PADDED COVER	◆
<b>LEGREST HANGER RH*</b>			
CEE0510	<input type="checkbox"/>	80° RH SWING IN/SWING OUT	◆
<b>LEGREST HANGER LH*</b>			
CEE0560	<input type="checkbox"/>	80° LH SWING IN/SWING OUT	◆
<b>FOOTREST RH*</b>			
CEE0600	<input type="checkbox"/>	FLIP UP RH FOOTPLATE	◆
<b>FOOTREST LH*</b>			
CEE0620	<input type="checkbox"/>	FLIP UP LH FOOTPLATE	◆
<b>ARMREST RH*</b>			
CEE0640	<input type="checkbox"/>	REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 0760, 0839)	○
CEE0650	<input type="checkbox"/>	REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 0760, 0839)	○
CEE0660	<input type="checkbox"/>	SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT (NC 0740, 0750, 0839)	◆
<b>ARMREST LH*</b>			
CEE0740	<input type="checkbox"/>	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEE0750	<input type="checkbox"/>	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEE0760	<input type="checkbox"/>	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT	◆
CEE0839	<input type="checkbox"/>	TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) (NC 1200, 1345)	○

<b>SEAT CUSHION</b>			
CEE0849	<input type="checkbox"/>	BASIC 50 MM FOAM CUSHION NYLON AND ANTISLIP COVER	◆
<b>SEAT CUSHION DEPTH</b>			
CEE0878	<input type="checkbox"/>	FOR 400 MM SEAT DEPTH	○
CEE0879	<input type="checkbox"/>	FOR 450 MM SEAT DEPTH	○
<b>COLOUR*</b>			
CEE0961	<input type="checkbox"/>	SLATE GREY (RAL7024 FINE TEXTURE)	○
CEE0964	<input type="checkbox"/>	PEARL WHITE	○
CEE0965	<input type="checkbox"/>	ELECTRIC GREEN	○
CEE0966	<input type="checkbox"/>	MATT BLACK	○
CEE0967	<input type="checkbox"/>	BLUE POP	○
CEE0968	<input type="checkbox"/>	HYPNOTIC PURPLE	○
CEE0969	<input type="checkbox"/>	RED MATT	○
CEE0972	<input type="checkbox"/>	WATER BLUE	○
CEE0980	<input type="checkbox"/>	ANTHRACITE BLACK	○
<b>CASTOR HOUSING*</b>			
CEE0995	<input type="checkbox"/>	STANDARD	◆
<b>FRONT CASTORS*</b>			
CEE1010	<input type="checkbox"/>	8" X 1.1/4" BLACK SOLID TYRE	◆
CEE1030	<input type="checkbox"/>	8" X 2" BLACK SOLID TYRE	○
<b>REAR WHEEL SIZE*</b>			
CEE1068	<input type="checkbox"/>	24"	◆
<b>REAR WHEELS*</b>			
CEE1071	<input type="checkbox"/>	STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATED SPOKES )	◆
<b>TYRES*</b>			
CEE1082	<input type="checkbox"/>	SOLID BLACK 1.3/8" / 35MM	◆
CEE1088	<input type="checkbox"/>	PNEUMATIC LOW PROFILE BLACK 1" / 25MM	○
<b>HANDRIM*</b>			
CEE1140	<input type="checkbox"/>	ALUMINIUM PREANODIZED HANDRIM	◆
<b>REAR WHEEL AXLE*</b>			
CEE1180	<input type="checkbox"/>	Q/R AXLE	◆
<b>REAR WHEEL POSITION*</b>			
CEE1190	<input type="checkbox"/>	STANDARD REAR WHEELS POSITION	◆
CEE1195	<input type="checkbox"/>	MID WHEEL PROPULSION REAR WHEEL POSITION	○
CEE1200	<input type="checkbox"/>	PASSIVE REAR WHEELS POSITION	○
<b>BRAKES*</b>			
CEE1270	<input type="checkbox"/>	PUSH TO LOCK BRAKES	◆
<b>TRANSPORTATION KIT*</b>			
CEE1345	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆
<b>SAFETY OPTIONS*</b>			
CEE1360	<input type="checkbox"/>	POSTURE BELT	◆

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

CEE0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEE0839	TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM)
CEE0878	FOR 400 MM SEAT DEPTH
CEE0879	FOR 450 MM SEAT DEPTH
CEE1195	MID WHEEL PROPULSION REAR WHEEL POSITION
CEE1200	PASSIVE REAR WHEELS POSITION
CEE1345	TRANSPORTATION BRACKETS