



Motability Retail Prescription Form

Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JAN - MAR 2024
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**TDX SP2 Low Rider M14 712116 - Weekly
Price £57.50**

Standard Features

- Max user weight 160kg (25.2 Stone)
- Crash tested ISO7176-Part 19
- Modulite seating
- Matrix Flo-tech Lite cushion
- Without push handles
- Flex 3 Adjustable backrest plates with laguna backrest cushion
- Elan Headrest
- Standard T-armrests, detachable with Long armpad
- Seat Lifter with powered Tilt (300mm, 30°) CoG Shift
- Seat height either 403mm
- Manual settable Recliner 0° - 30° in 3.75° steps
- Standard 80° Hangers with Angle and Depth Adjustable footplates(adapt)
- LiNX control REM211 with swing away bracket
- Solid Black tyres
- 6km/h (4mph) 4 Pole motors
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)



Standard Base price: £8,358

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
ARQ0072	TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height	<input type="radio"/>
User Weight*		
ARQ1918	User weight less 100 kg	<input type="radio"/>
ARQ1919	User weight 101 to 125 kg	<input type="radio"/>
ARQ1921	User weight 126 to 136 kg	<input type="radio"/>
ARQ1922	User weight 137 to 150 kg	<input type="radio"/>
ARQ1923	User weight 151 to 160 kg	<input type="radio"/>

SEATING TYPE*		
ARQ2100	<input type="checkbox"/>	Modulite Seating Unit £8,358
SEAT FRAME*		
ARQ2111	<input type="checkbox"/>	Telescopic Seat Frame with metal seat plate ◆
SEAT WIDTH*		
ARQ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19)) ○
ARQ2123	<input type="checkbox"/>	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21")) ○
ARQ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23")) ○
SEAT DEPTH RAIL LENGTH*		
ARQ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm) ○
ARQ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm) (NC 2134) ○
SEAT DEPTH*		
ARQ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm) (NC 2196) ○
ARQ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm) ○
ARQ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm) ○
ARQ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm) ○
SEAT CUSHION*		
ARQ2142	<input type="checkbox"/>	Matrx Flo-tech Lite ◆
BACKREST SUPPORT*		
ARQ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates ◆
BACKREST UPHOLSTERY / CUSHION*		
ARQ2173	<input type="checkbox"/>	Flex3 Adjustable Backrest Plates with Laguna Backrest cushion ◆
UPHOLSTERY TYPE*		
ARQ2152	<input type="checkbox"/>	Black multi stretch PU coated fabric ○
ARQ2155	<input type="checkbox"/>	Black Spacer fabric ○
ARMREST*		
ARQ2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment ◆
ARMPAD*		
ARQ2196	<input type="checkbox"/>	Long Armpad (375mm) ◆
HEADREST / NECKREST*		
ARQ2324	<input type="checkbox"/>	Elan Standard Pad 355mm (10") wide ◆
SEAT LIFTER / SEAT TILT*		
ARQ1094	<input type="checkbox"/>	Seat lifter with Electrical Tilt ◆
BACKREST RECLINER*		
ARQ0975	<input type="checkbox"/>	Manual Settable Recliner ◆
FOOTREST HANGER TYPE*		
ARQ0526	<input type="checkbox"/>	Side Mounted : Standard 80° hangers ◆
FOOTPLATE TYPE*		
ARQ0512	<input type="checkbox"/>	Fixed Footplate ◆
CONTROLS*		
ARQ2203	<input type="checkbox"/>	REM211 Drive & Actuators ◆
CONTROL BRACKETRY*		
ARQ1707	<input type="checkbox"/>	Swing Away RH Bracket ○
ARQ1708	<input type="checkbox"/>	Swing Away LH Bracket ○

SHROUD COLOUR*			
ARQ0228	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARQ2601	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARQ2602	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARQ2603	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARQ2607	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARQ2608	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARQ2609	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
RIM COLOUR*			
ARQ2611	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARQ2612	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARQ2613	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARQ2617	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARQ2618	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARQ2619	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
ARQ2620	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARQ2621	<input type="checkbox"/>	Without Rim insert	<input type="radio"/>
FORKS*			
ARQ2632	<input type="checkbox"/>	Double sided fork	<input type="checkbox"/>
CENTRE WHEELS (DRIVING WHEELS)*			
ARQ0733	<input type="checkbox"/>	355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	<input type="checkbox"/>
MOTOR TYPE*			
ARQ1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole)	<input type="checkbox"/>
CHARGER*			
ARQ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input type="checkbox"/>
LIGHTS*			
ARQ1254	<input type="checkbox"/>	Without Light system	<input type="checkbox"/>
TRANSPORTATION KIT*			
ARQ1494	<input type="checkbox"/>	4-Point tie-down	<input type="checkbox"/>
FIXATION BELT*			
ARQ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides)	<input type="checkbox"/>

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARQ2134 460mm (18") (410mm - 510mm)

ARQ2196 Long Armpad (375mm)