



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JAN - MAR 2024
Customer Ref:

Customer Ref*, Contact No*, Delivery Address*, Account No*, Order No*, Order Qty*, Order date*

TDX SP2 NB M3 735801 - Weekly Price £62.19

Standard Features

- Max User Weight 150kg (23.5 Stone)
Modulite Seating
Matrix Flo-tech cushion
Flex3 Adjustable Plates with Luguna back cushion
Standard T-armrests, detachable with Long armpad
Seat Lifter with electrical Tilt (300mm, 30°) with CoG shift
Manual settable backrest Recliner (0° - 30° in 3.75° steps)
Manual centre mounted hanger, with angle adjustable footplates
LiNX REM216 with swing away bracket
SolidBlack tyres
10km/h (6mph) 4 pole motors (including indicators, LED headlights and mirror)
4 Point tie-down system (Q'Straint QLK-150 compatible)
Lapbelt with autolock



Standard Base price: £8,927

Please make your selection with a cross in the box [] () Free of charge [] ♦ standard feature POA special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include BASE PRODUCT* (ARM2100), USER WEIGHT* (ARM1920, ARM1921, ARM1922), and BASE MODEL* (ARM0000).

SEATING TYPE*			
ARM2100	<input type="checkbox"/>	Modulite Seating Unit	◆
SEAT FRAME*			
ARM2111	<input type="checkbox"/>	Telescopic Seat Frame	◆
SEAT WIDTH*			
ARM2121	<input type="checkbox"/>	380mm (15")	○
ARM2122	<input type="checkbox"/>	430mm (17")	○
ARM2123	<input type="checkbox"/>	480mm (19")	○
ARM2124	<input type="checkbox"/>	530mm (21")	○
SEAT DEPTH RAIL LENGTH*			
ARM2129	<input type="checkbox"/>	Standard 410mm - 510mm	(NC 2138) ○
ARM2190	<input type="checkbox"/>	Extended 510mm - 560mm	(NC 2132, 2134) ○
SEAT DEPTH*			
ARM2132	<input type="checkbox"/>	410mm (16")	(NC 2196) ○
ARM2134	<input type="checkbox"/>	460mm (18")	○
ARM2136	<input type="checkbox"/>	510mm (20")	○
ARM2138	<input type="checkbox"/>	560mm (22")	○
SEAT CUSHION*			
ARM2142	<input type="checkbox"/>	Matrix Flo-Tech Lite	◆
BACKREST SUPPORT*			
ARM2162	<input type="checkbox"/>	Tall 540mm	◆
BACKREST UPHOLSTERY / CUSHION*			
ARM2173	<input type="checkbox"/>	Flex3 Backrest plates with Laguna cushion	◆
UPHOLSTERY TYPE*			
ARM2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	○
ARM2155	<input type="checkbox"/>	Black Spacer fabric	○
ARMREST*			
ARM2191	<input type="checkbox"/>	T-Armrest standard	◆
ARMPAD*			
ARM2196	<input type="checkbox"/>	Long 375mm Armpad	◆
HEADREST / NECKREST - Adequate headrest option must be selected*			
ARM2324	<input type="checkbox"/>	Elan standard pad 255mm (10")	◆
SEAT LIFTER / SEAT TILT*			
ARM1094	<input type="checkbox"/>	Seat Lifter with centre of gravity shift	◆
BACKREST RECLINER*			
ARM0975	<input type="checkbox"/>	Manual settable Recliner	◆
FOOTREST HANGER TYPE*			
ARM0556	<input type="checkbox"/>	Centre mounted : Manual	◆
FOOTPLATE TYPE*			
ARM0559	<input type="checkbox"/>	Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger only	◆
CONTROLS*			
ARM2204	<input type="checkbox"/>	REM216	◆
CONTROL BRACKETRY*			
ARM1707	<input type="checkbox"/>	Swing away RH Bracket	○
ARM1708	<input type="checkbox"/>	Swing away LH Bracket	○

SHROUD COLOUR*			
ARM0228	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARM2601	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARM2602	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARM2603	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARM2607	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARM2608	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARM2609	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
RIM COLOUR*			
ARM2611	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARM2612	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARM2613	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARM2617	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARM2618	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARM2619	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
ARM2620	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARM2621	<input type="checkbox"/>	Without rim insert	<input type="radio"/>
FORKS*			
ARM2632	<input type="checkbox"/>	Double sided	<input type="checkbox"/>
CENTRE WHEELS (DRIVING WHEELS)*			
ARM0733	<input type="checkbox"/>	Black tyres with gel foam inserts (solid)	<input type="checkbox"/>
MOTOR TYPE*			
ARM1908	<input type="checkbox"/>	10 km/h (6mph)	<input type="checkbox"/>
BATTERY*			
ARM1257	<input type="checkbox"/>	50 Ah Gel	<input type="checkbox"/>
CHARGER*			
ARM1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V)	<input type="checkbox"/>
TRANSPORTATION KIT*			
ARM1494	<input type="checkbox"/>	4-Point tie-down	<input type="checkbox"/>
FIXATION BELT*			
ARM1329	<input type="checkbox"/>	Lap Belt fitted to base of backrest	<input type="radio"/>
ARM2113	<input type="checkbox"/>	Lap Belt adjustable mounting on seat rail extrusion	<input type="checkbox"/>
REAR VIEW MIRROR*			
ARM1415	<input type="checkbox"/>	RH Mirror	<input type="radio"/>
ARM1339	<input type="checkbox"/>	LH Mirror	<input type="radio"/>

TOTAL PRICE	<input type="text"/>
--------------------	----------------------

Non Conforming Options:

ARM2132	410mm (16")
ARM2134	460mm (18")
ARM2138	560mm (22")
ARM2196	Long 375mm Armpad