



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Motability Retail Prescription Form

APRIL 2025  
Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## AVIVA RX40 Modulite M26 707017 - Weekly Price £60.49

### Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flex3 Tension Adjustable Upholstery
- Standard T-Armrests, detachable with long armpad
- Neckrest depth adjustable
- Powered "fixed pivot" Tilt 20°
- Manual Recliner
- Standard hanger 80° pair with Fixed Footplates
- LiNX control REM216 with swing away bracket
- 8" Castor, puncture proof solid tyre
- Solid inlay Dual density foam insert
- 12 km/h Ultra High Speed Package with LED Light system with Indicators, 30 Lux headlights & Mirror
- 73Ah batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £8,062**



Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA    special request

Customer weight =                      KG

| Code                  | Description  |        |
|-----------------------|--|--------|
| <b>BASE MODEL*</b>    |  |        |
| JBG0001               | AVIVA RX 40 Base Model Max user weight up to 136 kg<br>includes Performance Package: 120A Power Module + minimum 60 Ah Batteries | £8,062 |
| <b>CHASSIS WIDTH*</b> |  |        |
| JBL2810               | Chassis Narrow (560 mm/22")  | ◆      |
| <b>SEATING TYPE*</b>  |  |        |
| JBG2100               | Invacare® Modulite Seat Unit   | ◆      |

|   |                          |   |   |
|---|--------------------------|---|---|
| <b>SEAT FRAME*</b>  |                          |   |   |
| JBG2111   | <input type="checkbox"/> | Telescopic Seat Frame and Plate with Side Rails<br>width 380 to 530 mm Additional 50 mm by arms   | ◆ |
| <b>SEAT WIDTH*</b>  |                          |   |   |
| JBG2121   | <input type="checkbox"/> | 380 mm Seat Plate Can be reduced to 330 mm in<br>combination with inrigged armrests   | ○ |
| JBG2122   | <input type="checkbox"/> | 430 mm Seat Plate Can be reduced to 380 mm in<br>combination with inrigged armrests   | ○ |
| JBG2123   | <input type="checkbox"/> | 480 mm Seat Plate   | ○ |
| JBG2124   | <input type="checkbox"/> | 530 mm Seat Plate   | ○ |
| <b>SEAT DEPTH*</b>  |                          |   |   |
| JBG2132   | <input type="checkbox"/> | Seat Depth 410 mm   | ○ |
| JBG2134   | <input type="checkbox"/> | Seat Depth 460 mm   | ○ |
| JBG2136   | <input type="checkbox"/> | Seat Depth 510 mm   | ○ |
| <b>SEAT CUSHION*</b>  |                          |   |   |
| JBG2142   | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite  | ◆ |
| <b>BACKREST SUPPORT*</b>  |                          |   |   |
| JBG2162   | <input type="checkbox"/> | Backposts tall (540 mm) For Flex3 backrests/Tension<br>Adjustable Upholstery/Matrx Backs, measured from<br>seat plate to top of upholstery<br>Width adjustable cast bracket | ◆ |
| <b>PUSH HANDLES FOR BACKPOSTS</b>                                       |                          |   |   |
| JBG2108   | <input type="checkbox"/> | Backposts without push handles  | ◆ |
| <b>BACKREST UPHOLSTERY / CUSHION*</b>                                   |                          |   |   |
| JBG2173   | <input type="checkbox"/> | Flex3 adjustable Backrest Plates with Laguna Backrest<br>cushion  | ◆ |
| <b>UPHOLSTERY TYPE*</b>   |                          |   |   |
| JBG2152   | <input type="checkbox"/> | Black multi stretch PU coated fabric Upholstery   | ○ |
| JBG2155   | <input type="checkbox"/> | Black Spacer Fabric Upholstery  | ○ |
| <b>ARMREST (Pair)*</b>  |                          |   |   |
| JBG2191   | <input type="checkbox"/> | T-Arm detachable, Standard Adjustment range from<br>seat plate:<br>245 - 315 mm   | ◆ |
| <b>ARMPAD (Pair)*</b>   |                          |   |   |
| JBG2196   | <input type="checkbox"/> | Long Armpad (375 mm)  | ◆ |
| <b>HEADREST / NECKREST (Adequate headrest option must be selected)*</b> |                          |   |   |
| JBG1350   | <input type="checkbox"/> | Neckrest depth adjustable   | ◆ |
| <b>MANUAL / POWERED SEATING OPTIONS*</b>                                |                          |   |   |
| JBG1092   | <input type="checkbox"/> | Powered "fixed pivot" Tilt 20°  | ◆ |
| <b>SEAT TO GROUND HEIGHT*</b>   |                          |   |   |
| JBG2857   | <input type="checkbox"/> | Seat-to-Floor Height 435 mm   | ○ |
| JBG2858   | <input type="checkbox"/> | Seat-to-Floor Height 460 mm   | ○ |
| JBG2859   | <input type="checkbox"/> | Seat-to-Floor Height 485 mm   | ○ |
| <b>BACKREST RECLINER*</b>   |                          |   |   |
| JBG0975   | <input type="checkbox"/> | Manual Settable Recliner (in 3.75° steps) 90° - 120°  | ◆ |
| <b>LEG REST TYPE*</b>   |                          |   |   |
| JBG0526   | <input type="checkbox"/> | Standard hanger 80° pair Lower leg length: 350 - 450<br>mm  | ◆ |
| <b>FOOTPLATE TYPE*</b>  |                          |   |   |
| JBG0512   | <input type="checkbox"/> | Fixed Footplate pair  | ◆ |
| <b>LINX CONTROLS*</b>   |                          |   |   |
| JBG2204   | <input type="checkbox"/> | DLX-REM216 = Drive, Actuators & Lights 3 Drive<br>functions and power module incl. 2 Actuators & Lights   | ◆ |

| CONTROL BRACKETRY*  |                          |  |                                  |
|---------------------|--------------------------|--|----------------------------------|
| JBG1707             | <input type="checkbox"/> | Swing Away Bracket RH  | <input type="radio"/>            |
| JBG1708             | <input type="checkbox"/> | Swing Away Bracket LH  | <input type="radio"/>            |
| SHROUD COLOUR*      |                          |  |                                  |
| JBG2871             | <input type="checkbox"/> | Mid range - Black/Anthracite Shroud  | <input type="radio"/>            |
| JBG2872             | <input type="checkbox"/> | Mid range - Silver Shroud  | <input type="radio"/>            |
| JBG2873             | <input type="checkbox"/> | Mid range - Bengal Red Shroud  | <input type="radio"/>            |
| JBG2874             | <input type="checkbox"/> | Mid range - Mercury Blue Shroud  | <input type="radio"/>            |
| JBG2875             | <input type="checkbox"/> | Mid range - Midnight Purple Shroud   | <input type="radio"/>            |
| RIM COLOUR*         |                          |  |                                  |
| JBG2882             | <input type="checkbox"/> | Mid range - Black Rim Inserts  | <input type="radio"/>            |
| JBG2883             | <input type="checkbox"/> | Mid range - Silver Rim Inserts   | <input type="radio"/>            |
| JBG2884             | <input type="checkbox"/> | Mid range - Bengal Red Rim Inserts   | <input type="radio"/>            |
| JBG2885             | <input type="checkbox"/> | Mid range - Mercury Blue Rim Inserts   | <input type="radio"/>            |
| JBG2886             | <input type="checkbox"/> | Mid range - Midnight Purple Rim Inserts  | <input type="radio"/>            |
| FORKS*              |                          |  |                                  |
| JBG2632             | <input type="checkbox"/> | Double sided fork Without castor rim insert<br>(enlarges wide base to 656 mm)  | <input checked="" type="radio"/> |
| FRONT WHEEL*        |                          |  |                                  |
| JBG0836             | <input type="checkbox"/> | 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof<br>solid tyre  | <input checked="" type="radio"/> |
| DRIVE WHEEL*        |                          |  |                                  |
| JBG0705             | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual<br>density foam insert  | <input checked="" type="radio"/> |
| C.T.C SUSPENSION*   |                          |  |                                  |
| JBG2870             | <input type="checkbox"/> | CTC Suspension Control - Traction - Comfort  | <input checked="" type="radio"/> |
| SUSPENSION SETTING* |                          |  |                                  |
| JBG1910             | <input type="checkbox"/> | Soft (suggested for up to 70 kg User Weight)   | <input type="radio"/>            |
| JBG1912             | <input type="checkbox"/> | Medium (suggested for up to 90 kg User Weight)   | <input type="radio"/>            |
| JBG1911             | <input type="checkbox"/> | Firm (suggested for up to 110 kg User Weight)  | <input type="radio"/>            |
| JBG1915             | <input type="checkbox"/> | Very firm (suggested for up to 136 kg User Weight)   | <input type="radio"/>            |
| MOTOR TYPE*         |                          |  |                                  |
| JBG1909             | <input type="checkbox"/> | 12 km/h Ultra High Speed Package<br>Always combined with LED Light system with<br>Indicators, 30 Lux headlights                            | <input checked="" type="radio"/> |
| BATTERY TYPE*       |                          |  |                                  |
| JBG1260             | <input type="checkbox"/> | 73,5 Ah Gel  | <input checked="" type="radio"/> |
| CHARGER*            |                          |  |                                  |
| JBG1343             | <input type="checkbox"/> | Opti Charge Charger 8A (IP44)  | <input checked="" type="radio"/> |
| ANTI-TIPPERS*       |                          |  |                                  |
| JBG1205             | <input type="checkbox"/> | Anti-tipper  | <input checked="" type="radio"/> |
| TRANSPORTATION KIT* |                          |  |                                  |
| JBG1494             | <input type="checkbox"/> | 4-Point tie-down<br>(Tested version for using wheelchair as a vehicle seat<br>according to ISO 7176-19) Adequate headrest is<br>compulsory | <input checked="" type="radio"/> |
| LAP BELT*           |                          |  |                                  |
| JBG2113             | <input type="checkbox"/> | Lap belt with Auto Lock Flexible mount in side rails   | <input checked="" type="radio"/> |
| REAR VIEW MIRROR*   |                          |  |                                  |
| JBG1415             | <input type="checkbox"/> | Rear view mirror RH  | <input type="radio"/>            |
| JBG1339             | <input type="checkbox"/> | Rear view mirror LH  | <input type="radio"/>            |

|  |                    |
|--|--------------------|
|  |                    |
|  | <b>TOTAL PRICE</b> |

Non Conforming Options:



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):