



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Motability Retail Prescription Form

JULY 2025  
Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## AVIVA RX20 Standard Seat M20 707043 - Weekly Price £41.76

### Standard Features

- Max user weight 136kg
- Crash tested ISO 7176 -Part 19
- One piece seat plate
- Flo-tech lite cushion
- Tension adjustable backrest
- Neckrest
- Standard T-Armrests, detachable with Long armpad
- Powered "fixed pivot" Tilt 20°
- Manual settable backrest Recline (0° - 30°, in 3.75° steps)
- 80° hangers with fixed footplate
- Remote LiNX DLX-REM211 with swing away Bracket
- 8" Castor, puncture proof solid tyre
- Solid black tyres
- 6 km/h 4mph motors
- 60Ah batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock



**Standard Base price: £4,280**

Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA    special request

Customer weight =                      KG

| Code                  | Description   |        |
|-----------------------|---|--------|
| <b>BASE MODEL*</b>    |   |        |
| JBB0000               | AVIVA RX 20 Base Model Max user weight up to 136 kg<br>includes 75A Power Module + minimum 52Ah Batteries | £4,280 |
| <b>CHASSIS WIDTH*</b> |   |        |
| JBL2810               | Chassis Narrow (560 mm/22")   | ◆      |
| <b>SEATING TYPE*</b>  |   |        |
| JBB2105               | Modulite Standard Seat Standard Seat Plate and Backrest Frame in fixed Width                              | ◆      |
| <b>SEAT FRAME*</b>    |   |        |
| JBB2110               | One Piece Seat Plate  | ◆      |

|   |                          |  |                                  |
|---|--------------------------|--|----------------------------------|
| <b>SEAT WIDTH*</b>  |                          |  |                                  |
| JBB2121   | <input type="checkbox"/> | 380 mm Seat Plate  | <input type="radio"/>            |
| JBB2122   | <input type="checkbox"/> | 430 mm Seat Plate  | <input type="radio"/>            |
| <b>SEAT DEPTH*</b>  |                          |  |                                  |
| JBB2132   | <input type="checkbox"/> | Seat Depth 410 mm  | <input type="radio"/>            |
| JBB2134   | <input type="checkbox"/> | Seat Depth 460 mm  | <input type="radio"/>            |
| JBB2136   | <input type="checkbox"/> | Seat Depth 510 mm  | <input type="radio"/>            |
| <b>SEAT CUSHION*</b>  |                          |  |                                  |
| JBB2142   | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite   | <input checked="" type="radio"/> |
| <b>BACKREST SUPPORT*</b>  |                          |  |                                  |
| JBB2161   | <input type="checkbox"/> | Backrest frame, detachable, height 480 mm For Tension Adjustable Upholstery/Matrx Backs  | <input checked="" type="radio"/> |
| <b>BACKREST UPHOLSTERY / CUSHION*</b>                                   |                          |  |                                  |
| JBB2171   | <input type="checkbox"/> | Tension Adjustable Upholstery<br>nylon overlay short = 480 mm / 7 straps<br>measured from seat plate to top of upholstery 480 mm | <input checked="" type="radio"/> |
| <b>UPHOLSTERY TYPE*</b>   |                          |  |                                  |
| JBB2151   | <input type="checkbox"/> | Black Fabric Upholstery (Nylon w 2141; 2-way-stretch material w 2142, 2143)  | <input checked="" type="radio"/> |
| <b>ARMREST (Pair)*</b>  |                          |  |                                  |
| JBB2191   | <input type="checkbox"/> | T-Arm detachable, Standard Adjustment range from seat plate:<br>245 - 315 mm   | <input checked="" type="radio"/> |
| <b>ARMPAD (Pair)*</b>   |                          |  |                                  |
| JBB2196   | <input type="checkbox"/> | Long Armpad (375 mm)   | <input checked="" type="radio"/> |
| <b>HEADREST / NECKREST (Adequate headrest option must be selected)*</b> |                          |  |                                  |
| JBB1350   | <input type="checkbox"/> | Neckrest depth adjustable  | <input checked="" type="radio"/> |
| <b>MANUAL / POWERED SEATING OPTIONS*</b>                                |                          |  |                                  |
| JBB1092   | <input type="checkbox"/> | Powered "fixed pivot" Tilt 20°   | <input checked="" type="radio"/> |
| <b>SEAT TO GROUND HEIGHT*</b>   |                          |  |                                  |
| JBB2857   | <input type="checkbox"/> | Seat-to-Floor Height 435 mm  | <input type="radio"/>            |
| JBB2858   | <input type="checkbox"/> | Seat-to-Floor Height 460 mm  | <input type="radio"/>            |
| JBB2859   | <input type="checkbox"/> | Seat-to-Floor Height 485 mm  | <input type="radio"/>            |
| <b>BACKREST RECLINER*</b>   |                          |  |                                  |
| JBB0975   | <input type="checkbox"/> | Manual Settable Recliner (in 3.75° steps) 90° - 120°   | <input checked="" type="radio"/> |
| <b>LEG REST TYPE*</b>   |                          |  |                                  |
| JBB0526   | <input type="checkbox"/> | Standard hanger 80° pair Lower leg length: 350 - 450 mm  | <input checked="" type="radio"/> |
| <b>FOOTPLATE TYPE*</b>  |                          |  |                                  |
| JBB0512   | <input type="checkbox"/> | Fixed Footplate pair   | <input checked="" type="radio"/> |
| <b>LINX CONTROLS*</b>   |                          |  |                                  |
| JBB2203   | <input type="checkbox"/> | Remote LiNX<br>DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators                               | <input checked="" type="radio"/> |
| <b>CONTROL BRACKETRY*</b>   |                          |  |                                  |
| JBB1707   | <input type="checkbox"/> | Swing Away Bracket RH  | <input type="radio"/>            |
| JBB1708   | <input type="checkbox"/> | Swing Away Bracket LH  | <input type="radio"/>            |

| SHROUD COLOUR*      |                          |  |                                  |
|---------------------|--------------------------|--|----------------------------------|
| JBB2871             | <input type="checkbox"/> | Mid range - Black/Anthracite Shroud  | <input type="radio"/>            |
| JBB2872             | <input type="checkbox"/> | Mid range - Silver Shroud  | <input type="radio"/>            |
| JBB2873             | <input type="checkbox"/> | Mid range - Bengal Red Shroud  | <input type="radio"/>            |
| JBB2874             | <input type="checkbox"/> | Mid range - Mercury Blue Shroud  | <input type="radio"/>            |
| JBB2875             | <input type="checkbox"/> | Mid range - Midnight Purple Shroud   | <input type="radio"/>            |
| RIM COLOUR*         |                          |  |                                  |
| JBB2882             | <input type="checkbox"/> | Mid range - Black Rim Inserts  | <input type="radio"/>            |
| JBB2883             | <input type="checkbox"/> | Mid range - Silver Rim Inserts   | <input type="radio"/>            |
| JBB2884             | <input type="checkbox"/> | Mid range - Bengal Red Rim Inserts   | <input type="radio"/>            |
| JBB2885             | <input type="checkbox"/> | Mid range - Mercury Blue Rim Inserts   | <input type="radio"/>            |
| JBB2886             | <input type="checkbox"/> | Mid range - Midnight Purple Rim Inserts  | <input type="radio"/>            |
| FORKS*              |                          |  |                                  |
| JBB2632             | <input type="checkbox"/> | Double sided fork Without castor rim insert<br>(enlarges wide base to 656 mm)  | <input checked="" type="radio"/> |
| CASTORS*            |                          |  |                                  |
| JBB0828             | <input type="checkbox"/> | 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof<br>solid tyre  | <input checked="" type="radio"/> |
| DRIVE WHEELS*       |                          |  |                                  |
| JBB0705             | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual<br>density foam insert  | <input checked="" type="radio"/> |
| C.T.C SUSPENSION*   |                          |  |                                  |
| JBB2870             | <input type="checkbox"/> | CTC Suspension Control - Traction - Comfort  | <input checked="" type="radio"/> |
| SUSPENSION SETTING* |                          |  |                                  |
| JBB1910             | <input type="checkbox"/> | Soft (suggested for up to 70 kg User Weight)   | <input checked="" type="radio"/> |
| JBB1912             | <input type="checkbox"/> | Medium (suggested for up to 90 kg User Weight)   | <input type="radio"/>            |
| JBB1911             | <input type="checkbox"/> | Firm (suggested for up to 110 kg User Weight)  | <input type="radio"/>            |
| JBB1915             | <input type="checkbox"/> | Very firm (suggested for up to 136 kg User Weight)   | <input type="radio"/>            |
| MOTOR TYPE*         |                          |  |                                  |
| JBB1282             | <input type="checkbox"/> | 6 km/h   | <input checked="" type="radio"/> |
| BATTERY TYPE*       |                          |  |                                  |
| JBB1259             | <input type="checkbox"/> | 60 Ah Gel  | <input checked="" type="radio"/> |
| CHARGER*            |                          |  |                                  |
| JBB1343             | <input type="checkbox"/> | Opti Charge Charger 8A (IP44)  | <input checked="" type="radio"/> |
| ANTI-TIPPERS*       |                          |  |                                  |
| JBB1205             | <input type="checkbox"/> | Anti-tipper  | <input checked="" type="radio"/> |
| TRANSPORTATION KIT* |                          |  |                                  |
| JBB1494             | <input type="checkbox"/> | 4-Point tie-down<br>(Tested version for using wheelchair as a vehicle seat<br>according to ISO 7176-19) Adequate headrest is<br>compulsory | <input checked="" type="radio"/> |
| LAP BELT*           |                          |  |                                  |
| JBB1329             | <input type="checkbox"/> | Lap belt with Auto Lock Standard mount on backrest   | <input checked="" type="radio"/> |

|             |  |  |  |
|-------------|--|--|--|
| TOTAL PRICE |  |  |  |
|-------------|--|--|--|

Non Conforming Options:



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):