



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JULY 2025
Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

AVIVA RX40 Modulite M25 707020 - Weekly Price £60.14

Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flex3 Tension Adjustable Upholstery
- T-Armrests, detachable with Long armpad
- Neckrest depth adjustable
- Powered "fixed pivot" Tilt 20°
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- Standard hanger 80° pair with Fixed Footplates
- LiNX control REM216 with swing away bracket
- 8" Castor, puncture proof solid tyre
- Solid inlay Dual density foam insert
- 12 km/h Ultra High Speed Package with LED Light system with Indicators, 30 Lux headlights & Mirror
- 73Ah batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £8,062



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

| Code | Description | |
|-----------------------|--|-------------------------------------|
| BASE MODEL* | | |
| JBG0001 | AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries | £8,062 |
| CHASSIS WIDTH* | | |
| JBL2810 | Chassis Narrow (560 mm/22") | <input checked="" type="checkbox"/> |
| SEATING TYPE* | | |
| JBG2100 | Invacare® Modulite Seat Unit | <input checked="" type="checkbox"/> |

| | | | |
|---|--------------------------|---|---|
| SEAT FRAME* | | | |
| JBG2111 | <input type="checkbox"/> | Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms | ◆ |
| SEAT WIDTH* | | | |
| JBG2121 | <input type="checkbox"/> | 380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests | ○ |
| JBG2122 | <input type="checkbox"/> | 430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests | ○ |
| JBG2123 | <input type="checkbox"/> | 480 mm Seat Plate | ○ |
| JBG2124 | <input type="checkbox"/> | 530 mm Seat Plate | ○ |
| SEAT DEPTH* | | | |
| JBG2132 | <input type="checkbox"/> | Seat Depth 410 mm | ○ |
| JBG2134 | <input type="checkbox"/> | Seat Depth 460 mm | ○ |
| JBG2136 | <input type="checkbox"/> | Seat Depth 510 mm | ○ |
| SEAT CUSHION* | | | |
| JBG2142 | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite | ◆ |
| BACKREST SUPPORT* | | | |
| JBG2162 | <input type="checkbox"/> | Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket | ◆ |
| PUSH HANDLES FOR BACKPOSTS | | | |
| JBG2108 | <input type="checkbox"/> | Backposts without push handles | ◆ |
| BACKREST UPHOLSTERY / CUSHION* | | | |
| JBG2179 | <input type="checkbox"/> | Flex3 Tension Adjustable Upholstery height 540 mm (seat plate to top of upholstery) (NC 2155) | ◆ |
| UPHOLSTERY TYPE* | | | |
| JBG2152 | <input type="checkbox"/> | Black multi stretch PU coated fabric Upholstery | ○ |
| JBG2155 | <input type="checkbox"/> | Black Spacer Fabric Upholstery | ○ |
| ARMREST (Pair)* | | | |
| JBG2191 | <input type="checkbox"/> | T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm | ◆ |
| ARMPAD (Pair)* | | | |
| JBG2196 | <input type="checkbox"/> | Long Armpad (375 mm) | ◆ |
| HEADREST / NECKREST (Adequate headrest option must be selected)* | | | |
| JBG1350 | <input type="checkbox"/> | Neckrest depth adjustable | ◆ |
| MANUAL / POWERED SEATING OPTIONS* | | | |
| JBG1092 | <input type="checkbox"/> | Powered "fixed pivot" Tilt 20° | ◆ |
| SEAT TO GROUND HEIGHT* | | | |
| JBG2857 | <input type="checkbox"/> | Seat-to-Floor Height 435 mm | ○ |
| JBG2858 | <input type="checkbox"/> | Seat-to-Floor Height 460 mm | ○ |
| JBG2859 | <input type="checkbox"/> | Seat-to-Floor Height 485 mm | ○ |
| BACKREST RECLINER* | | | |
| JBG0975 | <input type="checkbox"/> | Manual Settable Recliner (in 3.75° steps) 90° - 120° | ◆ |
| LEG REST TYPE* | | | |
| JBG0526 | <input type="checkbox"/> | Standard hanger 80° pair Lower leg length: 350 - 450 mm | ◆ |
| FOOTPLATE TYPE* | | | |
| JBG0512 | <input type="checkbox"/> | Fixed Footplate pair | ◆ |
| LINX CONTROLS* | | | |
| JBG2204 | <input type="checkbox"/> | DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights | ◆ |

| | | | |
|----------------------------|--|--|---|
| CONTROL BRACKETRY* | | | |
| JBG1707 | | Swing Away Bracket RH | ○ |
| JBG1708 | | Swing Away Bracket LH | ○ |
| SHROUD COLOUR* | | | |
| JBG2871 | | Mid range - Black/Anthracite Shroud | ○ |
| JBG2872 | | Mid range - Silver Shroud | ○ |
| JBG2873 | | Mid range - Bengal Red Shroud | ○ |
| JBG2874 | | Mid range - Mercury Blue Shroud | ○ |
| JBG2875 | | Mid range - Midnight Purple Shroud | ○ |
| RIM COLOUR* | | | |
| JBG2882 | | Mid range - Black Rim Inserts | ○ |
| JBG2883 | | Mid range - Silver Rim Inserts | ○ |
| JBG2884 | | Mid range - Bengal Red Rim Inserts | ○ |
| JBG2885 | | Mid range - Mercury Blue Rim Inserts | ○ |
| JBG2886 | | Mid range - Midnight Purple Rim Inserts | ○ |
| FORKS* | | | |
| JBG2632 | | Double sided fork Without castor rim insert (enlarges wide base to 656 mm) | ◆ |
| CASTORS* | | | |
| JBG0836 | | 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre | ◆ |
| DRIVE WHEELS* | | | |
| JBG0705 | | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert | ◆ |
| C.T.C SUSPENSION* | | | |
| JBG2870 | | CTC Suspension Control - Traction - Comfort | ◆ |
| SUSPENSION SETTING* | | | |
| JBG1910 | | Soft (suggested for up to 70 kg User Weight) | ○ |
| JBG1912 | | Medium (suggested for up to 90 kg User Weight) | ○ |
| JBG1911 | | Firm (suggested for up to 110 kg User Weight) | ○ |
| JBG1915 | | Very firm (suggested for up to 136 kg User Weight) | ○ |
| MOTOR TYPE* | | | |
| JBG1909 | | 12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights | ◆ |
| BATTERY TYPE* | | | |
| JBG1260 | | 73,5 Ah Gel | ◆ |
| CHARGER* | | | |
| JBG1343 | | Opti Charge Charger 8A (IP44) | ◆ |
| ANTI-TIPPERS* | | | |
| JBG1205 | | Anti-tipper | ◆ |
| TRANSPORTATION KIT* | | | |
| JBG1494 | | 4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory | ◆ |
| LAP BELT* | | | |
| JBG2113 | | Lap belt with Auto Lock Flexible mount in side rails | ◆ |
| REAR VIEW MIRROR | | | |
| JBG1415 | | Rear view mirror RH | ○ |
| JBG1339 | | Rear view mirror LH | ○ |

Non Conforming Options:

| | |
|---------|--------------------------------|
| JBG2155 | Black Spacer Fabric Upholstery |
|---------|--------------------------------|



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):