



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Motability Retail Prescription Form

JULY 2025  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## AVIVA RX40 Modulite M26 707017 - Weekly Price £61.46

### Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flex3 Adjustable back plates with Laguna back
- Standard T-Armrests, detachable with long armpad
- Neckrest depth adjustable
- Powered "fixed pivot" Tilt 20°
- Manual Recliner
- Standard hanger 80° pair with Fixed Footplates
- LiNX control REM216 with swing away bracket
- 8" Castor, puncture proof solid tyre
- Solid inlay Dual density foam insert
- 12 km/h Ultra High Speed Package with LED Light system with Indicators, 30 Lux headlights & Mirror
- 73Ah batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



**Standard Base price: £8,315**

Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA    special request

Customer weight =                      KG

Code	Description	
<b>BASE MODEL*</b>		
JBG0001	AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries	£8,315
<b>CHASSIS WIDTH*</b>		
JBL2810	Chassis Narrow (560 mm/22")	<input checked="" type="checkbox"/>
<b>SEATING TYPE*</b>		
JBG2100	Invacare® Modulite Seat Unit	<input checked="" type="checkbox"/>

<b>SEAT FRAME*</b>			
JBG2111	<input type="checkbox"/>	Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	◆
<b>SEAT WIDTH*</b>			
JBG2121	<input type="checkbox"/>	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests	○
JBG2122	<input type="checkbox"/>	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	○
JBG2123	<input type="checkbox"/>	480 mm Seat Plate	○
JBG2124	<input type="checkbox"/>	530 mm Seat Plate	○
<b>SEAT DEPTH*</b>			
JBG2132	<input type="checkbox"/>	Seat Depth 410 mm	○
JBG2134	<input type="checkbox"/>	Seat Depth 460 mm	○
JBG2136	<input type="checkbox"/>	Seat Depth 510 mm	○
<b>SEAT CUSHION*</b>			
JBG2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	◆
<b>BACKREST SUPPORT*</b>			
JBG2162	<input type="checkbox"/>	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	◆
<b>PUSH HANDLES FOR BACKPOSTS</b>			
JBG2108	<input type="checkbox"/>	Backposts without push handles	◆
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
JBG2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	◆
<b>UPHOLSTERY TYPE*</b>			
JBG2152	<input type="checkbox"/>	Black multi stretch PU coated fabric Upholstery	○
JBG2155	<input type="checkbox"/>	Black Spacer Fabric Upholstery	○
<b>ARMREST (Pair)*</b>			
JBG2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
<b>ARMPAD (Pair)*</b>			
JBG2196	<input type="checkbox"/>	Long Armpad (375 mm)	◆
<b>HEADREST / NECKREST (Adequate headrest option must be selected)*</b>			
JBG1350	<input type="checkbox"/>	Neckrest depth adjustable	◆
<b>MANUAL / POWERED SEATING OPTIONS*</b>			
JBG1092	<input type="checkbox"/>	Powered "fixed pivot" Tilt 20°	◆
<b>SEAT TO GROUND HEIGHT*</b>			
JBG2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	○
JBG2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	○
JBG2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	○
<b>BACKREST RECLINER*</b>			
JBG0975	<input type="checkbox"/>	Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
<b>LEG REST TYPE*</b>			
JBG0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
<b>FOOTPLATE TYPE*</b>			
JBG0512	<input type="checkbox"/>	Fixed Footplate pair	◆
<b>LINX CONTROLS*</b>			
JBG2204	<input type="checkbox"/>	DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	◆

<b>CONTROL BRACKETRY*</b>			
JBG1707	<input type="checkbox"/>	Swing Away Bracket RH	<input type="radio"/>
JBG1708	<input type="checkbox"/>	Swing Away Bracket LH	<input type="radio"/>
<b>SHROUD COLOUR*</b>			
JBG2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	<input type="radio"/>
JBG2872	<input type="checkbox"/>	Mid range - Silver Shroud	<input type="radio"/>
JBG2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	<input type="radio"/>
JBG2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	<input type="radio"/>
JBG2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	<input type="radio"/>
<b>RIM COLOUR*</b>			
JBG2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	<input type="radio"/>
JBG2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	<input type="radio"/>
JBG2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	<input type="radio"/>
JBG2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	<input type="radio"/>
JBG2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	<input type="radio"/>
<b>FORKS*</b>			
JBG2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	<input checked="" type="radio"/>
<b>CASTORS*</b>			
JBG0836	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	<input checked="" type="radio"/>
<b>DRIVE WHEELS*</b>			
JBG0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	<input checked="" type="radio"/>
<b>C.T.C SUSPENSION*</b>			
JBG2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	<input checked="" type="radio"/>
<b>SUSPENSION SETTING*</b>			
JBG1910	<input type="checkbox"/>	Soft (suggested for up to 70 kg User Weight)	<input type="radio"/>
JBG1912	<input type="checkbox"/>	Medium (suggested for up to 90 kg User Weight)	<input type="radio"/>
JBG1911	<input type="checkbox"/>	Firm (suggested for up to 110 kg User Weight)	<input type="radio"/>
JBG1915	<input type="checkbox"/>	Very firm (suggested for up to 136 kg User Weight)	<input type="radio"/>
<b>MOTOR TYPE*</b>			
JBG1909	<input type="checkbox"/>	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	<input checked="" type="radio"/>
<b>BATTERY TYPE*</b>			
JBG1260	<input type="checkbox"/>	73,5 Ah Gel	<input checked="" type="radio"/>
<b>CHARGER*</b>			
JBG1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	<input checked="" type="radio"/>
<b>ANTI-TIPPERS*</b>			
JBG1205	<input type="checkbox"/>	Anti-tipper	<input checked="" type="radio"/>
<b>TRANSPORTATION KIT*</b>			
JBG1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	<input checked="" type="radio"/>
<b>LAP BELT*</b>			
JBG2113	<input type="checkbox"/>	Lap belt with Auto Lock Flexible mount in side rails	<input checked="" type="radio"/>
<b>REAR VIEW MIRROR*</b>			
JBG1415	<input type="checkbox"/>	Rear view mirror RH	<input type="radio"/>
JBG1339	<input type="checkbox"/>	Rear view mirror LH	<input type="radio"/>

	<b>TOTAL PRICE</b>

Non Conforming Options:



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):