



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JULY 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Erivo R11 Dual M43 1002906 - Weekly Price £39.45

Standard Features

- ERIVO R11 Dual
- REM 050 Dual Control included to be mounted in attendant position
- Bracket Kit for mounting REM 050 included
- 30Ah Battery
- 35KM range
- Lighting Package included
- Max user weight 150kg (23.5 Stone)
- 6kmh
- High quality foam seat cushion
- Solid seat plate
- Padded backrest with large built in pocket
- Fold-up armrests
- Fold up and height adjustable footplate (3 positions)
- LiNX REM060 remote with simple swivel arm
- 8" PU front castor
- Rear wheel suspension with 12.5" drive wheel, solid tyre



Standard Base price: £3,888

Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
ERIVO*		
1595188	<input type="checkbox"/> R11, 30Ah battery, with lighting package, 35KM range	£3,888
ACCESSORIES 1 (To be sent as spare parts for retrofitting. Additional lead times possible)*		
ERI5678	<input type="checkbox"/> 1602792 - REM 050 Dual Control (to work in conjunction with REM 060 - Bracket Kit for Attendant is required)	◆
ACCESSORIES 2 (To be sent as spare parts for retrofitting. Additional lead times possible)*		
ERI1234	<input type="checkbox"/> 1594974 - Bracket Kit for Attendant Mounting (includes 1000mm bus cable to fit REM 060 or REM 050 in the attendant position)	◆

ACCESSORIES SHIPMENT*

DNPS987		DO NOT PART SHIP - Accessories to be shipped with Erivo	◆
---------	--	---------------------------------------------------------	---

	TOTAL PRICE	
--	-------------	--

Non Conforming Options:



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):