



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JULY 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**TDX SP2 Low Rider M14 712116 - Weekly
Price £60.64**

Standard Features

- Max user weight 160kg (25.2 Stone)
- Crash tested ISO7176-Part 19
- Modulite seating
- Matrix Flo-tech Lite cushion
- Without push handles
- Flex 3 Adjustable backrest plates with laguna backrest cushion
- Elan Headrest
- Standard T-armrests, detachable with Long armpad
- Seat Lifter with powered Tilt (300mm, 30°) CoG Shift
- Seat height either 403mm
- Manual settable Recliner 0° - 30° in 3.75° steps
- Standard 80° Hangers with Angle and Depth Adjustable footplates(adapt)
- LiNX control REM211 with swing away bracket
- Solid Black tyres
- 6km/h (4mph) 4 Pole motors
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £8,216



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code Description

BASE MODEL*

ARQ0072	<input type="checkbox"/>	TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height	£8,216
---------	--------------------------	--	--------

USER WEIGHT*

ARQ1918	<input type="checkbox"/>	User weight less 100 kg	<input type="radio"/>
ARQ1919	<input type="checkbox"/>	User weight 101 to 125 kg	<input type="radio"/>
ARQ1921	<input type="checkbox"/>	User weight 126 to 136 kg	<input type="radio"/>
ARQ1922	<input type="checkbox"/>	User weight 137 to 150 kg	<input type="radio"/>
ARQ1923	<input type="checkbox"/>	User weight 151 to 160 kg	<input type="radio"/>

SEATING TYPE*			
ARQ2100	<input type="checkbox"/>	Modulite Seating Unit	◆
SEAT FRAME*			
ARQ2111	<input type="checkbox"/>	Telescopic Seat Frame with metal seat plate	◆
SEAT WIDTH*			
ARQ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19))	○
ARQ2123	<input type="checkbox"/>	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21"))	○
ARQ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23"))	○
SEAT DEPTH RAIL LENGTH*			
ARQ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm)	○
ARQ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm) (NC 2134)	○
SEAT DEPTH*			
ARQ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm) (NC 2196)	○
ARQ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm)	○
ARQ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm)	○
ARQ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm)	○
SEAT CUSHION*			
ARQ2142	<input type="checkbox"/>	Matrx Flo-tech Lite	◆
BACKREST SUPPORT*			
ARQ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	◆
BACKREST UPHOLSTERY / CUSHION*			
ARQ2173	<input type="checkbox"/>	Flex3 Adjustable Backrest Plates with Laguna Backrest cushion	◆
UPHOLSTERY TYPE*			
ARQ2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	○
ARQ2155	<input type="checkbox"/>	Black Spacer fabric	○
ARMREST*			
ARQ2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment	◆
ARMPAD*			
ARQ2196	<input type="checkbox"/>	Long Armpad (375mm)	◆
HEADREST / NECKREST*			
ARQ2324	<input type="checkbox"/>	Elan Standard Pad 355mm (10") wide	◆
MANUAL / POWERED SEATING OPTIONS*			
ARQ1094	<input type="checkbox"/>	Seat lifter with Electrical Tilt	◆
BACKREST RECLINER*			
ARQ0975	<input type="checkbox"/>	Manual Settable Recliner	◆
LEG REST TYPE*			
ARQ0526	<input type="checkbox"/>	Side Mounted : Standard 80° hangers	◆
FOOTPLATE TYPE*			
ARQ0512	<input type="checkbox"/>	Fixed Footplate	◆
LINX CONTROLS*			
ARQ2203	<input type="checkbox"/>	REM211 Drive & Actuators	◆
CONTROL BRACKETRY*			
ARQ1707	<input type="checkbox"/>	Swing Away RH Bracket	○
ARQ1708	<input type="checkbox"/>	Swing Away LH Bracket	○

SHROUD COLOUR*			
ARQ2640	<input type="checkbox"/>	Obsidian Black Frame Shrouds	<input type="radio"/>
ARQ2641	<input type="checkbox"/>	Ruby Red Frame Shrouds	<input type="radio"/>
ARQ2642	<input type="checkbox"/>	Purple Passion Frame Shrouds	<input type="radio"/>
ARQ2643	<input type="checkbox"/>	Peacock Blue Frame Shrouds	<input type="radio"/>
ARQ2644	<input type="checkbox"/>	Arctic Blue Frame Shrouds	<input type="radio"/>
RIM COLOUR*			
ARQ2650	<input type="checkbox"/>	Obsidian Black Rim Inserts	<input type="radio"/>
ARQ2651	<input type="checkbox"/>	Ruby Red Rim Inserts	<input type="radio"/>
ARQ2652	<input type="checkbox"/>	Purple Passion Rim Inserts	<input type="radio"/>
ARQ2653	<input type="checkbox"/>	Peacock Blue Rim Inserts	<input type="radio"/>
ARQ2654	<input type="checkbox"/>	Arctic Blue Rim Inserts	<input type="radio"/>
FORKS*			
ARQ2632	<input type="checkbox"/>	Double sided fork	<input type="checkbox"/>
CENTRE WHEELS (DRIVE WHEELS)*			
ARQ0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	<input type="checkbox"/>
MOTOR TYPE*			
ARQ1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole)	<input type="checkbox"/>
CHARGER*			
ARQ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input type="checkbox"/>
LIGHTS*			
ARQ1254	<input type="checkbox"/>	Without Light system	<input type="checkbox"/>
TRANSPORTATION KIT*			
ARQ1494	<input type="checkbox"/>	4-Point tie-down	<input type="checkbox"/>
FIXATION BELT*			
ARQ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides)	<input type="checkbox"/>
TOTAL PRICE			

Non Conforming Options:

ARQ2134	460mm (18") (410mm - 510mm)
ARQ2196	Long Armpad (375mm)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):