



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Motability Retail Prescription Form

JULY 2025  
Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

**TDX SP2 NB M4 1001336 - Weekly Price  
£50.41**

## Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash test ISO7176 Part-19
- Modulite Seating
- Matrix Flo-tech cushion
- Basic tension adjustable backrest
- Standard T-armrests, detachable, Long armpad (325mm)
- Manual settable seat Tilt (0° - 20° )
- Manual settable Recliner (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM216 with swing away bracket
- Solid Black tyres
- 10km/h (6mph) motors
- 4 Point tie-down brackets (Q'Straint QLK-150 compatible)
- Lap belt with auto lock (flexible mount in side rails)

**Standard Base price: £6,121**



Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA    special request

Customer weight =                      KG

| Code                 | Description                                          |
|----------------------|------------------------------------------------------|
| <b>BASE PRODUCT*</b> |                                                      |
| ARM2100              | TDX SP2 NARROW BASE £6,121                           |
| <b>USER WEIGHT*</b>  |                                                      |
| ARM1920              | User weight up to 125kg <input type="radio"/>        |
| ARM1921              | User weight 126kg to 136kg <input type="radio"/>     |
| ARM1922              | User weight 137kg to 150kg <input type="radio"/>     |
| <b>BASE MODEL*</b>   |                                                      |
| ARM0000              | TDX SP2 Narrow Base <input checked="" type="radio"/> |

|                                                                         |                          |                                        |   |
|-------------------------------------------------------------------------|--------------------------|----------------------------------------|---|
| <b>SEATING TYPE*</b>                                                    |                          |                                        |   |
| ARM2100                                                                 | <input type="checkbox"/> | Modulite Seating Unit                  | ◆ |
| <b>SEAT FRAME*</b>                                                      |                          |                                        |   |
| ARM2111                                                                 | <input type="checkbox"/> | Telescopic Seat Frame                  | ◆ |
| <b>SEAT WIDTH*</b>                                                      |                          |                                        |   |
| ARM2121                                                                 | <input type="checkbox"/> | 380mm (15")                            | ○ |
| ARM2122                                                                 | <input type="checkbox"/> | 430mm (17")                            | ○ |
| ARM2123                                                                 | <input type="checkbox"/> | 480mm (19")                            | ○ |
| ARM2124                                                                 | <input type="checkbox"/> | 530mm (21")                            | ○ |
| <b>SEAT DEPTH RAIL LENGTH*</b>                                          |                          |                                        |   |
| ARM2129                                                                 | <input type="checkbox"/> | Standard 410mm - 510mm (NC 2138)       | ○ |
| ARM2190                                                                 | <input type="checkbox"/> | Extended 510mm - 560mm (NC 2132, 2134) | ○ |
| <b>SEAT DEPTH*</b>                                                      |                          |                                        |   |
| ARM2132                                                                 | <input type="checkbox"/> | 410mm (16") (NC 2196)                  | ○ |
| ARM2134                                                                 | <input type="checkbox"/> | 460mm (18")                            | ○ |
| ARM2136                                                                 | <input type="checkbox"/> | 510mm (20")                            | ○ |
| ARM2138                                                                 | <input type="checkbox"/> | 560mm (22")                            | ○ |
| <b>SEAT CUSHION*</b>                                                    |                          |                                        |   |
| ARM2142                                                                 | <input type="checkbox"/> | Matrx Flo-Tech Lite                    | ◆ |
| <b>BACKREST SUPPORT*</b>                                                |                          |                                        |   |
| ARM2162                                                                 | <input type="checkbox"/> | Tall 540mm                             | ◆ |
| <b>BACKREST UPHOLSTERY / CUSHION*</b>                                   |                          |                                        |   |
| ARM2171                                                                 | <input type="checkbox"/> | Tension Adjustable Backrest            | ◆ |
| <b>UPHOLSTERY TYPE*</b>                                                 |                          |                                        |   |
| ARM2151                                                                 | <input type="checkbox"/> | Black fabric                           | ◆ |
| <b>ARMREST*</b>                                                         |                          |                                        |   |
| ARM2191                                                                 | <input type="checkbox"/> | T-Armrest standard                     | ◆ |
| <b>ARMPAD*</b>                                                          |                          |                                        |   |
| ARM2196                                                                 | <input type="checkbox"/> | Long 375mm Armpad                      | ◆ |
| <b>HEADREST / NECKREST (Adequate headrest option must be selected)*</b> |                          |                                        |   |
| ARM1495                                                                 | <input type="checkbox"/> | Push bar 22mm                          | ◆ |
| <b>MANUAL / POWERED SEATING OPTIONS*</b>                                |                          |                                        |   |
| ARM1091                                                                 | <input type="checkbox"/> | Manually settable                      | ◆ |
| <b>SEAT HEIGHT SETTING*</b>                                             |                          |                                        |   |
| ARM0354                                                                 | <input type="checkbox"/> | Seat height 450mm (17.75")             | ◆ |
| <b>BACKREST RECLINER*</b>                                               |                          |                                        |   |
| ARM0975                                                                 | <input type="checkbox"/> | Manual settable Recliner               | ◆ |
| <b>LEG REST TYPE*</b>                                                   |                          |                                        |   |
| ARM0526                                                                 | <input type="checkbox"/> | Side mounted : Standard 80° hanger     | ◆ |
| <b>FOOTPLATE TYPE*</b>                                                  |                          |                                        |   |
| ARM0512                                                                 | <input type="checkbox"/> | Fixed Footplate                        | ◆ |
| <b>LINX CONTROLS*</b>                                                   |                          |                                        |   |
| ARM2204                                                                 | <input type="checkbox"/> | REM216                                 | ◆ |
| <b>CONTROL BRACKETRY*</b>                                               |                          |                                        |   |
| ARM1707                                                                 | <input type="checkbox"/> | Swing away RH Bracket                  | ○ |
| ARM1708                                                                 | <input type="checkbox"/> | Swing away LH Bracket                  | ○ |

| SHROUD COLOUR*                |                          |                                                                     |                                  |
|-------------------------------|--------------------------|---------------------------------------------------------------------|----------------------------------|
| ARM2640                       | <input type="checkbox"/> | Obsidian Black Frame Shrouds                                        | <input type="radio"/>            |
| ARM2641                       | <input type="checkbox"/> | Ruby Red Frame Shrouds                                              | <input type="radio"/>            |
| ARM2642                       | <input type="checkbox"/> | Purple Passion Frame Shrouds                                        | <input type="radio"/>            |
| ARM2643                       | <input type="checkbox"/> | Peacock Blue Frame Shrouds                                          | <input type="radio"/>            |
| ARM2644                       | <input type="checkbox"/> | Arctic Blue Frame Shrouds                                           | <input type="radio"/>            |
| RIM COLOUR*                   |                          |                                                                     |                                  |
| ARM2650                       | <input type="checkbox"/> | Obsidian Black Rim Inserts                                          | <input type="radio"/>            |
| ARM2651                       | <input type="checkbox"/> | Ruby Red Rim Inserts                                                | <input type="radio"/>            |
| ARM2652                       | <input type="checkbox"/> | Purple Passion Rim Inserts                                          | <input type="radio"/>            |
| ARM2653                       | <input type="checkbox"/> | Peacock Blue Rim Inserts                                            | <input type="radio"/>            |
| ARM2654                       | <input type="checkbox"/> | Arctic Blue Rim Inserts                                             | <input type="radio"/>            |
| FORKS*                        |                          |                                                                     |                                  |
| ARM2632                       | <input type="checkbox"/> | Double sided                                                        | <input checked="" type="radio"/> |
| CENTRE WHEELS (DRIVE WHEELS)* |                          |                                                                     |                                  |
| ARM0705                       | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert | <input checked="" type="radio"/> |
| MOTOR TYPE*                   |                          |                                                                     |                                  |
| ARM1908                       | <input type="checkbox"/> | 10 km/h (6mph)                                                      | <input checked="" type="radio"/> |
| BATTERY*                      |                          |                                                                     |                                  |
| ARM1257                       | <input type="checkbox"/> | 50 Ah Gel                                                           | <input checked="" type="radio"/> |
| CHARGER*                      |                          |                                                                     |                                  |
| ARM1343                       | <input type="checkbox"/> | OptiCharge 8A (200 - 250V)                                          | <input checked="" type="radio"/> |
| TRANSPORTATION KIT*           |                          |                                                                     |                                  |
| ARM1494                       | <input type="checkbox"/> | 4-Point tie-down                                                    | <input checked="" type="radio"/> |
| FIXATION BELT*                |                          |                                                                     |                                  |
| ARM1329                       | <input type="checkbox"/> | Lap Belt fitted to base of backrest                                 | <input type="radio"/>            |
| ARM2113                       | <input type="checkbox"/> | Lap Belt adjustable mounting on seat rail extrusion                 | <input checked="" type="radio"/> |
| REAR VIEW MIRROR              |                          |                                                                     |                                  |
| ARM1415                       | <input type="checkbox"/> | RH Mirror                                                           | <input type="radio"/>            |
| ARM1339                       | <input type="checkbox"/> | LH Mirror                                                           | <input type="radio"/>            |

|             |  |  |  |
|-------------|--|--|--|
| TOTAL PRICE |  |  |  |
|-------------|--|--|--|

**Non Conforming Options:**

|         |                   |
|---------|-------------------|
| ARM2132 | 410mm (16")       |
| ARM2134 | 460mm (18")       |
| ARM2138 | 560mm (22")       |
| ARM2196 | Long 375mm Armpad |



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):