



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX20 Modulite M22 707038 - Weekly Price £51.71

Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flex3 adjustable Backrest Plates with Laguna Backrest cushion
- Standard T-Armrests, detachable with long armpad
- Neckrest depth adjustable
- Powered "fixed pivot" Tilt 20°
- Manual Settable Recliner (in 3.75° steps) 90° - 120°
- 80° hangers with fixed footplate
- Remote LiNX DLX-REM211 with swing away bracket
- 8" Castor, puncture proof solid tyre
- Solid inlay black drive wheels
- 6 km/h 4mph motors
- 60Ah batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,462



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBC0000	☐ AVIVA RX 20 Base Model Max user weight up to 136 kg	£6,462
CHASSIS WIDTH*		
JBL2810	☐ Chassis Narrow (560 mm/22")	◆
SEATING TYPE*		
JBC2100	☐ Invacare® Modulite Seat Unit	◆
SEAT FRAME*		
JBC2111	☐ Telescopic Seat Frame and Plate with Side Rails	◆

SEAT WIDTH*			
JBC2121	☐	380 mm Seat Plate	☐
JBC2122	☐	430 mm Seat Plate	☐
JBC2123	☐	480 mm Seat Plate	☐
JBC2124	☐	530 mm Seat Plate	☐
SEAT DEPTH*			
JBC2132	☐	Seat Depth 410 mm	☐
JBC2134	☐	Seat Depth 460 mm	☐
JBC2136	☐	Seat Depth 510 mm	☐
SEAT CUSHION*			
JBC2142	☐	Matrx Flo-tech Cushion Flo-tech Lite	☒
BACKREST SUPPORT*			
JBC2162	☐	Backposts tall (540 mm)	☒
PUSH HANDLES FOR BACKPOSTS			
JBC2108	☐	Backposts without push handles	☒
BACKREST UPHOLSTERY / CUSHION*			
JBC2173	☐	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	☒
UPHOLSTERY TYPE*			
JBC2152	☐	Black Black multi stretch PU coated fabric Upholstery	☐
JBC2155	☐	Black Spacer Fabric Upholstery	☐
ARMREST (Pair)*			
JBC2191	☐	T-Arm detachable, Standard	☒
ARMPAD (Pair)*			
JBC2196	☐	Long Armpad (375 mm)	☒
HEADREST / NECKREST (Adequate headrest option must be selected)*			
JBC1350	☐	Neckrest depth adjustable	☒
MANUAL / POWERED SEATING OPTIONS*			
JBC1092	☐	Powered "fixed pivot" Tilt 20°	☒
SEAT TO GROUND HEIGHT*			
JBC2857	☐	Seat-to-Floor Height 435 mm	☐
JBC2858	☐	Seat-to-Floor Height 460 mm	☐
JBC2859	☐	Seat-to-Floor Height 485 mm	☐
BACKREST RECLINER*			
JBC0975	☐	Manual Settable Recliner (in 3.75° steps) 90° - 120°	☒
LEG REST TYPE*			
JBC0526	☐	Standard hanger 80° pair Lower leg length: 350 - 450 mm	☒
FOOTPLATE TYPE*			
JBC0512	☐	Fixed Footplate pair	☒
LiNX CONTROLS*			
JBC2203	☐	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	☒
CONTROL BRACKETRY*			
JBC1707	☐	Swing Away Bracket RH	☐
JBC1708	☐	Swing Away Bracket LH	☐

SHROUD COLOUR*			
JBC2871	☐	Mid range - Black/Anthracite Shroud	☐
JBC2872	☐	Mid range - Silver Shroud	☐
JBC2873	☐	Mid range - Bengal Red Shroud	☐
JBC2874	☐	Mid range - Mercury Blue Shroud	☐
JBC2875	☐	Mid range - Midnight Purple Shroud	☐
RIM COLOUR*			
JBC2882	☐	Mid range - Black Rim Inserts	☐
JBC2883	☐	Mid range - Silver Rim Inserts	☐
JBC2884	☐	Mid range - Bengal Red Rim Inserts	☐
JBC2885	☐	Mid range - Mercury Blue Rim Inserts	☐
JBC2886	☐	Mid range - Midnight Purple Rim Inserts	☐
FORKS*			
JBC2632	☐	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	☒
CASTORS*			
JBC0828	☐	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	☒
DRIVE WHEELS*			
JBC0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	☒
C.T.C SUSPENSION*			
JBC2870	☐	CTC Suspension Control - Traction - Comfort	☒
SUSPENSION SETTING*			
JBC1910	☐	Soft (suggested for up to 70 kg User Weight)	☐
JBC1912	☐	Medium (suggested for up to 90 kg User Weight)	☐
JBC1911	☐	Firm (suggested for up to 110 kg User Weight)	☐
JBC1915	☐	Very firm (suggested for up to 136 kg User Weight)	☐
MOTOR TYPE*			
JBC1282	☐	6 km/h	☒
BATTERY TYPE*			
JBC1259	☐	60 Ah Gel	☒
CHARGER*			
JBC1343	☐	Opti Charge Charger 8A (IP44)	☒
ANTI-TIPPERS*			
JBC1205	☐	Anti-tipper	☒
TRANSPORTATION KIT*			
JBC1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	☒
LAP BELT*			
JBC2113	☐	Lap belt with Auto Lock Flexible mount in side rails	☒

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):