



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX20 Standard Seat M20 707043 - Weekly Price £41.76

Standard Features

- Max user weight 136kg
- Crash tested ISO 7176 -Part 19
- One piece seat plate
- Flo-tech lite cushion
- Tension adjustable backrest
- Neckrest
- Standard T-Armrests, detachable with Long armpad
- Powered "fixed pivot" Tilt 20°
- Manual settable backrest Recline (0° - 30°, in 3.75° steps)
- 80° hangers with fixed footplate
- Remote LiNX DLX-REM211 with swing away Bracket
- 8" Castor, puncture proof solid tyre
- Solid black tyres
- 6 km/h 4mph motors
- 60Ah batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock



Standard Base price: £4,280

Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBB0000	AVIVA RX 20 Base Model Max user weight up to 136 kg includes 75A Power Module + minimum 52Ah Batteries	£4,280
CHASSIS WIDTH*		
JBL2810	Chassis Narrow (560 mm/22")	◆
SEATING TYPE*		
JBB2105	Modulite Standard Seat Standard Seat Plate and Backrest Frame in fixed Width	◆
SEAT FRAME*		
JBB2110	One Piece Seat Plate	◆

SEAT WIDTH*			
JBB2121	☐	380 mm Seat Plate	☐
JBB2122	☐	430 mm Seat Plate	☐
SEAT DEPTH*			
JBB2132	☐	Seat Depth 410 mm	☐
JBB2134	☐	Seat Depth 460 mm	☐
JBB2136	☐	Seat Depth 510 mm	☐
SEAT CUSHION*			
JBB2142	☐	Matrx Flo-tech Cushion Flo-tech Lite	☒
BACKREST SUPPORT*			
JBB2161	☐	Backrest frame, detachable, height 480 mm For Tension Adjustable Upholstery/Matrx Backs	☒
BACKREST UPHOLSTERY / CUSHION*			
JBB2171	☐	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm	☒
UPHOLSTERY TYPE*			
JBB2151	☐	Black Fabric Upholstery (Nylon w 2141; 2-way-stretch material w 2142, 2143)	☒
ARMREST (Pair)*			
JBB2191	☐	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	☒
ARMPAD (Pair)*			
JBB2196	☐	Long Armpad (375 mm)	☒
HEADREST / NECKREST (Adequate headrest option must be selected)*			
JBB1350	☐	Neckrest depth adjustable	☒
MANUAL / POWERED SEATING OPTIONS*			
JBB1092	☐	Powered "fixed pivot" Tilt 20°	☒
SEAT TO GROUND HEIGHT*			
JBB2857	☐	Seat-to-Floor Height 435 mm	☐
JBB2858	☐	Seat-to-Floor Height 460 mm	☐
JBB2859	☐	Seat-to-Floor Height 485 mm	☐
BACKREST RECLINER*			
JBB0975	☐	Manual Settable Recliner (in 3.75° steps) 90° - 120°	☒
LEG REST TYPE*			
JBB0526	☐	Standard hanger 80° pair Lower leg length: 350 - 450 mm	☒
FOOTPLATE TYPE*			
JBB0512	☐	Fixed Footplate pair	☒
LINX CONTROLS*			
JBB2203	☐	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	☒
CONTROL BRACKETRY*			
JBB1707	☐	Swing Away Bracket RH	☐
JBB1708	☐	Swing Away Bracket LH	☐

SHROUD COLOUR*			
JBB2871	☐	Mid range - Black/Anthracite Shroud	☐
JBB2872	☐	Mid range - Silver Shroud	☐
JBB2873	☐	Mid range - Bengal Red Shroud	☐
JBB2874	☐	Mid range - Mercury Blue Shroud	☐
JBB2875	☐	Mid range - Midnight Purple Shroud	☐
RIM COLOUR*			
JBB2882	☐	Mid range - Black Rim Inserts	☐
JBB2883	☐	Mid range - Silver Rim Inserts	☐
JBB2884	☐	Mid range - Bengal Red Rim Inserts	☐
JBB2885	☐	Mid range - Mercury Blue Rim Inserts	☐
JBB2886	☐	Mid range - Midnight Purple Rim Inserts	☐
FORKS*			
JBB2632	☐	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	☒
CASTORS*			
JBB0828	☐	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	☒
DRIVE WHEELS*			
JBB0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	☒
C.T.C SUSPENSION*			
JBB2870	☐	CTC Suspension Control - Traction - Comfort	☒
SUSPENSION SETTING*			
JBB1910	☐	Soft (suggested for up to 70 kg User Weight)	☒
JBB1912	☐	Medium (suggested for up to 90 kg User Weight)	☐
JBB1911	☐	Firm (suggested for up to 110 kg User Weight)	☐
JBB1915	☐	Very firm (suggested for up to 136 kg User Weight)	☐
MOTOR TYPE*			
JBB1282	☐	6 km/h	☒
BATTERY TYPE*			
JBB1259	☐	60 Ah Gel	☒
CHARGER*			
JBB1343	☐	Opti Charge Charger 8A (IP44)	☒
ANTI-TIPPERS*			
JBB1205	☐	Anti-tipper	☒
TRANSPORTATION KIT*			
JBB1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	☒
LAP BELT*			
JBB1329	☐	Lap belt with Auto Lock Standard mount on backrest	☒

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):