



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 Modulite M27 707015 - Weekly Price £62.33

Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flo-tech lite cushion
- Flex3 Laguna backrest
- Standard T-Armrests, detachable with Long armpad
- Neckrest depth adjustable
- Seat lifter with powered tilt (300 mm, 30°) with CoG shift
- Manual settable backrest Recline
- Standard hanger 80° pair with Fixed Footplates
- LiNX control REM216 with swing away bracket
- 9" Castor, puncture proof solid tyre
- Solid inlay Dual density foam insert
- 12 km/h Ultra High Speed Package with LED Light system with Indicators, 30 Lux headlights & Mirror
- 73Ah batteries
- Anti tippers
- 6-Point tie-down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £9,804



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBG0001	AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries	£9,804
CHASSIS WIDTH*		
JBL2810	Chassis Narrow (560 mm/22")	◆
SEATING TYPE*		
JBG2100	Invacare® Modulite Seat Unit	◆

SEAT FRAME*			
JBG2111	☐	Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	◆
SEAT WIDTH*			
JBG2121	☐	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests	○
JBG2122	☐	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	○
JBG2123	☐	480 mm Seat Plate	○
JBG2124	☐	530 mm Seat Plate	○
SEAT DEPTH*			
JBG2132	☐	Seat Depth 410 mm	○
JBG2134	☐	Seat Depth 460 mm	○
JBG2136	☐	Seat Depth 510 mm	○
SEAT CUSHION*			
JBG2142	☐	Matrx Flo-tech Cushion Flo-tech Lite	◆
BACKREST SUPPORT*			
JBG2162	☐	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	◆
PUSH HANDLES FOR BACKPOSTS			
JBG2108	☐	Backposts without push handles	◆
BACKREST UPHOLSTERY / CUSHION*			
JBG2173	☐	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	◆
UPHOLSTERY TYPE*			
JBG2152	☐	Black multi stretch PU coated fabric Upholstery	○
JBG2155	☐	Black Spacer Fabric Upholstery	○
ARMREST (Pair)*			
JBG2191	☐	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
ARMPAD (Pair)*			
JBG2196	☐	Long Armpad (375 mm)	◆
HEADREST / NECKREST (Adequate headrest option must be selected)*			
JBG1350	☐	Neckrest depth adjustable	◆
MANUAL / POWERED SEATING OPTIONS*			
JBG1094	☐	Seat lifter with powered tilt (300 mm, 30°) with CoG shift	◆
SEAT TO GROUND HEIGHT*			
JBG2857	☐	Seat-to-Floor Height 435 mm	○
JBG2858	☐	Seat-to-Floor Height 460 mm	○
JBG2859	☐	Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBG0975	☐	Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
LEG REST TYPE*			
JBG0526	☐	Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
FOOTPLATE TYPE*			
JBG0512	☐	Fixed Footplate pair	◆
LINX CONTROLS*			
JBG2204	☐	DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	◆

CONTROL BRACKETRY*			
JBG1707	☐	Swing Away Bracket RH	☐
JBG1708	☐	Swing Away Bracket LH	☐
SHROUD COLOUR*			
JBG2871	☐	Mid range - Black/Anthracite Shroud	☐
JBG2872	☐	Mid range - Silver Shroud	☐
JBG2873	☐	Mid range - Bengal Red Shroud	☐
JBG2874	☐	Mid range - Mercury Blue Shroud	☐
JBG2875	☐	Mid range - Midnight Purple Shroud	☐
RIM COLOUR*			
JBG2882	☐	Mid range - Black Rim Inserts	☐
JBG2883	☐	Mid range - Silver Rim Inserts	☐
JBG2884	☐	Mid range - Bengal Red Rim Inserts	☐
JBG2885	☐	Mid range - Mercury Blue Rim Inserts	☐
JBG2886	☐	Mid range - Midnight Purple Rim Inserts	☐
FORKS*			
JBG2632	☐	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	☒
CASTORS*			
JBG0836	☐	9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	☒
DRIVE WHEELS*			
JBG0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	☒
C.T.C SUSPENSION*			
JBG2870	☐	CTC Suspension Control - Traction - Comfort	☒
SUSPENSION SETTING*			
JBG1910	☐	Soft (suggested for up to 70 kg User Weight) (NC 1094)	☐
JBG1912	☐	Medium (suggested for up to 90 kg User Weight) (NC 1094)	☐
JBG1911	☐	Firm (suggested for up to 110 kg User Weight)	☐
JBG1915	☐	Very firm (suggested for up to 136 kg User Weight)	☐
MOTOR TYPE*			
JBG1909	☐	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	☒
BATTERY TYPE*			
JBG1260	☐	73,5 Ah Gel	☒
CHARGER*			
JBG1343	☐	Opti Charge Charger 8A (IP44)	☒
ANTI-TIPPERS*			
JBG1205	☐	Anti-tipper	☒
TRANSPORTATION KIT*			
JBG1401	☐	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory	☒
LAP BELT*			
JBG2113	☐	Lap belt with Auto Lock Flexible mount in side rails	☒
REAR VIEW MIRROR			
JBG1415	☐	Rear view mirror RH	☐
JBG1339	☐	Rear view mirror LH	☐

	TOTAL PRICE

Non Conforming Options:

JBG1094	Seat lifter with powered tilt (300 mm, 30°) with CoG shift
---------	--



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):