



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm4 Modulite M31 1001116 - Weekly Price £61.94

Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flex 3 Laguna Backrest
- Standard T-Armrests, detachable with Long armpad
- Electrical Tilt 0° to 25°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F with Fixed Footplate
- LiNX Control REM216 with swing away bracket
- Pneumatic Castors
- Pneumatic Drive Wheels
- 12 km/h (8mph) Motors, LED light system with indicators
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers



Standard Base price: £8,334

Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ATB0000	Storm 4 Modulite	£8,334
SEATING*		
ATB2100	Modulite Seating Unit (Flex3 seating unit)	◆
SEATING OPTION*		
ATB2111	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆

SEAT WIDTH*			
ATBXXXX	☐	Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	☐
ATB2121	☐	Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest)	☐
ATB2122	☐	Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)	☐
ATB2123	☐	Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest)	☐
ATB2124	☐	Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	☐
SEAT DEPTH*			
ATB2132	☐	410mm (16") (adjustable 410mm - 510mm stepless adjustment) (NC 2196)	☐
ATB2134	☐	460mm (18") (adjustable 410mm - 510mm stepless adjustment)	☐
ATB2136	☐	510mm (20") (adjustable 410mm - 510mm stepless adjustment)	☐
SEAT CUSHION*			
ATB2142	☐	Matrx Flo-tech Lite	☒
BACKREST SUPPORT*			
ATB2162	☐	Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	☒
PUSH HANDLES FOR BACKPOSTS*			
ATB2107	☐	Backposts WITH push handles	☒
BACKREST UPHOLSTERY / CUSHION*			
ATB2173	☐	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	☒
UPHOLSTERY TYPE*			
ATB2152	☐	Black multi stretch PU coated fabric	☐
ATB2155	☐	Black Spacer fabric	☐
ARMREST (PAIR)*			
ATB2191	☐	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	☒
ARMPAD (PAIR)*			
ATB2196	☐	Long Armpad (375mm) (NC 2132)	☒
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ATB1350	☐	neckrest: height and depth adjustable	☒
SEAT HEIGHT AND SEAT RISER*			
ATB1023	☐	Fixed Seat Height 450mm (17.75")	☐
ATB1025	☐	Fixed Seat Height 480mm (19")	☐
SEAT TILT*			
ATB1014	☐	Electrical Tilt 0° to 25°	☒
BACKREST RECLINER*			
ATB0975	☐	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	☒
LEG REST TYPE*			
ATB0503	☐	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	☐
FOOTPLATE TYPE*			
ATB0512	☐	Fixed Footplate	☐

LINX CONTROLS*			
ATB2204	☐	REM216 Drive, Actuators & Lights	◆
CONTROL BRACKETRY*			
ATB1707	☐	Swing Away RH Bracket	○
ATB1708	☐	Swing Away LH Bracket	○
SHROUD COLOUR*			
ATB0202	☐	Blue	○
ATB0228	☐	Black	○
ATB0204	☐	Transparent Cherry Red	○
CASTORS*			
ATB0851	☐	3.00 - 6 Grey Pneumatic	◆
DRIVE WHEELS*			
ATB0771	☐	3.00 - 8 Grey Pneumatic	◆
MOTOR TYPE*			
ATB1909	☐	12 km/h (8mph) (includes LED light system with indicators) (NC 1437)	◆
BATTERY*			
ATB1364	☐	70 Ah Gel	◆
SUSPENSION*			
ATB1240	☐	Suspension low (up to 80kg user weight)	○
ATB1241	☐	Suspension middle (from 70kg up to 100kg user weight)	○
ATB1242	☐	Suspension hard (from 90kg up to 150kg user weight)	○
CHARGER*			
ATB1343	☐	OptiCharge Charger 8A (200 - 250V)	◆
LIGHTS*			
ATB1437	☐	LED Light system with indicators (NC 1909)	◆
ANTI TIPPER*			
ATB1205	☐	Anti-tipper (standard version)	◆
TRANSPORTATION KIT*			
ATB1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
FIXATION BELT*			
ATB2113	☐	Lap belt with auto lock (flexible mount in side rails)	◆
REAR VIEW MIRROR			
ATB1415	☐	Rear view Mirror RH - Must be selected with 10 km/h Motor	○
ATB1339	☐	Rear view Mirror LH - Must be selected with 10 km/h Motor	○

TOTAL PRICE			
--------------------	--	--	--

Non Conforming Options:

ATB2132	410mm (16") (adjustable 410mm - 510mm stepless adjustment)
ATB2196	Long Armpad (375mm)
ATB1909	12 km/h (8mph) (includes LED light system with indicators)
ATB1437	LED Light system with indicators



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):