



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm4 X-Plore Modulite M34 707458 - Weekly Price £58.72

Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO7176-Part 19
- Modulite seating system
- Flex3 Laguna backrest
- Matrix Flo-tech Lite cushion
- Backposts 540mm (Tall)
- Standard T-armrests, detachable with Long armpads
- Manual settable seat Tilt
- Manual settable backrest Recline
- Angle adjustable legrest hangers (VARI F) with Fixed footplates
- LINX REM216 control with swing away bracket
- Puncture proof (sealed) drive wheels and front castors
- 12 km/h (8mph) includes LED light system with indicators, Mirror
- Lapbelt (adjustable both sides)

Standard Base price: £7,885



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ATE0017	Storm 4 X-plore	£7,885
ATZ0017	Special Storm 4 X-plore	POA
SEATING TYPE*		
ATE2100	Invacare Modulite 'Flex3'	◆
SEAT FRAME*		
ATE2111	Telescopic Seat Frame with metal seat plate	◆

SEAT WIDTH*			
ATEXXXX	☐	Seat width 330mm (13")	☐
ATE2121	☐	Set at 380mm	☐
ATE2122	☐	Set at 430mm	☐
ATE2123	☐	Set at 480mm	☐
ATE2124	☐	Set at 530mm	☐
SEAT DEPTH*			
ATE2132	☐	410mm (16") (NC 2196)	☐
ATE2134	☐	460mm (18")	☐
ATE2136	☐	510mm (20")	☐
SEAT CUSHION*			
ATE2142	☐	Matrx Flo-tech Lite	☒
BACKREST SUPPORT*			
ATE2162	☐	Backposts 540mm (Tall)	☒
PUSH HANDLES FOR BACKPOSTS*			
ATE2107	☐	Backposts WITH push handles	☒
BACKREST UPHOLSTREY / CUSHION*			
ATE2173	☐	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	☒
UPHOLSTERY COLOUR*			
ATE2152	☐	Black multi stretch PU coated fabric	☐
ATE2155	☐	Black Spacer fabric	☐
ARMREST*			
ATE2191	☐	T-Armrest Standard	☒
ARMPAD*			
ATE2196	☐	Long Armpad (375mm)	☒
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ATE1495	☐	Push bar	☒
SEAT HEIGHT SETTING & SEAT RISER*			
ATE1023	☐	Fixed Seat Height 450mm	☐
ATE1025	☐	Fixed Seat Height 480mm	☐
SEAT TILT*			
ATE1002	☐	Manual Seat Tilt by spindle	☒
BACKREST RECLINER*			
ATE0975	☐	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	☒
LEG REST TYPE*			
ATE0503	☐	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	☒
FOOTPLATE TYPE*			
ATE0512	☐	Fixed Footplate	☒
LINX CONTROLS*			
ATE2204	☐	REM216 Drive, Actuators & Lights	☒
CONTROL BRACKETRY*			
ATE1707	☐	Swing Away RH Bracket	☐
ATE1708	☐	Swing Away LH Bracket	☐

SHROUD COLOUR*			
ATE0202	☐	Blue	☐
ATE0228	☐	Black	☐
ATE0204	☐	Transparent Cherry Red	☐
CASTORS*			
ATE0853	☐	3.00 - 6 Grey Puncture Proof (sealed)	☒
DRIVE WHEELS*			
ATE0773	☐	3.00 - 8 Grey Puncture Proof (sealed)	☒
MOTOR TYPE*			
ATE1909	☐	12 km/h (8mph) includes LED light system with indicators	☒
BATTERY*			
ATE1364	☐	73.5 Ah GEL	☒
REAR SUSPENSION*			
ATE1233	☐	Hitec spring X-plore	☐
FRONT SUSPENSION*			
ATE1238	☐	Front suspension (suggested for up to 100kg)	☐
ATE1239	☐	Front suspension (suggested for above 100kg)	☐
CHARGER*			
ATE1343	☐	OptiCharge Charger 8A (200 - 250V)	☒
ANTI-TIPPERS*			
ATE1205	☐	Anti-tipper (standard version)	☒
TRANSPORTATION*			
ATE1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory).	☒
FIXATION BELT*			
ATE1329	☐	Lapbelt (adjustable both sides) fitted to base of backrest	☒
REAR VIEW MIRROR*			
ATE1415	☐	Rear view Mirror RH - Must be selected with 10 km/h Motor	☐
ATE1339	☐	Rear view Mirror LH - Must be selected with 10 km/h Motor	☐

			TOTAL PRICE	
--	--	--	--------------------	--

Non Conforming Options:

ATE2196	Long Armpad (375mm)
---------	---------------------



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):