



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm4 X-Plore Modulite M37 1001120 - Weekly Price £62.23

Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO7176-Part 19
- Modulite seating system
- Flex3 Laguna Backrest
- Standard T-Armrests, detachable with Long Armpads
- Neckrest
- Electrical Seat Riser, Electrical Tilt from 0° to 25°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F hangers with Fixed Footplate
- LiNX control REM216 with swing away bracket
- Puncture proof castors
- Puncture proof Drive Wheels
- 12 km/h (8mph) Motors, LED light system with indicators
- 4 Point-tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers
- X-plore suspension 50mm all round

Standard Base price: £9,653



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ATE0017	☐ Storm 4 X-plore	£9,653
SEATING TYPE*		
ATE2100	☐ Invacare Modulite 'Flex3'	◆
SEAT FRAME*		
ATE2111	☐ Telescopic Seat Frame with metal seat plate	◆

SEAT WIDTH*			
ATEXXXX	☐	Seat width 330mm (13")	☐
ATE2121	☐	Set at 380mm	☐
ATE2122	☐	Set at 430mm	☐
ATE2123	☐	Set at 480mm	☐
ATE2124	☐	Set at 530mm	☐
SEAT DEPTH*			
ATE2132	☐	410mm (16") (NC 2196)	☐
ATE2134	☐	460mm (18")	☐
ATE2136	☐	510mm (20")	☐
SEAT CUSHION*			
ATE2142	☐	Matrx Flo-tech Lite	☒
BACKREST SUPPORT*			
ATE2162	☐	Backposts 540mm (Tall)	☒
PUSH HANDLES FOR BACKPOSTS*			
ATE2107	☐	Backposts WITH push handles	☒
BACKREST UPHOLSTREY / CUSHION*			
ATE2173	☐	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	☒
UPHOLSTERY COLOUR*			
ATE2152	☐	Black multi stretch PU coated fabric	☐
ATE2155	☐	Black Spacer fabric	☐
ARMREST*			
ATE2191	☐	T-Armrest Standard	☒
ARMPAD*			
ATE2196	☐	Long Armpad (375mm)	☒
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ATE1350	☐	Neckrest: height and depth Adjustable	☒
SEAT HEIGHT SETTING & SEAT RISER*			
ATE1040	☐	Electrical Seat Riser (440mm - 690mm)	☒
SEAT TILT*			
ATE1014	☐	Electrical Tilt from 0° to 25°	☒
BACKREST RECLINER*			
ATE0975	☐	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	☒
LEG REST TYPE*			
ATE0503	☐	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	☒
FOOTPLATE TYPE*			
ATE0512	☐	Fixed Footplate	☒
LINX CONTROLS*			
ATE2204	☐	REM216 Drive, Actuators & Lights	☒
CONTROL BRACKETRY*			
ATE1707	☐	Swing Away RH Bracket	☐
ATE1708	☐	Swing Away LH Bracket	☐
SHROUD COLOUR*			
ATE0202	☐	Blue	☐
ATE0228	☐	Black	☐
ATE0204	☐	Transparent Cherry Red	☐

CASTORS*			
ATE0853	☐	3.00 - 6 Grey Puncture Proof (sealed)	◆
DRIVE WHEELS*			
ATE0773	☐	3.00 - 8 Grey Puncture Proof (sealed)	◆
MOTOR TYPE*			
ATE1909	☐	12 km/h (8mph) includes LED light system with indicators	◆
BATTERY*			
ATE1364	☐	73.5 Ah GEL	◆
REAR SUSPENSION*			
ATE1233	☐	Hitec spring X-plore	○
FRONT SUSPENSION*			
ATE1238	☐	Front suspension (suggested for up to 100kg)	○
ATE1239	☐	Front suspension (suggested for above 100kg)	○
CHARGER*			
ATE1343	☐	OptiCharge Charger 8A (200 - 250V)	◆
LIGHTS*			
ATZ1006	☐	Lights Included with Motor Package 10 & 12kmh	◆
ANTI-TIPPERS*			
ATE1205	☐	Anti-tipper (standard version)	◆
TRANSPORTATION*			
ATE1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory).	◆
FIXATION BELT*			
ATE2113	☐	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	◆
REAR VIEW MIRROR			
ATE1415	☐	Rear view Mirror RH - Must be selected with 10 km/h Motor	○
ATE1339	☐	Rear view Mirror LH - Must be selected with 10 km/h Motor	○

			TOTAL PRICE	
--	--	--	--------------------	----------------------

Non Conforming Options:

ATE2196	Long Armpad (375mm)
---------	---------------------



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):