



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**TDX SP2 HD M16 1001587 - Weekly Price
£51.84**

Standard Features

- Max user weight 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with long armpad
- Neckrest
- Powereed tilt 20°
- Manual settable backrest Recline 0° - 30° in 3.75° steps
- 80° hangers with fixed footplate
- LiNX control REM211 swing away bracket
- Solid Black tyres
- 6km/h (4mph) 4 pole motors
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,410



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARJ0000	☐ TDX SP2 Heavy Duty	£6,410
USER WEIGHT*		
ARJ1918	☐ User weight less 100 kg	☐
ARJ1919	☐ User weight 101 to 125 kg	☐
ARJ1921	☐ User weight 126 to 136 kg	☐
ARJ1922	☐ User weight 137 to 150 kg (NC 1092)	☐
ARJ1923	☐ User weight 151 to 160 kg (NC 1092)	☐
ARJ1924	☐ User weight 161 to 180 kg (NC 1092)	☐

SEATING TYPE*			
ARJ2100	☐	Modulite Seating unit	◆
SEAT WIDTH*			
ARJ2121	☐	Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2122	☐	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2123	☐	Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2124	☐	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	○
SEAT DEPTH RAIL LENGTH*			
ARJ2129	☐	Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138, 2196)	○
ARJ2130	☐	Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134)	○
SEAT DEPTH SETTING*			
ARJ2132	☐	410mm (16") (410mm - 510mm simple stepless adjustment from the side) (NC 2196)	○
ARJ2134	☐	460mm (18") (410mm - 510mm simple stepless adjustment from the side)	○
ARJ2136	☐	510mm (20") (410mm - 510mm simple stepless adjustment from the side)	○
ARJ2138	☐	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)	○
SEAT CUSHION*			
ARJ2142	☐	Matrx Flo-tech Lite	◆
BACKREST SUPPORT*			
ARJ2162	☐	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	◆
PUSH HANDLES FOR BACKPOSTS*			
ARJ2107	☐	Backposts WITH push handles	◆
BACKREST UPHOLSTERY / CUSHION*			
ARJ2173	☐	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	◆
UPHOLSTERY COLOUR*			
ARJ2152	☐	Black Black multi stretch PU coated fabric Upholstery	○
ARJ2155	☐	Black Spacer fabric	○
ARMREST*			
ARJ2191	☐	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	◆
ARMPAD*			
ARJ2196	☐	Long Armpad (375mm) (NC 2132)	◆
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ARJ1350	☐	Neckrest height and depth adjustable	◆
MANUAL / POWERED SEATING OPTIONS*			
ARJ1092	☐	Electric Tilt (0° - 20°) Seat Height 450mm	◆
ARJ1095	☐	Electric HD Tilt (0° - 20°) Seat Height 450mm (upto 180kg Max User Weight)	○

SEAT HEIGHT SETTING			
ARJ0354	☐	Seat height set at 450mm (17.75") (from floor to front of seat plate at 0° tilt)	☐
ARJ0364	☐	Seat height set at 475mm (18.75") (from floor to front of seat plate at 0° tilt)	☐
BACKREST RECLINER*			
ARJ0975	☐	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	☒
LEG REST TYPE*			
ARJ0526	☐	Side Mounted : Standard 80° hangers	☒
FOOTPLATE TYPE*			
ARJ0512	☐	Fixed Footplate	☒
LINX CONTROLS*			
ARJ2203	☐	REM211 Drive & Actuators	☒
CONTROL BRACKETRY*			
ARJ1707	☐	Swing Away RH Bracket	☐
ARJ1708	☐	Swing Away LH Bracket	☐
SHROUD COLOUR*			
ARJ2640	☐	Obsidian Black Frame Shrouds	☐
ARJ2641	☐	Ruby Red Frame Shrouds	☐
ARJ2642	☐	Purple Passion Frame Shrouds	☐
ARJ2643	☐	Peacock Blue Frame Shrouds	☐
ARJ2644	☐	Arctic Blue Frame Shrouds	☐
RIM COLOUR*			
ARJ2650	☐	Obsidian Black Frame Shrouds	☐
ARJ2651	☐	Ruby Red Frame Shrouds	☐
ARJ2652	☐	Purple Passion Frame Shrouds	☐
ARJ2653	☐	Peacock Blue Frame Shrouds	☐
ARJ2654	☐	Arctic Blue Frame Shrouds	☐
FORKS*			
ARJ2632	☐	Double sided fork	☒
CENTRE WHEELS (DRIVE WHEELS)*			
ARJ0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	☒
MOTOR TYPE*			
ARJ1282	☐	6 km/h (4mph) (4 pole)	☒
BATTERY*			
ARJ1260	☐	73.5 Ah (Gel)	☒
CHARGER*			
ARJ1343	☐	OptiCharge Charger 8A (200 - 250V)	☒
TRANSPORTATION KIT*			
ARJ1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	☒
FIXATION BELT*			
ARJ2113	☐	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	☒

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:

ARJ2132	410mm (16") (410mm - 510mm simple stepless adjustment from the side)
ARJ2134	460mm (18") (410mm - 510mm simple stepless adjustment from the side)
ARJ2138	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)
ARJ2196	Long Armpad (375mm)
ARJ1092	Electric Tilt (0° - 20°) Seat Height 450mm



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):