



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**TDX SP2 HD M17 1001588 - Weekly Price
£63.28**

Standard Features

- Max user weight 160kg (25 stone)
- Crash tested ISO 7176 -Part 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with long armpad
- Elan Headrest
- Seat lifter with electrical tilt (300 mm, 30°) with CoG shift
- Manual settable backrest Recline 0° - 30° in 3.75° steps
- Centre mounted Manual legrest
- LiNX control REM216 swing away bracket
- Solid Black tyres
- 12km/h (8mph) 4 pole motors
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- indicators, LED headlights and Mirror
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £8,939



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARJ0000	☐ TDX SP2 Heavy Duty	£8,939
USER WEIGHT*		
ARJ1918	☐ User weight less 100 kg	☐
ARJ1919	☐ User weight 101 to 125 kg	☐
ARJ1921	☐ User weight 126 to 136 kg	☐
ARJ1922	☐ User weight 137 to 150 kg	☐
SEATING TYPE*		
ARJ2100	☐ Modulite Seating unit	☒

TDX SP2 HD M17 1001588 - Weekly Price £63.28

SEAT WIDTH*			
ARJ2121		Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2122		Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2123		Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2124		Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	○
SEAT DEPTH RAIL LENGTH*			
ARJ2129		Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138, 2196)	○
ARJ2130		Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134)	○
SEAT DEPTH SETTING*			
ARJ2132		410mm (16") (410mm - 510mm simple stepless adjustment from the side) (NC 2196)	○
ARJ2134		460mm (18") (410mm - 510mm simple stepless adjustment from the side)	○
ARJ2136		510mm (20") (410mm - 510mm simple stepless adjustment from the side)	○
ARJ2138		560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)	○
SEAT CUSHION*			
ARJ2142		Matrx Flo-tech Lite	◆
BACKREST SUPPORT*			
ARJ2162		Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	◆
PUSH HANDLES FOR BACKPOSTS*			
ARJ2107		Backposts WITH push handles	◆
BACKREST UPHOLSTERY / CUSHION*			
ARJ2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion	◆
UPHOLSTERY COLOUR*			
ARJ2152		Black Black multi stretch PU coated fabric Upholstery	○
ARJ2155		Black Spacer fabric	○
ARMREST*			
ARJ2191		T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	◆
ARMPAD*			
ARJ2196		Long Armpad (375mm) (NC 2132)	◆
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ARJ2324		Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware	◆
MANUAL / POWERED SEATING OPTIONS*			
ARJ1094		Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift	◆
BACKREST RECLINER*			
ARJ0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
LEG REST TYPE*			
ARJ0556		Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate)	◆

FOOTPLATE TYPE*		
ARJ0559	☐	Angle Adjustable Footplate for Centre Mounted Hanger ◆
LINX CONTROLS*		
ARJ2204	☐	REM216 Drive, Actuators & Lights ◆
CONTROL BRACKETRY*		
ARJ1707	☐	Swing Away RH Bracket ○
ARJ1708	☐	Swing Away LH Bracket ○
SHROUD COLOUR*		
ARJ2640	☐	Obsidian Black Frame Shrouds ○
ARJ2641	☐	Ruby Red Frame Shrouds ○
ARJ2642	☐	Purple Passion Frame Shrouds ○
ARJ2643	☐	Peacock Blue Frame Shrouds ○
ARJ2644	☐	Arctic Blue Frame Shrouds ○
RIM COLOUR*		
ARJ2650	☐	Obsidian Black Frame Shrouds ○
ARJ2651	☐	Ruby Red Frame Shrouds ○
ARJ2652	☐	Purple Passion Frame Shrouds ○
ARJ2653	☐	Peacock Blue Frame Shrouds ○
ARJ2654	☐	Arctic Blue Frame Shrouds ○
FORKS*		
ARJ2632	☐	Double sided fork ◆
CENTRE WHEELS (DRIVE WHEELS)*		
ARJ0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert ◆
MOTOR TYPE*		
ARJ1224	☐	12 km/h (8mph) (includes light system and indicators) G-Trac recommended ◆
BATTERY*		
ARJ1260	☐	73.5 Ah (Gel) ◆
CHARGER*		
ARJ1343	☐	OptiCharge Charger 8A (200 - 250V) ◆
TRANSPORTATION KIT*		
ARJ1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) ◆
FIXATION BELT*		
ARJ2113	☐	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion ◆
REAR VIEW MIRROR*		
ARJ1415	☐	Rear view Mirror RH - Must be selected with 10 km/h Motor ○
ARJ1339	☐	Rear view Mirror LH - Must be selected with 10 km/h Motor ○

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARJ2132	410mm (16") (410mm - 510mm simple stepless adjustment from the side)
ARJ2134	460mm (18") (410mm - 510mm simple stepless adjustment from the side)
ARJ2138	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)
ARJ2196	Long Armpad (375mm)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):