



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 Low Rider M14 712116 - Weekly Price £60.64

Standard Features

- Max user weight 160kg (25.2 Stone)
- Crash tested ISO7176-Part 19
- Modulite seating
- Matrix Flo-tech Lite cushion
- Without push handles
- Flex 3 Adjustable backrest plates with laguna backrest cushion
- Elan Headrest
- Standard T-armrests, detachable with Long armpad
- Seat Lifter with powered Tilt (300mm, 30°) CoG Shift
- Seat height either 403mm
- Manual settable Recliner 0° - 30° in 3.75° steps
- Standard 80° Hangers with Angle and Depth Adjustable footplates(adapt)
- LiNX control REM211 with swing away bracket
- Solid Black tyres
- 6km/h (4mph) 4 Pole motors
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £8,216



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code Description

BASE MODEL*

ARQ0072	☐	TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height	£8,216
---------	--------------------------	--	--------

USER WEIGHT*

ARQ1918	☐	User weight less 100 kg	☐
ARQ1919	☐	User weight 101 to 125 kg	☐
ARQ1921	☐	User weight 126 to 136 kg	☐
ARQ1922	☐	User weight 137 to 150 kg	☐
ARQ1923	☐	User weight 151 to 160 kg	☐

SEATING TYPE*			
ARQ2100	☐	Modulite Seating Unit	◆
SEAT FRAME*			
ARQ2111	☐	Telescopic Seat Frame with metal seat plate	◆
SEAT WIDTH*			
ARQ2122	☐	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19))	○
ARQ2123	☐	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21"))	○
ARQ2124	☐	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23"))	○
SEAT DEPTH RAIL LENGTH*			
ARQ2129	☐	Standard Depth Package (16" - 20", 410mm - 510mm)	○
ARQ2130	☐	Extended Depth Package (20" - 22", 510mm - 560mm) (NC 2134)	○
SEAT DEPTH*			
ARQ2132	☐	410mm (16") (410mm - 510mm) (NC 2196)	○
ARQ2134	☐	460mm (18") (410mm - 510mm)	○
ARQ2136	☐	510mm (20") (410mm - 510mm)	○
ARQ2138	☐	560mm (22") (510mm - 560mm)	○
SEAT CUSHION*			
ARQ2142	☐	Matrx Flo-tech Lite	◆
BACKREST SUPPORT*			
ARQ2162	☐	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	◆
BACKREST UPHOLSTERY / CUSHION*			
ARQ2173	☐	Flex3 Adjustable Backrest Plates with Laguna Backrest cushion	◆
UPHOLSTERY TYPE*			
ARQ2152	☐	Black multi stretch PU coated fabric	○
ARQ2155	☐	Black Spacer fabric	○
ARMREST*			
ARQ2191	☐	T-Armrest Standard, Range of adjustment	◆
ARMPAD*			
ARQ2196	☐	Long Armpad (375mm)	◆
HEADREST / NECKREST*			
ARQ2324	☐	Elan Standard Pad 355mm (10") wide	◆
MANUAL / POWERED SEATING OPTIONS*			
ARQ1094	☐	Seat lifter with Electrical Tilt	◆
BACKREST RECLINER*			
ARQ0975	☐	Manual Settable Recliner	◆
LEG REST TYPE*			
ARQ0526	☐	Side Mounted : Standard 80° hangers	◆
FOOTPLATE TYPE*			
ARQ0512	☐	Fixed Footplate	◆
LINX CONTROLS*			
ARQ2203	☐	REM211 Drive & Actuators	◆
CONTROL BRACKETRY*			
ARQ1707	☐	Swing Away RH Bracket	○
ARQ1708	☐	Swing Away LH Bracket	○

SHROUD COLOUR*			
ARQ2640	☐	Obsidian Black Frame Shrouds	☐
ARQ2641	☐	Ruby Red Frame Shrouds	☐
ARQ2642	☐	Purple Passion Frame Shrouds	☐
ARQ2643	☐	Peacock Blue Frame Shrouds	☐
ARQ2644	☐	Arctic Blue Frame Shrouds	☐
RIM COLOUR*			
ARQ2650	☐	Obsidian Black Rim Inserts	☐
ARQ2651	☐	Ruby Red Rim Inserts	☐
ARQ2652	☐	Purple Passion Rim Inserts	☐
ARQ2653	☐	Peacock Blue Rim Inserts	☐
ARQ2654	☐	Arctic Blue Rim Inserts	☐
FORKS*			
ARQ2632	☐	Double sided fork	☒
CENTRE WHEELS (DRIVE WHEELS)*			
ARQ0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	☒
MOTOR TYPE*			
ARQ1282	☐	6 km/h (4mph) (4 pole)	☒
CHARGER*			
ARQ1343	☐	OptiCharge Charger 8A (200 - 250V)	☒
LIGHTS*			
ARQ1254	☐	Without Light system	☒
TRANSPORTATION KIT*			
ARQ1494	☐	4-Point tie-down	☒
FIXATION BELT*			
ARQ2113	☐	Lapbelt (adjustable both sides)	☒
TOTAL PRICE			

Non Conforming Options:

ARQ2134	460mm (18") (410mm - 510mm)
ARQ2196	Long Armpad (375mm)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):