



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**TDX SP2 Max M15 1001086 - Weekly Price
£49.14**

Standard Features

- Max user weight 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Tension Adjustable backrest
- Armrest with reinforced depth adjustable bracket with Short armpad
- Powered 20 degree seat Tilt
- Manual settable backrest Recline
- VARI F hangers with fixed footplate
- LiNX Control REM211 with swing away bracket
- Solid Black Tyres
- 6 km/h (4mph) 4 Pole Motors
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £5,882



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARH0000	TDX SP2 Max	£5,882
ARZ0000	Special TDX SP2 Max (Please provide details in special section)	POA
SEATING TYPE*		
ARH1145	MAX Seating system	☒
SEAT WIDTH SETTING*		
ARH0378	Seat Width of Seat Frame 550mm (+100mm adjustment in arms)	☐
ARH0380	Seat Width of Seat Frame 610mm (+100mm adjustment in arms)	☐

SEAT DEPTH SETTING SETTABLE IN 3 STEPS*			
ARH0396	☐	Seat Depth 560mm	☐
ARH0397	☐	Seat Depth 595mm	☐
ARH0398	☐	Seat Depth 630mm	☐
SEAT CUSHION*			
ARH1484	☐	MAX Seat Cushion	☐
ADAPTOR TO REDUCE THE SEAT DEPTH BY 130MM			
ARH0399	☐	Adaptor to shorten the seat depth by 130mm (Required if seat width 490mm is selected in combination with seat depth 595mm or 630mm)	☐
BACKREST UPHOLSTERY*			
ARH1106	☐	Tension Adjustable	☒
UPHOLSTERY COLOUR*			
ARH0278	☐	Black multi stretch PU coated fabric	☒
ARMREST INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*			
ARH0635	☐	Fixed Armrest, 250mm - 340mm tool adjustable in height, detachable	☒
ARMPADS*			
ARH0649	☐	Long armpad (375mm)	☒
PUSH BAR*			
ARH1307	☐	Push bar	☒
HEADREST (Adequate headrest option must be selected)*			
ARH1432	☐	Headrest: height, depth and sideways adjustable (Azalea type)	☒
MANUAL / POWERED SEATING OPTIONS*			
ARH1095	☐	Electric heavy duty basic Tilt (0° - 20°) select either 470mm or 495mm seat to ground height. (Max user weight 180kg)	☒
SEAT HEIGHT SETTING*			
ARH1024	☐	Fixed seat height 470mm (18.5")	☐
ARH1046	☐	Fixed seat height 495mm (19.5")	☐
BACKREST RECLINER*			
ARH0958	☐	Manually settable recliner by spindle (0° - 25°)	☒
LEG REST TYPE*			
ARH0503	☐	Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	☒
FOOTPLATE TYPE*			
ARH0512	☐	Fixed Footplate	☒
LINX CONTROLS*			
ARH2203	☐	REM211 Drive & Actuators	☒
CONTROL BRACKETRY*			
ARH1707	☐	Swing Away RH Bracket	☐
ARH1708	☐	Swing Away LH Bracket	☐
SHROUD COLOUR*			
ARH2640	☐	Obsidian Black Frame Shrouds	☐
ARH2641	☐	Ruby Red Frame Shrouds	☐
ARH2642	☐	Purple Passion Frame Shrouds	☐
ARH2643	☐	Peacock Blue Frame Shrouds	☐
ARH2644	☐	Arctic Blue Frame Shrouds	☐

RIM COLOUR*			
ARH2650	☐	Obsidian Black Rim Inserts	☐
ARH2651	☐	Ruby Red Rim Inserts	☐
ARH2652	☐	Purple Passion Rim Inserts	☐
ARH2653	☐	Peacock Blue Rim Inserts	☐
ARH2654	☐	Arctic Blue Rim Inserts	☐
FORKS*			
ARH2632	☐	Double sided fork	☒
CENTRE WHEELS (DRIVE WHEELS)*			
ARH0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam	☒
MOTOR TYPE*			
ARH1282	☐	6 km/h (4mph) (4 pole)	☒
BATTERY*			
ARH1260	☐	73.5 Ah	☒
CHARGER*			
ARH1343	☐	OptiCharge Charger 8A (200 - 250V)	☒
LIGHTS*			
ARH1254	☐	Without light system	☒
TRANSPORTATION KIT			
ARH1494	☐	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	☒
FIXATION BELT*			
ARH1329	☐	MAX Lapbelt (adjustable both sides) fitted to base of backrest	☒

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):