



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**TDX SP2 NB M3 735801 - Weekly Price
£65.66**

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Modulte Seating
- Matrix Flo-tech cushion
- Flex3 Adjustable Plates with Luguna back cushion
- Standard T-armrests, detachable with Long armpad
- Seat Lifter with electrical Tilt (300mm, 30°) with CoG shift
- Manual settable backrest Recliner (0° - 30° in 3.75° steps)
- Manual centre mounted hanger, with angle adjustable footplates
- LiNX REM216 with swing away bracket
- SolidBlack tyres
- 10km/h (6mph) 4 pole motors (including indicators, LED headlights and mirror)
- 4 Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock

Standard Base price: £8,927



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARM0000	☐ TDX SP2 Narrow Base	£8,927
USER WEIGHT*		
ARM1920	☐ User weight up to 125kg	☐
ARM1921	☐ User weight 126kg to 136kg	☐
ARM1922	☐ User weight 137kg to 150kg	☐
SEATING TYPE*		
ARM2100	☐ Modulte Seating Unit	☒

SEAT FRAME*			
ARM2111	☐	Telescopic Seat Frame	◆
SEAT WIDTH*			
ARM2121	☐	380mm (15")	○
ARM2122	☐	430mm (17")	○
ARM2123	☐	480mm (19")	○
ARM2124	☐	530mm (21")	○
SEAT DEPTH RAIL LENGTH*			
ARM2129	☐	Standard 410mm - 510mm (NC 2138)	○
ARM2190	☐	Extended 510mm - 560mm (NC 2132, 2134)	○
SEAT DEPTH*			
ARM2132	☐	410mm (16") (NC 2196)	○
ARM2134	☐	460mm (18")	○
ARM2136	☐	510mm (20")	○
ARM2138	☐	560mm (22")	○
SEAT CUSHION*			
ARM2142	☐	Matrx Flo-Tech Lite	◆
BACKREST SUPPORT*			
ARM2162	☐	Tall 540mm	◆
BACKREST UPHOLSTERY / CUSHION*			
ARM2173	☐	Flex3 Backrest plates with Laguna cushion	◆
UPHOLSTERY TYPE*			
ARM2152	☐	Black Black multi stretch PU coated fabric Upholstery	○
ARM2155	☐	Black Spacer fabric	○
ARMREST*			
ARM2191	☐	T-Armrest standard	◆
ARMPAD*			
ARM2196	☐	Long 375mm Armpad	◆
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ARM2324	☐	Elan standard pad 255mm (10")	◆
MANUAL / POWERED SEATING OPTIONS*			
ARM1094	☐	Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift	◆
BACKREST RECLINER*			
ARM0975	☐	Manual settable Recliner	◆
LEG REST TYPE*			
ARM0556	☐	Centre mounted : Manual	◆
FOOTPLATE TYPE*			
ARM0559	☐	Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger only	◆
LINX CONTROLS*			
ARM2204	☐	REM216	◆
CONTROL BRACKETRY*			
ARM1707	☐	Swing away RH Bracket	○
ARM1708	☐	Swing away LH Bracket	○

SHROUD COLOUR*			
ARM2640	☐	Obsidian Black Frame Shrouds	☐
ARM2641	☐	Ruby Red Frame Shrouds	☐
ARM2642	☐	Purple Passion Frame Shrouds	☐
ARM2643	☐	Peacock Blue Frame Shrouds	☐
ARM2644	☐	Arctic Blue Frame Shrouds	☐
RIM COLOUR*			
ARM2650	☐	Obsidian Black Rim Inserts	☐
ARM2651	☐	Ruby Red Rim Inserts	☐
ARM2652	☐	Purple Passion Rim Inserts	☐
ARM2653	☐	Peacock Blue Rim Inserts	☐
ARM2654	☐	Arctic Blue Rim Inserts	☐
FORKS*			
ARM2632	☐	Double sided	☒
CENTRE WHEELS (DRIVE WHEELS)*			
ARM0705	☐	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	☒
MOTOR TYPE*			
ARM1908	☐	10 km/h (6mph)	☒
BATTERY*			
ARM1257	☐	50 Ah Gel	☒
CHARGER*			
ARM1343	☐	OptiCharge 8A (200 - 250V)	☒
TRANSPORTATION KIT*			
ARM1494	☐	4-Point tie-down	☒
FIXATION BELT*			
ARM1329	☐	Lap Belt fitted to base of backrest	☐
ARM2113	☐	Lap Belt adjustable mounting on seat rail extrusion	☒
REAR VIEW MIRROR*			
ARM1415	☐	RH Mirror	☐
ARM1339	☐	LH Mirror	☐

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:

ARM2132	410mm (16")
ARM2134	460mm (18")
ARM2138	560mm (22")
ARM2196	Long 375mm Armpad



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):