



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Motability Retail Prescription Form

February 2026  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## AVIVA RX20 Standard Seat M20 707043 - Weekly Price £43.77

### Standard Features

- Max user weight 136kg
- Crash tested ISO 7176 -Part 19
- One piece seat plate
- Flo-tech lite cushion
- Tension adjustable backrest
- Neckrest
- Standard T-Armrests, detachable with Long armpad
- Powered "fixed pivot" Tilt 20°
- Manual settable backrest Recline (0° - 30°, in 3.75° steps)
- 80° hangers with fixed footplate
- Remote LiNX DLX-REM211 with swing away Bracket
- 8" Castor, puncture proof solid tyre
- Solid black tyres
- 6 km/h 4mph motors
- 60Ah batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock



**Standard Base price: £4,280**

Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA    special request

Customer weight =                      KG

Code	Description	
<b>BASE MODEL*</b>		
JBB0000	AVIVA RX 20 Base Model Max user weight up to 136 kg includes 75A Power Module + minimum 52Ah Batteries	£4,280
<b>CHASSIS WIDTH*</b>		
JBL2810	Chassis Narrow (560 mm/22")	◆
<b>SEATING TYPE*</b>		
JBB2105	Modulite Standard Seat Standard Seat Plate and Backrest Frame in fixed Width	◆
<b>SEAT FRAME*</b>		
JBB2110	One Piece Seat Plate	◆

<b>SEAT WIDTH*</b>			
JBB2121	<input type="checkbox"/>	380 mm Seat Plate	<input type="radio"/>
JBB2122	<input type="checkbox"/>	430 mm Seat Plate	<input type="radio"/>
<b>SEAT DEPTH*</b>			
JBB2132	<input type="checkbox"/>	Seat Depth 410 mm	<input type="radio"/>
JBB2134	<input type="checkbox"/>	Seat Depth 460 mm	<input type="radio"/>
JBB2136	<input type="checkbox"/>	Seat Depth 510 mm	<input type="radio"/>
<b>SEAT CUSHION*</b>			
JBB2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	<input checked="" type="radio"/>
<b>BACKREST SUPPORT*</b>			
JBB2161	<input type="checkbox"/>	Backrest frame, detachable, height 480 mm For Tension Adjustable Upholstery/Matrx Backs	<input checked="" type="radio"/>
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
JBB2171	<input type="checkbox"/>	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm	<input checked="" type="radio"/>
<b>UPHOLSTERY TYPE*</b>			
JBB2151	<input type="checkbox"/>	Black Fabric Upholstery (Nylon w 2141; 2-way-stretch material w 2142, 2143)	<input checked="" type="radio"/>
<b>ARMREST (Pair)*</b>			
JBB2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	<input checked="" type="radio"/>
<b>ARMPAD (Pair)*</b>			
JBB2196	<input type="checkbox"/>	Long Armpad (375 mm)	<input checked="" type="radio"/>
<b>HEADREST / NECKREST (Adequate headrest option must be selected)*</b>			
JBB1350	<input type="checkbox"/>	Neckrest depth adjustable	<input checked="" type="radio"/>
<b>MANUAL / POWERED SEATING OPTIONS*</b>			
JBB1092	<input type="checkbox"/>	Powered "fixed pivot" Tilt 20°	<input checked="" type="radio"/>
<b>SEAT TO GROUND HEIGHT*</b>			
JBB2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	<input type="radio"/>
JBB2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	<input type="radio"/>
JBB2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	<input type="radio"/>
<b>BACKREST RECLINER*</b>			
JBB0975	<input type="checkbox"/>	Manual Settable Recliner (in 3.75° steps) 90° - 120°	<input checked="" type="radio"/>
<b>LEG REST TYPE*</b>			
JBB0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	<input checked="" type="radio"/>
<b>FOOTPLATE TYPE*</b>			
JBB0512	<input type="checkbox"/>	Fixed Footplate pair	<input checked="" type="radio"/>
<b>LINX CONTROLS*</b>			
JBB2203	<input type="checkbox"/>	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	<input checked="" type="radio"/>
<b>CONTROL BRACKETRY*</b>			
JBB1707	<input type="checkbox"/>	Swing Away Bracket RH	<input type="radio"/>
JBB1708	<input type="checkbox"/>	Swing Away Bracket LH	<input type="radio"/>

SHROUD COLOUR*			
JBB2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	<input type="radio"/>
JBB2872	<input type="checkbox"/>	Mid range - Silver Shroud	<input type="radio"/>
JBB2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	<input type="radio"/>
JBB2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	<input type="radio"/>
JBB2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	<input type="radio"/>
RIM COLOUR*			
JBB2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	<input type="radio"/>
JBB2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	<input type="radio"/>
JBB2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	<input type="radio"/>
JBB2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	<input type="radio"/>
JBB2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	<input type="radio"/>
FORKS*			
JBB2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	<input checked="" type="radio"/>
CASTORS*			
JBB0828	<input type="checkbox"/>	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	<input checked="" type="radio"/>
DRIVE WHEELS*			
JBB0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	<input checked="" type="radio"/>
C.T.C SUSPENSION*			
JBB2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	<input checked="" type="radio"/>
SUSPENSION SETTING*			
JBB1910	<input type="checkbox"/>	Soft (suggested for up to 70 kg User Weight)	<input checked="" type="radio"/>
JBB1912	<input type="checkbox"/>	Medium (suggested for up to 90 kg User Weight)	<input type="radio"/>
JBB1911	<input type="checkbox"/>	Firm (suggested for up to 110 kg User Weight)	<input type="radio"/>
JBB1915	<input type="checkbox"/>	Very firm (suggested for up to 136 kg User Weight)	<input type="radio"/>
MOTOR TYPE*			
JBB1282	<input type="checkbox"/>	6 km/h	<input checked="" type="radio"/>
BATTERY TYPE*			
JBB1259	<input type="checkbox"/>	60 Ah Gel	<input checked="" type="radio"/>
CHARGER*			
JBB1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	<input checked="" type="radio"/>
ANTI-TIPPERS*			
JBB1205	<input type="checkbox"/>	Anti-tipper	<input checked="" type="radio"/>
TRANSPORTATION KIT*			
JBB1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	<input checked="" type="radio"/>
LAP BELT*			
JBB1329	<input type="checkbox"/>	Lap belt with Auto Lock Standard mount on backrest	<input checked="" type="radio"/>

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):