



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

February 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 Modulite M24 707026 - weekly Price £58.85

Standard Features

- Max user weight 136kg
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Tension adjustable backrest
- Standard T-Armrests, detachable with long armpad
- Manually settable Seat Tilt 0 - 15°
- Manual settable backrest Recline (0° - 30°, in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX Control REM216 with swing away Bracket
- 9" Castor, puncture proof solid tyre
- Solid inlay Dual density foam insert
- 12 km/h Ultra High Speed Package
- LED Light system with Indicators, 30 Lux headlights
- 73,5 Ah Gel Batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



Standard Base price: £7,424

Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBF0001	AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries	£7,424
CHASSIS WIDTH*		
JBL2810	Chassis Narrow (560 mm/22")	<input checked="" type="checkbox"/>
SEATING TYPE*		
JBF2100	Invacare® Modulite Seat Unit	<input checked="" type="checkbox"/>

SEAT FRAME*			
JBF2111		Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	○
JBF2112		Telescopic Seat Frame + Sling Seat with Side Rails width 380 to 530 mm; + 30 mm seat height Additional 50 mm by arms	○
SEAT WIDTH*			
JBF2121		380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests	○
JBF2122		430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	○
JBF2123		480 mm Seat Plate	○
JBF2124		530 mm Seat Plate	○
SEAT DEPTH*			
JBF2132		Seat Depth 410 mm	○
JBF2134		Seat Depth 460 mm	○
JBF2136		Seat Depth 510 mm	○
SEAT CUSHION*			
JBF2142		Matrx Flo-tech Cushion Flo-tech Lite	◆
BACKREST SUPPORT*			
JBF2163		Backposts short (480 mm) For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	◆
JBF2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	○
BACKREST UPHOLSTERY / CUSHION*			
JBF2171		Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm	◆
UPHOLSTERY TYPE*			
JBF2151		Black Fabric Upholstery (Nylon w 2141; 2-way-stretch material w 2142, 2143)	◆
ARMREST (Pair)*			
JBF2191		T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
ARMPAD (Pair)*			
JBF2196		Long Armpad (375 mm)	◆
MANUAL / POWERED SEATING OPTIONS*			
JBF1091		Manually settable Seat Tilt 0 - 15°	◆
SEAT TO GROUND HEIGHT*			
JBF2857		Seat-to-Floor Height 435 mm	○
JBF2858		Seat-to-Floor Height 460 mm	○
JBF2859		Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBF0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
LEG REST TYPE*			
JBF0526		Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
FOOTPLATE TYPE*			
JBF0512		Fixed Footplate pair	◆

LINX CONTROLS*			
JBF2204		Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	◆
CONTROL BRACKETRY*			
JBF1707		Swing Away Bracket RH	○
JBF1708		Swing Away Bracket LH	○
SHROUD COLOUR*			
JBF2871		Mid range - Black/Anthracite Shroud	○
JBF2872		Mid range - Silver Shroud	○
JBF2873		Mid range - Bengal Red Shroud	○
JBF2874		Mid range - Mercury Blue Shroud	○
JBF2875		Mid range - Midnight Purple Shroud	○
RIM COLOUR*			
JBF2882		Mid range - Black Rim Inserts	○
JBF2883		Mid range - Silver Rim Inserts	○
JBF2884		Mid range - Bengal Red Rim Inserts	○
JBF2885		Mid range - Mercury Blue Rim Inserts	○
JBF2886		Mid range - Midnight Purple Rim Inserts	○
FORKS*			
JBF2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
CASTORS*			
JBF0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
DRIVE WHEELS*			
JBF0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
C.T.C SUSPENSION*			
JBF2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBF1910		Soft (suggested for up to 70 kg User Weight)	○
JBF1912		Medium (suggested for up to 90 kg User Weight)	○
JBF1911		Firm (suggested for up to 110 kg User Weight)	○
JBF1915		Very firm (suggested for up to 136 kg User Weight)	○
MOTOR TYPE*			
JBF1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	◆
BATTERY TYPE*			
JBF1260		73,5 Ah Gel	◆
CHARGER*			
JBF1343		Opti Charge Charger 8A (IP44)	◆
ANTI-TIPPERS*			
JBF1205		Anti-tipper	◆
TRANSPORTATION KIT*			
JBF1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	◆
LAP BELT*			
JBF1329		Lap belt with Auto Lock Standard mount on backrest	◆

REAR VIEW MIRROR*			
JBF1415		Rear view mirror RH	○
JBF1339		Rear view mirror LH	○

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):