



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

February 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 Modulite M25 707020 - Weekly Price £69.34

Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flex3 Tension Adjustable Upholstery
- T-Armrests, detachable with Long armpad
- Neckrest depth adjustable
- Powered "fixed pivot" Tilt 20°
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- Standard hanger 80° pair with Fixed Footplates
- LiNX control REM216 with swing away bracket
- 8" Castor, puncture proof solid tyre
- Solid inlay Dual density foam insert
- 12 km/h Ultra High Speed Package with LED Light system with Indicators, 30 Lux headlights & Mirror
- 73Ah batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £8,062



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code **Description**

BASE MODEL*

JBG0001	<input type="checkbox"/>	AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries	£8,062
---------	--------------------------	--	--------

CHASSIS WIDTH*

JBL2810	<input type="checkbox"/>	Chassis Narrow (560 mm/22")	<input checked="" type="checkbox"/>
---------	--------------------------	-----------------------------	-------------------------------------

SEATING TYPE*

JBG2100	<input type="checkbox"/>	Invacare® Modulite Seat Unit	<input checked="" type="checkbox"/>
---------	--------------------------	------------------------------	-------------------------------------

SEAT FRAME*			
JBG2111	<input type="checkbox"/>	Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	◆
SEAT WIDTH*			
JBG2121	<input type="checkbox"/>	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests	○
JBG2122	<input type="checkbox"/>	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	○
JBG2123	<input type="checkbox"/>	480 mm Seat Plate	○
JBG2124	<input type="checkbox"/>	530 mm Seat Plate	○
SEAT DEPTH*			
JBG2132	<input type="checkbox"/>	Seat Depth 410 mm	○
JBG2134	<input type="checkbox"/>	Seat Depth 460 mm	○
JBG2136	<input type="checkbox"/>	Seat Depth 510 mm	○
SEAT CUSHION*			
JBG2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	◆
BACKREST SUPPORT*			
JBG2162	<input type="checkbox"/>	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	◆
PUSH HANDLES FOR BACKPOSTS			
JBG2108	<input type="checkbox"/>	Backposts without push handles	◆
BACKREST UPHOLSTERY / CUSHION*			
JBG2179	<input type="checkbox"/>	Flex3 Tension Adjustable Upholstery height 540 mm (seat plate to top of upholstery) (NC 2155)	◆
UPHOLSTERY TYPE*			
JBG2152	<input type="checkbox"/>	Black multi stretch PU coated fabric Upholstery	○
JBG2155	<input type="checkbox"/>	Black Spacer Fabric Upholstery	○
ARMREST (Pair)*			
JBG2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
ARMPAD (Pair)*			
JBG2196	<input type="checkbox"/>	Long Armpad (375 mm)	◆
HEADREST / NECKREST (Adequate headrest option must be selected)*			
JBG1350	<input type="checkbox"/>	Neckrest depth adjustable	◆
MANUAL / POWERED SEATING OPTIONS*			
JBG1092	<input type="checkbox"/>	Powered "fixed pivot" Tilt 20°	◆
SEAT TO GROUND HEIGHT*			
JBG2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	○
JBG2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	○
JBG2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBG0975	<input type="checkbox"/>	Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
LEG REST TYPE*			
JBG0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
FOOTPLATE TYPE*			
JBG0512	<input type="checkbox"/>	Fixed Footplate pair	◆
LINX CONTROLS*			
JBG2204	<input type="checkbox"/>	DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	◆

CONTROL BRACKETRY*			
JBG1707	<input type="checkbox"/>	Swing Away Bracket RH	<input type="radio"/>
JBG1708	<input type="checkbox"/>	Swing Away Bracket LH	<input type="radio"/>
SHROUD COLOUR*			
JBG2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	<input type="radio"/>
JBG2872	<input type="checkbox"/>	Mid range - Silver Shroud	<input type="radio"/>
JBG2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	<input type="radio"/>
JBG2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	<input type="radio"/>
JBG2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	<input type="radio"/>
RIM COLOUR*			
JBG2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	<input type="radio"/>
JBG2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	<input type="radio"/>
JBG2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	<input type="radio"/>
JBG2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	<input type="radio"/>
JBG2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	<input type="radio"/>
FORKS*			
JBG2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	<input checked="" type="radio"/>
CASTORS*			
JBG0836	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	<input checked="" type="radio"/>
DRIVE WHEELS*			
JBG0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	<input checked="" type="radio"/>
C.T.C SUSPENSION*			
JBG2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	<input checked="" type="radio"/>
SUSPENSION SETTING*			
JBG1910	<input type="checkbox"/>	Soft (suggested for up to 70 kg User Weight)	<input type="radio"/>
JBG1912	<input type="checkbox"/>	Medium (suggested for up to 90 kg User Weight)	<input type="radio"/>
JBG1911	<input type="checkbox"/>	Firm (suggested for up to 110 kg User Weight)	<input type="radio"/>
JBG1915	<input type="checkbox"/>	Very firm (suggested for up to 136 kg User Weight)	<input type="radio"/>
MOTOR TYPE*			
JBG1909	<input type="checkbox"/>	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	<input checked="" type="radio"/>
BATTERY TYPE*			
JBG1260	<input type="checkbox"/>	73,5 Ah Gel	<input checked="" type="radio"/>
CHARGER*			
JBG1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	<input checked="" type="radio"/>
ANTI-TIPPERS*			
JBG1205	<input type="checkbox"/>	Anti-tipper	<input checked="" type="radio"/>
TRANSPORTATION KIT*			
JBG1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	<input checked="" type="radio"/>
LAP BELT*			
JBG2113	<input type="checkbox"/>	Lap belt with Auto Lock Flexible mount in side rails	<input checked="" type="radio"/>
REAR VIEW MIRROR			
JBG1415	<input type="checkbox"/>	Rear view mirror RH	<input type="radio"/>
JBG1339	<input type="checkbox"/>	Rear view mirror LH	<input type="radio"/>

	TOTAL PRICE

Non Conforming Options:

JBG2155	Black Spacer Fabric Upholstery
---------	--------------------------------



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):