



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Motability Retail Prescription Form

January 2026  
Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## E-Fix & Action 3NG M82 707362 - Weekly Price £51.22

### Standard Features

- Action 3NG + E-Fix
- Max user weight 125kg (19.7 stone)
- Crash tested ISO- 7176-Part 19
- Seat widths 380mm to 505mm (15" -20")
- Seat depths 400mm or 450mm (16" or 18")
- Seat heights 460mm to 510mm (18" -20")
- Nylon with hook and loop strips seat upholstery
- Fixed 7°Backrest 2 heights 400mm/430mm
- Lumbar support padded backrest with shaped cover
- 80°swing in/swing out hangers
- Flip up footplates
- Removable full arms
- Basic 50mm foam cushion
- Adjustable calf strap
- Choice of 8 standard colours
- 8" Black solid front wheels, with 24" Black solid tyre spoked rim
- Aluminium Pre-anodized handrim
- Quick release axle
- Push to lock brakes
- Transportation kit (including buckle belt)
- Alber Efix E35 1566293 (24" wheels)
- Brackets
- Alber anti-tippers
- Adapter Kit
- Mobility Plus package
- Fitting at Invacare



Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA special request

Customer weight =                      KG

| Code                 | Description                                     |             |
|----------------------|---|-------------|
| <b>BASE PRODUCT*</b> |   |             |
| CEE0003              | <input type="checkbox"/> ACTION 3NG SELF PROPEL | (NC 1195) ◆ |
| <b>PROPULSION*</b>   |   |             |
| CEE9259              | <input type="checkbox"/> E-Fix E36              | £6,086      |

| SEAT WIDTH*          |  |   |                       |   |
|----------------------|--|---|-----------------------|---|
| CEE0180              |  | 15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)   |                       | ○ |
| CEE0190              |  | 16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)   |                       | ○ |
| CEE0200              |  | 17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)   |                       | ◆ |
| CEE0210              |  | 18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)   |                       | ○ |
| CEE0220              |  | 19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)   |                       | ○ |
| CEE0230              |  | 20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)   |                       | ○ |
| SEAT DEPTH*          |  |   |                       |   |
| CEE0240              |  | 16" / 400 MM SEAT DEPTH   | (NC 0879)             | ○ |
| CEE0250              |  | 18" / 450 MM SEAT DEPTH   | (NC 0878)             | ◆ |
| SEAT HEIGHT*         |  |   |                       |   |
| CEE0260              |  | 20" / 510 MM SEAT HEIGHT  |                       | ○ |
| CEE0270              |  | 19" / 485 MM SEAT HEIGHT  |                       | ○ |
| CEE0280              |  | 18" / 460 MM SEAT HEIGHT  |                       | ◆ |
| SEAT UPHOLSTERY*     |  |   |                       |   |
| CEE0330              |  | NYLON W/O HOOK AND LOOP STRIPS  |                       | ◆ |
| BACKREST*            |  |   |                       |   |
| CEE0410              |  | FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES  |                       | ◆ |
| BACKREST UPHOLSTERY* |  |   |                       |   |
| CEE0474              |  | LUMBAR SUPPORT BI MATERIAL PADDED COVER   |                       | ◆ |
| LEGREST HANGER RH*   |  |   |                       |   |
| CEE0510              |  | 80° RH SWING IN/SWING OUT   |                       | ◆ |
| LEGREST HANGER LH*   |  |   |                       |   |
| CEE0560              |  | 80° LH SWING IN/SWING OUT   |                       | ◆ |
| FOOTREST RH*         |  |   |                       |   |
| CEE0600              |  | FLIP UP RH FOOTPLATE  |                       | ◆ |
| FOOTREST LH*         |  |   |                       |   |
| CEE0620              |  | FLIP UP LH FOOTPLATE  |                       | ◆ |
| ARMREST RH*          |  |   |                       |   |
| CEE0640              |  | REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED   | (NC 0760, 0839)       | ○ |
| CEE0650              |  | REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED  | (NC 0760, 0839)       | ○ |
| CEE0660              |  | SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT   | (NC 0740, 0750, 0839) | ◆ |
| ARMREST LH*          |  |   |                       |   |
| CEE0740              |  | REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED   |                       | ○ |
| CEE0750              |  | REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED  |                       | ○ |
| CEE0760              |  | SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT   |                       | ◆ |
| CEE0839              |  | TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) | (NC 1200, 1345)       | ○ |

|                             |                          |  |   |
|-----------------------------|--------------------------|--|---|
| <b>SEAT CUSHION</b>         |                          |  |   |
| CEE0849                     | <input type="checkbox"/> | BASIC 50 MM FOAM CUSHION NYLON AND ANTISLIP COVER                    | ◆ |
| <b>SEAT CUSHION DEPTH</b>   |                          |  |   |
| CEE0878                     | <input type="checkbox"/> | FOR 400 MM SEAT DEPTH  | ○ |
| CEE0879                     | <input type="checkbox"/> | FOR 450 MM SEAT DEPTH  | ○ |
| <b>COLOUR*</b>              |                          |  |   |
| CEE0961                     | <input type="checkbox"/> | SLATE GREY (RAL7024 FINE TEXTURE)                                    | ○ |
| CEE0964                     | <input type="checkbox"/> | PEARL WHITE  | ○ |
| CEE0965                     | <input type="checkbox"/> | ELECTRIC GREEN   | ○ |
| CEE0966                     | <input type="checkbox"/> | MATT BLACK   | ○ |
| CEE0974                     | <input type="checkbox"/> | BLUE   | ○ |
| CEE0976                     | <input type="checkbox"/> | DARK PURPLE  | ○ |
| CEE0969                     | <input type="checkbox"/> | RED MATT   | ○ |
| CEE0972                     | <input type="checkbox"/> | WATER BLUE   | ○ |
| CEE0980                     | <input type="checkbox"/> | ANTHRACITE BLACK   | ○ |
| <b>CASTOR HOUSING*</b>      |                          |  |   |
| CEE0995                     | <input type="checkbox"/> | STANDARD   | ◆ |
| <b>FRONT CASTORS*</b>       |                          |  |   |
| CEE1010                     | <input type="checkbox"/> | 8" X 1.1/4" BLACK SOLID TYRE   | ◆ |
| CEE1030                     | <input type="checkbox"/> | 8" X 2" BLACK SOLID TYRE   | ○ |
| <b>REAR WHEEL SIZE*</b>     |                          |  |   |
| CEE1068                     | <input type="checkbox"/> | 24"  | ◆ |
| <b>REAR WHEELS*</b>         |                          |  |   |
| CEE1071                     | <input type="checkbox"/> | STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES ) | ◆ |
| <b>TYRES*</b>               |                          |  |   |
| CEE1082                     | <input type="checkbox"/> | SOLID BLACK 1.3/8" / 35MM  | ◆ |
| CEE1088                     | <input type="checkbox"/> | PNEUMATIC LOW PROFILE BLACK 1" / 25MM                                | ○ |
| <b>HANDRIM*</b>             |                          |  |   |
| CEE1140                     | <input type="checkbox"/> | ALUMINIUM PREANODIZED HANDRIM  | ◆ |
| <b>REAR WHEEL AXLE*</b>     |                          |  |   |
| CEE1180                     | <input type="checkbox"/> | Q/R AXLE   | ◆ |
| <b>REAR WHEEL POSITION*</b> |                          |  |   |
| CEE1190                     | <input type="checkbox"/> | STANDARD REAR WHEELS POSITION  | ◆ |
| CEE1195                     | <input type="checkbox"/> | MID WHEEL PROPULSION REAR WHEEL POSITION                             | ○ |
| CEE1200                     | <input type="checkbox"/> | PASSIVE REAR WHEELS POSITION   | ○ |
| <b>BRAKES*</b>              |                          |  |   |
| CEE1270                     | <input type="checkbox"/> | PUSH TO LOCK BRAKES  | ◆ |
| <b>TRANSPORTATION KIT*</b>  |                          |  |   |
| CEE1345                     | <input type="checkbox"/> | TRANSPORTATION BRACKETS  | ◆ |
| <b>SAFETY OPTIONS*</b>      |                          |  |   |
| CEE1360                     | <input type="checkbox"/> | POSTURE BELT   | ◆ |

|  |  |  |                    |                      |
|--|--|--|--------------------|----------------------|
|  |  |  | <b>TOTAL PRICE</b> | <input type="text"/> |
|--|--|--|--------------------|----------------------|

**Non Conforming Options:**

|         |   |
|---------|---|
| CEE0740 | REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED   |
| CEE0750 | REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED  |
| CEE0760 | SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT   |
| CEE0839 | TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) |
| CEE0878 | FOR 400 MM SEAT DEPTH   |
| CEE0879 | FOR 450 MM SEAT DEPTH   |
| CEE1195 | MID WHEEL PROPULSION REAR WHEEL POSITION  |
| CEE1200 | PASSIVE REAR WHEELS POSITION  |
| CEE1345 | TRANSPORTATION BRACKETS   |



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):