



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

January 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

E-Motion / DuoDrive with Action 3NG M123 & M124 706257 / 706256

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO- 7176-Part 19
- Self Propulsion with 24" Quick release rear wheels
- Choice of 6 seat widths 380mm to 505mm (15" to 20")
- Choice of 2 seat depths 400mm or 450mm (16" or 18")
- Choice of 6 seat heights 385mm to 510mm (15" to 20")
- Black nylon upholstery
- Fixed 7° 430mm backrest
- Lumbar support padded backrest with shaped cover
- 80° Swing in / Swing out legrest hangers
- Flip up footplates
- Swing away and detachable desk armrests
- Choice of 8 colours
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and buckle belt (adjustable both sides)
- Adjustable heel straps
- User tool kit
-
- Alber E-Motion M25/ DuoDrive
- 24" Stainless steel handrims
- Schwalbe Rightrun Black
- Anti-tippers with jack-up function (pair)
- Fitted at Invacare Included



Please make your selection with a cross in the box	☐ Free of charge	☒ standard feature	☐ POA	☐ special request
--	--------------------------------------	--	------------------------------	--

Customer weight = KG

Code **Description**

BASE PRODUCT*	
CEE0003	☐ ACTION 3NG SELF PROPEL ◆

PROPULSION*			
CEE0100		SELF PROPULSION	◆
CEE9250		E-Motion M25 Weekly Price £55.34 Including Axle adaption kit, Anti tippers with Jack up function M123	£6,478
CEE9251		E-Motion DuoDrive Weekly Price £58.27 Including Axle adaption kit, Anti tippers with Jack up function M124	£6,937
SEAT WIDTH*			
CEE0180		15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0190		16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0200		17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	◆
CEE0210		18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0220		19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0230		20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
SEAT DEPTH*			
CEE0240		16" / 400 MM SEAT DEPTH (NC 0879)	○
CEE0250		18" / 450 MM SEAT DEPTH (NC 0878)	◆
SEAT HEIGHT*			
CEE0260		20" / 510 MM SEAT HEIGHT (NC 0511, 0561)	○
CEE0270		19" / 485 MM SEAT HEIGHT (NC 0511, 0561)	○
CEE0280		18" / 460 MM SEAT HEIGHT	◆
CEE0290		17" / 435 MM SEAT HEIGHT (NC 1010)	○
CEE0300		16" / 410 MM SEAT HEIGHT (NC 1010)	○
CEE0310		15" / 385 MM SEAT HEIGHT , (NC 0511, 0561, 0660, 0670, 0760, 0770, 0839, 1010, 1068)	○
SEAT UPHOLSTERY*			
CEE0330		NYLON W/O HOOK AND LOOP STRIPS	◆
CEE0340		NYLON WITH HOOK AND LOOP STRIPS	○
CEE0345		TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM	○
BACKREST*			
CEE0380		FIXED 0° BACKREST 2 HEIGHTS 16" / 400 MM or 17" / 430 MM	○
CEE0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT	○
CEE0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES	○
CEE0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES	○
CEE0428		FIXED 7° BACKREST ACTIVE 14" / 350 MM HT ADJ TO 16" / 400 MM FIX PUSH HANDLES	○
CEE0455		FOLDING 7° BACKREST 17"	○
CEE0460		FOLDING 7° BACKREST 20"	○
BACKREST UPHOLSTERY*			
CEE0465		WITHOUT BACK UPHOLSTERY	○
CEE0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER	○
LEGREST HANGER RH*			
CEE0510		80° RH SWING IN/SWING OUT (NC 66666)	◆
CEE0511		90° RH SWING IN/SWING OUT (NC 0600, 1010, 66666)	○

LEGREST HANGER LH*				
CEE0560		80° LH SWING IN/SWING OUT		◆
CEE0561		90° LH SWING IN/SWING OUT	(NC 0620, 1010)	○
FOOTREST RH*				
CEE0600		FLIP UP RH FOOTPLATE		◆
CEE66666		NO FOOTREST		○
FOOTREST LH*				
CEE0620		FLIP UP LH FOOTPLATE		◆
CEE55555		NO FOOTREST		○
ARMREST RH*				
CEE0640		REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0839)	○
CEE0650		REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0839)	○
CEE0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT	(NC 0740, 0750, 0839, 0955)	◆
CEE0670		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0740, 0750, 0839, 0955)	○
ARMREST LH*				
CEE0740		REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED		○
CEE0750		REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED		○
CEE0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT	(NC 0955)	◆
CEE0770		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0955)	○
CEE0839		TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM)	(NC 0955, 1200, 1345)	○
SEAT CUSHION DEPTH				
CEE0878		FOR 400 MM SEAT DEPTH		○
CEE0879		FOR 450 MM SEAT DEPTH		○
OTHER SEATING OPTIONS				
CEE0930		CALF STRAP		○
CEE0955		W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY		○
COLOUR*				
CEE0961		SLATE GREY (RAL7024 FINE TEXTURE)		○
CEE0964		PEARL WHITE		○
CEE0965		ELECTRIC GREEN		○
CEE0966		MATT BLACK		○
CEE0974		BLUE		○
CEE0976		DARK PURPLE		○
CEE0969		RED MATT		○
CEE0972		WATER BLUE		○
CEE0980		ANTHRACITE BLACK		○
CASTOR HOUSING*				
CEE0995		STANDARD		◆
FRONT CASTORS*				
CEE1010		8" X 1.1/4" BLACK SOLID TYRE		◆

REAR WHEEL SIZE*		
CEE1068	☐	24" ◆
REAR WHEELS*		
CEE1071	☐	STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES) ◆
TYRES*		
CEE1082	☐	SOLID BLACK 1.3/8" / 35MM ◆
CEE1088	☐	PNEUMATIC LOW PROFILE BLACK 1" / 25MM ○
HANDRIM*		
CEE1140	☐	ALUMINIUM PREANODIZED HANDRIM ◆
REAR WHEEL AXLE*		
CEE1180	☐	Q/R AXLE ◆
REAR WHEEL POSITION*		
CEE1190	☐	STANDARD REAR WHEELS POSITION ◆
CEE1200	☐	PASSIVE REAR WHEELS POSITION ○
BRAKES*		
CEE1270	☐	PUSH TO LOCK BRAKES ◆
CEE1275	☐	PULL TO LOCK BRAKES ○
TRANSPORTATION KIT*		
CEE1345	☐	TRANSPORTATION BRACKETS ◆
SAFETY OPTIONS*		
CEE1360	☐	POSTURE BELT ◆

TOTAL PRICE		
-------------	--	--

Non Conforming Options:	
CEE0511	90° RH SWING IN/SWING OUT
CEE0561	90° LH SWING IN/SWING OUT
CEE0600	FLIP UP RH FOOTPLATE
CEE66666	NO FOOTREST
CEE0620	FLIP UP LH FOOTPLATE
CEE0660	SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT
CEE0670	SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT
CEE0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEE0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CEE0839	TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM)
CEE0878	FOR 400 MM SEAT DEPTH
CEE0879	FOR 450 MM SEAT DEPTH
CEE0955	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY
CEE1010	8" X 1.1/4" BLACK SOLID TYRE
CEE1068	24"
CEE1200	PASSIVE REAR WHEELS POSITION
CEE1345	TRANSPORTATION BRACKETS



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):