



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

January 2026
Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

TDX SP2 Low Rider M14 712116 - Weekly Price £67.41

Standard Features

- Max user weight 160kg (25.2 Stone)
- Crash tested ISO7176-Part 19
- Modulite seating
- Matrix Flo-tech Lite cushion
- Without push handles
- Flex 3 Adjustable backrest plates with laguna backrest cushion
- Elan Headrest
- Standard T-armrests, detachable with Long armpad
- Seat Lifter with powered Tilt (300mm, 30°) CoG Shift
- Seat height either 403mm
- Manual settable Recliner 0° - 30° in 3.75° steps
- Standard 80° Hangers with Angle and Depth Adjustable footplates(adapt)
- LiNX control REM211 with swing away bracket
- Solid Black tyres
- 6km/h (4mph) 4 Pole motors
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £8,216



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code Description

BASE MODEL*

| | | | |
|---------|--------------------------|--|--------|
| ARQ0072 | <input type="checkbox"/> | TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height | £8,216 |
|---------|--------------------------|--|--------|

USER WEIGHT*

| | | | |
|---------|--------------------------|---------------------------|-----------------------|
| ARQ1918 | <input type="checkbox"/> | User weight less 100 kg | <input type="radio"/> |
| ARQ1919 | <input type="checkbox"/> | User weight 101 to 125 kg | <input type="radio"/> |
| ARQ1921 | <input type="checkbox"/> | User weight 126 to 136 kg | <input type="radio"/> |
| ARQ1922 | <input type="checkbox"/> | User weight 137 to 150 kg | <input type="radio"/> |
| ARQ1923 | <input type="checkbox"/> | User weight 151 to 160 kg | <input type="radio"/> |

| | | | |
|--|--------------------------|--|---|
| SEATING TYPE* | | | |
| ARQ2100 | <input type="checkbox"/> | Modulite Seating Unit | ◆ |
| SEAT FRAME* | | | |
| ARQ2111 | <input type="checkbox"/> | Telescopic Seat Frame with metal seat plate | ◆ |
| SEAT WIDTH* | | | |
| ARQ2122 | <input type="checkbox"/> | Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19)) | ○ |
| ARQ2123 | <input type="checkbox"/> | Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21")) | ○ |
| ARQ2124 | <input type="checkbox"/> | Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23")) | ○ |
| SEAT DEPTH RAIL LENGTH* | | | |
| ARQ2129 | <input type="checkbox"/> | Standard Depth Package (16" - 20", 410mm - 510mm) | ○ |
| ARQ2130 | <input type="checkbox"/> | Extended Depth Package (20" - 22", 510mm - 560mm) (NC 2134) | ○ |
| SEAT DEPTH* | | | |
| ARQ2132 | <input type="checkbox"/> | 410mm (16") (410mm - 510mm) (NC 2196) | ○ |
| ARQ2134 | <input type="checkbox"/> | 460mm (18") (410mm - 510mm) | ○ |
| ARQ2136 | <input type="checkbox"/> | 510mm (20") (410mm - 510mm) | ○ |
| ARQ2138 | <input type="checkbox"/> | 560mm (22") (510mm - 560mm) | ○ |
| SEAT CUSHION* | | | |
| ARQ2142 | <input type="checkbox"/> | Matrx Flo-tech Lite | ◆ |
| BACKREST SUPPORT* | | | |
| ARQ2162 | <input type="checkbox"/> | Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates | ◆ |
| BACKREST UPHOLSTERY / CUSHION* | | | |
| ARQ2173 | <input type="checkbox"/> | Flex3 Adjustable Backrest Plates with Laguna Backrest cushion | ◆ |
| UPHOLSTERY TYPE* | | | |
| ARQ2152 | <input type="checkbox"/> | Black multi stretch PU coated fabric | ○ |
| ARQ2155 | <input type="checkbox"/> | Black Spacer fabric | ○ |
| ARMREST* | | | |
| ARQ2191 | <input type="checkbox"/> | T-Armrest Standard, Range of adjustment | ◆ |
| ARMPAD* | | | |
| ARQ2196 | <input type="checkbox"/> | Long Armpad (375mm) | ◆ |
| HEADREST / NECKREST* | | | |
| ARQ2324 | <input type="checkbox"/> | Elan Standard Pad 355mm (10") wide | ◆ |
| MANUAL / POWERED SEATING OPTIONS* | | | |
| ARQ1094 | <input type="checkbox"/> | Seat lifter with Electrical Tilt | ◆ |
| BACKREST RECLINER* | | | |
| ARQ0975 | <input type="checkbox"/> | Manual Settable Recliner | ◆ |
| LEG REST TYPE* | | | |
| ARQ0526 | <input type="checkbox"/> | Side Mounted : Standard 80° hangers | ◆ |
| FOOTPLATE TYPE* | | | |
| ARQ0512 | <input type="checkbox"/> | Fixed Footplate | ◆ |
| LINX CONTROLS* | | | |
| ARQ2203 | <input type="checkbox"/> | REM211 Drive & Actuators | ◆ |
| CONTROL BRACKETRY* | | | |
| ARQ1707 | <input type="checkbox"/> | Swing Away RH Bracket | ○ |
| ARQ1708 | <input type="checkbox"/> | Swing Away LH Bracket | ○ |

| | | |
|--------------------------------------|--------------------------|--|
| SHROUD COLOUR* | | |
| ARQ2640 | <input type="checkbox"/> | Obsidian Black Frame Shrouds ○ |
| ARQ2641 | <input type="checkbox"/> | Ruby Red Frame Shrouds ○ |
| ARQ2642 | <input type="checkbox"/> | Purple Passion Frame Shrouds ○ |
| ARQ2643 | <input type="checkbox"/> | Peacock Blue Frame Shrouds ○ |
| ARQ2644 | <input type="checkbox"/> | Arctic Blue Frame Shrouds ○ |
| RIM COLOUR* | | |
| ARQ2650 | <input type="checkbox"/> | Obsidian Black Rim Inserts ○ |
| ARQ2651 | <input type="checkbox"/> | Ruby Red Rim Inserts ○ |
| ARQ2652 | <input type="checkbox"/> | Purple Passion Rim Inserts ○ |
| ARQ2653 | <input type="checkbox"/> | Peacock Blue Rim Inserts ○ |
| ARQ2654 | <input type="checkbox"/> | Arctic Blue Rim Inserts ○ |
| FORKS* | | |
| ARQ2632 | <input type="checkbox"/> | Double sided fork ◆ |
| CENTRE WHEELS (DRIVE WHEELS)* | | |
| ARQ0705 | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert ◆ |
| MOTOR TYPE* | | |
| ARQ1282 | <input type="checkbox"/> | 6 km/h (4mph) (4 pole) ◆ |
| CHARGER* | | |
| ARQ1343 | <input type="checkbox"/> | OptiCharge Charger 8A (200 - 250V) ◆ |
| LIGHTS* | | |
| ARQ1254 | <input type="checkbox"/> | Without Light system ◆ |
| TRANSPORTATION KIT* | | |
| ARQ1494 | <input type="checkbox"/> | 4-Point tie-down ◆ |
| FIXATION BELT* | | |
| ARQ2113 | <input type="checkbox"/> | Lapbelt (adjustable both sides) ◆ |
| TOTAL PRICE | | |

Non Conforming Options:

| | |
|---------|-----------------------------|
| ARQ2134 | 460mm (18") (410mm - 510mm) |
| ARQ2196 | Long Armpad (375mm) |



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):