



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Motability Retail Prescription Form

January 2026  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**TDX SP2 HD M16 1001587 - Weekly Price  
£56.27**

## Standard Features

- Max user weight 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with long armpad
- Neckrest
- Powereed tilt 20°
- Manual settable backrest Recline 0° - 30° in 3.75° steps
- 80° hangers with fixed footplate
- LiNX control REM211 swing away bracket
- Solid Black tyres
- 6km/h (4mph) 4 pole motors
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £6,410**



Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA    special request

Customer weight =                      KG

Code	Description	
<b>BASE MODEL*</b>		
ARJ0000	<input type="checkbox"/> TDX SP2 Heavy Duty	£6,410
<b>USER WEIGHT*</b>		
ARJ1918	<input type="checkbox"/> User weight less 100 kg	<input type="radio"/>
ARJ1919	<input type="checkbox"/> User weight 101 to 125 kg	<input type="radio"/>
ARJ1921	<input type="checkbox"/> User weight 126 to 136 kg	<input type="radio"/>
ARJ1922	<input type="checkbox"/> User weight 137 to 150 kg (NC 1092)	<input type="radio"/>
ARJ1923	<input type="checkbox"/> User weight 151 to 160 kg (NC 1092)	<input type="radio"/>
ARJ1924	<input type="checkbox"/> User weight 161 to 180 kg (NC 1092)	<input type="radio"/>

SEATING TYPE*			
ARJ2100	<input type="checkbox"/>	Modulite Seating unit	◆
SEAT WIDTH*			
ARJ2121	<input type="checkbox"/>	Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2123	<input type="checkbox"/>	Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	○
SEAT DEPTH RAIL LENGTH*			
ARJ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138, 2196)	○
ARJ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134)	○
SEAT DEPTH SETTING*			
ARJ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm simple stepless adjustment from the side) (NC 2196)	○
ARJ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm simple stepless adjustment from the side)	○
ARJ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm simple stepless adjustment from the side)	○
ARJ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)	○
SEAT CUSHION*			
ARJ2142	<input type="checkbox"/>	Matrx Flo-tech Lite	◆
BACKREST SUPPORT*			
ARJ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	◆
PUSH HANDLES FOR BACKPOSTS*			
ARJ2107	<input type="checkbox"/>	Backposts WITH push handles	◆
BACKREST UPHOLSTERY / CUSHION*			
ARJ2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	◆
UPHOLSTERY COLOUR*			
ARJ2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	○
ARJ2155	<input type="checkbox"/>	Black Spacer fabric	○
ARMREST*			
ARJ2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	◆
ARMPAD*			
ARJ2196	<input type="checkbox"/>	Long Armpad (375mm) (NC 2132)	◆
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ARJ1350	<input type="checkbox"/>	Neckrest height and depth adjustable	◆
MANUAL / POWERED SEATING OPTIONS*			
ARJ1092	<input type="checkbox"/>	Electric Tilt (0° - 20°) Seat Height 450mm	◆
ARJ1095	<input type="checkbox"/>	Electric HD Tilt (0° - 20°) Seat Height 450mm (upto 180kg Max User Weight)	○

SEAT HEIGHT SETTING			
ARJ0354	<input type="checkbox"/>	Seat height set at 450mm (17.75") (from floor to front of seat plate at 0° tilt)	<input type="radio"/>
ARJ0364	<input type="checkbox"/>	Seat height set at 475mm (18.75") (from floor to front of seat plate at 0° tilt)	<input type="radio"/>
BACKREST RECLINER*			
ARJ0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="radio"/>
LEG REST TYPE*			
ARJ0526	<input type="checkbox"/>	Side Mounted : Standard 80° hangers	<input checked="" type="radio"/>
FOOTPLATE TYPE*			
ARJ0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="radio"/>
LINX CONTROLS*			
ARJ2203	<input type="checkbox"/>	REM211 Drive & Actuators	<input checked="" type="radio"/>
CONTROL BRACKETRY*			
ARJ1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ARJ1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>
SHROUD COLOUR*			
ARJ2640	<input type="checkbox"/>	Obsidian Black Frame Shrouds	<input type="radio"/>
ARJ2641	<input type="checkbox"/>	Ruby Red Frame Shrouds	<input type="radio"/>
ARJ2642	<input type="checkbox"/>	Purple Passion Frame Shrouds	<input type="radio"/>
ARJ2643	<input type="checkbox"/>	Peacock Blue Frame Shrouds	<input type="radio"/>
ARJ2644	<input type="checkbox"/>	Arctic Blue Frame Shrouds	<input type="radio"/>
RIM COLOUR*			
ARJ2650	<input type="checkbox"/>	Obsidian Black Frame Shrouds	<input type="radio"/>
ARJ2651	<input type="checkbox"/>	Ruby Red Frame Shrouds	<input type="radio"/>
ARJ2652	<input type="checkbox"/>	Purple Passion Frame Shrouds	<input type="radio"/>
ARJ2653	<input type="checkbox"/>	Peacock Blue Frame Shrouds	<input type="radio"/>
ARJ2654	<input type="checkbox"/>	Arctic Blue Frame Shrouds	<input type="radio"/>
FORKS*			
ARJ2632	<input type="checkbox"/>	Double sided fork	<input checked="" type="radio"/>
CENTRE WHEELS (DRIVE WHEELS)*			
ARJ0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	<input checked="" type="radio"/>
MOTOR TYPE*			
ARJ1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole)	<input checked="" type="radio"/>
BATTERY*			
ARJ1260	<input type="checkbox"/>	73.5 Ah (Gel)	<input checked="" type="radio"/>
CHARGER*			
ARJ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input checked="" type="radio"/>
TRANSPORTATION KIT*			
ARJ1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	<input checked="" type="radio"/>
FIXATION BELT*			
ARJ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	<input checked="" type="radio"/>

TOTAL PRICE			
-------------	--	--	--

**Non Conforming Options:**

ARJ2132	410mm (16") (410mm - 510mm simple stepless adjustment from the side)
ARJ2134	460mm (18") (410mm - 510mm simple stepless adjustment from the side)
ARJ2138	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)
ARJ2196	Long Armpad (375mm)
ARJ1092	Electric Tilt (0° - 20°) Seat Height 450mm



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):