



DESIGNED FOR

Life

Küschall Newsletter Winter Edition

Welcome to the winter issue of our

Küschall newsletter!

Hello and welcome to the winter edition of the Küschall newsletter. I have some exciting updates I'd like to share with you. The K-Series, our most popular rigid wheelchair has

undergone quite a transformation. It's lighter than ever, more configuration possibilities with exciting new options and colours. I would also like to introduce the new Champion SK, the rigid chair that folds even more! Check it out, because with folding front frames it's now even smaller. Finally, our R33 has been re-designed. Look out for the new and improved R33, the

unique shock absorbing wheelchair. More adjustment possibilities and lighter weight are just some of the enhancements made to this unique wheelchair.

From our regular conversations and feedback, I know you will welcome these improvements and new options and enjoy catching up on the latest Küschall has on offer. As always, please

feel free to get in touch with me if you have any questions or queries, I'm always happy to hear from you!

Sharon Jones
PRODUCT SPECIALIST KÜSCHALL

Direct line: 01656 776216,
Mobile 07792 620394
Email: sjones@invacare.com

Küschall Champion SK SMALL BUT SMART

The active rigid wheelchair that folds, is now even smaller.

Evolved from our iconic Küschall Champion, it's designed to be the most compact, practical and lightweight folding wheelchair. The Champion SK with its folding front frames makes it effortless to transport in the smallest of spaces.

Highly configurable with a wide range of adjustments, options and accessories; this smart, durable and practical wheelchair is perfect for everyday life and an essential travel companion for any holiday destination.



NEW

Design, Durability & Driving performance



Kuschall K-Series

ACTIVE INDEPENDENCE

With its minimal and open frame design the Kuschall K-Series makes it the wheelchair of choice for the most active and independent person.

With the new design, the K-Series is approximately 1/2kg lighter. Combined with an increased maximum user weight and configuration possibilities, the K-Series can be built to accommodate an even wider range of clients.

20 years on, the iconic K-Series is still the original rigid and fully adjustable wheelchair making it a versatile and practical choice. Over the years, Kuschall have built and improved on the strong identity of this popular wheelchair and today brings us the latest in the series.

What's new?

Increased maximum user weight to 130kg

NEW seat width range 320mm - 500mm

NEW seat depth range 350mm - 525mm

Exciting new frame colours

In addition, there are new options which are available across the whole of the Kuschall range! See next page for details.

Kuschall R33

DISTINCTIVE STYLE - UNRIVALLED COMFORT

The unique style of the R33 has been improved and enhanced.

Nothing else compares to the R33 with its most unique feature the LIS® Low Impact System. Designed to take the bumps out of everyday life the shock absorbing performance of the R33 makes this the most comfortable rigid active wheelchair around.

The seat module now boasts more adjustment possibilities and is over 1/2kg lighter, making this an ideal choice for the independent, active person looking for something unique.


A fixed seat option means the individual style of the R33 can be enjoyed even if there is no requirement for the shock absorber.

Combining this with a wide range of options, accessories and colours, the R33 can be matched to suit all personalities and styles.



NEW OPTIONS

Lighter weight Comfort & Everyday Castors. Available in a range of sizes.



NEW Aluminium Softroll Castors with black or silver hubs.

Available in 4" and 5".




Curve L Handrim

The unique shape of these handrims is wider to increase the grip of the hand on the handrim.

The benefits include:

- Reduced fatigue
- Cramp reduction
- Relaxed and lasting ride experience



Gekko Handrim

The Gekko Handrim features an integrated rubber strip to provide extra grip and comfort.






NEW Tension Backrest Upholstery

Black, nylon premium tension adjustable upholstery with reflective trim.



NEW colours



For full option details please refer to 2017 Küschall prescription forms.

Get in touch:

Marketing Communications Team - 01656 776200
email: uk@invacare.com

