

Essential Range

Essential Basic



Essential Care





Essential Basic

The **Essential Basic** mattress offers the user a comfortable and supportive surface. This is an ideal entry level mattress for individuals who do not require therapy for pressure ulcer prevention.

- Multi-stretch vapor-permeable, water resistant polyurethane fabric cover
- Zipped around 3 sides for ease of inspection and cleaning
- Fire retardant

Essential Care

The **Essential Care** mattress is a pressure redistributing support surface, developed for home care and residential care environments. A castellated foam surface provides a cost effective solution for comfort and pressure redistribution to meet the needs of individuals 'At Risk/High Risk'^{1/2} of developing pressure ulcers.

- Cut foam castellations
- Multi-stretch vapor-permeable, water resistant polyurethane fabric cover
- Zipped around 3 sides for ease of inspection and cleaning
- Fire retardant

Size	Max User weight	Order Code
200 x 88 x 15cm	140kg (22 stone)	1586660

Size	Max User weight	Order Code
200 x 88 x 15cm	140kg (22 stone)	1586659



1. Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.
2. The risk category in this instance is based on the Waterlow Scale

Mattresses is manufactured to: BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

Tel: +44 1 656 776 222
 Fax: +44 1656 776 220
 E-mail: uk@invacare.com
 Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Limited
 Pencoed Technology Park
 Pencoed - Bridgend CF35
 5AQ - United Kingdom

Tel. +353 1 8107084
 Fax +353 1 8107085
 Email: ireland@invacare.com
www.invacare.ie

Invacare Ireland Ltd
 Unit 5 - Seatown
 Business Campus
 Seatown Road
 SWORDS - County
 Dublin - Ireland

