

matrix[®]
SEATING SERIES


Flo-tech[™] Lite Visco



Low profile visco elastic low memory
foam cushion offering exceptional comfort



The *Invacare* **Matrx Flo-tech Lite Visco** cushion is made of high-density visco elastic foam, allowing it to conform to the user's shape. This affords exceptional levels of comfort and increases pressure redistribution. This low profile cushion format is ideal for both chair and wheelchair use, facilitating side transfers.



Features and Options



Ingress Resistant Cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle.

Available in 61 sizes

The *Invacare* **Matrx Flo-tech Lite Visco** is available in a range of 61 sizes, from 250 x 250 to 480 x 530 mm. (See size chart for details).



Low memory foam designed to maximise body contact area

The visco-elastic foam conforms to the user's body, maximising the contact area, delivering exceptional levels of comfort and optimising pressure redistribution.

Three foam densities available

Available in three densities to accommodate a range of user weights – Soft (up to 32 kg), Medium (up to 51 kg) and Firm (up to 100 kg).

Contoured seating surface with pommel and side support

Offers mild lateral support and pelvic stability to enhance the user's posture.

Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local *Invacare* website.

Various sizes	Height	Total product weight	Max. user weight ¹	Max. cleaning temperature
A range of 61 sizes	78 mm	1.26 kg	100 kg	80°

Manufactured to European Standards:
Interior Foam: BS EN 1021-1, BS EN 1021-2,
BS 5852: 1982 Crib 5

1. Guide maximum user weight limits are set assuming that the cushion selected is of appropriate width and depth for the intended user/s.

Colours



QR
Code

Invacare Ltd
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel +44 1 656 776 222
E-mail: uk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS, County Dublin - Ireland
Tel. +353 1 8107084
Email: Ireland@invacare.com
www.invacare.ie



© 2018 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Lite Visco - 1628867 UK - 10/2018